



# Amines, C<sub>1</sub> – C<sub>2</sub>

## Analysis of impurities in dimethylamine

### Application Note

Materials Testing & Research

#### Authors

Agilent Technologies, Inc.

#### Introduction

Amines are difficult to analyze due to their strong basic nature. Capillary columns must be base-modified to elute amines with acceptable recovery. For highly volatile amines, including ammonia, the siloxane-based phases do not provide enough retention. The Agilent PoraPLOT for Amines porous polymer provides a high retention combined with a high inertness for amines.

Volatile amines elute at low levels which makes impurity analysis possible as shown in this application. Volatile amines can be analyzed at relatively high temperatures due to the high retention of the PoraPLOT Amines. If, besides these amines, alcohols and/or water must be measured, a 5 µm film Agilent CP-Sil 5 CB is recommended, operated at temperatures around 30 °C.




**Agilent Technologies**

## Conditions

Technique : GC-capillary  
Column : Agilent PoraPLOT for Amines, 0.32 mm x 25 m, fused silica PLOT (df = 10  $\mu$ m) (Part no. CP7591)  
Temperature : 110 °C (8 min)  $\rightarrow$  250 °C, 15 °C/min  
Carrier Gas : H<sub>2</sub>, 95 kPa (0.95 bar, 13 psi)  
Injector : Split,  
T = 250 °C  
Detector : FID  
T = 250 °C  
Sample Size : 0.2  $\mu$ L  
Concentration Range : 10 - 100 ppm impurities  
Solvent Sample : dimethylamine balance

## Peak identification

1. methylamine
2. unknown
3. dimethylamine
4. trimethylamine
5. ethylamine



[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01398



**Agilent Technologies**