



## PAHs

# Analysis of EPA 610 polycyclic aromatic hydrocarbons

## Application Note

Environmental

### Authors

Agilent Technologies, Inc.

### Introduction


For the analysis of trace PAH often the on-column injection technique is used because in hot injection ports discrimination, adsorption or decomposition can occur. Best columns for on-column injection are 0.53 mm and 0.32 mm id capillary columns because the inside diameter of the capillary is easily accessed by the syringe needles used for on-column injection. Also, high flow rates are easily obtained resulting in low elution temperatures.



**Agilent Technologies**

## Conditions

Technique : GC  
Column : Agilent VF-5ms, 0.53 mm x 30 m fused silica  
(df = 0.5  $\mu$ m) (Part no. CP8974)  
Temperature : 50 °C, 6 min  $\rightarrow$  300 °C, 10 °C/min  
Carrier Gas : Helium, 50 cm/s  
Detector : FID  
T = 300 °C  
Injector : Direct, 2  $\mu$ L  
Concentration Range : ca. 10 pg/ $\mu$ L



45 min

## Analysis of standard mixture

Column : Agilent VF-5ms, 0.53 mm x 30 m fused silica  
(df = 0.5  $\mu$ m) (Part no. CP8974)

Temperature : 40 °C, 4 min  $\rightarrow$  270 °C, 10 °C/min

Carrier Gas : H<sub>2</sub>, 50 kPa


Detector : FID  
T = 300 °C

Injector : Split

Concentration range : ca. 10 ng of each component on the column

## Peak identification

1. naphthalene
2. acenaphthylene
3. acenaphthene
4. fluorene
5. phenanthrene
6. anthracene
7. fluoranthene
8. pyrene
9. chrysene
10. benzo(a)anthracene
11. benzo(k)fluoranthene
12. benzo(b)fluoranthene
13. benzo(a)pyrene
14. indeno(1,2,3-c,d)pyrene
15. dibenzo(a,h)anthracene
16. benzo(g,h,i)perylene



**[www.agilent.com/chem](http://www.agilent.com/chem)**

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A02253



**Agilent Technologies**