

# Basic aromatic compounds

## Application Note

Environmental

### Authors

Agilent Technologies, Inc.

### Introduction

Eight basic aromatics are separated by GC using the stabilized 50% phenyl PDMS phase of Agilent VF-17ms in less than 20 minutes.



**Agilent Technologies**


## Conditions

Technique : GC  
Column : Agilent VF-17ms, 0.25 mm x 30 m fused silica  
(df = 0.25 µm) (Part No. CP8982)  
Temperature : 50 °C + 10 °C/min → 300 °C  
Carrier Gas : Helium, 70 kPa  
Injector : Splitter, 1:100  
Detector : FID  
Sample Size : 1 µL  
Concentration Range : 200 µg/mL

Courtesy : J. Peene, Agilent application laboratory,  
Middelburg, The Netherlands

## Peak identification

1. aniline
2. benzyl alcohol
3. 4-chloroaniline
4. 2-methylnaphthalene
5. 2-nitroaniline
6. dibenzofuran
7. 3-nitroaniline
8. 4-nitroaniline



[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A02281



**Agilent Technologies**