

# Analysis of Petroleum Fractions by ASTM D2887

Peter Morgan, Thermo Fisher Scientific, Runcorn, Cheshire, UK

## Key Words

Simulated distillation, D2887, TRACE TR-SimDist

## Abstract

ASTM method D2887 describes a simulated distillation protocol for the separation of petroleum fractions over the boiling points range of 55.5 to 538 °C. The Thermo Scientific™ TRACE™ TR-SimDist GC column has been shown to achieve this separation in less than 20 minutes.

## Introduction

ASTM method D2887 describes the gas chromatography (GC) conditions for simulated distillation of petroleum fractions in the boiling point range of 55.5 to 538 °C.[1]

Simulated distillation by GC is used to determine the boiling point range of a petroleum sample. Column heating causes the sequential elution of sample components based on their boiling point. A calibration curve is constructed that relates carbon number and boiling point to retention time and allows for the determination of boiling points for the components in an unknown test sample.

ASTM method D2887 presents a chromatographic challenge as the boiling point range is wide. To resolve hydrocarbons in the C5–C12 range, a particularly thick film is used. To elute heavier hydrocarbons, a high temperature is required.



## Experimental Details

### Chemicals and Reagents

1. A commercially available standard containing all components in Table 1 was prepared in carbon disulfide at a concentration of 800 µg/mL. This concentration was found to give the optimum sensitivity and peak shape.

2. A commercially available gas oil reference mix, diluted twenty fold in carbon disulfide, was used to demonstrate the performance of the TRACE TR-SimDist GC column in a real sample.

Peak Number	Compound Name	Boiling Point (°C)	Peak Number	Compound Name	Boiling Point (°C)
1	n-pentane (C5)	36	10	n-hexadecane (C16)	287
2	n-hexane (C6)	69	11	n-octadecane (C18)	316
3	n-heptane (C7)	98	12	n-eicosane (C20)	344
4	n-octane (C8)	126	13	n-tetracosane (C24)	391
5	n-nonane (C9)	151	14	n-octacosane (C28)	431
6	n-decane (C10)	174	15	n-dotriacontane (C32)	466
7	n-undecane (C11)	196	16	n-hexatriacontane (C36)	496
8	n-dodecane (C12)	216	17	n-tetracontane (C40)	522
9	n-tetradecane (C14)	254	18	n-tetratetracontane (C44)	545

Table 1: ASTM D2887 Calibration standard containing C5–C44

### Sample Handling Equipment

Vials and closures: 9 mm screw cap vial with Blue Silicone/PTFE Closure 60180-599

### Separation Conditions

		Part Number
Instrumentation:	Thermo Scientific TRACE GC	
Column:	TRACE TR-SimDist 10 m x 0.53 mm x 2.65 µm	260S348P
Septum:	Thermo Scientific BTO Septa, 12.7 mm	31303228
Liner:	Thermo Scientific PTV Silcosteel Liner for PTV simulated on-column injection	45322052
Column ferrules:	Graphite ferrules to fit 0.53 mm ID columns	29053488
Injection syringe:	10 µL fixed needle 80 mm cone tip	36502019
Carrier gas:	Helium	
Column flow:	12 mL/min (constant flow)	
Oven temperature:	40 °C hold 1 minute, 20 °C/min to 370 °C, hold 5 minutes	
Injector type:	Programmable temperature vaporizer (PTV)	
Injector mode:	Simulated on-column	
Injection volume:	1 µL	
Injector temperature:	40 °C, 10 °C/s to 370 °C, hold 20 minutes	
Detector details:	Flame ionization detector (FID), 370 °C	

### Data Processing

Software: Thermo Scientific Xcalibur™

## Results

The lower boiling point components in the C5–C10 region exhibited good retention and resolution (Figure 1) because the 2.65  $\mu\text{m}$  film gave greater retention than a thinner film column would.




Figure 1: Compounds C5–C10

The resolution between the C16 and C18 peaks (10 and 11) was calculated to be 6.12, exceeding the minimum specification of 3 as set out in the ASTM method (Figure 2).




Figure 2: Resolution between the C16 and C18 (peaks 10 and 11)

Peak shapes across the boiling point range were symmetrical and showed no evidence of overload at a concentration of 800  $\mu\text{g}/\text{mL}$ . This allowed a calibration curve to be constructed of retention time against boiling point (boiling point data taken from ASTM method) showing a linearity of greater than 0.99 (Figure 3).




Figure 3: Boiling point calibration curve

The calibration curve was used to determine the boiling point range of a reference gas oil sample. This sample was found to contain compounds with boiling points between 36 and 475  $^{\circ}\text{C}$  (Figure 4).




Figure 4: Calibration standard at 800  $\mu\text{g}/\text{mL}$  (top chromatogram) overlaid with reference gas oil sample diluted 20 times in  $\text{CS}_2$  (bottom chromatogram)

## Conclusion

- The TRACE TR-SimDist GC column meets the criteria for the analysis of hydrocarbons with boiling points between 55.5 and 538 °C as required by ASTM D2887.
- A combination of thick film and high temperature allows the separation of these n-alkanes with a broad boiling point range with good peak shape.

## References

1. ASTM Standard D2887-08, 2008, “Standard Test Method for Boiling Range Distribution of Petroleum Fractions by Gas Chromatography,” ASTM International, West Conshohocken, PA, 2008, DOI: 10.1520/D2887-08, [www.astm.org](http://www.astm.org).

[thermoscientific.com/gccolumns](http://thermoscientific.com/gccolumns)

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**USA and Canada** +1 800 332 3331  
**France** +33 (0)1 60 92 48 34  
**Germany** +49 (0) 2423 9431 20 or 21  
**United Kingdom** +44 (0)1928 534110  
**Japan** +81 3 5826 1615

**China** +86 21 68654588 +86 10 84193588  
+86 20 83145199 800 810 5118  
**India** +91 22 6742 9494 +91 27 1766 2352  
**Australia** 1 300 735 292 (free call domestic)  
**New Zealand** 0800 933 966 (free call domestic)  
**All Other Enquiries** +44 (0) 1928 534 050

**Technical Support**  
**North America** +1 800 332 3331  
**Outside North America** +44 (0) 1928 534 440

**Thermo**  
SCIENTIFIC

Part of Thermo Fisher Scientific