


# Reproducibility of EGA Curves In Evolved Gas Analysis (EGA)

EGA is a thermal analysis technique in which a sample is continuously heated and gases evolved are not separated but directly detected. Data obtained by this technique provides useful information such as thermal desorption conditions of volatile components released from sample and start/end temperatures of pyrolysis<sup>1)</sup>. The reproducibility of EGA curves utilizing Double-Shot Pyrolyzer® is described here. Fig. 1. shows EGA curves of polystyrene (PS) obtained by heating a sample from 100°C to 600°C at 20°C/min. The process was repeated 5 times. A peak arising from thermal decomposition of PS is observed at 340°C~470°C in each run. The average temperature of 5 peaks was 441.6°C+/-1°C (top of peak). This clearly shows that excellent reproducibility has been obtained.

1) *Double-Shot Pyrolyzer® Technical Note, PYT-004*



**Fig. 1 EGA Curves of PS**

Pyrolysis temp.: 100–600°C (20°C/min), Carrier gas : He 50kPa, Split ratio: ca. 1/50  
 EGA capillary tube : Id 0.15mm, Length 2.5m (UADTM-2.5N)  
 GC oven temp.: 300°C, Injection temp.:320°C, Sample : 30µg, Detector : FID

**Keywords :** Evolved Gas Analysis, Reproducibility

**Products used :** Multi-functional pyrolyzer, UADTM-2.0N

**Applications :** General Polymer Analysis

**Related technical notes :**

Please forward your inquiries via our web page or send us a fax message.

**R&D and manufactured by :**  
**Frontier Laboratories Ltd.**

Phone: (81)24-935-5100 Fax: (81)24-935-5102  
<http://www.frontier-lab.com/>

