


APPLICATION OF LARGE VOLUME SPME FIBERS

Jason S. Herrington, Colton Myers, Chris Rattray, Linx Wacłaski, Jaap de Zeeuw, Gary Stidsen, Erin Herrold


Restek Corporation, Bellefonte, PA, USA

The SPME Arrow was designed to overcome 2 major short-comings associated with traditional SPME fibers: small phase volumes and limited mechanical robustness.


The following figure provides a side-by-side comparison:




With over 4x the phase surface area and 6x the volume, SPME Arrows were compared head-to-head with traditional SPME fibers for the following applications:



CTC RTC Parameters	Agilent 7890B/5977B GC-MS Parameters
HS-SPME	Inlet
SPME: DVB/Carbon/PDMS	250 °C
Agitator Speed: 250 rpm	Split: 50:1
Agitator Temperature: 30 °C	Topaz 1.8 mm ID Straight/SPME Liner (cat.# 23280)
Incubation Time: 120 s	Column
Heater Stirrer Speed: 1000 rpm	Rtx-VMS, 40 m, 0.18 mm ID, 1.00 µm (cat.# 49915)
Heater Stirrer Temperature: 30 °C	Oven
Extraction Time: 450 s	Insert: GC Accelerator (cat.# 23849)
Vial Penetration Depth: 35 mm	Ramp: 30 °C/hold 2 min to 180 °C at 30 °C/min
Injector Penetration Depth: 50 mm	Carrier Gas
Desorption Time: 5 s	Type: Helium
Pre Conditioning: True	Mode: Constant Flow
Post Conditioning: False	Flow Rate: 1.44 mL/min
Conditioning Time: 60 s	Detector
Conditioning Temperature: 280 °C	Type: Single Quadrupole MS
Pre Conditioning: True	Mode: Scan 12.1 scans/sec (35 to 226 amu)
Post Conditioning: False	Transfer Line Temp: 250 °C
Conditioning Time: 60 s	Source Temp: 220 °C
Conditioning Temperature: 280 °C	Quad Temp: 150 °C
	Electron Energy: 70 eV
	Tune Type: BFB
	Ionization Code: EI



CTC RTC Parameters	Agilent 7890B/5977B GC-MS Parameters
HS-SPME	Inlet
SPME: Carbon	250 °C
Agitator Speed: 250 rpm	Split: 100:1
Agitator Temperature: 30 °C	Topaz 1.8 mm ID Straight/SPME Liner (cat.# 23280)
Incubation Time: 120 s	Column
Heater Stirrer Speed: 1000 rpm	Rtx-VMS, 40 m, 0.18 mm ID, 1.00 µm (cat.# 49915)
Heater Stirrer Temperature: 30 °C	Oven
Extraction Time: 240 s	Insert: GC Accelerator (cat.# 23849)
Vial Penetration Depth: 35 mm	Ramp: 30 °C/hold 2 min to 180 °C at 30 °C/min
Injector Penetration Depth: 50 mm	Carrier Gas
Desorption Time: 5 s	Type: Helium
Pre Conditioning: True	Mode: Constant Flow
Post Conditioning: True	Flow: 1.00 mL/min
Conditioning Time: 60 s	Detector
Conditioning Temperature: 280 °C	Type: Single Quadrupole MS-MS
	Mode: Scan 6.1 scans/sec (29 to 250 amu)
	Transfer: 250 °C
	Source: 220 °C
	Quad: 150 °C
	Electron: 70 eV
	Tune: BFB
	Ionization: EI



CTC RTC Parameters	Agilent 7890B/5977B GC-MS Parameters
HS-SPME	Inlet
SPME: 100 µm PDMS	280 °C
Agitator Speed: 250 rpm	Split: 10:1
Agitator Temperature: 30 °C	Topaz 1.8 mm ID Straight/SPME Liner (cat.# 23280)
Incubation Time: 120 s	Column
Heater Stirrer Speed: 1000 rpm	Sola
Heater Stirrer Temperature: 30 °C	Poli-Sil MS, 20 m, 0.15 mm ID, 0.15 µm (Restek cat.# 43816)
Extraction Time: 450 s	Oven
Vial Penetration Depth: 35 mm	Insert: GC Accelerator (cat.# 23849)
Injector Penetration Depth: 50 mm	Ramp: 70 °C to 115 °C at 95 °C/min, to 175 °C at 65 °C/min
Desorption Time: 60 s	to 225 °C at 45 °C/min to 300 °C at 10 °C/min (hold 1.0 min)
Pre Conditioning: True	Carrier Gas
Post Conditioning: False	Type: Helium
Conditioning Time: 60 s	Mode: Constant Flow
Conditioning Temperature: 280 °C	Flow: 1.20 mL/min
	Detector
	Type: Single Quadrupole MS-MS
	Mode: Scan 12.7 scans/sec (100 to 280 amu)
	Transfer: 280 °C
	Source: 280 °C
	Quad: 200 °C
	Electron: 70 eV
	Tune: D51P0
	Ionization: EI

The SPME Arrow generated 4 – 6x the response of the traditional SPME fiber, which is consistent with the increase in phase surface area and volume. In addition, after hundreds of samples the SPME Arrow does not suffer the same fate as traditional SPME fibers:



For all 3 applications, the SPME Arrow generally demonstrated an increase in response over the traditional SPME fiber. There was a consistent trend that as molecular weight increased, the SPME Arrow's advantage diminished. Regardless, the Arrow's mechanical robustness never lessened.