


# Basics of SPE Technology & Mechanisms

Pieter Grobler,  
Sigma-Aldrich RSA




# Agenda

- The Importance of Sample Prep
- Overview of SPE Technology
- SPE Strategies
- Understanding Retention Mechanisms



# Analytical Chromatography Heaven

- Short run times
- Baseline resolution
- Symmetric peak shape
- Good S/N ratio
- No misleading peaks
- High precision/accuracy



# The Importance of Sample Preparation




 **SUPELCO**

  
**SIGMA-ALDRICH**


# Real World & Real Samples

## The Importance of Sample Preparation

Urine Sample without SPE



Urine Sample with SPE




# Why is sample preparation required?

Collected Sample




GC, HPLC, or LC-MS/MS Analysis



Current Sample = Unsuitable for further analysis!!!... Why?


- **Too dirty**- contains other sample matrix components that interfere with the analysis
- **Too dilute**- analyte(s) not concentrated enough for quantitative detection
- Present **sample matrix not compatible** with or harmful to the chromatographic column/system

# Sources of Chromatographic Errors



(R.E. Majors, LCGC Magazine, 1991, 1997, 2002)

# Time Spend on Analytical Process



(R.E. Majors, LC/GC Magazine, 1991, 1997, 2002)



# Many Tools/Technology for Sample Prep

- Dilute and Shoot
- Filtration
- Protein Precipitation
- Equilibrium dialysis/  
ultrafiltration
- Liquid Liquid Extraction
- Solid Phase Microextraction  
(SPME)
- Solid Phase Extraction  
(SPE)
- Turbulent Flow  
Chromatography
- Monolithic Chromatography
- Immunoaffinity

**Simpler**  
**Generic**  
**Methodology**



**More**  
**Complicated**

**Requires Method Dev**

**Less Selective**

**Minimal Sample Cleanup &  
Concentration**



**Greater Selectivity**

**Optimal Sample Cleanup &  
Concentration**



**SIGMA-ALDRICH**

# Separatory Funnel/LLE = Old Technology




- Works for many samples
- Large solvent consumption
  - Disposal of solvent
- Vigorous shaking/mixing
- Waiting for layers to separate
- Phase emulsions
- Longer Rotovap Times
- Separatory funnel is spacey equipment
  - (sample throughput, Automatisation?)

# Purpose of Solid Phase Extraction (SPE)

Prior to the actual analysis, SPE is most commonly used to...


1. **Clean Up** - Strip the analyte(s) away from endogenous interferences.
2. **Concentrate** analytes(s) for better sensitivity.
3. **Exchange** sample environments for better chromatography  
-e.g., analytes in serum => analytes in mobile phase.



# Overview of Solid Phase Extraction (SPE)




# Basic SPE Concept



- Another **form of chromatography**
- Hardware = plastic (polypropylene) or glass
- Sorbent held in place by two PE frits
- Packing material is very similar to HPLC
  - Often irregular shape vs. spherical (HPLC)
  - Much larger particle size ( $>50\mu\text{m}$ ) vs. HPLC ( $\leq 5\mu\text{m}$ )
  - SPE particle size distribution much broader than HPLC
- Use it only once

# SPE Vacuum Manifold



**Sample introduction**

**Indiv. Port Valves**

**Sample collection tubes (volumetric flasks)**


**Waste reservoir**

**SPE tubes**

**Vacuum manifold**

**Vacuum line and gauge**

# SPE Tube Device Processing Equipment



995-0134

**SUPELCO**

sigma-aldrich.com



**SIGMA-ALDRICH**

# Most Common SPE Robots for Automated SPE



Zymark RapidTrace System



TomTec Quadra System



Gilson SPE 215 System



Code 802 & 803 "Tab-less" 1 & 3mL racks



SIGMA-ALDRICH

 SUPELCO



# Types of SPE Tubes/Cartridges

SPE tubes are available in two materials:

- **Polypropylene (serological grade)**

- Most common
- Suitable for most SPE applications
- Inexpensive



An assortment of Supelco SPE tubes. Second tubes in from either side are glass.

- **Glass (serological grade)**

- Greater solvent resistance than plastic
- No phthalates or plasticizers to leach into sample
- Can be heated
- More expensive than plastic
- Common in environmental analysis




# SPE Bed Weight/Tube Size Selection

- **Smaller tube dimensions (1 mL) contain smaller bed weights.**
  - **reduced elution volumes** which can be beneficial
- **3 mL SPE tubes are most common size**
- **6 mL SPE tubes when one or more steps require volumes greater than 3 mL.**
- **12, 20, and 60 mL tubes contain larger bed weights allow to use SPE as a prep purification or modified LPLC/Flash technique.**

Bed Weight	Tube Volume	Minimum Elution Vol.	Bed Capacity*
50-100 mg	1 mL	100-200 $\mu$ L	2.5-10 mg
500 mg	3 mL	1-3 mL	25-100 mg
0.5-1 g	6 mL	2-6 mL	25-100 mg
2 g	12 mL	10-20 mL	0.1-0.2 g
5 g	20 mL	20-40 mL	1.25-2.5 g
10 g	60 mL	40-100 mL	0.5-1 g

\* This value depends on the analyte and sample matrix. As a rule of thumb, the bed capacity can be estimated with ~5% of the bed weight.



# Common SPE Hardware


## Funnels

Büchner format  
ideal for large  
sample volumes



## Tubes

Glass or plastic, tubes  
are the most common  
SPE format




## 96-well plates



## SPE Disks

# Disk & 96-Well Plates Manifold

## ENVI-Disk™



SPE packing embedded  
in glass fiber matrix



## 96-well plates

**SUPELCO**

  
**SIGMA-ALDRICH**

# SPE Advantages & Disadvantages

## Disadvantages

- **Perceived difficulty** to master its usage (method development)
  - Wide range of chemistries, many choices for manipulating solvent and pH conditions make it **difficult to grasp**
- More steps and MD time required
- Greater cost per sample (really?)

## Advantages

- **Greater selectivity**- paramount importance (e.g. bioanalysis (pg/mL))
- Wide variety of sample matrices
- High recoveries & good reproducibility
- Amenable to **automation**
- Low solvent volumes

# Three different SPE Strategies

Which one to choose depends on the goal of the extraction.

## 1. Bind & Elute Strategy

- Most common
- Bind: **Analytes bind** to tube, unwanted **matrix comp.** are **washed off**
- Elute: **Eluant** changed to **remove analytes** from tube
- Analytes are concentrated via evaporation prior to HPLC or GC analysis

## 2. Interference Removal Strategy

- **Bind** all unwanted **matrix components** and allow **analytes** to **pass through** during the sample loading stage
- Like chemical filtration

## 3. Fractionation Strategy (Form of Bind Elute)

- Retain and **sequentially elute** different **classes of compounds** by modifying eluant pH or % organic

# General Steps of an SPE Procedure (Bind & Elute)

1. Sample Pre-treatment
2. Conditioning & Equilibration
3. Sample Load
4. Washing
5. Elution
6. Evaporation

## 1) Sample Pre-treatment:

Dependent on analyte, sample matrix, and nature of retention chemistry; involves **pH adjustment**, **centrifugation**, **filtration**, **dilution**, **buffer addition**, etc..

## 2a) Conditioning:

Solvent is passed through the SPE material to **wet** the bonded **functional groups** => ensures consistent interaction.

## 2b) Equilibration:

Sorbent/ phase is treated with a solution that is similar (in polarity, pH, etc.) to the sample matrix => **maximizes retention**.

# General Steps of an SPE Procedure (Bind & Elute)

## 3) Sample Load:

Introduction of the sample = **analytes** of interest are **bound/** extracted **onto the phase/** sorbent.

## 4) Washing:

Selectively **remove** unwanted **interferences** co-extracted with the analyte **without** prematurely **eluting analytes** of interest.

## 5) Elution:

**Removing analytes** of interest with a solvent that overcomes the primary and secondary retention interactions b/w sorbent and analytes of interest.


## 6) Evaporation

of eluent/ **reconstitution** with **mobile phase** (optional).




# Bind-elute strategy diagram

(Filtered) sample with internal standard (IS) → Analytes of interest in suitable matrix




# Interference removal strategy diagram

Sample with Internal Standard in Matrix → Matrix adsorbed → Analytes & IS pass



# Fraction strategy diagram

## Form of Bind and Elute Strategy with multiple elution steps



# Understanding Retention Mechanisms



# Reversed-Phase SPE

**Retention Mechanism:** Non-polar or hydrophobic interactions

- Van der Waals or dispersion forces

**Sample Matrix:** Aqueous samples

- Biological fluids (serum, plasma, urine)
- Aqueous extracts of tissues
- Environmental water samples
- Wine, beer and other aqueous samples

**Analyte Characteristics:** Analytes exhibiting non-polar functionalities

- Most organic analytes
- Alkyl, aromatic, alicyclic functional groups


**Elution Scheme:** Disrupt reversed-phase interaction with solvent or solvent mixtures of adequate non-polar character

- Methanol, acetonitrile, dichloromethane
- Buffer/solvent mixtures

**Common Applications:** Drugs and metabolites in biological fluids

- Environmental pollutants in water
- Aqueous extracts of tissues and solids

Aqueous Sample Matrix/Mobile Phase Environment



# Example RP SPE Protocol

## 1. Sample Pre-Treatment

- Dilute samples 1:1 with buffer (10mM ammonium acetate)
- pH manipulation important for ionizable analytes
- Filter or centrifuge out particulates

## 2. Condition & Equilibrate

- Condition with 1-2 tube volumes MeOH or MeCN
- Equilibrate with 1-2 tube volumes buffer

## 3. Load sample (consistent rate; 1-2 drops per second)

# Example RP SPE Protocol

## 4. Wash sorbent (elutes co-retained interferences)

- Critical for improving selectivity
- 5-20% MeOH common
- Dilute MeOH in buffer used during sample load


## 5. Elute analytes of interest

- MeOH or MeCN most common
- pH manipulation can improve recovery  
(adjust pH opposite to load conditions)


## 6. Evaporate/reconstitute as necessary

# C18 vs. C8 vs. Ph vs. CN

More polar RP sorbents (e.g. CN, Ph) can offer better selectivity



Blank urine extract on conventional **C18**



Blank urine extract on **DSC-CN**



# C18 vs. C8 vs. Ph vs. CN

- More **polar RP sorbents**
  - can offer **better selectivity**
  - Often allow for weaker & smaller elution volumes
  - Greater **risk of premature** analyte **elution** during wash step
    - Requires weaker wash solvents
  - Less risk of sorbent over drying
  
- More **non-polar RP sorbents**
  - Have **broader** analyte **retention range**
  - Greater risk insufficient clean-up
  - Allows for **stronger wash** solvents
  - May require increased elution volume

# Useful RP SPE Tips

- Drugs in biological fluids risk **drug-protein binding** effect
  - Disrupt during sample pre-treatment using 40uL 2% disodium EDTA or 2% formic acid per 100uL plasma
- **Sorbent over drying** only a concern during first conditioning step
  - Only critical with C18 & only critical in first conditioning step
  - Phase just needs to be moist during sample addition
  - All other steps non-critical
- If eluate evaporation necessary, **dry SPE tube with vacuum** for 10-15 min. prior to elution to **remove residual moisture**
- Pass **DCM** through SPE before conditioning to remove **SPE tube impurities** for highly sensitive analyses
- **Reduce bed weight** to minimize elution volume
- **Increase bed weight** to retain more polar compounds

# Normal-Phase SPE

## Retention Mechanism: Polar Interactions

- Hydrogen bonding, pi-pi, dipole-dipole, and induced dipole-dipole

## Sample Matrix: Non-polar samples

- Organic extracts of solids
- Very non-polar solvents
- Fatty oils, hydrocarbons

## Analyte Characteristics: Analytes exhibiting polar functionalities


- Hydroxyl groups, carbonyls, amines, double bonds
- Hetero atoms (O, N, S, P)
- Functional groups with resonance properties

## Elution Scheme: Polar interactions disrupted with a more polar solvent or solution

- Acetonitrile, methanol, isopropanol
- Combinations of buffer/solvent or solvent/solvent mixtures

## Common Applications:

- Cleanup of organic extracts of soils and sludge
- Fractionation of petroleum hydrocarbons
- PCBs in transformer oil
- Isolation of compounds in cosmetics



# Example NP SPE Protocol

## 1. Sample Pre-Treatment

- Liq samples **extracted/diluted with non-polar solvent** (e.g. hexane, DCM)
- Solid samples (soil, sediment, etc.) **extracted** (soxhlet, sonication, etc.) **with non-polar solvent, and concentrated**
- **Dry solvent** extract with Na-sulfate or Mg-sulfate
  - Residual moisture can greatly affect analyte retention

## 2. Condition & Equilibrate

- w/ 1-2 tube volumes non-polar solvent

## 3. Load sample (consistent rate; 1-2 drops per second)

- Sample should **not** be in MeCN or MeOH

# Example NP SPE Protocol

## 4. Wash sorbent (elutes co-retained interferences)

- Use a **more polar solvent**, but not so polar as to elute analytes of interest
- Fractionation common in NP SPE

## 5. Elute analytes of interest with polar solvent

- MeOH, MeCN, Acetone, IPA are common

## 6. Evaporate/reconstitute as necessary

# Common Normal Phase Solvents

Solvent	Elutropic (e°) or elution strength on silica		
Hexane	0.00	Promotes Normal-Phase Retention	
Isooctane	0.00		
Carbon tetrachloride	0.14	↓	
Toluene	0.22		
Benzene	0.27		
Tert-butyl methyl ether	0.29		
Chloroform	0.31		
Methylene chloride (dichloromethane)	0.32		
Diethyl ether	0.29		
Ethyl acetate	0.43		
Tetrahydrofuran	0.35		
Acetone	0.45		
Acetonitrile	0.50		
40% methanol in acetonitrile	0.67		
20% methanol in diethyl ether	0.65		
20% methanol in methylene chloride	0.63		
Isopropanol	0.63		
Methanol	0.73		Promotes Normal-Phase Elution
Water	>0.73		
Acetic acid	>0.73		



# Ion-Exchange SPE

**Retention Mechanism:** Electrostatic attraction of charged functional groups of the analyte(s) to oppositely charged functional groups on the sorbent. Combination of reversed-phase and ion-exchange for mixed-mode

**Sample Matrix:** Aqueous or organic samples of low salt concentration (< 0.1M)

- Biological fluids
- Solution phase synthesis reactions

**Analyte Characteristics:**


- Use cation-exchange for isolating basic compounds: primary, secondary, tertiary, and quarternary amines
- Use anion-exchange for isolating acidic compounds: carboxylic acids, sulphonic acids, and phosphates

**Elution Scheme:** Electrostatic interactions disrupted via:

- pH modification to neutralize compound and/or sorbent functional groups
- Increase salt concentration (> 1M); or use a more selective counter-ion to compete for ion-exchange binding sites

**Common Applications:**

- Drugs of abuse and pharmaceutical compounds in biological fluids
- Fatty acids removal in food/agricultural samples
- Cleanup of synthetic reactions
- Organic acids from urine
- Herbicides in soil



# Example IOX SPE Protocol

## 1. Sample Pre-Treatment:

- Basic compounds: dilute w/ 10-25mM buffer (e.g., potassium phosphate, ammonium acetate), pH 3-6
- Acidic compounds: dilute with 10-25mM buffer (e.g. acetate), pH 7-9
- BOTH sorbent functional group & analyte must be ionized

## 2. Condition & Equilibrate

- Condition with 1-2 tube volumes MeOH or MeCN
- Equilibrate with 1-2 tube volumes buffer (used during sample pre-treatment)

## 3. Load sample (consistent rate; 1-2 drops per second)



# Example IOX SPE Protocol

## 4. Wash sorbent (elutes co-retained interferences)

- Wash interferences with buffer
- Wash with 100% MeOH to remove hydrophobic interferences

## 5. Elute analytes of interest

- Adjust **pH opposite to load conditions** (e.g. 2-5% ammonium hydroxide for basic compounds)
- May require organic modifier (50-100% MeOH)

## 6. Evaporate/reconstitute as necessary

# Useful IOX SPE Tips

- IOX kinetics slower than RP & NP => reduce flow rate
- Strong vs. weak ion-exchangers
  - Strong = sorbent functional group always ionized regardless of pH
  - Weak = sorbent functional group has controllable pKa; commonly used for extracting strong analytes
- Counter-Ion Selectivity in IOX

## For Cation Exchangers:

- $\text{Ca}^{2+} > \text{Mg}^{2+} > \text{K}^+ > \text{Mn}^{2+} > \text{RNH}_3^{2+} > \text{NH}_4^+ > \text{Na}^+ > \text{H}^+ > \text{Li}^+$

## For Anion Exchangers:

- Benzene Sulphonate > Citrate >  $\text{HSO}_4^- > \text{NO}_3^- > \text{HSO}_3^- > \text{NO}_2^- > \text{Cl}^- > \text{HCO}_3^- > \text{HPO}_4^- > \text{Formate} > \text{Acetate} > \text{Propionate} > \text{F}^- > \text{OH}^-$


# The Critical Role of pH in SPE

Neutral State (Blue) = promotes hydrophobic (RP) interaction

Ionized State (Green) = promotes electrostatic (IOX) interaction

## Ionization of Acidic & Basic Molecules-



**Acids** (e.g., carboxylic acids): (e.g.,  $R-COOH \leftrightarrow R-COO^-$ )



pKa of most acids (e.g. -COOH) is 3-5

- Presence of halogen atom near a carboxy group strengthens acid effect (electron sink)
- e.g., acetic acid (pKa 4.75), monochloroacetic acid (pKa 2.85), dichloroacetic acid (pKa 1.48)

**Bases** (e.g., amines): (e.g.,  $R-NH_3^+ \leftrightarrow R-NH_2$ )



pKa of most amines is 8-11

- Aromatic (electron sink) amines have a lower pKa than aliphatic amines
- e.g., Aromatic amines- aniline (pKa 4.6), pyridine (pKa 5.2); Aliphatic amines- dimethylamine (pKa 10.7)

# SPE Phase Selection

Your Sample Matrix is:

**Aqueous**  
(biological fluids, water, aqueous extracts of tissues, etc.)

**Organic**  
(organic extracts of tissues, hexane, dichloromethane, etc.)

Recommended Retention Mechanisms:

**Reversed-Phase**  
See page 27 for more details

**Ion-Exchange**  
See page 28 for more details

**Normal-Phase**  
See page 29 for more details

Analyte Characteristics:

**Moderately polar to non-polar compounds**

**Weak cations/ anions**

**Strong cations/ anions**

**Polar to moderately polar compounds**

Application:

Pharma.	Environ.
DSC-18	ENVI™-18
DSC-18Lt	ENVI-8
DSC-8	ENVI-Chrom P
DSC-Ph	LC-18
DSC-CN	LC-8
DPA-6S	LC-Ph
	LC-CN

Pharma.	Environ.
DSC-SAX	LC-SAX
DSC-SCX	LC-SCX

Pharma.	Environ.
DSC-WCX	LC-WCX
DSC-NH <sub>2</sub>	LC-NH <sub>2</sub>
	PSA

Pharma.	Environ.
DSC-Si	ENVI-Florisil®
DSC-CN	LC-Alumina
DSC-Diol	LC-Florisil
DSC-NH <sub>2</sub>	LC-Si
	LC-NH <sub>2</sub>
	LC-Diol
	LC-CN
	PSA

Recommended SPE Phases:

# New SPE Brochure 2007

- T402150 (FEB)
- 28 pages
- Complete list of SPE products and accessories

