

Analysis of Artificial Colorants in Haw Jelly Using SampliQ WAX Coupled with High Performance Liquid Chromatography (HPLC)

Authors

Congcong Zhang and
Xia Yang
Agilent Technologies, Inc.

Abstract

This study developed and validated a method for the quantitative analysis of six artificial colorants in haw jelly candies using the Agilent SampliQ WAX coupled with HPLC analysis. The method delivered a reliable solution with excellent recoveries and reproducibility that can be extended to colorant analysis in other fruit jelly candies.

Experimental

Target analytes

The six target analytes in this application note include Citrine, Amaranth, Carmine, Sunset Yellow, Allura Red, and Brilliant Blue artificial colorants. Figure 1 shows the chemical structures.









Figure 1. Chemical structures of target analytes.

HPLC conditions

Parameter	Value																		
Column	Agilent InfinityLab Poroshell 120 EC-C18, 100 × 4.6 mm, 2.7 μm (p/n 695975-902)																		
Flow Rate	1 mL/min																		
Column Temperature	35 °C																		
Injection Volume	10 μL																		
Detection Wavelength	425 nm (citrine), 500 nm (amaranth, carmine, sunset yellow, allura red), 630 nm (brilliant blue)																		
Mobile Phase	A) 20 mM ammonium acetate pH = 5 B) ACN																		
Gradient	<table border="1"> <thead> <tr> <th>Time (min)</th> <th>%A</th> <th>%B</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>100</td> <td>2</td> </tr> <tr> <td>6.0</td> <td>100</td> <td>30</td> </tr> <tr> <td>11.0</td> <td>88</td> <td>90</td> </tr> <tr> <td>11.5</td> <td>88</td> <td>100</td> </tr> <tr> <td>12.0</td> <td>2</td> <td>2</td> </tr> </tbody> </table>	Time (min)	%A	%B	0	100	2	6.0	100	30	11.0	88	90	11.5	88	100	12.0	2	2
Time (min)	%A	%B																	
0	100	2																	
6.0	100	30																	
11.0	88	90																	
11.5	88	100																	
12.0	2	2																	
Post Time	3 minutes																		

Instrument method

The samples were run on an Agilent 1260 Infinity II LC system with a diode array detector (DAD). Agilent ChemStation software was used for data acquisition and processing.

Sample extraction

Weigh 2 g of haw jelly piece, then add 20 mg of pectinase and 5 mL water followed by 15 minutes of

sonication mixing for enzyme digestion. Add 10 mL of extraction solution (ammonia:water:acetonitrile = 2:48:50), vortex for one minute, then centrifuge for 10 minutes to get the supernatant. Compensate the supernatant to 30 mL with water. Adjust the pH of the sample solution to pH = 6 with formic acid; the sample is now ready for SPE cleanup, following the procedure shown in Figure 2.




Figure 2. Sample preparation workflow chart.

Results and discussion

Table 1. Method recovery (REC) and RSD.

Analytes	0.5 µg/g Spiking		1 µg/g Spiking		5 µg/g Spiking	
	Rec (%)	RSD (%)	Rec (%)	RSD (%)	Rec (%)	RSD (%)
Citrine	99.4	7.9	96.6	2.5	92.2	2.1
Amaranth	95.5	7.4	109.1	3.4	90.7	2.6
Carmines	111.3	6.3	89.3	2.8	88.9	2.0
Sunset Yellow	88.0	8.5	94.0	3.4	90.5	2.1
Allura Red	102.1	7.9	107.3	4.5	90.8	4.4
Brilliant Blue	98.9	8.9	79.7	4.0	88.3	3.3




Figure 3. Chromatograms of six artificial pigments in neat standards solution (A) and haw jelly sample (other three figures: blue—neat standard; red—matrix blank; green—pre-spiked sample with 5 µg/g).

Conclusion

This method using Agilent SampliQ WAX SPE cartridges coupled with HPLC for analysis of six artificial colorants in haw jelly candies provided excellent recoveries of approximately 79% to 112%. This method can be a good reference to detect these six artificial colorants in other fruit jelly candies.

www.agilent.com/chem

DE.3471875

This information is subject to change without notice.

© Agilent Technologies, Inc. 2020
Printed in the USA, March 26, 2020
5994-1843EN