

Shimadzu Electronic Balances General Catalog



Supporting laboratories in the future, utilizing and building on 100 plus years of experience and knowledge

Shimadzu began manufacturing balances in 1918. For more than 100 years, we have been at the forefront of providing precision, quality solutions for the most challenging R&D and QA/QC requirements. Our steadfast customer-focused commitment and unwavering dedication to technical excellence are both hallmarks of our history and the principles that guide us into the future.

SHIMADZU ELECTRONIC BALANCES

A Tradition of Weighing Expertise

Established in 1875 in Kyoto, Japan, Shimadzu Corporation is one of the pioneers of scientific precision instruments.

Top-pan and torsion balance production started in 1918, and equal-beam analytical balances were introduced in 1925. Since their release, the continuous improvement of Shimadzu balances has contributed to research and development across all industries.


Around the turn of the 20th century, precision weighing was a time-consuming practice performed only by experienced operators. Placing the sample and small masses on pans hung from a beam scale with a moving indicator was a tedious process. Shimadzu strove continuously to streamline weighing procedures. The introduction of the direct reading analytical balance (patented in Japan in 1948) signified a new era in weighing technology. In the Type L balance, the sensitive mass-loading work was replaced by

convenient dial operations. This reduced weighing time by 66% and, subsequently, reduced demand for conventional balances.

Shimadzu then added the top-loading direct reading balance with Roberval's mechanism in 1959. Until recently, many of these instruments were still utilized in modern laboratories. Shimadzu continued to pioneer technologies, releasing its first electronic balance in 1971—the Digibalance. This release marked a milestone in precision weighing, introducing simplicity and ease of use to analytical weighing.

Six years later (1977), the application of microprocessors in electronic balances further enhanced weighing performance. The compact ED Series provided substantial improvements in sensitivity, resolution, and stability.





More recently, Shimadzu has introduced user-friendly instruments and features to the market, such as:

- Temperature-based fully-automatic calibration in 1985, the first one-piece force cell (OPF, later renamed UniBloc™) in 1989
- The high-sensitivity AEM-5200 Micro Balance in 1993
- The unique Windows® Direct feature perfectly suited for the computerized laboratory of the 21st Century.

Its most recent achievement is the AP Series, advanced performance balances featuring UniBloc and a high response speed, and which are applicable for a wide range of applications.

Moving forward, Shimadzu is committed to providing innovative products for the analytical marketplace.

Contents


P 04	- UniBloc
P 06	- Functions Menu
P 08	- Product Lineup
P 12	- AP Series
P 18	- STABLO-AP (Ionizer)
P 20	- AU Series AT-R Series
P 24	- UP Series UW/UX Series TW/TX/TXB Series
P 30	- BL Series
P 31	- ELB Series
P 32	- BW-K/BX-K Series
P 33	- EP-100/EP-110
P 34	- MOC63u MOC-120H
P 37	- Specific Gravity Analyzers
P 42	- Animal Balances
P 44	- Analytical Network Data System
P 46	- Optional Accessories
P 48	- Physical Dimensions



UniBloc Goes Beyond the Concept of Precision Balances


The aluminum UniBloc integrated sensor mechanism is precisely machined from a single block of aluminum alloy.

Because it integrates all the components used in the previous sensor unit without the need for any fastening screws, it is fastened without deflecting the elastic fulcrum, resulting in an exceptionally impact-resistant sensor that is unlikely to ever fail.



UniBloc

Illustrations of Balance Mechanism, Assembled Balance, and UniBloc Balance



High responsiveness


Control speed has been improved and the fastest high-speed response is realized by reducing the mass of the lever system. (Comparison with previous model)

High temperature resistance

As a result of the integral block structure, the humidity in the whole cell is more uniform, which improves the temperature characteristics.

High impact resistance

The shock resistance is significantly improved by the compact aluminum integrally formed cell.



One Million-Cycle Endurance Testing

Engineered with the impact-resistant UniBloc sensor, they pass endurance testing with one million cycles.*



Appearance of Endurance Testing

* Endurance Testing

Method: A 1-kg weight is moved on and off the pan at 1.5-second cycles.

Pass/fail criterion: Satisfies Shimadzu standards for routine inspections after the on-off cycles.

Note: Japan's Measurement Act specifies endurance testing to 100,000 cycles.


UP series balances are tested to over ten times that level.

UniBloc family of balances

Experience for yourself the performance of UniBloc.



AP series



AU series



AT-R series



UP series



UW/UX series



TW/TX series



BW-K/BX-K series



MOC-120H



MOC63u

Shimadzu Balances Offer a Diverse Range of Functions

01

High Level Functionality


Easy Setting

During operation, if you want to make the display slightly more stable, or alternatively, want to improve the response speed, you can make one-touch adjustments without interrupting measurement. A special indicator is provided that instantly shows the adjustment status.



Menu Operation Key


Keys exclusively for menu operations are arranged separately from the measurement keys. Menu can be operated intuitively using the cross-shaped key layout.



Perfect Self Calibration (PSC)

Electronic balances are precision instruments very susceptible to changes in room temperature. Sensitivity must be calibrated every time the balance is used since changes in room temperature influence mass measurement values, which are not supposed to change. The balance detects changes in room temperature that affect sensitivity, and automatically starts calibration using built-in weights. As a result, sensitivity errors are always kept within a constant range.

This allows the operator to concentrate on measurement tasks without having to worry about sensitivity calibration.




The perfect self calibration (PSC) function keeps the sensitivity error within a constant range at all times.

Clock-CAL

The balance starts calibration using built-in weights at preset times. If you set calibration times before important measurements (e.g. before starting work in the morning or during the lunch or evening break), the balance will automatically start calibration when the preset time is reached.

This lets you take stable, reliable measurements without worrying about sensitivity calibration.



02

Durability

UniBloc

Developed by Shimadzu, UniBloc is created by high-precision electric discharge wire processing applied to a block of aluminum alloy and replaces the conventional sensor block assembly. As such, it uses no springs or screws. This uniform structure dramatically improves response and temperature characteristics, and the simple yet compact design enhances impact resistance. Balances equipped with UniBloc provide highly reliable mass measurement even during prolonged use.

03

Convenient Functions

Internal Calibration

The balance has built-in motor-driven calibration weights. Sensitivity can be calibrated as needed with a single key press.

Dry Battery Operation

The balance can also run on dry cell batteries, enabling use when no power is available.



Formulation Mode

This is convenient when formulating (preparing) multiple substances.

Interval Timer

Data can be automatically output at time intervals set in the range from 1 second to 99 minutes 59 seconds. This function can be combined with Balance Keys.

Checkweighing

When upper and lower thresholds are set, the balance indicates if the sample weight is within the range (GO), over (HI) or under (LO).

Backlight LCD

The backlight LCD display can be clearly read in the darkest of environments.



Analog Bar Display

Remaining weighing capacity can be seen by the analog bar graph at a glance.

04

ISO / GLP / GMP




Built-in Clock

With the optional printer connected, data can be recorded with date and time stamps. Calibration reports can also be date- and time-stamped, which is ideal for establishing measurement and traceability requirements specified by GLP, GMP and ISO 9001.



ISO Calibration Report

Simply connect an optional printer to automatically print out which balance was calibrated when, and the calibration results. No difficult settings are required. Furthermore, the current date and time can be printed at any time during measurement.



(AUW Series Printout Sample)

05

PC Connection



RS-232C Interface

Equipped with an RS-232C interface for easy integration with other devices and computers.



Computer Connection Function

Systems can be connected to a computer via an optional cable/adaptor kit. For more details, visit the Shimadzu website.



Balance Keys



Multi-Balance Collect

Multi-Balance Collect and Balance Keys software run on Windows®. Directly importing balance data into a PC using these free wares eliminates transcription errors, improves work efficiency, and increases data reliability. Multi-Balance Collect can also collect data continuously at fixed intervals. It is recommended when you want to record changes over time, such as weight changes due to evaporation, on multiple balances.

06

Network



LabSolutions™ Balance

Connecting with LabSolutions lets you save data from balances, HPLC and other analytical instruments to a database and create reports automatically. Uniform data management means no transcription errors and is perfect for security.

07

Applications



Standard Below-weigh Hook

Also compatible with suspension measurements.



Piece Counting

A built-in piece counting function enables balances to be used as parts counters (piece scales).



Specific Gravity Measurement

A specific gravity calculation function based on the immersion method is built in.

Just attach the optional Specific Gravity Measurement Kit to use a balance as a specific gravity meter.




Built-in Animal Measurement Mode

The weight of mice, rats, rabbits, and other small animals can be measured. Stable measurements are obtained even if the animal moves.



Carat Measurement

Results can be displayed in carats when measuring precious stones. Change among many weighing units and functions with a single touch.



08

Printer




Printer

Using an optional printer can enhance support for ISO/GLP/GMP.



Product Lineup







































	Up to 100 g	More than 100 g	More than 300 g	More than 500 g	More than 1 kg	More than 3 kg
0.01 mg	<p>AP125WD-AD (120 g/0.1 mg) (52 g/0.01 mg)</p> <p>AP125WD (120 g/0.1 mg) (52 g/0.01 mg)</p> <p>AUW220D (220 g/0.1 mg) (82 g/0.01 mg)</p> <p>AUW120D (120 g/0.1 mg) (42 g/0.01 mg)</p>	<p>AP225W-AD (220 g/0.01 mg)</p> <p>AP135W-AD (135 g/0.01 mg)</p> <p>AP225WD-AD (220 g/0.1 mg) (102 g/0.01 mg)</p> <p>AP225W (220 g/0.01 mg)</p> <p>AP135W (135 g/0.01 mg)</p> <p>AP225WD (220 g/0.1 mg) (102 g/0.01 mg)</p>				
0.1 mg	<p>ATX84R (82 g/0.1 mg)</p> <p>ATY64R (62 g/0.1 mg)</p>	<p>AP124W</p> <p>AP124X</p> <p>AP124Y (120 g/0.1 mg)</p> <p>AP224W-AD</p> <p>AP224W</p> <p>AP224X</p> <p>AP224Y (220 g/0.1 mg)</p> <p>AUW220</p> <p>AUX220</p> <p>AUY220 (220 g/0.1 mg)</p> <p>ATX124R</p> <p>ATY124R (120 g/0.1 mg)</p> <p>ATX224R</p> <p>ATY224R (220 g/0.1 mg)</p> <p>AUW120</p> <p>AUX120</p> <p>AUY120 (120 g/0.1 mg)</p> <p>TWC623L</p> <p>TXC623L (620 ct/0.001 ct)</p>	<p>AP324W-AD</p> <p>AP324W</p> <p>AP324X</p> <p>AP324Y (320 g/0.1 mg)</p> <p>AUX320</p> <p>AUW320 (320 g/0.1 mg)</p> <p>ATX324R</p> <p>ATY324R (320 g/0.1 mg)</p>			
0.001 g (1 mg)		<p>UP223X</p> <p>UP223Y (220 g/0.001 g)</p> <p>UW220H</p> <p>UX220H (220 g/0.001 g)</p> <p>TW223L</p> <p>TX223L (220 g/0.001 g)</p> <p>BL-220H (220 g/0.001 g)</p>	<p>UP423X</p> <p>UP423Y (420 g/0.001 g)</p> <p>UW420H</p> <p>UX420H (420 g/0.001 g)</p> <p>TW423L</p> <p>TX423L (420 g/0.001 g)</p> <p>TW323L</p> <p>TX323L (320 g/0.001 g)</p> <p>BL-320H (320 g/0.001 g)</p>	<p>UP623X</p> <p>UP623Y (620 g/0.001 g)</p> <p>UP823X</p> <p>UP823Y (820 g/0.001 g)</p> <p>UW820H</p> <p>UX820H (820 g/0.001 g)</p> <p>UW620H</p> <p>UX620H (620 g/0.001 g)</p>	<p>UP1023X</p> <p>UP1023Y (1020 g/0.001 g)</p> <p>UW1020H</p> <p>UX1020H (1020 g/0.001 g)</p>	










* The UW/UX series, TXB series and ELB series will be discontinued when the stocks are sold out.

	Up to 100 g	More than 100 g	More than 300 g	More than 500 g	More than 1 kg	More than 3 kg
0.01 g (10 mg)		ELB200 (200 g/0.01 g) ELB120 (120 g/0.01 g) UW220HV (220 g/0.01 g)	UP422X UP422Y (420 g/0.01 g) UW420S UX420S (420 g/0.01 g) TXB422L (420 g/0.01 g) BL-320S (320 g/0.01 g) ELB300 (300 g/0.01 g) UW420HV (420 g/0.01 g)	UP822X UP822Y (820 g/0.01 g) UW820S UX820S (820 g/0.01 g) TXB622L (620 g/0.01 g) BL-620S (620 g/0.01 g) ELB600 (600 g/0.05 g) UW620HV (620 g/0.01 g)	UP2202X UP2202Y (2200 g/0.01 g) UW2200H UX2200H (2200 g/0.01 g) TX2202L (2200 g/0.01 g) BL-2200L (2200 g/0.01 g)	UP4202X UP4202Y (4200 g/0.01 g) UP6202X UP6202Y (6200 g/0.01 g) UW4200H UX4200H (4200 g/0.01 g) UW6200H UX6200H (6200 g/0.01 g) TX4202L (4200 g/0.01 g) TX3202L (3200 g/0.01 g) BL-3200H (3200 g/0.01 g)
0.1 g (100 mg)				ELB600S (600 g/0.1 g) UW820SV (820 g/0.1 g)	ELB2000 (2000 g/0.1 g) ELB1200 (1200 g/0.1 g) UW2200HV (2200 g/0.1 g)	UP4201X UP4201Y (4200 g/0.1 g) UW4200S UX4200S (4200 g/0.1 g) TXB4201L (4200 g/0.1 g) BL-3200S (3200 g/0.1 g) ELB3000 (3000 g/0.1 g) UW4200HV (4200 g/0.1 g)
0.1 g (100 mg)	UP8201X UP8201Y (8200 g/0.1 g) UW8200S UX8200S (8200 g/0.1 g) TXB6201L (6200 g/0.1 g) UW6200HV (6200 g/0.1 g)	BW22KH BX22KH (22 kg/0.1 g) BW12KH BX12KH (12 kg/0.1 g)	BW32KH BX32KH (32 kg/0.1 g)			
1 g	ELB6000S (6000 g/1 g) UW8200SV (8200 g/1 g)	ELB12K (12 kg/1 g)	BW32KS BX32KS (32 kg/1 g)	BW52KS BX52KS (52 kg/1 g)		

Product Lineup

									
 UniBloc	UniBloc	✓	✓	✓	✓	✓	✓	✓	✓
 PSC	Perfect Self Calibration	✓	✓	✓	✓	✓	✓		✓ (ATX-R only)
 Clock-CAL	Clock-CAL	✓	✓	✓	✓	✓			
 Internal CAL	Internal Calibration	✓	✓	✓	✓	✓	✓		✓ (ATX-R only)
 Single-Lever CAL	Single-Lever CAL								
 Built-in Clock	Built-in Clock	✓	✓	✓	✓	✓	✓		
 ISO Calibration Report	ISO Calibration Report	✓	✓	✓	✓	✓	✓		
 Printer*1	Printer*1	✓	✓	✓	✓	✓	✓	✓	✓
 PC	Computer Connection Function	✓	✓	✓	✓	✓	✓	✓	✓
 Multi-Balance Collect	Multi-Balance Collect	✓	✓	✓	✓	✓	✓	✓	✓
 Balance Keys	Balance Keys	✓	✓	✓	✓	✓	✓	✓	✓
 RS-232C INTERFACE	Built-in RS-232C Interface	✓	✓	✓	✓	✓	✓	✓	✓
 USB INTERFACE	Built-in USB Interface	✓	✓	✓	✓				✓
 Touchless Sensor	Touchless Sensor	✓							
 Smart AutoDoor	Smart Auto Door	✓							
 Menu Operation Key	Menu Operation Key	✓	✓	✓	✓				
 Easy Setting	Easy Setting	✓	✓	✓	✓				✓
 Back Light	Backlight					✓ (AUW only)			
 Organic EL Display	Organic EL Display	✓	✓	✓	✓				
 Analog Bar Display	Analog Bar Display	✓	✓	✓	✓	✓	✓	✓	
 Checkweighing	Checkweighing	✓	✓	✓	✓				✓
 Comparator Output	Comparator Output								
 PCS	Piece Counting	✓	✓	✓	✓	✓	✓	✓	✓
 CARAT	Carat Measurement	✓	✓	✓	✓	✓	✓	✓	✓
 Specific Gravity	Specific Gravity Measurement	✓	✓	✓	✓	✓	✓	✓	✓
 DRY Icon	Dry Battery Operation								
 Standard Below-weigh Hook	Standard Below-weigh Hook	✓	✓	✓	✓	✓	✓	✓	
 Built-in Animal Measurement Mode	Built-in Animal Measurement Mode								
 Formulation Mode	Formulation Mode	✓	✓	✓	✓	✓	✓	✓	✓
 Interval Timer	Interval Timer	✓	✓	✓	✓	✓	✓		

	 UP-X UW*2 *3	 UP-Y UX*3 *4	 TW TX	 TXB*3	 BW-K BX-K	 ELB*3	 BL	 UW-V	 MOC63u MOC-120H
	✓	✓	✓		✓			✓	✓
	✓							✓	
	✓							✓	
	✓		✓ (TW-N only)					✓	
					✓ (BW-K only)				
	✓	✓			✓			✓	✓
	✓	✓			✓			✓	✓
	✓	✓	✓	✓	✓	✓	✓	✓	✓*5
	✓	✓	✓	✓	✓	✓*6	✓*6	✓	✓
	✓	✓	✓	✓	✓	✓	✓	✓	✓ (MOC63u only)
	✓	✓	✓	✓	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓	✓*6	✓*6	✓	✓
									✓ (MOC63u only)
									
			✓	✓					
	✓	✓	✓	✓				✓	✓
									
	✓	✓			✓		✓	✓	
	✓	✓	✓	✓	✓			✓	
	✓	✓						✓	
	✓	✓	✓	✓	✓	✓	✓	✓	
	✓	✓	✓	✓	✓		✓	✓	
	✓	✓			✓	✓		✓	
				✓		✓			
	✓	✓			✓*7	✓*7		✓	
	✓	✓			✓				
	✓	✓	✓	✓	✓			✓	
	✓	✓			✓				✓

*1 Option *2 The pictures show the UP-X. *3 The UW/UX series, TXB series and ELB series will be discontinued when the stocks are sold out. *4 The pictures show the UP-Y. *5 A dedicated printer is available for the MOC-120H. *6 Requires optional I/O-RS conversion cable or interface IFB-102A. *7 Requires an optional below-weigh hook.



AP W-AD

AP with Automatic Door

AP W-AD Series

Advanced Performance UniBloc Balances

Provides high-speed response and high stability

New automatic door functionality makes weighing operations even more convenient

01

Touchless sensors and a Smart Auto Door improve hygiene and lower contamination risk to provide a superior operating environment.

02

An ionizer and adjustable windbreak plate reduce static electricity and convection effects to provide highly stable and reliable measurements.

03

LabSolutions Balance supports weighing data integrity.



Visit our website for more information.

Product



Advanced Performance UniBloc Balances

AP Series

Product 

Provides High-Speed Response and High Stability

New automatic door functionality makes weighing operations even more convenient



Minimum display 0.1 mg model

AP324W	AP324X	AP324Y
AP224W	AP224X	AP224Y
AP124W	AP124X	AP124Y

Minimum display 0.01 mg model

AP225W	AP135W
--------	--------

Minimum display 0.01 mg/0.1 mg model

AP225WD	AP125WD
---------	---------

Minimum display 0.1 mg model















































































AP324W-AD	AP224W-AD
-----------	-----------

Minimum display 0.01 mg model

AP225W-AD	AP135W-AD
-----------	-----------

Minimum display 0.01 mg/0.1 mg model

AP225WD-AD	AP125WD-AD
------------	------------

AP W-AD Series	                    
AP W Series	                   
AP X Series	                   
AP Y Series	                



Aluminum integrated mass sensor

Excellent response time and stable temperature characteristics.



Built-in Clock

If a printer (optional) is connected, then data can be marked with the date and time. Calibration results can also be marked with date and time, which is perfect for managing measurements or establishing traceability as required by GLP/GMP/ISO 9001 standards.



USB Interface

A USB connector is built in for connecting to a PC.



ISO Calibration Report

Simply connect an optional printer to automatically print out which balance was calibrated when, and the calibration results. No difficult settings are required. Furthermore, the current date and time can be printed at any time during measurement.



Organic EL Display

Adopts an "Organic EL Display" that is clear even in dark places.

AP W-AD Series

Model Name	W-AD Series					
	AP225W-AD	AP135W-AD	AP225WD-AD	AP125WD-AD	AP324W-AD	AP224W-AD
Capacity	220 g	135 g	220 g / 102 g	120 g / 52 g	320 g	220 g
Minimum Display	0.01 mg		0.1 mg / 0.01 mg		0.1 mg	
Pan Size	Approx. 91 mm dia.					
Dimensions	Approx. W215 × D411 × H346 mm (power supply unit included)				Approx. W215 × D367 × H346 mm	
Weight	Approx. 9.7 kg			Approx. 8.6 kg		
Calibration Weight	Built-in					
External Calibration Weight Range for Span Calibration (recommended weight value)	95 to 220.00090 g (200 g)	45 to 135.00090 g (100 g)	95 to 220.00090 g (200 g)	45 to 120.00090 g (100 g)	95 to 320.0090 g (300 g)	95 to 220.0090 g (200 g)
Repeatability (at weighing capacity) *1	0.015 mg (to 20 g) 0.03 mg (to 100 g) 0.05 mg (to weighing capacity)	0.05 mg	0.1 mg / 0.05 mg	0.1 mg / 0.02 mg	0.15 mg	0.1 mg
Repeatability (for Low Loads)	0.01 mg (for 5 g load) *1				0.1 mg (for 20 g load)	0.1 mg (for 10 g load)
Minimum Weight	20 mg *1				200 mg	
Linearity *1	±0.1 mg	±0.2 mg / ±0.1 mg		±0.2 mg / ±0.05 mg	±0.3 mg	±0.2 mg
Response Time for Trace Measurements *2,3	2 sec.					
Response Time *3	8 sec.	2 sec. / 8 sec.			2 sec.	
Operational Temperature/ Humidity Range	5 to 40 °C at 20 to 85 % RH *4					
Sensitivity Stability Against Temperature Range	±2 ppm/°C (10 to 30 °C)					
Display	OEL display (dot matrix)					
Rated Electric Power Supply	DC 12 V 1.5 A					
Power Supply Input (AC Adapter) *5	AC 100 to 240V 480 mA 50/60 Hz					
Input/Output Terminal	RS-232C (D-sub 9P plug), USB Host (Type A), USB device (Type B), Ionizer					
Functions, Options	USB Host (Type A), USB Device (Type B), Recipe Compounding, HPLC Buffer Solution Preparation, mol Conversion Function, Sample (Concentration) Preparation, Inspection Support Function, Clock-CAL, Automatic Doors, Touchless Sensors, Ionizer *6					Adjustable Internal Windbreak Plate

*1 Measurement conditions of W-AD series (0.01 mg models only) are as follows:

- Set the adjustable windbreak plate in the lowest position
- With a shield plate configured around the pan

*2 The response time for displaying 90 % of the added sample amount value in trace measurements (from 1 mg).

*3 The response time value is typical.

*4 Non-condensing.

*5 Depending on the attached AC adapter.

*6 See more details of ionizer on page 18.



Touchless Sensors (AP W-AD series only)

Doors open/close by waving a hand over the left and right infrared sensors that enables door operation without touching the balance.



Smart Auto Door (AP W-AD series only)

The automatic doors include automatic learning functionality that enables freely setting how far to open/close each glass door. That minimizes external air effects and increases operational efficiency.



Computer Connection Function

Systems can be connected to a computer via an RS-232C cable or using a USB–serial adapter kit. For more details, visit the Shimadzu website.

AP W/X/Y Series

Model name	W Series						
	AP225W	AP135W	AP225WD	AP125WD	AP324W	AP224W	AP124W
Capacity	220 g	135 g	220 g / 102 g	120 g / 52 g	320 g	220 g	120 g
Minimum Display	0.01 mg		0.1 mg / 0.01 mg		0.1 mg		
Pan Size	Approx. 91 mm dia.						
Dimensions	Approx. W213 × D411 × H345 mm (power supply unit included)				Approx. W213 × D367 × H345 mm		
Weight	Approx. 7.9 kg				Approx. 7.0kg		
Calibration Weight	Built-in						
External Calibration Weight Range for Span Calibration (recommended weight value)	95 to 220.00090 g (200 g)	45 to 135.00090 g (100 g)	95 to 220.00090 g (200 g)	45 to 120.00090 g (100 g)	95 to 320.0090 g (300 g)	95 to 220.0090 g (200 g)	45 to 120.0090g (100 g)
Repeatability (at weighing capacity)	0.015 mg (to 20 g) 0.03 mg (to 100 g) 0.05 mg (to weighing capacity)	0.05 mg	0.1 mg / 0.05 mg	0.1 mg / 0.02 mg	0.15 mg	0.1 mg	
Repeatability (for Low Loads)	0.01 mg (for 5 g load)				0.1 mg (for 20 g load)	0.1 mg (for 10 g load)	0.1 mg (for 5 g load)
Minimum Weight	20 mg				200 mg		
Linearity	±0.1 mg		±0.2 mg / ±0.1 mg		±0.2 mg / ±0.05 mg	±0.3 mg	±0.2 mg
Response Time for Trace Measurements *1,2	2 sec.						
Response Time *2	8 sec.		2 sec./8 sec.		2 sec.		
Operational Temperature/Humidity Range	5 to 40 °C at 20 to 85 % RH *3						
Sensitivity Stability Against Temperature Range	±2 ppm/°C (10 to 30 °C)						
Display	OEL display (dot matrix)						
Rated Electric Power Supply	DC 12 V 1.0A						
Power Supply Input (AC Adapter) *4	AC100 to 240V 320 mA 50/60 Hz						
Input/Output Terminal	RS-232C (D-sub 9P plug)		USB host (Type A)		USB device (Type B)		Ionizer
Functions, Options	USB Host (Type A)		✓				
	USB Device (Type B)		✓				
	Recipe Compounding		✓				
	HPLC Buffer Solution Preparation		✓				
	mol Conversion Function		✓				
	Sample (Concentration) Preparation		✓				
	Inspection Support Function		✓				
	Clock-CAL		✓				
	Internal Windbreak Plate		✓		✓ (optional)		
	Ionizer				✓ (optional)		

Model name	X Series			Y Series		
	AP324X	AP224X	AP124X	AP324Y	AP224Y	AP124Y
Capacity	320 g	220 g	120 g	320 g	220 g	120 g
Minimum Display	0.1mg					
Pan Size	Approx. 91 mm dia.					
Dimensions	Approx. W213 × D367 × H345 mm					
Weight	Approx. 7.0 kg			Approx. 6.5 kg		
Calibration Weight	Built-in			None		
External Calibration Weight Range for Span Calibration (recommended weight value)	95 to 320.0090 g (300 g)	95 to 220.0090 g (200 g)	45 to 120.0090 g (100 g)	95 to 320.0090 g (300 g)	95 to 220.0090 g (200 g)	45 to 120.0090 g (100 g)
Repeatability (at weighing capacity)	0.15 mg	0.1 mg		0.15 mg	0.1 mg	
Repeatability (for Low Loads)	0.1 mg (for 20 g load)	0.1 mg (for 10 g load)	0.1 mg (for 5 g load)	0.1 mg (for 20 g load)	0.1 mg (for 10 g load)	0.1 mg (for 5 g load)
Minimum Weight	200 mg					
Linearity	±0.3 mg	±0.2 mg		±0.3 mg	±0.2 mg	
Response Time for Trace Measurements *1,2	2 sec.					
Response Time *2	2 sec.					
Operational Temperature/Humidity Range	5 to 40°C at 20 to 85% RH *3					
Sensitivity Stability Against Temperature Range	±2 ppm/°C (10 to 30°C)					
Display	OEL display (dot matrix)					
Rated Electric Power Supply	DC 12 V 1.0 A					
Power Supply Input (AC Adapter) *4	AC100 to 240V 320 mA 50/60 Hz					
Input/Output Terminal	RS-232C (D-sub 9P plug)		USB host (Type A)		USB device (Type B) Ionizer	
Functions, Options	USB Host (Type A)		✓			
	USB Device (Type B)		✓			
	Recipe Compounding		✓			
	HPLC Buffer Solution Preparation		✓			
	mol Conversion Function		✓			
	Sample (Concentration) Preparation		✓			
	Inspection Support Function		✓			
	Clock-CAL		✓			
	Internal Windbreak Plate		✓		✓ (optional)	
	Ionizer		✓ (optional)		✓ (optional)	

*1 The response time for displaying 90 % of the added sample amount value in trace measurements (from 1 mg).

*3 Non-condensing.

*2 The response time value is typical.

*4 Depends on the attached AC adapter.

Various accessories and options suitable for semi-micro measurement are available.

Multi-Stand

(included standard with 0.01 mg models of W-AD series only)

If placing weighing paper, microtubes, or other containers that exceed the pan diameter, or when weighing long rod-like samples, attach a specialized multi-stand to easily weigh samples.

- Example Using a Multi-Stand



Internal Windbreak Plate

(included standard with 0.01 mg models of W series only)



The plate suppresses the influence of convection and airflow within the weighing chamber to improve the stability and response of measurement values.



Internal Windbreak Plate (for W/X/Y series models)

Static Electricity Remover (Ionizer)

STABLO^{AP}

Use as an external stand configuration or install it inside the balance unit.



When using a stand



Built-in

Other Optional Products



SMK-601
Specific Gravity Measurement Kit



AP Holder



EP-100



EP-110

Optional Accessories

Description
STABLO ^{AP} Ionizer Static Electricity Remover
EP-100 Electronic Printer
EP-110 Electronic Printer
Label Roll Paper for EP-100/110 (10 Rolls)
Internal Windbreak Plate (for W/X/Y Series) *1
SMK-601 Specific Gravity Measurement Kit
AP Holder *2
Multi-Stand *3

Description
Shield plate
AC Adapter (for W/X/Y Series)
AC Adapter (for W-AD Series Balances)
AC Adapter (for W-AD Series STABLO-AP Ionizers)
Display Protective Cover (Set of 5)
USB Cable Assembly (2 m) with Core
RS-IO Adapter Cable (for Connecting EP-80/90)

*1 Included standard with 0.01 mg models of W series only

*2 Included standard with AP225W-AD/AP225W models

*3 Included standard with 0.01 mg models of W-AD series only

Static Remover (Ionizer)



An excellent solution against static electricity.

STABLO-AP provides reliable measurement by removing static electricity.

- 1 On Stand
- 2 Hand-held
- 3 Built-in Balance (AP series)




Easily install in the balance

Features of STABLO-AP

Static Electricity Removal by Ion Irradiation

With the high-frequency AC corona discharge method, Shimadzu's STABLO-AP ionizer provides a stable ion balance and excellent static removal performance on samples and containers. Precision weighing work becomes remarkably efficient. Electrodes are safely housed inside the unit.


Concept of static removal



AC Method Produces Excellent Ion Balance

AC method: AC voltage is applied on the discharge needle and a well-balanced mixture of positive/negative ions is emitted in rapid alternation from one electrode.


DC method: DC voltage is applied to a couple of electrodes. One is positive and the other is negative. Each electrode emits ions of one polarity only. An effective static removal angle is limited if the two electrodes are distanced. As electrodes deteriorate, the initial ion balance is lost.



What is "ion balance"?


Ion balance is the balance of positive and negative ions that are supplied by an ionizer. If ion balance is poor, static electricity is not removed or inverse charging may result.

If ion balance is poor.....



Applications


Static electricity keeps the sample out of the ampoule



The sample is hard to handle because it adheres to the ampoule inlet and sides.



STABLO-AP removes the charge from the ampoule.



The static charge is gone in seconds. This improves productivity.


Plastic wrap sticks to rubber gloves



Plastic wrap adheres to rubber gloves, making it difficult to work with them.



Fasten STABLO-AP to the stand, and remove the static from the gloves.



The static is removed in about 10 seconds, and the plastic wrap no longer sticks.

STABLO-AP is convenient when using an electronic balance



When the powder in the Petri dish is electrically charged, and the numerical value fluctuates



When the powdered medicine paper is electrically charged, and the numerical value is unstable



When the measurement spoon is electrically charged, and bringing it near the pan affects the numerical value

Specifications

Ion Generation Method	AC corona discharge method
Ion Balance	±10 V
Effective Static Removal Range	Approx. 400 mm from the outlet
Static Elimination Time	Approx. 1 second (typical value) (from ±1000 V to ±100 V)
Ozone Concentration	0.06 ppm
Electrode Probes	Tungsten (durability: 30,000 hours)
Weight	Approx. 710 g (Main unit: 395 g, Stand: 315 g)
Operating Temperature and Humidity	0 °C to + 40 °C, 25 % RH to 85 % RH (non-condensing)
Rated Electric Power Supply	DC 24 V, 1.0 A

Multi Functional Analytical Balances

AU Series

Product 



































































UniBloc Analytical Balances

AUW-D series dual-range semi-micro balances
 AUW/AUX/AUY series analytical balances



Excellent Weighing Performance


- Compact UniBloc mechanism and digital processing technology produce fast response and stability at the same time.
- Measurement results can be transmitted to Excel® or other Windows® applications.
- Piece counting, various mass units, below-weigh hook, and specific gravity measurement software are all standard features.

AUW-D Series																			
AUW Series																			
AUX Series																			
AUY Series																			

Model Name	AUW-D Series		AUW Series			AUX Series			AUY Series	
	AUW220D	AUW120D	AUW320	AUW220	AUW120	AUX320	AUX220	AUX120	AUY220	AUY120
Capacity	220 g / 82 g	120 g / 42 g	320 g	220 g	120 g	320 g	220 g	120 g	220 g	120 g
Minimum Display	0.1 mg / 0.01 mg		0.1 mg							
Pan Size	Approx. 80 mm dia.									
Dimensions	Approx. W217 × D356 × H338 mm									
Weight	Approx. 7.0 kg									


PSC (Perfect Self-Calibration)
(AUW-D/AUW/AUX series only)

Automatically calibrates the balance if an ambient temperature change occurs that could affect sensitivity.



Clock-CAL

This automatically calibrates the balance at pre-specified times (such as before starting work, during lunch, or after work hours).



Internal Calibration
(AUW-D/AUW/AUX series only)

The balance has built-in motor-driven calibration weights. Sensitivity can be calibrated as needed with a single key press.

ISO Calibration Report

Simply connect an optional printer to automatically print out which balance was calibrated when, and the calibration results. No difficult settings are required. Furthermore, the current date and time can be printed at any time during measurement.

Interval Timer
(AUW-D/AUW/AUX series only)

Data can be automatically output at time intervals set in the range from 1 second to 99 minutes 59 seconds. This function can be combined with Balance Keys.

Specific Gravity Measurement

Installing the optional SMK-401 specific gravity measurement kit transforms the balance into a dedicated instrument for measuring specific gravity or density. Specific gravity measurement software is already installed in the Shimadzu balance (refer to page 38).

Piece Counting

A built-in piece counting function enables balances to be used as parts counters (piece scales).

Carat Measurement

Results can be displayed in carats when measuring precious stones. Change among many weighing units and functions with a single touch.

RS-232C Interface

Equipped with an RS-232C interface for easy integration with other devices and computers.

Printer

By adding an optional printer, the measurement results can be automatically printed (refer to page 33).

Optional Accessories

Description
STABLO-AP Ionizer Static Electricity Remover
EP-100 Electronic Printer
EP-110 Electronic Printer
SMK-401 Specific Gravity Measurement Kit
AP Holder
In-use Protective Cover
RS-232C Cable
USB Conversion Cable
AKB-301 Application Keyboard
FSB-102PK Foot Switch
FSB-102TK Foot Switch
AC Adapter (Provided with the main unit)




AKB-301
Application Keyboard



AP Holder



STABLO™-AP



SMK-401
Specific Gravity Measurement Kit

Basic Models of Analytical Balances


AT-R Series

Product 

Basic Model with Improved Convenience

- Freely set responsiveness or stability
- Perform stable measurements in combination with an ionizer
- Connect with a computer via USB or RS-232C




ATX-R Series															
ATY-R Series															

Model Name	ATX-R Series				ATY-R Series			
	ATX324R	ATX224R	ATX124R	ATX84R	ATY324R	ATY224R	ATY124R	ATY64R
Capacity	320 g	220 g	120 g	82 g	320 g	220 g	120 g	62 g
Minimum Display	0.1 mg							
Pan Size	Approx. 91 mm dia.							
Dimensions	Approx. W213 × D356 × H338 mm							
Weight	Approx. 6.2 kg				Approx. 6.0 kg			

PSC (Perfect Self-Calibration)
(ATX-R only)

Automatically calibrates the balance if an ambient temperature change occurs that could affect sensitivity.



Internal Calibration
(ATX-R only)

The balance has built-in motor-driven calibration weights. Sensitivity can be calibrated as needed by a single key press.

Checkweighing

When upper and lower thresholds are set, the balance indicates if the sample weight is within the range (GO), over (HI) or under (LO).

Computer Connection Function

Systems can be connected to a computer via an optional cable/adaptor kit. For more details, visit the Shimadzu website.

RS-232C Interface

Equipped with an RS-232C interface for easy integration with other devices and computers.

USB Interface

A USB connector is built in for connecting to a PC.

Piece Counting

A built-in piece counting function enables balances to be used as parts counters (piece scales).

Carat Measurement

Results can be displayed in carats when measuring precious stones. Change among many weighing units and functions with a single touch.

Percentage Measurement


Measures a percentage value with respect to a preset reference.

Formulation Mode

This is convenient when formulating (preparing) multiple substances.

Easy Setting

During operation, if you want to make the display slightly more stable, or alternatively, want to improve the response speed, you can make one-touch adjustments without interrupting measurement. A special indicator is provided that instantly shows the adjustment status.



Optional Accessories

Description
STABLO-AP Ionizer Static Electricity Remover
EP-100 Electronic Printer
EP-110 Electronic Printer
USB Cable Set
SMK-501 Specific Gravity Measurement Kit
AP Holder
Protective Cover (5 pc set)
AC Adapter (Provided with the main unit)



STABLO™-AP



SMK-501
Specific Gravity Measurement Kit



AP Holder

Multi Functional Top-Loading Balances

UP Series

Product 

Top-loading Balances with the Fastest Response Performance in Their Class



Large Pan with 0.1 g Minimum Display Value

UP4201X UP4201Y
UP8201X UP8201Y



Small Pan with 0.01 g Minimum Display Value

UP422X UP422Y
UP822X UP822Y



Large Pan with 0.01 g Minimum Display Value

UP2202X UP2202Y UP4202X UP4202Y UP6202X UP6202Y



Small Pan with 0.001 g Minimum Display Value

UP223X UP223Y UP423X UP423Y
UP623X UP623Y UP823X UP823Y
UP1023X UP1023Y

* Windbreak can be removed.

UP-X Series																					
UP-Y Series																					

UP-X Series

Model name	UP223X	UP423X	UP623X	UP823X	UP1023X	UP2202X	UP4202X	UP6202X	UP422X	UP822X	UP4201X	UP8201X
Capacity	220 g	420 g	620 g	820 g	1020 g	2200 g	4200 g	6200 g	420 g	820 g	4200 g	8200 g
Minimum Display	0.001 g					0.01 g				0.1 g		
Pan Size	Approx. W113 × D110 mm					Approx. W176 × D189 mm			Approx. W113 × D110 mm		Approx. W176 × D189 mm	
Dimensions	Approx. W191 × D317 × H79 mm					Approx. W191 × D317 × H82 mm			Approx. W191 × D317 × H79 mm			
Weight	Approx. 3.3 kg					Approx. 4.7 kg			Approx. 3.3 kg		Approx. 4.7 kg	

UP-Y Series

Model name	UP223Y	UP423Y	UP623Y	UP823Y	UP1023Y	UP2202Y	UP4202Y	UP6202Y	UP422Y	UP822Y	UP4201Y	UP8201Y
Capacity	220 g	420 g	620 g	820 g	1020 g	2200 g	4200 g	6200 g	420 g	820 g	4200 g	8200 g
Minimum Display	0.001 g					0.01 g				0.1 g		
Pan Size	Approx. W113 × D110 mm					Approx. W176 × D189 mm			Approx. W113 × D110 mm		Approx. W176 × D189 mm	
Dimensions	Approx. W191 × D317 × H79 mm					Approx. W191 × D317 × H82 mm			Approx. W191 × D317 × H79 mm			
Weight	Approx. 2.5 kg					Approx. 2.9 kg			Approx. 2.5 kg		Approx. 2.9 kg	

PSC **PSC (Perfect Self-Calibration)**
(UP-X series only)

Automatically calibrates the balance if an ambient temperature change occurs that could affect sensitivity.

Clock-CAL
(UP-X series only)

This automatically calibrates the balance at pre-specified times (such as before starting work, during lunch, or after work hours).

Internal Calibration
(UP-X series only)

The balance has built-in motor-driven calibration weights. Sensitivity can be calibrated as needed by a single key press.

Built-in Clock

If a printer (optional) is connected, then data can be marked with the date and time. Calibration results can also be marked with date and time, which is perfect for measurement and traceability requirements specified by GLP/GMP/ISO 9001 standards.

ISO Calibration Report

Simply connect an optional printer to automatically print out which balance was calibrated when, and the calibration results. No difficult settings are required. Furthermore, the current date and time can be printed at any time during measurement.

Backlight LCD

The backlight LCD display can be clearly read in the darkest of environments.

Checkweighing

When upper and lower thresholds are set, the balance indicates if the sample weight is within the range (GO), over (HI) or under (LO).

Comparator Output

Acceptable, high, or low weight results or other pass/fail results can be indicated by an audible sound or output externally as a contact signal. (Optional comparator audible indicator or relay output interface required.)

Carat Measurement

Results can be displayed in carats when measuring precious stones. Change among many weighing units and functions with a single touch.

Piece Counting

A built-in piece counting function enables the balance to be used as a parts counter (counting balance).

Computer Connection Function

Systems can be connected to a computer via an RS-232C cable or using a USB-serial adapter kit. For more details, visit the Shimadzu website.

Specific Gravity Measurement

The built-in specific gravity measurement function is based on the liquid immersion method. By installing the specific gravity measurement kit (optional), the balance can be used as a hydrometer (refer to page 40).

Response for Weighing Minute Quantities


Large-pan model with 0.01 g minimum display value

Small-pan model with 0.001 g minimum display value

Measurement Conditions:

Large-pan model with 0.01 g minimum display value and 0.1 g load

Small-pan model with 0.001 g minimum display value and 0.01 g load



Optional Accessories

Description
STABLO-AP Ionizer Static Electricity Remover
EP-100 Electronic Printer
EP-110 Electronic Printer
RS-232C Cable 25P-9P (1.5 m)
USB-Serial Replacement Kit
AKB-301 Application Keyboard
Windbreak Set *1
Windbreak Set (For large pan) *2
WBC-102 Glass Windbreak
WBC-502 Large Size Windbreak
In-use Protective Cover (3 pcs) (For large-pan models with 0.01 g min. display value)
In-use Protective Cover (5 pcs) (For large-pan models with 0.1 g min. display value)

Description
In-use Protective Cover (5 pcs) (For small-pan model)
In-use Protective Cover (5 pcs) (For display and key part)
Animal Bucket Set (For large-pan models with 0.01 g minimum display value)
Animal Bucket Set (For large-pan models with 0.1 g minimum display value)
SMK-101A Specific Gravity Measurement Kit (For large-pan model) (Up to approximately 2 kg for weighing capacity)
SMK-102 Specific Gravity Measurement Kit (For small-pan model)
IFB-RY1 Relay Output Interface
RY1 Connection Cable
AC Adapter *3

*1 Included standard with small pan models with 0.001 g minimum display value.

*2 Included standard with large pan models with 0.01 g minimum display value.

*3 Included standard with main unit.

Multi Functional Top-Loading Balances

UW/UX Series

Product 

Top-loading Balances with UniBloc Technology



UW420S UX420S
UW820S UX820S

UW220H UW620H UW1020H
UW420H UW820H
UX220H UX620H UX1020H
UX420H UX820H

UW4200S UW8200S UW2200H UW4200H UW6200H
UX4200S UX8200S UX2200H UX4200H UX6200H

* Windbreak can be removed.

UW Series																				
UX Series																				

UW Series


Model Name	UW220H	UW420H	UW620H	UW820H	UW1020H	UW420S	UW820S	UW2200H	UW4200H	UW6200H	UW4200S	UW8200S
Capacity	220 g	420 g	620 g	820 g	1020 g	420 g	820 g	2200 g	4200 g	6200 g	4200 g	8200 g
Minimum Display	0.001 g					0.01 g					0.1 g	
Pan Size	Approx. W113 × D110 mm						Approx. W176 × D189 mm					
Dimensions	Approx. W191 × D317 × H79 mm											
Weight	Approx. 3.4 kg						Approx. 4.6 kg					

UX Series

Model Name	UX220H	UX420H	UX620H	UX820H	UX1020H	UX420S	UX820S	UX2200H	UX4200H	UX6200H	UX4200S	UX8200S
Capacity	220 g	420 g	620 g	820 g	1020 g	420 g	820 g	2200 g	4200 g	6200 g	4200 g	8200 g
Minimum Display	0.001 g					0.01 g					0.1 g	
Pan Size	Approx. W113 × D110 mm						Approx. W176 × D189 mm					
Dimensions	Approx. W191 × D317 × H79 mm											
Weight	Approx. 2.7 kg						Approx. 2.9 kg					


PSC (Perfect Self-Calibration)
(UW only)

Automatically calibrates the balance if an ambient temperature change occurs that could affect sensitivity.



Clock-CAL (UW only)

This automatically calibrates the balance at pre-specified times (such as before starting work, during lunch, or after work hours).



Internal Calibration
(UW only)

The balance has built-in motor-driven calibration weights. Sensitivity can be calibrated as needed by a single key press.

ISO Calibration Report

Simply connect an optional printer to automatically print out which balance was calibrated when, and the calibration results. No difficult settings are required. Furthermore, the current date and time can be printed any time during measurement.

Carat Measurement

Results can be displayed in carats when measuring precious stones. Change among many weighing units and functions with a single touch.

Backlight LCD

The backlight LCD display can be clearly read in the darkest of environments.

Checkweighing

When upper and lower thresholds are set, the balance indicates if the sample weight is within the range (GO), over (HI) or under (LO).

Piece Counting

A built-in piece counting function enables balances to be used as parts counters (piece scales).

Specific Gravity Measurement

The built-in specific gravity measurement function is based on the liquid immersion method. By installing the specific gravity measurement kit (optional), the balance can be used as a hydrometer (refer to page 41).

Computer Connection Function

Systems can be connected to a computer via an RS-232C cable or using a USB-serial adapter kit. For more details, visit the Shimadzu website.



Data transfer port of UW/UX Series

Optional Accessories

Description
STABLO-AP Ionizer Static Electricity Remover
EP-100 Electronic Printer
EP-110 Electronic Printer
IFB-102A RS-232C Interface (Needed only for multiple connection)
Small Size Windbreak (For models with capacity of 220 g to 1020 g only) (Std. acc. for models with readability of 0.001 g)
Glass Windbreak (For models with capacity of 220 g to 1020 g only)
Large Size Windbreak (For all models)
SMK-101 Specific Gravity Measurement Kit (For large pan 170 × 180 mm)
SMK-102 Specific Gravity Measurement Kit (For small pan 108 × 105 mm)
Protective In-use Cover for Key Panel and Display (5 pcs)
Small Animal Bucket Set (For large-pan models only)
FSB-102PK Foot Switch (For printing)
FSB-102PK Foot Switch (For taring)
RS-232C Cable, for IBM PC/AT Compatibles (25P-9P, Null modem, 1.5 m)
RS-232C Cable, for Multiple Connections (25P-25P, Null modem, 1.5 m)
AKB-301 Application Keyboard



WBC-102
Glass Windbreak



WBC-502
Large Size Windbreak



STABLO™-AP



AKB-301
Application Keyboard




























TW/TX/TXB Series

Product 

The beginning of the standard.
Extremely capable,
but easy to operate.



TX3202L

TW Series														
TX Series														
TXB Series														

TW Series

Model Name	TW223L	TW323L	TW423L	TWC623L
Capacity	220 g	320 g	420 g	620 ct (124 g)
Minimum Display	0.001 g			0.001 ct (0.0002 g)
Pan Size	Approx. 110 mm dia.			Approx. 80 mm dia.
Dimensions	Approx. W206 × D291 × H241 mm			
Weight	Approx. 4.2 kg			Approx. 4.1 kg



TW223L TW323L TW423L
TX223L TX323L TX423L



TX2202L TX3202L
TX4202L



TWC623L TXC623L

TX Series

Model Name	TX223L	TX323L	TX423L	TX2202L	TX3202L	TX4202L	TXC623L
Capacity	220 g	320 g	420 g	2200 g	3200 g	4200 g	620 ct (124 g)
Minimum Display	0.001 g			0.01 g			0.001 ct (0.0002 g)
Pan Size	Approx. 110 mm dia.			Approx. W167 × D181			Approx. 80 mm dia.
Dimensions	Approx. W206 × D291 × H241 mm			Approx. W200 × D291 × H80 mm			Approx. W206 × D291 × H241 mm
Weight	Approx. 3.8 kg			Approx. 2.8 kg			Approx. 3.8 kg

TXB Series (This series will be discontinued when the stocks are sold out.)

Model Name	TXB422L	TXB622L	TXB4201L	TXB6201L
Capacity	420 g	620 g	4200 g	6200 g
Minimum Display	0.01 g		0.1 g	
Pan Size	Approx. 110 mm dia.		Approx. 160 mm dia.	
Dimensions	Approx. W199 × D260 × H77 mm			
Weight	Approx. 1.5 kg			



TXB6201L TXB4201L



TXB622L TXB422L



Internal Calibration

The balance has built-in motor-driven calibration weights. Sensitivity can be calibrated as needed by a single key press. (TW only)



Easy Setting

During operation, if you want to make the display slightly more stable, or alternatively, want to improve the response speed, you can make one-touch adjustments without interrupting measurement. A special indicator is provided that instantly shows the adjustment status.



Piece Counting

A built-in piece counting function enables balances to be used as parts counters (piece scales).



Carat Measurement

Results can be displayed in carats when measuring precious stones. Change among many weighing units and functions with a single touch.



Checkweighing

When upper and lower thresholds are set, the balance indicates if the sample weight is within the range (GO), over (HI) or under (LO).




Dry Battery Operation

The balance can also run on dry cell batteries, enabling use when no power is available.



Menu Operation Key

Keys exclusively for menu operations are arranged separately from the measurement keys. Menus can be operated intuitively using the cross-shaped key layout.



Optional Accessories

Description
EP-100 Electronic Printer
EP-110 Electronic Printer
RS-232C Cable

Description
In-use Protective Cover (5 pcs)
In-use Protective Cover for Display
USB Conversion Kit
AC Adapter (Included standard with main unit)

Basic Top-Loading Balances

BL Series

Product 

High-Resolution Balances Made Affordable



Piece Counting

A built-in piece counting function enables balances to be used as parts counters (piece scales).



Analog Bar Display

Remaining weighing capacity can be seen by the analog bar graph at a glance.

BL Series



Model Name	BL-220H	BL-320H	BL-320S	BL-620S	BL-2200H	BL-3200H	BL-3200S
Capacity	220 g	320 g		620 g	2200 g	3200 g	
Minimum Display	0.001 g		0.01 g			0.1 g	
Pan Size	Approx. W100 × D100 mm			Approx. W159 × D126 mm			
Dimensions	Approx. W171 × D243 × H77 mm			Approx. W171 × D243 × H78 mm			
Weight	Approx. 2.2 kg						

Optional Accessories

Description
EP-100 Electronic Printer
EP-110 Electronic Printer
In-use Protective Cover
I/O-RS Cable
USB Serial Adapter
Simplified Windshield (Main) *1
Simplified Windshield (Lid) *1
AC Adapter (Included standard with main unit)



Data transfer port of BL Series

*1: Simplified windshield can be used for BL-220H, 320H, 320S. Also, it is included as standard for BL-220H and 320H.

Basic Portable Electronic Balances

ELB Series

Product 

| Precision without Compromise



ELB120
ELB200
ELB300



ELB600S
ELB600S5
ELB600
ELB1200
ELB2000
ELB3000
ELB12K



Dry Battery Operation

The balance can also run on dry cell batteries, enabling use when no power is available.



Specific Gravity Measurement

The built-in specific gravity measurement function is based on the liquid immersion method. By installing the specific gravity measurement kit (optional), the balance can be used as a hydrometer.



Piece Counting

A built-in piece counting function enables balances to be used as parts counters (piece scales).

ELB Series



Model Name	ELB120	ELB200	ELB300	ELB600	ELB600S	ELB1200	ELB2000	ELB3000	ELB6000S	ELB12K
Capacity	120 g	200 g	300 g	600 g		1200 g	2000 g	3000 g	6000 g	12 kg
Minimum Display	0.01 g			0.05 g	0.1 g				1 g	
Pan Size	Approx. 110 mm dia.			Approx. W170 × D130 mm						
Dimensions	Approx. W188 × D216 × H58 mm									
Weight	Approx. 1.3 kg									

Optional Accessories

Description
EP-100 Electronic Printer
EP-110 Electronic Printer
SMK-201 Specific Gravity Measurement Kit (Except for ELB120, 200, 300 for rectangular-pan models only)
Carrying Case
In-use Protective Cover
Below-weigh Hook (Except for ELB12K)
I/O-RS Cable
USB Serial Adapter
AC Adapter (Included standard with main unit)



Data transfer port of ELB Series

Precision Balances for Heavy Samples

BW-K/BX-K Series

Product 

































The Shimadzu precision platform balances have been engineered with the innovative UniBloc mechanism since 1989. Powerful features support any imaginable weighing application. The BW-K Series includes internal calibration weight.



BW12KH

 Single-Lever CAL

The balance has built-in calibration weights. Sensitivity is calibrated with a simple lever operation. Sensitivity can be calibrated easily, whenever needed.

BW-K Series																
BX-K Series																


Model Name	BW-K Series					BX-K Series				
	BW12KH	BW22KH	BW32KH	BW32KS	BW52KS	BX12KH	BX22KH	BX32KH	BX32KS	BX52KS
Capacity	12 kg	22 kg	32 kg		52 kg	12 kg	22 kg	32 kg		52 kg
Minimum Display	0.1 g		1 g			0.1 g		1 g		
Pan Size	Approx. W347 × D248 mm									
Dimensions	Approx. W382 × D366 × H125 mm					Approx. W360 × D366 × H115 mm				
Weight	Approx. 16.5 kg					Approx. 10.5 kg				

Optional Accessories

Description
EP-100 Electronic Printer
EP-110 Electronic Printer
FSB-102PK Foot Switch (For printing)
AKB-301 Application Keyboard
USB-Serial Conversion Kit
RS-232C Cable
Below-weigh Hook
AC Adapter (Included standard with main unit)



AKB-301 Application Keyboard



Hook assy for below weighing

Enhanced Support for ISO/GLP/GMP

EP-100/EP-110

Product 

- Built-In clock
- Customized printing
- Easy communication settings



Compatible Balance Models

AP, AU, AT-R, AT, UP, UW/UX, TX, TXB, BX/BW-K, BL, and ELB series, and MOC63u moisture analyzers.

Note: The automatic setting function cannot be used with models that do not include the PRINT key, such as ELB series balances and MOC63u moisture analyzers.

Specifications

Model Name	EP-100	EP-110
Display		OLED 128 x 64 Dot Matrix Display Easy-to-understand fluorescent dot matrix display
Protected Date Setting		Password protectable (six-character)
Printing	Paper for printing: Regular paper (does not fade with age) Method: 8-pin reciprocating impact dot matrix Speed: Approx. 1.7 lines/sec. Printer head life: 1 million lines Character size: Approx. W1.7 x H2.6 mm	
Interface	USB B-Type female, RS-232 (D-sub 9-pin male)	
Power Supply	AC adapter: Input 100 to 240 V AC, 50/60 Hz; Output 12 V DC/1500 mA Power consumption: 8 W (while printing) Standby power: 0.5 W (when not printing)	
Battery		1500 to 2500 mAh capacity rechargeable nickel-metal hydride (NiMH) batteries can be used (four AA cells). Note: Dry cell batteries cannot be used.
Installation Environment	Temperature: 5 to 45 °C; Humidity: 10 to 80 % No condensation	

		Printout Samples	
		Normal Mode	Statistical Calculation Mode
Manufacturer Information	Shimadzu Corporation	Shimadzu Corporation	Shimadzu Corporation
Device Name	Model: AP225W	Model: AP225W	Model: AP225W
Serial No.	S/N: 000000000	S/N: 000000000	S/N: 000000000
Sample Name (ID)	Device ID: 0000	Device ID: 0000	Device ID: 0000
Date	Date: 2022/12/28	Date: 2022/12/28	Date: 2022/12/28
Measurement Start Time	Start Time: 13:52:52	Start Time: 13:52:52	Start Time: 15:33:35
Measurement Values	001: 100.00029g 002: 100.00036g 003: 100.00037g	001: 100.00029g 002: 100.00036g 003: 100.00037g	001: 100.00035g 002: 100.00032g 003: 100.00032g
Measurement End Time	End Time: 13:53:04	End Time: 13:53:04	End Time: 15:33:45
Signature Field	Signature:	Signature:	Signature:

Output Items

Item	Symbol	Remarks
Title (Header)		Manufacturer information, device name, serial number (S/N), date, and measurement start time
Number of samples	N	
Total value	T	
Maximum value	MAX	
Minimum value	MIN	
Range	RNG	= MAX - MIN
Mean value	MEAN	= T / N
Standard deviation	SD	$\sqrt{\sum(X_i - \text{MEAN})^2 / (N-1)}$
Coefficient of variation	CV	(SD / MEAN x 100)%
Data suffix (footer)		Measurement end time and signature field

Maintenance Parts

Description
Recording Paper
Labeling Paper Rolls
Ink Ribbon
AC Adapter
Connection Cable

Features unique to EP-110

- Supports GLP/GMP using password protection-based date/time alternation prevention
- Enhanced visibility for OLED display
- Powered by rechargeable batteries

Moisture Analyzer

MOC63u

Product

Makes Moisture Content Measurements Quick and Easy

- The moisture ratio is found by heating the sample with the built-in halogen heater to drive out the moisture.
- The measurement procedure is simple. Just close the heater cover to start the measurement (automatic starting mode).
- Measurements are faster than the loss on drying method using a dryer.
- A USB connector makes connecting to a PC easy (built-in Balance Keys function).
- Equipped with the UniBloc aluminum block mass sensor.



MOC63u



Model Name		MOC63u
Capacity	Max. Sample Quantity	60 g
	Min. Sample Quantity	0.02 g
Minimum Display	Mass	0.001 g
	Moisture Ratio	0.01 %
Repeatability *1		0.15 % (2 g), 0.05 % (5 g), 0.02 % (10 g)
Heat Source	Method	Halogen (straight tube)
	Power	Rated at 400 W
Temperature Settings		50 to 200 °C (1 °C interval) (up to 1 hour for settings over 180 °C)
Display		Backlit LCD
Pan Size		Approx. 95 mm dia.
Dimensions		Approx. W202 × D336 × H157 mm
Weight		4.2 kg
Rated Power		430 VA
Ambient Temperature		5 to 40 °C, relative humidity of 85 % max.
Measurement Modes		Standard drying mode (Automatic ending/timed ending) Rapid drying mode (Automatic ending/timed ending) Slow drying mode (Automatic ending/timed ending) Step (3-stage) drying mode (Automatic ending/timed ending)
Time Settings		1 to 240 min, or continuous (up to 12 hours)
External Output		USB Data I/O printer (EP-100/EP-110) output RS-232C (D-sub 9P)
Storage of Measurement Conditions		10 sets
Data Memory		100 items
Standard Accessories		Sample pans (3 aluminum pans), pan holder, windbreak, board, aluminum sheets (50), pan handler, power cable, spare fuses (2), protective display cover, hexagonal wrench

*1 The repeatability (standard deviation) value is from a standard measurement (sample: sodium tartrate dihydrate). This value is not guaranteed for all samples, environments, and measurement conditions.

Optional Accessories

Description
EP-100 Electronic Printer
EP-110 Electronic Printer
Protective Display Cover (5 pcs)
Aluminum Pans (Disposable) (50 pcs)
Fiberglass Sheets (For liquid sample measurements) (100 pcs)
Temperature Calibration Kit
Sample Pan (Stainless steel) (5 pcs)
Sample Pan (Aluminum) (5 pcs)
RS-232C Cable
USB Cable Set
Sample Pan Handler (Stainless steel)
Halogen Heater (For replacement) *2
Power Cable

*2 The halogen heater can be removed and replaced by the user.

Note: For delivery related matters, contact your Shimadzu representative.

Warning

- Use this balance to heat samples to evaporate moisture for measurement.
- The built-in heater will be hotter than the set temperature.
- Samples must not be measured if there is a risk of an explosion or fire, or a dangerous chemical reaction from heating.




This product is certified as Shimadzu's Eco-Products Plus. Energy savings: Power consumption reduced by 30% as compared to a conventional Shimadzu product.

A Total of Five Modes Makes This Balance Compatible with a Variety of Sample Measurements

Ending Modes


Automatic Ending Mode

This automatically ends measurement when the moisture change (% margin) over 30 seconds drops below a set value.



Timed Ending Mode

This automatically ends measurement after a preset amount of time (t1).




The sample is easy to see!
Wide observation window

Alternate Drying Modes


Rapid Drying Mode

The sample is dried at the highest temperature for the initial drying stage, and when the moisture has been reduced, it returns to the set temperature, shortening the measurement time.




Slow Drying Mode

This gently heats samples that might form a surface film or are prone to degrading at high temperatures.



Step Drying Mode


Drying conditions are changed step by step for samples that contain a lot of moisture, such as surface water or crystallization water.



USB INTERFACE

A Wealth of PC Connection Functions

A built-in USB connector allows connecting to a PC. It can be used in conjunction with the Balance Keys function. For USB port connections, check the Shimadzu website or contact your Shimadzu representative.




Compatible with printers (EP-100/EP-110). Equipped with I/O port.

Built-in RS-232C interface is standard.

Equipped with USB interface. Importing data to a PC is easy.

Printout Sample

• Sample Measurement Results Output



Using the EP-100/EP-110

AP Series

Product 

When combined with an optional specific gravity measurement kit, balances can be used to measure specific gravity. Operations are simplified by a text-based navigation function. By using sinkers, the specific gravity of liquid can be measured as well. This allows easily measuring the specific gravity of metals, rubbers, plastics, and other materials.



SMK-601 Specific Gravity Measurement Kit



Model Name	W-AD Series					
	AP225W-AD	AP135W-AD	AP225WD-AD	AP125WD-AD	AP324W-AD	AP224W-AD
Capacity	220 g	135 g	220 g / 102 g	120 g / 52 g	320 g	220 g
Minimum Display	0.01 mg		0.1 mg / 0.01 mg		0.1 mg	
Pan Size	Approx. 91 mm dia.					
Dimensions	Approx. W215 × D411 × H346 mm (incl. power supply unit)				Approx. W215 × D367 × H346 mm	
Weight	Approx. 9.7 kg				Approx. 8.6 kg	

Model Name	W Series						X Series			Y Series			
	AP225W	AP135W	AP125WD	AP225WD	AP124W	AP224W	AP324W	AP124X	AP224X	AP324X	AP124Y	AP224Y	AP324Y
Capacity	220 g	135 g	120 g / 52 g	220 g / 102 g	120 g	220 g	320 g	120 g	220 g	320 g	120 g	220 g	320 g
Minimum Display	0.01 mg		0.1 mg / 0.01 mg					0.1 mg					
Pan Size							Approx. 91 mm dia.						
Dimensions	Approx. W213 × D411 × H345 mm (incl. power supply unit)				Approx. W213 × D367 × H345 mm								
Weight	Approx. 7.9 kg				Approx. 7.0 kg				Approx. 6.5 kg				

Specific Gravity Measurement Kit

Model	SMK-601
-------	---------

Note: The optional liquid density sinker is required for liquid density measurements.

Optional Accessories


Description
Liquid Density Sinker
Petri Dish, Square

P. 37-41 common


What's a Specific Gravity Analyzer?



1 First measure the in-air weight.



2 Then place it in the container filled with water, as instructed on the screen.



3 The specific gravity value is displayed using simple steps.

AU Series

Product 

Attach the optional SMK-401 Specific Gravity Measurement Kit to a balance in the AU series, and set the balance to specific gravity measurement mode. You can then use the balance as a specific gravity analyzer, capable of automatically calculating and displaying specific gravity values.

Liquid density can also be measured by using an optional sinker.

Various balances are available, including a semi-micro (0.01 mg) model. Choose the model best suited to the sample amount and required precision for your application.



AUW Series + SMK-401

Two kinds of weighing pans as standard.



For standard sample



For floating sample

Model Name	AUW-D Series		AUW Series			AUX Series			AUY Series	
	AUW220D	AUW120D	AUW320	AUW220	AUW120	AUX320	AUX220	AUX120	AUY220	AUY120
Capacity	220 g / 82 g	120 g / 42 g	320 g	220 g	120 g	320 g	220 g	120 g	220 g	120 g
Minimum Display	0.1 mg / 0.01 mg		0.1 mg							
Pan Size	Approx. 80 mm dia.									
Dimensions	Approx. W217 × D356 × H338 mm									
Weight	Approx. 7 kg									

Specific Gravity Measurement Kit

Model	SMK-401
-------	---------

Optional Accessories

Description
Liquid Density Sinker
Petri Dish, Square

AT-R Series

Product 

Attach the optional SMK-501 Specific Gravity Measurement Kit to a balance in the AT-R series, and set the balance to specific gravity measurement mode. You can then use the balance as a specific gravity analyzer, capable of automatically calculating and displaying specific gravity values.

Liquid density can also be measured by using an optional sinker.

Various balances are available. Choose the model best suited to the sample amount and required precision for your application.



SMK-501 Specific Gravity Measurement Kit



AT-R Series + SMK-501

Model Name	ATX-R Series				ATY-R Series			
	ATX324R	ATX224R	ATX124R	ATX84R	ATY324R	ATY224R	ATY124R	ATY64R
Capacity	320 g	220 g	120 g	82 g	320 g	220 g	120 g	62 g
Minimum Display	0.1 mg							
Pan Size	Approx. 91 mm dia.							
Dimensions	Approx. W213 × D356 × H338 mm							
Weight	Approx. 6.2 kg				Approx. 6.0 kg			

Specific Gravity Measurement Kit

Model	SMK-501
-------	---------

Optional Accessories

Description
Liquid Density Sinker
Petri Dish, Square

UP Series

Product 

Attach the optional SMK-101A/102 Specific Gravity Measurement Kit to a balance in the UP series, and set the balance to specific gravity measurement mode. You can then use the balance as a specific gravity analyzer, capable of automatically calculating and displaying specific gravity values. Liquid density can also be measured by using an optional sinker.

Various balances are available. Choose the model best suited to the sample amount and required precision for your application.

The large submersible pan makes it easy to measure bulky samples.



UP Series + SMK-101A

Large-pan Models

Model Name						Models with built-in calibration weights				
	UP2202Y	UP4202Y	UP6202Y	UP4201Y	UP8201Y	UP2202X	UP4202X	UP6202X	UP4201X	UP8201X
Capacity	2200 g	4200 g	6200 g	4200 g	8200 g	2200 g	4200 g	6200 g	4200 g	8200 g
Minimum Display	0.01 g			0.1 g		0.01 g			0.1 g	
Pan Size	Approx. 176 × 189 mm					Approx. 176 × 189 mm				

Small-pan Models

Model Name								Models with built-in calibration weights						
	UP223Y	UP423Y	UP623Y	UP823Y	UP1023Y	UP422Y	UP822Y	UP223X	UP423X	UP623X	UP823X	UP1023X	UP422X	UP822X
Capacity		420 g	620 g	820 g	1020 g	420 g	820 g		420 g	620 g	820 g	1020 g	420 g	820 g
Minimum Display		0.001 g				0.01 g			0.001 g				0.01 g	
Pan Size	Approx. 113 × 110 mm							Approx. 113 × 110 mm						

Specific Gravity Measurement Kit

Model	SMK-101A *1*2	SMK-102 *1*3
	for large pan type	for small pan type


Optional Accessories

Description
Liquid Density Sinker for SMK-101A/102

*1 The optional liquid density sinker is required for liquid density measurements.

*2 For UP-X/UP-Y series large-pan (176 × 189 mm) models. The actual capacity is 100 g less than the capacity of the balance.

*3 For UP-X/UP-Y series small-pan (113 × 110 mm) models. The actual capacity is 290 g less than the capacity of the balance. Cannot be attached to the UP223X/UP223Y.

In addition, balances with the  mark are equipped with a specific gravity calculation function, so they can be used for specific gravity measurement.

UW/UX Series

Product 

Attach the optional SMK-101/102 Specific Gravity Measurement Kit to a balance in the UW/UX series, and set the balance to specific gravity measurement mode. You can then use the balance as a specific gravity analyzer, capable of automatically calculating and displaying specific gravity values.

Liquid density can also be measured by using an optional sinker.

Various balances are available. Choose the model best suited to the sample amount and required precision for your application.

The large submersible pan makes it easy to measure bulky samples.



UW/UX Series + SMK-101

UW Series

Model Name	UW220H	UW420H	UW620H	UW820H	UW1020H	UW420S	UW820S	Models with built-in calibration weights				
								UW2200H	UW4200H	UW6200H	UW4200S	UW8200S
Capacity		420 g	620 g	820 g	1020 g	420 g	820 g	2200 g	4200 g	6200 g	4200 g	8200 g
Minimum Display		0.001 g	0.001 g	0.001 g	0.001 g	0.01 g	0.01 g	0.01 g	0.01 g	0.01 g	0.1 g	0.1 g
Pan Size		Approx. 113 × 110 mm						Approx. 176 × 189 mm				

UX Series

Model Name	UX220H	UX420H	UX620H	UX820H	UX1020H	UX420S	UX820S	Models with built-in calibration weights				
								UX2200H	UX4200H	UX6200H	UX4200S	UX8200S
Capacity		420 g	620 g	820 g	1020 g	420 g	820 g	2200 g	4200 g	6200 g	4200 g	8200 g
Minimum Display		0.001 g	0.001 g	0.001 g	0.001 g	0.01 g	0.01 g	0.01 g	0.01 g	0.01 g	0.1 g	0.1 g
Pan Size		Approx. 113 × 110 mm						Approx. 176 × 189 mm				

Specific Gravity Measurement Kit

Model	SMK-101 ^{*1*2} for large pan type	SMK-102 ^{*1*3} for small pan type


Optional Accessories

Description
Liquid Density Sinker for SMK-101/102

*1 The optional liquid density sinker is required for liquid density measurements.

*2 For UW/UX series large-pan (176 × 189 mm) models. The actual capacity is 100 g less than the capacity of the balance.

*3 For UW/UX series small-pan (113 × 110 mm) models. The actual capacity is 290 g less than the capacity of the balance. Cannot be attached to the UW/UX 220H.

In addition, balances with the  mark are equipped with a specific gravity calculation function, so they can be used for specific gravity measurement.

Quick, Stable Measurements of Animal Weight

UP Series

Product 

Attach an optional animal bucket to a UP series balance and set the unit to animal mode. The balance can now be used as a user-friendly animal balance.



Bucket for small animals

Deep round bucket

Rectangular bucket



Clock-CAL

This automatically calibrates the balance at pre-specified times (such as before starting work, during lunch, or after work hours).



Internal Calibration

The balance has built-in motor-driven calibration weights. Sensitivity can be calibrated as needed with a single key press.



Built-in Animal Measurement Mode

- When the animal is unloaded, residual weight from excretions and other materials is automatically subtracted and the display is set to zero. The next animal can be loaded without pressing the TARE button, which increases efficiency.
- Specially developed software allows quick, stable measurement of the weight of moving animals.

UP-X Series																			
UP-Y Series																			

UP-X/UP-Y Series

Model Name	Models with built-in calibration weights					Standard models				
	UP2202X	UP4202X	UP6202X	UP4201X	UP8201X	UP2202Y	UP4202Y	UP6202Y	UP4201Y	UP8201Y
Capacity	2200 g	4200 g	6200 g	4200 g	8200 g	2200 g	4200 g	6200 g	4200 g	8200 g
Minimum Display	0.01 g			0.1 g		0.01 g			0.1 g	
Dimensions	Approx. W191 × D317 × H82 mm			Approx. W191 × D317 × H79 mm		Approx. W191 × D317 × H82 mm			Approx. W191 × D317 × H79 mm	
Weight	Approx. 4.6 kg					Approx. 2.9 kg			Approx. 2.9 kg	

Bucket

Small Animal Bucket	Shape: round / Size: bottom 110 dia. × top 200 dia. × height 130 mm
Deep Round Bucket	Shape: round / Size: bottom 155 dia. × top 195 dia. × height 200 mm
Rectangular Bucket *1	Shape: rectangular / Size: bottom 250 × 210 mm; top 290 × 250 mm; height 150 mm

*1 The rectangular bucket can only be attached to the UP8201X and UP8201Y.
Note: The rectangular bucket can also be used with the UW/UX series.

Three movement levels can be selected corresponding to the animal movement.

Animals can be measured whether they are docile or extremely active.

When the animal is loaded and the stability mark is displayed, the weight is output automatically.

Increases efficiency by eliminating needless operations.

When the animal is unloaded, residual weight from excretions and other materials is automatically subtracted and the display is set to zero.

The next animal can be loaded without pressing the TARE button, which increases efficiency.

Quick, Stable Measurements of Animal Weight

BW-K/BX-K Series

Product 

Attach an optional animal bucket to a BW-K/BX-K series balance and set the unit to animal mode. The balance can now be used as a user-friendly animal balance.



Medium bucket set







Small bucket set



Built-in Animal Measurement Mode

- Models with a range of capacities are available. Ideal for medium weight measurements of rabbits and small dogs.
- Specially developed software allows quick, stable measurement of the weight of moving animals.
- When the animal is unloaded, residual weight from excretions and other materials is automatically subtracted and the display is set to zero. The next animal can be loaded without pressing the TARE button, which increases efficiency.

BW-K Series															
BX-K Series															

BW-K/BX-K Series

Model Name	Models with built-in calibration weights					Standard models				
	BW12KH	BW22KH	BW32KH	BW32KS	BW52KS	BX12KH	BX22KH	BX32KH	BX32KS	BX52KS
Capacity*1 *2	12 kg	22 kg	32 kg		52 kg	12 kg	22 kg	32 kg		52 kg
Minimum Display	0.1 g		1 g			0.1 g		1 g		
Dimensions	Approx. W382 × D366 × H125 mm					Approx. W382 × D366 × H115 mm				
Weight	Approx. 16.5 kg					Approx. 10.5 kg				

*1 When an animal small bucket is attached, the capacity will be reduced about 2 kg from the value indicated.
 *2 When an animal medium bucket is attached, the capacity will be reduced about 6 kg from the value indicated.

Bucket

Small Bucket (mainly suited to rabbits)	Shape: rectangular / Size: bottom 305 × 215 mm; top 335 × 245 mm; height 215 mm
Medium Bucket (mainly suited to small dogs)*3	Shape: rectangular / Size: bottom 335 × 245 mm; top 445 × 295 mm; height 345 mm

*3 The bucket cannot be attached to the BW12KH or BX12KH.

Optional Accessories for the UP-X/UP-Y & UW/UX & BW-K/BX-K

Description
EP-100 Electronic Printer
EP-110 Electronic Printer
RS-232C Cable
USB-Serial Conversion Kit




For Customers at Pharmaceutical Companies — ER/ES Regulatory Compliance — LabSolutions Balance



In recent years, data tampering has caused a decline in the reliability of measurement data. To ensure the reliability of measurement data, as data integrity, it is important to retain not only numeric measurement results, but also other measurement information, such as who measured the data, when, using which instruments, and under what conditions. Information about the operations involved is also important, including information about transcribing measurement values. Such information about measurements is referred to as metadata, such that measurement results are considered reliable (with data integrity ensured) only if they include corresponding metadata. The same applies to data measured using an analytical balance. LabSolutions Balance is software designed for customers that need to ensure the integrity of analytical balance data in the same manner as it pertains to data from LC, GC, and other analytical instruments.

LabSolutions Balance Functionality


- LabSolutions Balance eliminates the need to enter weighing data manually and the risk of transcription errors. All weighing data is saved in a safe database.
- A spreadsheet report of tamper-proof weighing data and analytical data is automatically created.
- Spreadsheet reports can be customized to customer requirements, such as by combining weighing data with HPLC or other analytical results for system suitability tests, content uniformity tests, or elution tests.



Integrated Report Creation Function*3 Combines Analysis Results from HPLC and Weighing Results from a Balance


Creation of Report Template

Enables creating the report by reading sample data and confirming the sample report at the same time.



*1 Metadata refers to information about corresponding data, such as measurement date/time and sample information.
 *2 LabSolutions Balance Ver. 1.0.5 or later
 *3 Multi-data report creation (optional) is necessary to use this function.

Weighing Process Flow Using LabSolutions Balance




Wireless Networking Capability and Tablet Computer Support Enable Convenient Operation in Confined Spaces

Tablet computer compatibility (with wireless networking^{*4}) is convenient for weighing rooms or other locations with limited space.

Weighing data can be transmitted or saved via the wireless network. Of course, it also supports desktop computers.

Key Specifications

OS	Windows® 10 Pro / 11 Pro
Compatible Analytical Balance Models	Shimadzu AP, AU, AT-R, AT, UP, UW/UX, and BW-K/BX-K series
Other Functionality	Simultaneous connection of up to two analytical balances, PDF file creation, and optional LIMS interface supported



Windows® tablet computer (functionality verified using Surface Pro7+)

^{*4} A wireless router and serial device server are required for using wireless networking functionality.

Balance Data Collection Software

Multi-Balance Collect / Balance Keys

Features

Multi-Balance Collect and Balance Keys software run on Windows®. Directly importing balance data into a PC using these free wares eliminates transcription errors, improves work efficiency, and increases data reliability.

Multi-Balance Collect can also collect data continuously at fixed intervals. It is recommended when you want to record changes over time, such as weight changes due to evaporation, on multiple balances.

Accessories for Shimadzu Balances

		Analytical Balances			Electronic Balances			Basic Top-Loading Balances		Precision Platform Balances	Certified Scales and Balances	Moisture Analyzer		Specific Gravity Analyzer				Animal Balances	
		AP ATX-R ATY-R	AU	ATX ATY	UP-X UP-Y UW UX	TW-N TX-N TXB	BL ELB	BW-K BX-K	UW-V	MOC-120H	MOC 63u	AP ATX-R ATY-R	AU UP-X UP-Y UW UX	ELB	UP-X UP-Y UW UX	BW-K BX-K			
Printer	EP-100	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
	EP-110	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
	Printer for MOC-120H AC adapter is included.											✓							
RS-232C Interface	IFB-102A-UNC		✓	✓ ^{*1}	✓	✓	✓ ^{*1}	✓	✓	✓	✓	✓	✓	✓ ^{*1}	✓	✓			
	I/O-RS cable			✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
Application keyboard	AKB-301 For single-point calibration, quantity measurement, threshold setting for pass/fail judgment, and other numerical configurations			✓		✓			✓	✓				✓	✓	✓			
Windbreak	WBC-102 Windbreak Designed for UP-X/UP-Y and UX/UW series models with a maximum capacity of 1020g					✓				✓									
	WBC-502 Large windbreak For UP-X/UP-Y and UX/UW Series					✓				✓									
USB conversion kit with RS-232C cable			✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			

		Analytical Balances				Electronic Balances			Basic Top-Loading Balances		Precision Platform Balances	Certified Scales and Balances	Moisture Analyzer		Specific Gravity Analyzer					Animal Balances			
		AP	AU	ATX-R ATY-R	ATX ATY	UP-X UP-Y	UW UX	TW-N TX-N TXB	BL	ELB	BW-K BX-K	UW-V	MOC-120H	MOC 63u	AP	AU	ATX-R ATY-R	UP-X UP-Y	UW UX	ELB	UP-X UP-Y	UW UX	BW-K BX-K
Foot switch	for TARE	FSB-102PK		✓		✓	✓			✓					✓	✓	✓	✓	✓	✓	✓	✓	
	for print	FSB-102TK		✓		✓	✓			✓					✓	✓	✓	✓	✓	✓	✓	✓	
Specific gravity measurement kit	SMK-101A					✓											✓		✓				
	SMK-101						✓				✓							✓		✓			
	SMK-102						✓	✓				✓						✓	✓	✓	✓	✓	
	SMK-201									✓										✓			
	SMK-401			✓												✓							
	SMK-501				✓												✓						
SMK-601		✓											✓										

*1 To connect to a computer, a separate RS-232C cable is required. *2 USB serial adapter and RS-232C cable for MOC are needed.



IFB-102A-UNC



I/O-RS cable

Optional Accessories List

Balances	Optional accessories
AP W-AD / AP W / AP X / AP Y	STABLO-AP Ionizer Static Electricity Remover
	EP-100 / EP-110 Electronic Printer
	Label Roll Paper for EP-100/110 (10 Rolls)
	Internal Windbreak Plate (for W/X/Y Series) *1
	SMK-601 Specific Gravity Measurement Kit
	AP Holder *2
	Multi-Stand *3
	Shield plate
	AC Adapter (for W/X/Y Series)
	AC Adapter (for W-AD Series Balances)
	AC Adapter (for W-AD Series STABLO-AP Ionizers)
	Display Protective Cover (Set of 5)
	USB Cable Assembly (2 m) with Core
I/O-RS Adapter Cable (for Connecting EP-80/90)	
AUW-D / AUW / AUX / AUY Series	STABLO-AP Ionizer Static Electricity Remover
	EP-100 / EP-110 Electronic Printer
	SMK-401 Specific Gravity Measurement Kit
	AP Holder
	In-use Protective Cover
	RS-232C Cable
	USB Conversion Cable
	AKB-301 Application Keyboard
	FSB-102PK Foot Switch
	FSB-102TK Foot Switch
	AC Adapter *4
ATX-R / ATY-R Series	STABLO-AP Ionizer Static Electricity Remover
	EP-100 / EP-110 Electronic Printer
	USB Cable Set
	SMK-501 Specific Gravity Measurement Kit
	AP Holder
ATX / ATY Series	EP-100 / EP-110 Electronic Printer
	IFB-102A-UNC
	USB Conversion Kit
	In-use Protective Cover (5 pcs)
	I/O-RS Cable
TX / TW / TXB / TXC / TWC Series	EP-100 / EP-110 Electronic Printer
	RS-232C Cable
	In-use Protective Cover (5 pcs)
	In-use Protective Cover for Display
	USB Conversion Kit
UP Series	STABLO-AP Ionizer Static Electricity Remover
	EP-100 / EP-110 Electronic Printer
	RS-232C Cable 25P-9P (1.5 m)
	USB-Serial Replacement Kit
	AKB-301 Application Keyboard
	Windbreak Set *5
	Windbreak Set (for large pan) *6
	WBC-102 Glass Windbreak
	WBC-502 Large Size Windbreak
	In-use Protective Cover (3 pcs) (For Large-Pan Models with 0.01 g Min. Display Value)
	In-use Protective Cover (5 pcs) (For Large-Pan Models with 0.1 g Min. Display Value)
	In-use Protective Cover (5 pcs) (For Small-Pan Model)
	In-use Protective Cover (5 pcs) (For Display and Key Part)
	Animal Bucket Set (for Large-Pan Models with 0.01 g Minimum Display Value)
	Animal Bucket Set (for Large-Pan Models with 0.1 g Minimum Display Value)
	SMK-101A Specific Gravity Measurement Kit (For Large-Pan Model) (Up to approximately 2 kg for weighing capacity)
	SMK-102 Specific Gravity Measurement Kit (For Small-Pan Model)
	AC Adapter *4

*1 Included standard with 0.01 mg models of W series only.

*2 Included standard with AP225W-AD/AP225W models.

*3 Included standard with 0.01 mg models of W-AD series only.

*4 Included standard with main unit.

*5 Included standard with small-pan models with 0.001 g minimum display value.

*6 Included standard with large-pan models with 0.01 g minimum display value.

*7 Std Acc. for models with readability of 1 mg.

*8 Simplified windshield can be used for BL-220H, 320H, 320S. Also, it is included as standard for BL-220H and 320H.

Balances	Optional accessories
UW / UX Series	STABLO-AP Ionizer Static Electricity Remover
	EP-100 / EP-110 Electronic Printer
	RS-232C Interface IFB-102A (for multiple connections)
	Small Size Windbreak (for models with capacity of 220 to 1020 g only) *7
	Glass Windbreak (for models with capacity of 220 to 1020 g only)
	Large Size Windbreak (for all models)
	SMK-101 Specific Gravity Measurement Kit (For large pan 170 × 180 mm)
	SMK-102 Specific Gravity Measurement Kit (For small pan 108 × 105 mm)
	Protective In-use Cover for Key Panel and Display (5 pcs)
	Small Animal Bucket Set (For large-pan models only)
	Angle Adjuster and Wall Hook for Remote Display
	Stand for Remote Display (1 m high)
	FSB-102PK Foot Switch (For printing)
	FSB-102PK Foot Switch (For taring)
	RS-232C Cable, for IBM PC/AT Compatibles (25P-9P, Null modem, 1.5 m)
	RS-232C Cable, for Multiple Connections (25P-25P, Null modem, 1.5 m)
	AKB-301 Application Keyboard
BL Series	EP-100 / EP-110 Electronic Printer
	In-use Protective Cover (5 pcs)
	I/O-RS Cable
	USB Serial Adapter
	Simplified Windshield (Main) *8
ELB Series	EP-100 / EP-110 Electronic Printer
	In-use Protective Cover (5 pcs)
	SMK-201 Specific Gravity Measurement Kit (Cannot be used with small-pan models)
	Carrying Case
	Below-weigh Hook (Except for ELB12K)
BW-K / BX-K Series	EP-100 / EP-110 Electronic Printer
	FSB-102PK Foot Switch (For printing)
	AKB-301 Application Keyboard
	USB-Serial Conversion Kit
	RS-232C Cable
	Below-weigh Hook
MOC63u	EP-100 / EP-110 Electronic Printer
	Protective Display Cover (5 pcs)
	Aluminum Pans (Disposable) (50 pcs)
	Fiberglass Sheets (For liquid sample measurements) (100 pcs)
	Temperature Calibration Kit
	Sample Pan (Stainless steel) (5 pcs)
	Sample Pan (Aluminum) (5 pcs)
	RS-232C Cable
	USB Cable Set
	Sample Pan Handler (Stainless steel)
	Halogen Heater (For replacement) *9
Power Cable	
MOC-120H	Printer (Includes a connection cord and 1 roll of thermal printer paper.) *10
	Printer Paper (10 rolls)
	RS-232C Cable
	Sample Pan
	Aluminum Sheets (500 pcs)
	Temperature Calibration Kit *11

*9 The halogen heater can be removed and replaced by the user. For delivery related matters, contact your Shimadzu representative.

*10 The dedicated printer will be discontinued when the stocks are sold out.

*11 Temperature calibration may be necessary depending on the measurement sample and the measurement conditions. Temperature calibration makes it possible to more accurately control the drying temperature of the measurement sample.


Physical Dimensions

Dimensions (unit: mm) and weight are approximate.
 Appearance and specifications are subject to change without prior notice.


AP Series

Note: For functions, features, and specifications, please refer to pages 12 to 17.


AP W-AD Series (Minimum display 0.01 mg model)




AP W-AD Series (Minimum display 0.1 mg model)



AP W Series (Minimum display 0.01 mg model)



AP W Series (Minimum display 0.1 mg model)
 AP X Series, AP-Y Series




Dimensions (unit: mm) and weight are approximate.
 Appearance and specifications are subject to change without prior notice.


UP and UW/UX Series

Note: For functions, features, and specifications, please refer to pages 24 to 27.


Small-pan model with glass windbreak (optional accessory) (Types with a capacity of 1,020 g or less.)



Large-pan model with large windbreak (optional accessory) (Types with a capacity of 2,200 g or less.)



Small-pan model (Types with a capacity of 1,020 g or less.)




* The 0.1 g minimum display model is 79, while the 0.01 g minimum display model is 82.


* Figure shows combination with simple windbreak (standard only for models with minimum display of 0.001 g). The delivered windbreak may differ slightly in size and shape.

Below-weigh Hook


AP Series



AU Series




UP Series UW/UX Series



STABLO-AP


Note: For functions, features, and specifications, please refer to pages 18, 19.



Dimensions (unit: mm) and weight are approximate.
 Appearance and specifications are subject to change without prior notice.


AU Series

Note: For functions, features, and specifications, please refer to pages 20, 21.




AT-R Series

Note: For functions, features, and specifications, please refer to pages 22, 23.




AP Series + Specific gravity measurement kit

Note: For functions, features, and specifications, please refer to page 37.




AU Series + Specific gravity measurement kit

Note: For functions, features, and specifications, please refer to page 38.




AT-R Series + Specific gravity measurement kit

Note: For functions, features, and specifications, please refer to page 39.



UP Series + Specific gravity measurement kit

Note: For functions, features, and specifications, please refer to page 40.




Dimensions (unit: mm) and weight are approximate.
 Appearance and specifications are subject to change without prior notice.

TX-N Series

For functions, features, and specifications, please refer to pages 28, 29.


Large-pan model
 (Types with a capacity of 2,200 g or less.)



TW-N/TX-N Series


For functions, features, and specifications, please refer to pages 28, 29.

Small-pan model
 (Types with a capacity of 420 g or less.)




TWC-N/TXC-N Series

For functions, features, and specifications, please refer to pages 28, 29.





TXB Series

For functions, features, and specifications, please refer to pages 28, 29.





BL Series

Note: For functions, features, and specifications, please refer to page 30.




ELB Series

Note: For functions, features, and specifications, please refer to page 31.




BW-K/BX-K Series


Note: For functions, features, and specifications, please refer to page 32.

**EP-100/EP-110**


Note: For functions, features, and specifications, please refer to page 33.

**MOC63u**

Note: For functions, features, and specifications, please refer to pages 34, 35.

**MOC-120H**

Note: For functions, features, and specifications, please refer to page 36.



Shimadzu Electronic Balances Demonstration Movies

[Movie](#)


UniBloc, LabSolutions, STABLO and eco logo are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.
Windows, Excel and Surface are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.



Shimadzu Corporation

www.shimadzu.com/an/**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.