

Identification of Phospholipid Molecular Species Using Neutral Loss Survey and MS³ Analysis

Technical Report vol.2



1. Introduction

To elucidate the functions of phospholipids, it is important to conduct analysis not only of their classes and sub-classes, but of their molecular species as well. At Shimadzu, we found that electrospray ionization (ESI) MS³ analysis is effective for more detailed and accurate annotation of each molecular species. We established a system for phospholipid molecular species analysis using a neutral loss (NL) survey experiment of the head-group-related mass values, and succeeding MS³ analyses by selecting the resulting product ions as precursor ions for MS³ analyses (Fig. 1). Using this new method, 34 different molecular species of phosphatidylcholine (PC) were identified even without pre-separation by LC. We then turned our attention to establishing a system for analysis of phosphatidylethanolamine (PE) and phosphatidylserine (PS) in a mixture of lipids.




Fig. 1 System for Analysis of Phospholipid Molecular Species - Neutral Loss (NL) Survey and MS³ Analysis -

In the NL survey + MS³ analysis, information for class and fatty acyl chains were obtained to make detailed analysis of diradyl phospholipids

2. Method

All of the phospholipids were extracted from rat brain (approx. 2 g), liver (approx. 5 g) and calf serum (100 μL) according to Bligh and Dyer's method. The ESI-MS analysis was conducted using a Shimadzu LCMS-IT-TOF (ion trap-time-of-flight mass spectrometer). The extracted phospholipids were directly subjected to ESI MS² and MS³ analysis, using a Si60 column (1 x 100 mm, Nomura Chem., Japan). The mobile phase consisted of acetonitrile/ methanol (spiked with 0.1 % ammonium and 0.3 % acetate).

3. Results

The [M-phosphorylethanolamine] ([M-(Pi-EthN)]⁺) ion corresponding to the [diglyceride-OH]⁺ ion is generated by MS² analysis of the [M+H]⁺ ion of PE in the positive mode (Fig. 2). By conducting MS³ analysis, selecting this [M-(Pi-EthN)]⁺ as the second precursor ion, [fatty acid (FA)-OH]⁺ ions are generated through neutral loss of monoglyceride (MG) moieties, and [MG-H]⁺ ions are generated by neutral loss of FA-related moieties, allowing effective identification of

PE species fatty acyl chains. By conducting MS² of the [M-H]⁻ ion in the negative mode, the [M-serine]⁻ ion, corresponding to the [phosphatidic acid-H]⁻ ion, was generated as a product ion (Fig. 3). By selecting this [M-serine]⁻ ion as the second precursor ion, the fatty acyl chains of the PS species were effectively identified by MS³ analysis, which generated the corresponding [FA-H]⁻ ion.

Sample : diacyl 14:0-14:0 PE (standard)




Fig. 2 Neutral Loss Survey and MS³ Analysis of PE
Two fatty acyl chains of PE were identified using a combination of NL survey of 141 u (C₂H₃NH₂+H) and MS³ of DG (MS² product ion).

Sample : diacyl 16:0-16:0 PS (standard)




Fig. 3 Neutral Loss Survey and MS³ Analysis of PS
Two fatty acyl chains of PS were identified using a combination of NL survey of 87 u [C₂H₃(NH₂)COOH+H] and MS³ of PA (MS² product ion).

Using this new method in conjunction with the MS² method in rapid analysis allowed identification of 7 molecular species of PS (Table 1, Figs. 4, 5). Using MS³ in combination with NL survey allowed highly accurate identification of two fatty acyl chains of the phospholipids. In addition to the NL survey information obtained in MS² analysis, the LCMS-IT-TOF proved very useful in assuring reliable identification of the two fatty acyl chains by providing excellent mass accuracy for the MS³ product ions (Table 2). 132 (one hundred thirty-two) phospholipids (including PC, sphingomyelins,

lysophosphatidylcholine, PE, lysoPE, PS, phosphatidylinositol, phosphatidylglycerol, and triglyceride) were identified in the lipid mixture derived from the rat liver, as summarized in Table 3. The possibility of using this method for quantitative analysis of the metabolome was also investigated. A 1.3-fold increase or 20% decrease in detection was obtained using this method. It is expected that this new method will be effective for comprehensive lipid metabolome analysis.




Fig. 4 Retention Times of Respective Lipid Classes Using a Si Column

	Molecular Species	Theoretical m/z	Difference (Da)	Mass Accuracy (ppm)
1	diacyl 3 8 : 4 PS (1 8 : 0 – 2 0 : 4)	810.5285	0.0017	2.1
2	diacyl 4 0 : 6 PS (1 8 : 0 – 2 2 : 6)	834.5285	0.0020	2.4
3	diacyl 4 0 : 5 PS (1 8 : 0 – 2 2 : 5)	836.5442	0.0034	4.0
4	diacyl 3 6 : 1 PS (1 8 : 0 – 1 8 : 1)	788.5442	0.0019	2.4
5	diacyl 3 6 : 2 PS (1 8 : 0 – 1 8 : 2)	786.5285	0.0069	8.8
6	diacyl 3 6 : 4 PS (1 6 : 0 – 2 0 : 4)	782.4972	0.0037	4.7
7	diacyl 3 8 : 6 PS (1 6 : 0 – 2 2 : 6)	806.4972	0.0015	1.8

Table 1 Identification of PS in Rat Liver




Fig. 5 Mass Spectrum of PS in Rat Liver

Number of possible candidates : m/z of MS

	Possible candidates for peak at m/z 810		theoretical m/z	difference (Da)	mass accuracy (ppm)
1	diacyl 38 : 4 PS	-H	810.5285	0.0017	2.1 ↑10ppm
2	alkacyl 41 : 10 DiMePE	-H	810.5438	0.0170	20.9
3	alkacyl 43 : 10 PE	-H	810.5438	0.0170	20.9 ↑20ppm
4	diacyl 40 : 10 DiMePE	-H	810.5074	0.0194	24.0
5	diacyl 42 : 10 PE	-H	810.5074	0.0194	24.0
6	alkacyl 39 : 4 PS	-H	810.5649	0.0381	47.0
7	alkacyl 40 : 11 PS	-H	810.4710	0.0558	68.9
8	diacyl 39 : 3 DiMePE	-H	810.6013	0.0745	91.9
9	diacyl 41 : 3 PE	-H	810.6013	0.0745	91.9
10	diacyl 39 : 11 PS	-H	810.4346	0.0922	113.8 ↑130ppm
11	alkacyl 40 : 3 DiMePE	-H	810.6377	0.1109	136.8
12	alkacyl 42 : 3 PE	-H	810.6377	0.1109	136.8

Reduction in number of possible candidates

10ppm : 1
20ppm : 3
130ppm : 10

percentage of reduction

10ppm : 100% (=1/1*100)
20ppm : 33% (=1/3*100)
130ppm : 10% (=1/10*100)

Number of possible candidates : m/z of MS³

	Possible candidates for peak at m/z 283		theoretical m/z	difference (Da)	mass accuracy (ppm)
1	acyl 18 : 0 FA	-H	283.2637	0.0000	0.0 ↑10ppm
2	alk 19 : 0 FA	-H	283.3001	0.0364	128.5
3	alk 20 : 7 FA	-H	283.2062	0.0575	203.1
4	acyl 19 : 7 FA	-H	283.1698	0.0939	331.6

Reduction in number of possible candidates

10ppm : 1

	Possible candidates for peak at m/z 303		theoretical m/z	difference (Da)	mass accuracy (ppm)
1	acyl 20 : 4 FA	-H	303.2324	0.0028	9.3 ↑10ppm
2	alk 21 : 4 FA	-H	303.2688	0.0392	129.2

Reduction in number of possible candidates

10ppm : 1

Table 2 Possible Peak Identifications for MS and MS³ Ions

Class of phospholipid	Identification summary	Class of phospholipid	Identification summary	Class of phospholipid	Identification summary
PC	10	LPE	9	PG	11
SM	3	PS	7	LPG	5
LPC	8	PI	2	TG	67
PE	7	LPI	3	Total	132

Table 3 Summary of Identified Phospholipids in a Lipid Mixture from Rat Liver





Fig. 6 A possibility of this method for use in quantitative analysis of a metabolome.

4. Conclusion

1) By selecting the proper conditions for neutral loss scanning of 141 u or 87 u, PE or PS species were identified separately from other phospholipids (Figs. 2, 3).
 2) The new systematic analysis of individual classes of phospholipids by conditional NL survey (MS & MS²) combined with subsequent MS³ has been shown to be a very effective method (Fig. 6). This method will be useful for lipidome (lipid metabolome) analysis.

3) When using the IT-TOF, very high mass accuracy was obtained in MS, MS² and MS³ without using internal standards (Table 2), demonstrating that NL survey in combination with MS³ provides high mass accuracy identification of a set of two FA of phospholipids.



LCMS-IT-TOF LIQUID CHROMATOGRAPH MASS SPECTROMETER

Mayuko Ishida^{1,2}, Shinichi Yamaguchi³, Junichi Taniguchi³, Junko Iida³, Kozo Miseki³, Osamu Nishimura¹, Takao Shimizu², Ryo Taguchi^{2,4}

¹Life Science Laboratory, Life Science Business Unit, Analytical and Measuring Instruments Division, Shimadzu Corporation, Japan; ²Graduate School of Medicine, The University of Tokyo; ³MS/GC Business Unit, Analytical and Measuring Instruments Division, Shimadzu Corporation, Japan; ⁴Crest JST

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. We maintain a global network of sales, service, technical support and applications centers on six continents, and have established long-term relationships with a host of highly trained distributors located in over 100 countries. For information about Shimadzu, and to contact your local office, please visit our Web site at www.shimadzu.com



SHIMADZU CORPORATION. International Marketing Division

3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan

Phone: 81(3)3219-5641 Fax: 81(3)3219-5710

URL <http://www.shimadzu.com>

The contents of this brochure are subject to change without notice.