

High-Speed Video Camera

# HyperVision HPV-X3



# HyperVision HPV-X3

Extreme Speed with High Resolution

## Visualization Technology that Drives Science and Technology

Visualization technology has been responsible for dramatic advances in medical care and the industrial sector. For example, the invention of the microscope allowed humans to observe microscopic objects otherwise too small for the naked eye, and radiography systems and infrared cameras are able to create images from wavelengths of light outside of the visible spectrum. Similarly, high-speed cameras allow humans to capture images of phenomena that are otherwise too quick for human perception. As an established tool in the field of ultra high-speed visualization, the HyperVision series high-speed video camera helps improve our understanding of ultra high-speed phenomena in a variety of fields.

**01 High Image Resolution Across an Impressive Range of Recording Speeds**

— 20 Mfps at 300,000 pixels

**02 External Input/Output Functionality for Synchronized Image Capture**

— Synchronization accurate to 5 nsec

**03 Flexible Design and Superior Ease-of-Use**

— Meets varied application requirements




# High Image Resolution Across an Impressive Range of Recording Speeds


## Burst Method for Ultra High-Speed Recording

High-speed cameras typically store images in memory located separate from the image sensor. Image data is transferred serially from each sensor pixel to the memory in a sequential arrangement through output connections that number relatively few compared to the number of pixels in the image. This configuration makes ultra high-speed recording at 1 Mfps or above hard to achieve.

The burst method used by Shimadzu places enough memory to record the frame capture capacity of the image sensor directly on the sensor and links every pixel to this memory through individual connections. This configuration allows signals to be transferred from the pixels to memory in a fully parallel arrangement for ultra high-speed recording at 20 Mfps. This design removes the limitations imposed by sequential signal transfer via a restricted number of connections, thereby allowing for ultra high-speed, high-resolution image recording.



## Three Times Higher Image Sensor Resolution (300,000 pixels)



The new image sensor has three times the number of pixels and six times the amount of memory thanks to a reduction in memory size. This technology is used to achieve an improved image resolution of 300,000 pixels with no negative effect on

frame rate. The higher image sensor resolution also provides more accurate measurements in applications that use Digital Image Correlation (DIC) technology.


Note: FTCMOS and FTCMOS3 sensors were developed through joint research with Professor Shigetoshi Sugawa of Tohoku University.  
Patents: 04931160, 04844853, 04844854



## External Input/Output Functionality for Synchronized Image Capture


Synchronizing image capture with the subject and an illumination source is an extremely important aspect of high-speed visualization. In addition to having an existing external output feature that synchronizes subject illumination by sending the image capture timing signal to an external illumination device, the HPV-X3 comes with external input and frame synchronization features that allow

256 individual frames to be synchronized to an external signal. With this technology, the camera can begin image exposure in response to a timing signal with an accuracy of 5 nsec. The timing accuracy of the HPV-X3 has also been improved from 10 nsec to 5 nsec. These improvements provide the user with highly reliable synchronized image capture.



## Flexible Design and Superior Ease-of-Use

The HPV-X3 is designed to combine flexibility with excellent ease-of-use. The camera comes with a variety of features that ensure it meets the needs of a wide range of users.




## Applications in a Wide Range of Fields

### Material

### DIC Analysis for High-Speed Tensile Testing of CFRPs


Both static and dynamic material properties, such as impact characteristics, are important for understanding the behavior of materials. Carbon fiber reinforced plastics (CFRPs) exhibit brittle fracture behavior with fracture progression that occurs instantaneously upon damage, and observing this phenomenon requires high-speed video cameras with excellent recording speeds and resolution. The improved resolution of the image sensor in the HPV-X3 improves camera performance for DIC analysis.



Recording speed: 20 Mfps, Test speed: 10 m/s, Specimen width: 12 mm

### HITS-TX Series High-Speed Impact Testing Machine

Utilizes hydraulic controls to test materials at any speed from 0.0001 m/s up to 20 m/s. The HITS-TX series is ideal for high-speed tensile testing of plastics and composite materials.





## Software Development Kit (SDK) for Improved System Development

For seamless integration of the HPV-X3 with commercial DIC software and user-developed analytical software, a software development kit (SDK)<sup>1)</sup> has been published.


1) Purchase of a separate SDK license certification kit is required drop down to next line operate the HPV-X3 with an SDK.

## Camera Synchronization and Dual-Camera Control

The control software can control two cameras. This allows the user to capture images with two cameras simultaneously and then play the recorded images on a single PC.



## Improved Image Sensor Resolution for Improved DIC Analysis Performance

The resolution of the HPV-X3 image sensor is three times that of its predecessor. The resulting improvement in DIC analysis performance was verified by simultaneously capturing images of a single specimen with the new and old sensor and comparing the results. The resulting DIC image clearly shows a build-up of strain in the material just before crack initiation.




## Observing Blast and Shock Waves During Detonation of Micro-Explosives

A silver azide pellet was detonated with a laser and the resulting blast and shock wave propagation was visualized in Schlieren images. The shock wave propagated around the blast wave and its reflection was visualized clearly in an aluminum alloy plate.



Recording speed: 1 Mfps, Field of view width: approx. 250 mm

Images of the area around a silver azide pellet during detonation were captured at 20 Mfps. The images captured a blast wave that appeared approx. 450 ns after laser irradiation of the pellet, which was followed by the progression of a shock wave around the blast wave.





Recording speed: 20 Mfps, Field of view width: approx. 5 mm

Images captured by: Specially Appointed Associate Professor Kiyonobu Otani, Institute of Fluid Science, Tohoku University




## Expansion and Contraction of Bubbles in Polyvinyl Alcohol (PVA) Gel

Images were captured of bubbles being formed while PVA gel was irradiated with a laser. Bubbles were observed to repeatedly expand and contract inside the gel. Images captured by the camera show the progression of shock waves produced when the bubbles were formed and collapsed.




### Bubble Formation



Recording speed: 20 Mfps, Field of view width: approx. 75 mm

### Bubble Collapse






Recording speed: 20 Mfps, Field of view width: approx. 20 mm

Images captured by: Associate Professor Tokitada Hashimoto, Department of Mechanical Engineering, Faculty of Science and Engineering, Saga University

## Observing the High-Frequency Oscillation of Microbubbles

Images were captured of microbubbles that formed in water when the water was irradiated and heated locally with a laser. The microbubbles first expanded and then contracted. The images show a jet flow that occurs during contraction as the bubble disappears.






Recording speed: 20 Mfps, Field of view width: approx. 110  $\mu\text{m}$

Images captured by: Associate Professor Kyoko Namura, Department of Micro Engineering, Graduate School of Engineering, Kyoto University

## Observing Crack Progression During Ring-on-Ring Testing of Glass

Ring-on-ring testing was performed on reinforced glass and images were captured of the cracks that occurred during failure. (Reference standard: ASTM C1499)



Recording speed: 10 Mfps, Field of view width: approx. 45 mm

# Specifications

P/N: S348-00021-41  
Model Name: HPV-X3

Camera Head		
Lens Mount <sup>1)</sup>	Nikon F-mount	
Image Sensor <sup>2)</sup>	FTCMOS3 image sensor (approx. 30 × 23 mm)	
Pixel Size	48 μm × 48 μm	
Recording Speed <sup>3)</sup> (frame rate)	20 Mfps Variable recording speed between 60 fps and 10 Mfps in 5 ns steps	
Recording Capacity	256 frames	
Resolution <sup>4)</sup>	300,000 pixels; 628 horizontal x 480 vertical	
Color/Gradations <sup>5)</sup>	Monochrome, 10 bits	
Exposure Time <sup>6)</sup>	20 Mfps, fixed at approx. 25 ns Variable in a 5 ns interval starting from 50 ns in a range from 60 fps to 10 Mfps	
External Trigger Input	Two channels (TRIG, STANDBY) TTL level (5 V), capable of either positive or negative polarity Or contact input Time resolution 5 ns	
Recording Mode	Internal trigger, external trigger, continuous trigger	
Camera Synchronization Function <sup>7)</sup>	Capable of synchronized recording with 2 cameras connected	
Frame Synchronization Function <sup>7)</sup>	Signal/level	One channel (F.SYNC) TTL level (5 V), capable of either positive or negative polarity Or contact input Time resolution 5 ns
	Recording speed	Range between 60 fps and 10 Mfps
Optional Outputs	Two channels (exposure start timing, trigger detection timing, or other outputs depending on settings)	
Trigger Point Setting	Can be set to any frame from the second frame onwards	
Interface <sup>8)</sup>	1000Base-T 1 port	
External Monitor Output <sup>9)</sup>	HDMI output	
Data Memory Format	10-bit dedicated format, BMP, AVI, JPEG, TIFF (8-bit and 16-bit formats supported)	
Power Supply Unit		
Power Rating	Single phase 100-240 V AC ± 10 %, 200 VA, 50/60 Hz	
Required Specifications for the Control PC		
OS	Windows® 11 Pro 64-bit (Ver. 22H2 or later)	
CPU	Intel® Core™ i5-1245U Max. 4.2 GHz or faster	
Memory	8 GB or more	
HDD	500 GB or more (SSD is recommended)	
Screen Size	1920 × 1080 dots or more (15.6-inch wide or larger is recommended)	
Interface	1000Base-T	
Optic Drive	Super multi drive or higher	
Other Peripherals	Mouse and keyboard	
Environmental Conditions		
Operating Temperature Range	5° to 40° C	
Operating Humidity Range	35 % to 75 % RH with no condensation	
Storage Temperature Range	0° to 50° C	
Storage Humidity Range	20 % to 80 % RH with no condensation	
Size and Weight		
Camera Head	WHD 200 × 270 × 402 mm, approx. 9.2 kg	
Power Supply Unit	WHD 154 × 199 × 398 mm, approx. 6.0 kg	
Length of Interface Cable between Camera and Control PC	Approx. 2 m	
Length of Cable between Camera and Power Supply Unit	Approx. 2.8 m	

Note 1) Shimadzu does not guarantee that all F-mount lenses can be attached.

Note 2) The FTCMOS3 image sensor used in this instrument is manufactured using high-accuracy technology, but defective pixels may exist. Note that this is not a defect or failure of the product.

Please note that some image sensor attributes are not published. When purchasing an HPV-X3, please request an image capture demonstration to verify product operation in your operating environment.

Note 3) The recording speed is a reference value. It is not guaranteed to be an accurate value for the time interval between recording frames.

Note 4) Stored images will be 628 pixels (horizontal) × 480 pixels (vertical).

Note 5) 10-bit refers to the data format. It does not indicate a guarantee of data precision.

Note 6) These exposure times are rough indications and are not guaranteed as exact exposure time ratios for all recording speeds.


Note 7) The synchronization function value is for reference only. It is not guaranteed to be an accurate value for the time interval between recording frames.

Note 8) For Ethernet, only 1 Gbps is supported. 100 Mbps/10 Mbps does not work properly.


Note 9) The output signal is 640 (horizontal) × 480 (vertical) in VGA.

## External Dimensions

### Camera Head



### Power Supply Unit



#### Export Control Regulations

Any export of a Shimadzu High-Speed Camera, HPV-X3/HPV-X2/HPV-X/HPV-2, is subject to export control regulations of the nation, based on Part 2 of the NSG guideline, 5.B.3. Please contact a sales agent or representative of Shimadzu should you have any questions.

Intel and intel Core are trademarks of Intel Corporation or its subsidiaries.

Windows is either a registered trademark or a trademark of Microsoft Corporation in the United States and/or other countries.



Shimadzu Corporation

[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

#### For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.