

Exploring the application of a universal method for pesticide screening in foods using a high data acquisition speed MS/MS

ISCC 2012 P. 54

Anja Grüning¹, Susanne Kräher¹, Ichiro Hirano²,
Yuka Fujito³, Kiyomi Arakawa², Yusuke Inohana²

¹Shimadzu Europa GmbH, Albert-Hahn-Straße 6-10,
Germany

²Shimadzu Corporation, 1, Nishinokyo-Kuwabaracho,
Nakagyo-ku, Japan

³Shimadzu Analytical & Measuring Center Inc.,
1, Shimoaicho, Nishinokyo, Nakagyo-ku, Japan

Exploring the application of a universal method for pesticide screening in foods using a high data acquisition speed MS/MS


Introduction

Effective management, use, and disposal of agrochemicals, particularly pesticides, is an increasingly important health and environment issue in developing countries where economies may be heavily reliant on agriculture. The conventional approach to monitor these pesticides is to develop highly optimized triple quadrupole MRM methods to achieve the required levels of sensitivity, selectivity and

speed of analysis whilst still providing confidence in pesticide identification. In this study LCMS technology, developed for ultra-fast scanning MRM analysis, allows the possibility of a single generic 'universal' method. High speed MRM analysis and a generic parameters were used for screening 172 pesticides (344 MRM transitions) with 5 msec dwell and 1 msec pause times in food matrices.

Materials and Methods

Sample Preparation



Analytical Conditions

HPLC : Nexera UHPLC system	
Column	: Shim-pack XR-ODSII (75 mm x 2 mmI.D., 2.2 um)
Mobile phase	: A ; 5 mM ammonium acetate – water B ; 5 mM ammonium acetate – methanol
Gradient program	: 30% B (0 min.) → 80% B (4 min.) → 95% B (10-15 min.) → 30% B (15.01-20 min.)
Flow rate	: 0.2 mL / min.
Column temperature	: 40°C
MS : LCMS-8040 Triple quadrupole mass spectrometer	
Ionization	: ESI (Positive / Negative)
Ion spray voltage	: +4.5 kV / -3.5 kV
MRM	: 344 MRM transitions (2 MRMs / compound)

Features of LCMS-8040

- **5 times** higher sensitivity compared to LCMS-8030
- An ultra fast scan speed of **15000 u / sec.**
- An ultra fast polarity switching of **15 msec.**
- An ultra fast MRM transition speed of **555 ch./ sec.**



Fig. 1 LCMS-8040 Triple Quadrupole Mass Spectrometer

Setting of MRM analysis & integration parameter

- In this study, no scheduling of MRM transitions was applied; thereby creating a universal generic method.

The figure shows two software interface screenshots:

- Instrument Parameters View (Realtime analysis)**: Shows a list of 177 compounds. A red box highlights 167 compounds from the 'Residual pesticide method package' and 5 additional compounds. The 'Positive' and 'Negative' ionization modes are indicated.
- MS Data Processing View (Postrun analysis)**: Shows a peak list with 'Largest Peak' selected for identification. Parameters include Dwell 5 msec, Pause 1 msec, and Loop time 2.058 sec.

Additional annotations include: 'Acquisition time set to the entire chromatographic run (0-20 min.)' and 'Largest peak' selected for identification.

Exploring the application of a universal method for pesticide screening in foods using a high data acquisition speed MS/MS

Results

Screening of 10 pesticides in food matrices




Fig. 2 MRM chromatogram of 10 pesticides

Table 2 Result of 10 pesticides screening (10 ppb spiked in each matrices)

Compounds	Leek	Paprika	Green tea leaves
Carbofuran	✓	✓	✓
Chlorfluazuron	✓	✓	✓
Fosthiazate	✓	✓	✓
Hexythiazox	✓	✓	✓
Indanofan	✓	✓	✓
Lufenuron	✓	✓	✓
Mevinphos	✓	✓	✓
Propoxur	✓	✓	✓
Pyrimidifen	✓	✓	✓
Tricyclazole	✓	✓	✓
False positives*	8	7	10

- All peaks were automatically selected as the target compound to permit automatic identification of target analytes without retention time information. (* Number of false positives out of 172 screened pesticides.)

Exploring the application of a universal method for pesticide screening in foods using a high data acquisition speed MS/MS




Fig. 3 Result of automatic identification (10 ppb spiked in the green tea leaves)

Conclusion

- Pesticides spiked in all matrices at 10 ppb (10 compounds) could be automatically detected using fast 5 msec MRM with 15 msec polarity switching without retention time information.



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.
The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.