

# From Instrument to Column Tracking Down the Problem

Melissa Goodlad  
Applications Engineer  
Columns and Supplies Technical Support



# Troubleshooting Topics

## System pressure

- Increased pressure
- Low pressure
- Pressure fluctuations

## Separation

- Changing retention time
- Loss of resolution

## Peak shape

- Tailing
- Broadening
- Fronting
- Peak splitting and doubling

## Detection

- Noisy baseline
- Reduced intensity or sensitivity
- Drifting baseline

# Agilent Lab Advisor

- Tools for calibration, diagnosis, and maintenance
- Daily instrument tests
- General calibration and maintenance procedures
- Advanced version also available for expert level troubleshooting
- EMF (Early Maintenance Feedback) shows the number of valve switches or pumped solvent

**Agilent Lab Advisor** System Overview Advanced Version

System Name	System Information	EMF	Status
1100	Agilent LC / 134.40.26.21	EMF ?	Offline <input type="button" value="Connect"/>
1290 Training	Agilent LC/CE / 192.168.254.11	EMF ?	Offline <input type="button" value="Connect"/>
Afraco	Agilent LC / 0030D32F48C4.gemany.agilent.com	EMF ?	Offline <input type="button" value="Connect"/>
App - Prime	Agilent LC / 10.68.8.30	EMF ?	Offline <input type="button" value="Connect"/>
<b>App-Bio</b>	<b>Agilent LC / 10.68.8.36</b>	<b>EMF</b>	<b>Not Ready <input type="button" value="Disconnect"/></b>

Component	Serial #	Firmware	Details	EMF
<b>G5654A 1260 Bio-inert Pump</b>	DEAGH00169	D.07.22 [0001]	- Seal Wash Pump - Degasser - AIV - LAN Settings (IP: 192.168.254.11, SM: 255.255.255.0, GW: not specified, MAC: 0030D330FF46, Init Mode: ...)	EMF ?
<b>G5668A 1260 Bio Multisampler</b>	DEAGJ00123	D.07.20 [0007]	- Sample Thermostat (Product# 20448, Serial# DEBAT07529) - Multi-wash - Multisampler Parameter Right (Bio Needle Seat 0.17mm, Bio Sample Loop Flex 100uL right, Bio Analytical He... - Injection Valve (5067-4263 - "Bio Injection Valve 600bar")	EMF ?
<b>G7116A 1260 MCT</b>	DEAEM00622	C.07.20 [0002]	- Hosted by 'G5668A:DEAGJ00123' with Firmware Revision 'D.07.20 [0007]' - Left Column Tag Reader - Column Selection Valve - Valve Head (10 Ports, 2 Positions, Product# 5067-4132, Serial# 0003064327, Max Pressure 600 bar)	EMF ? - Thermo off
<b>G7115A 1260 DAD WR</b>	DEAC600102	D.07.20 [0007]	- Flow Cell (Product# G5615-60022, Serial# DE311R1091, Path Length 10.00 mm, Volume 13.00 u) - UV Lamp (Product# 2140-0820, Serial# F93997) - LAN Settings (IP: 10.68.8.36 'ap-07', SM: 255.255.254.0, GW: 10.68.8.1, MAC: 0030D32F3303, Init Mode: ... - Current LAN Controller: 10.68.9.102 'AL-27', 10.68.9.102 'AL-27', 10.68.10.162 'CND64912BC'	EMF ? - UV lamp not ready


  

Bin Pump	Agilent LC / Localhost	EMF ?	Offline <input type="button" value="Connect"/>
CompactLC	Agilent LC / localhost	EMF ?	Offline <input type="button" value="Connect"/>

Connection Address: 10.68.8.36 Version B.02.11 [292] - Advanced | Licenses 54/35 ..


# System Pressure


# Understand Your LC System and Follow the Flow Path



# Changes in System Pressure


## Increased Pressure/Overpressure and Blockages

Potential cause		Recommended action
●	Clogging of filter frits in the high-pressure flow path	<ul style="list-style-type: none"><li>• Identify the culprit by logical elimination process and replace affected part.</li><li>• Use clean, prefiltered solvent</li><li>• Prevent algae growth in water</li></ul> 
●	Plugging of capillaries, needles, and needle seats	
●	Wrong solvent	<ul style="list-style-type: none"><li>• Check for correct mobile phase</li><li>• <b>Check</b> solvent reservoir and tube connections</li></ul>




# Blockages and Clogging

Characteristics	
Parts affected	<p>Blockages:</p> <ul style="list-style-type: none"> <li>• Capillaries, needle, and needle seat</li> <li>• Detector flow cells</li> </ul> <p>Clogging:</p> <ul style="list-style-type: none"> <li>• Filter frits (inline filter, column filter)</li> </ul>
Characteristic	●
Identification	<ul style="list-style-type: none"> <li>• Check easily accessible points: needle seat, purge valve, column</li> <li>• Disconnect capillaries one-by-one, starting at detector and moving back toward the pump</li> </ul>
Possible root cause	<ul style="list-style-type: none"> <li>• Debris from mechanically worn parts (needle seat material, rotor seal at injection valve)</li> <li>• Coring of vial septa material</li> </ul>
Instant action/First aid	<ul style="list-style-type: none"> <li>• Replace part</li> <li>• Backflush affected part</li> </ul>
Preventive measures	<ul style="list-style-type: none"> <li>• Proper preventive maintenance schedules, replace worn parts regularly</li> <li>• Use high quality septa</li> <li>• Install inline filters</li> </ul>





**Blockages:** Instant pressure increase step



**Clogging:** Constant pressure increase over time

# Checking for Blocked Needle Seat





# Checking for Blocked Needle Seat

Track the typical operating pressure for a given application

To troubleshoot:

- Create high pressure on the system by turning on flow
- While the pressure is climbing, move the sampler to the “Bypass” position
- Watch the pressure when the valve switches to “Bypass”
- If the pressure drops immediately, then the source of the high pressure is in the portion of the flow path specific to the bypass position
  - Needle Seat
    - Where the sample first meets the mobile phase
    - **Most commonly clogged piece of tubing in an LC**
  - Needle
    - Less commonly clogged
    - Watch for issues with septa

The screenshot displays the 'Control' panel of an Agilent instrument. Under the 'Method Parameters' section, the 'Injection valve (Single Needle)' is set to 'Mainpass' (indicated by a selected radio button). Below this, the 'Needle Wash' section is expanded to show 'Set Needle Wash Multi Mode Parameters'. The parameters are: Solvent 1 (dropdown), Channel: A (dropdown), Time [s]: 30 (text input), Seat back flush: Off (dropdown), Needle wash at flush port: On (dropdown), and Solvent Name: (text input). The 'Configuration' and 'Special Commands' sections are partially visible at the bottom.

# Checking for Blocked Needle Seat

## Loop


- Not commonly clogged
- Watch for issues with sample

## Metering head

- Never exposed to sample
- Consider solvent issues

## Injection valve

- Most common issue is with rotor seal
- Look for scratches on stator face




# Checking for Blocked Needle Seat

Driving the needle into debris may result in a clogged needle or seat

- In well plate samplers, bottom sensing gives the most consistent position
- But this isn't recommended if there's debris in the bottom of the sample vial
- For vial samplers, a zero offset is approximately 2 mm from the bottom of a 2 mL vial


Needle height position without bottom sensing:  
G1367E/G4226A 54 vial tray = 4 mm  
G1367E/G4226A 100 vial tray = 2.5 mm  
G7167X 54 vial tray = 5 mm



# Locating a Clog

If the pressure doesn't drop when the valve switches to "Bypass", the issue is likely outside the sampler.

- Purge valves:
  - With 1260 Infinity II model pumps, open the manual purge valve. The pressure should drop to between 0 and 5 bar. **If the pressure is higher than this, the PTFE filter may be clogged.**
  - With 1290 Infinity II model pumps, purging is done through an automated valve, activated using software. 1290 binary pumps have the same PTFE filter, 1290 quaternary pumps have a 5 µm filter frit.



**G7120A 1290 High Speed Pump**

Serial #	DEA0000000
Firmware:	B.07.20 [0007]

▼ **Controls**

▼ **Control**

Pump:  On  Off  Standby  Initializing

Purge + Prime

Purge Process:  On  Off

Prime:  On  Off

# Locating a Clog



PTFE replacement on a 1260 pump:

1. Remove pump outlet and purge waste tubing
2. Unscrew the purge valve using a 14 mm wrench
3. Remove the gold seal cap
4. Remove the frit
5. Install the new frit, slot side up
6. Replace the gold seal cap
7. Reinstall the valve

Re-align the waste tubing in the correct orientation during installation.




# Locating a Clog

If the pressure doesn't drop when the valve switches to "Bypass", the issue is likely outside the sampler.

## Column:

- Open the fitting at the inlet of the column.
- Pumping 1 mL/min of water through an Agilent LC with 0.17 mm id tubing typically shows a pressure of 40 bar.
- If the pressure is much higher than this, a capillary may be clogged. If the pressure appears "normal" the issue may be with the column.




# Locating a Clog

If the pressure doesn't drop when the valve switches to "Bypass", the issue is likely outside the sampler.

Working backwards from detectors:


- Clogs are located by opening a fitting, typically at most a half turn.
- If the pressure drops, the clog is downstream from the fitting or towards the detector. If pressure remains high, the clog is upstream or towards the pump.



# Changes in System Pressure

## Low pressure



Potential cause	Recommended action
● Leak in high-pressure flow path	<ul style="list-style-type: none"><li>• Visual inspection of flow path</li><li>• Instrument diagnostic tests</li></ul>
● Wrong mobile phase	<ul style="list-style-type: none"><li>• Check for correct mobile phase</li><li>• <b>Check</b> solvent reservoir and tube connections</li></ul>





# Leaks



Characteristics	
Parts affected	<ul style="list-style-type: none"><li>• Potentially all parts in the flow path</li><li>• High potential at frequently operated fitting connections (such as the column inlet) and parts with high mechanical stress (rotor seal, needle, and needle seat)</li></ul>
Characteristic	<ul style="list-style-type: none"><li>• Lower pressure</li><li>• Potentially impacting retention times and peak shape</li></ul>
Identification	<ul style="list-style-type: none"><li>• Drops of solvent or residues of salt</li><li>• System diagnostic tests</li></ul>
Possible root cause	<ul style="list-style-type: none"><li>• Loose or bad fitting connections</li><li>• Cracked capillaries</li><li>• Worn needle and needle seat</li></ul>
Instant action/first aid	<ul style="list-style-type: none"><li>• Replace affected parts</li><li>• Renew or redo fitting connection</li></ul>
Preventive measures	<ul style="list-style-type: none"><li>• Use proper fitting connections</li><li>• Replace fittings and wear parts in time</li></ul>




If dimension X is too long, leaks will occur

# How Do I Locate a Leak?

- Each Agilent LC module is equipped with a leak sensor
- If liquid is detected, the entire LC stack will shut down
- The LC will not start up again until the sensor has been dried and returned to temperature



# Overtightened Fittings



# Changes in System Pressure


## Pressure fluctuations

Potential Cause		Recommended Action
●	Air in the system	<ul style="list-style-type: none"><li>• Prime and flush instrument</li><li>• Check for sufficient solvent supply</li><li>• Check for correct plumbing (SSV/MCGV)</li><li>• Check for correct degassing</li></ul>
●	Malfunctions at pump head	<ul style="list-style-type: none"><li>• Perform pump head diagnostic tests LA</li><li>• Replace defective parts</li><li>• Implement proper maintenance schedule</li></ul>
●	Cavitation effects	<ul style="list-style-type: none"><li>• Check for flow restrictions (solvent bottle to pump head inlet)</li><li>• Clean or replace parts</li><li>• Verify that solvent supply is positioned above pump inlet</li></ul>

### Important to know



Pressure fluctuations typically also impact the UV signal due to refractive index effects.




# Peak Shape




# Changes in Peak Shape

## Peak tailing

If applicable to some peaks	Recommended Action
Secondary interactions	<ul style="list-style-type: none"> <li>• Check pH of mobile phase (most likely)</li> </ul>
Small peak eluting on tail of larger peak	<ul style="list-style-type: none"> <li>• Pump malfunction</li> </ul>



If applicable to all peaks	Recommended Action
Poor tubing connections; high dispersion volume	<ul style="list-style-type: none"> <li>• Minimize number of connections</li> <li>• Check connections/fitting condition and proper seat of fittings</li> <li>• Use fittings with spring-loaded function</li> </ul>
Column damage	<ul style="list-style-type: none"> <li>• Use specialty, polymeric or sterically protected column</li> <li>• Column cleaning</li> </ul>



# InfinityLab Quick Connect and Quick Turn Fittings


- Spring-loaded design
- Easy to use
- Works for all column types
- Reusable
- Consistent ZDV connection

## Quick Connect Fitting

- Finger tight up to 1300 bar
- Hand tighten the nut, then depress the lever

## Quick Turn Fitting


- Finger tight up to 400 bar
- Up to 1300 bar with a wrench
- Compact design




# Peak Tailing: Column Contamination

Column: StableBond SB-C8, 4.6 x 250 mm, 5 mm      Mobile phase: 20% H2O : 80% MeOH      Flow rate: 1.0 mL/min  
Temperature: R.T.      Detection: UV 254 nm      Sample: 1. Uracil    2. Phenol    3. 4-Chloronitrobenzene    4. Toluene


QC test forward direction



QC test reverse direction



QC test after cleaning  
100% IPA, 35 °C






# Changes in Peak Shape

## Peak splitting/doubling



Potential Cause	Recommended Action
Partially plugged column frit	<ul style="list-style-type: none"><li>• Backflush column (if applicable)</li><li>• Use inline filter</li><li>• Use guard column</li></ul>
Column void	<ul style="list-style-type: none"><li>• Replace column</li><li>• Use guard column</li><li>• Use less aggressive mobile phase conditions</li></ul>
Sample volume overload	<ul style="list-style-type: none"><li>• Use smaller injection volume</li></ul>
Sample solvent incompatibility with mobile phase	<ul style="list-style-type: none"><li>• Use mobile phase or weaker miscible solvent as injection solvent</li></ul>
Issues with injection valve	<ul style="list-style-type: none"><li>• Check injector valve parts</li><li>• Replace worn parts</li></ul>



# Changes in Peak Shape

## Fronting


Potential Cause	Recommended Action
Channeling in column	<ul style="list-style-type: none"> <li>• Replace column</li> <li>• Use guard columns</li> </ul>
Column overload	<ul style="list-style-type: none"> <li>• Use higher capacity column (increase length, diameter, or change to high-capacity material)</li> <li>• Decrease sample amount</li> </ul>




# Changes in Peak Shape

## Peak broadening


Potential Cause	Recommended Action
Injection volume too large	<ul style="list-style-type: none"><li>Decrease injection volume</li></ul>
Long retention times	<ul style="list-style-type: none"><li>Use gradient elution or stronger mobile phase</li></ul>
System settings	<ul style="list-style-type: none"><li>Check data collection rate</li><li>Adjust the detector setting or time constant to the fastest possible value without compromising signal-to-noise.</li></ul>
Viscosity of mobile phase too high	<ul style="list-style-type: none"><li>Increase column temperature</li></ul>
Detector cell volume too large	<ul style="list-style-type: none"><li>Use smallest possible cell volume</li></ul>
Improper fittings and connections	<ul style="list-style-type: none"><li>Ensure that your fitting connections are correct</li></ul>
Extra tubing volume on system	<ul style="list-style-type: none"><li>Ensure that the tubing is narrow and as short as possible to avoid extra volume</li></ul>
Sample diluent too strong	<ul style="list-style-type: none"><li>Reduce diluent strength</li></ul>




# Strong Diluents Can Disrupt Equilibration – Isocratic Method



# Strong Diluents Can Disrupt Equilibration – Gradient Analysis



# Comparison of Peak Shape at Low and High Loads Broadening and Tailing




High sample loads give broad or broad and tailing peaks

Dextromethorphan is 35% broader at high load

0.005 mg/mL dextromethorphan (4.1 uL injection volume)

Low sample loads provide symmetrical, nontailing peaks with narrow peak widths




# Changes in Separation

# Changes in Separation


## Retention time changing

Potential Cause	Recommended Action
Inconsistent online mobile phase mixing	Ensure gradient system <b>is</b> delivering constant composition check vs. manual <b>preparation</b> of mobile phase
Flow rate changing	Check 'pressure fluctuation'
Column temperature varying	Thermostat column and ensure constant lab temperature
Equilibration time insufficient with gradient run or change in isocratic mobile phase	Flush with at least 10 column volumes after solvent change or gradient conclusion
Selective evaporation of mobile phase component	Keep solvent reservoirs covered prepare fresh mobile phase
Buffer capacity insufficient	Use > 20 mM concentration of buffer
Contamination buildup	Occasionally flush column with strong solvent to remove contaminants
First few injections – adsorption on active sites	Condition column by initial injection of concentrated sample
Column overloaded with sample	Decrease injection volume or concentration
Mobile phase composition changing	Follow 'best practices'






# Mobile Phase Preparation




Agilent 1100 with quaternary pump  
ZORBAX Eclipse XDB-C8 Rapid-Resolution (3.5 $\mu$ m), 4.6 x 50 mm  
Agilent Part No. 935967-906  
Dial-a-Mix= A: water B: MeOH, pump 50% B  
Vol/Vol=250 mL water + 250 mL MeOH, pump 100%  
Mix-to-Mark = 250 mL MeOH, fill to 500 mL with water, pump 100%  
Premixed (w/w) = 200 g MeOH + 200 g water, pump 100%  
UV 254 nm  
1 mL/ min.

# Retention Time Shift – Selectivity Differences Due to Incorrect pH

pH 4.5 shows selectivity change from lot-to-lot for basic compounds




**pH 3.0 - Lot 1**




pH 3.0 shows no selectivity change from lot-to-lot

**pH 4.5 - Lot 2**



**pH 3.0 - Lot 2**




- For method ruggedness
  - Test three different column lots
  - Compare  $R_s$  for the three lots
    - If  $\Delta R_s$  is too large, modify method

# Changes in Separation

## Ghost peaks, carry over

	Potential Cause	Recommended Action
●	Peaks from previous injection	<ul style="list-style-type: none"><li>• Flush column to remove contaminants</li><li>• Check with blank injection</li></ul>
●	Specific interaction with metal surfaces	<ul style="list-style-type: none"><li>• Passivate instrument</li><li>• Use InfinityLab Deactivator Additive</li><li>• Use bio-inert LC equipment</li></ul>
●	Contamination or unknown interferences in samples	<ul style="list-style-type: none"><li>• Proper sample cleanup</li></ul>

**BIO  
INERT**



# Mobile Phase Hygiene


Contaminated mobile phases can cause:

- Lower sensitivity
- Rising/drifted baselines
- Higher noise
- Ghost peaks with gradient separations

Often the issue is confused with autosampler carryover.

It can be identified by repeating the gradient run without sample injection or increasing the pre-run equilibration.


Always run multiple blanks before standards or samples to distinguish gradient artifacts from possible carryover.




# Mobile Phase Hygiene: Glassware

## Improper cleaning of solvent bottles can cause contamination of mobile phases and result in gradient artifacts

- Wash solvent bottles with hot water, deionized water, and organic solvent (IPA or acetonitrile).
- Leave glassware inverted on paper towels on a bench or clean pegboard dowels to dry.
- **Avoid using detergents.** If it is necessary to use detergents to get glassware clean, rewash with plenty of hot and cold water so that all detergent residues are removed. Follow with deionized water and organic (IPA or acetonitrile) rinses.
- Store glassware inverted on shelves or in drawers, or cover openings.



System blank injection,  
water/ACN gradient on a C18 column,  
PEG contamination



# Mobile Phase Hygiene: Solvent Purity and Buffer Preparation

- Use HPLC grade organic mobile phases
- Use HPLC grade water or Milli-Q DI water
- Use HPLC grade reagents, including salts, ion pair reagents, and base and acid modifiers
- Always rinse pH electrodes thoroughly when measuring/adjusting **the** pH of **the** mobile phase
- Prepare fresh buffers to avoid contaminants from the growth of bacteria or algae
- Filter your mobile phase buffer with 0.45  $\mu\text{m}$  filter before use
- Solvent filters installed at the end of solvent lines should be replaced periodically





100% ACN

90%ACN+10% buffer  
(10 mM phosphate)

# Autosampler Carryover


## Common sources:

- Exterior of needle (use needle wash)
- Worn needle seat
- Worn rotor seal
- Poorly made fitting





# Autosampler Carryover






# Changes in Detection

# Changes in Detection

## Noisy baseline

	Potential Cause	Recommended Action
●	Gas bubbles in mobile phase	<ul style="list-style-type: none"> <li>• Apply degassing</li> <li>• Check degasser performance</li> </ul>
●	Low difference between sample and mobile phase absorbance	<ul style="list-style-type: none"> <li>• Check absorbance values of sample vs. mobile phase</li> </ul>
●	Contamination	<ul style="list-style-type: none"> <li>• Use degassed HPLC-grade solvents</li> <li>• Flush system</li> <li>• Clean up the sample</li> </ul>
●	Detector optics	<ul style="list-style-type: none"> <li>• Perform intensity test</li> <li>• Check signal with flow cell removed if possible</li> <li>• Replace lamp</li> </ul>
	Pressure instability	<ul style="list-style-type: none"> <li>• Check 'pressure fluctuation'</li> </ul>



# UV Lamp Tests

How do I know if my UV lamp is good?

- Visual inspection of an equilibrated baseline
- Accumulated UV lamp on-time from RFI tag or Lab Advisor
- Lab Advisor intensity test
- Lab Advisor ASTM drift and noise test
- Lab Advisor cell test


The screenshot shows the Lab Advisor software interface with a table of instrument counters. The table has columns for Title, Value, Unit, Limit, and Progress. The progress column includes a progress bar and a percentage. The table is filtered to show 'All Counters'.

		Title	Value	Unit	Limit	Progress		
 <b>G4220A 1290 Bin Pump</b> Serial # DEBAA00157			0	Hour	3000	0%	★	↺
			3.45	Liter	50	6%	★	↺
			8.64	Liter	50	17%	★	↺
			631	Count	15000	4%	★	↺
		Liquimeter (A+B)	12.09	Liter	0	0%	★	↺
 <b>G4226A 1290 ALS</b> Serial # DE93000560			0	Count	1000	0%	★	↺
		Needle into seat counter	1191	Count	1500	79%	★	↺
			1.53	Hour	3000	0%	★	↺
		Valve switching counter	2418	Count	60000	4%	★	↺
 <b>G4212A 1290 DAD</b> Serial # DEBAF00163		Accumulated UV lamp on-time	2519.65	Hour	2000	100%	★	↺
		UV lamp ignition counter	28	Count	1500	1%	★	↺
		UV lamp on-time	360.65	Hour	0	0%	★	↺
 <b>G4208A 1200 Instant Pilot</b> Serial # PP55055002								

# UV Lamp Tests

Diode array and multiple wavelength

Counters and hours

		Title	Value	Unit	Limit	Progress	
G7117B							
	<b>G7117B</b> Serial #	<b>1290 DAD</b> PPBAW00058	Accumulated UV Lamp On-Time	3.17	h	0	0%
			Number of UV Lamp Ignitions	2	Count	0	0%

The useable lifetime of a deuterium lamp will depend on it's use:

- How many hours has it been on?
- How many times has it been ignited?
- What wavelength is being used?

# UV Lamp Tests

Diode array and multiple wavelength

Intensity test

Intensity Test 1260 HPLC » 1260 HPLC » G7115A:DEAC600377

General Limits Signals

<b>Test Name</b>	Intensity Test	<b>Description</b>	The test scans the Intensity spectrum generated by the UV and VIS Lamp.
<b>Module</b>	G7115A:DEAC600377 (1260 DAD WR)		
<b>Status</b>	<b>Passed</b>		
<b>Start Time</b>	2/21/2019 3:21:31 PM		
<b>Stop Time</b>	2/21/2019 3:22:15 PM		

Test Procedure

- ✓ 1. Check Prerequisites...
- ✓ 2. Remove Flow Cell.
- ✓ 3. Scan Intensity Spectrum...
- ✓ 4. Evaluate Data...


Result

Name	Value
Accumulated UV Lamp Burn Time	126.35 h
UV Lamp On-Time	0.02 h
Accumulated Vis Lamp Burn Time	262.45 h
Vis Lamp On-Time	0.02 h
Lowest Intensity in Range 190 - 220 nm	36289 Counts
Lowest Intensity in Range 190 - 220 nm	2000 Counts
Lowest Intensity in Range 221 - 350 nm	21963 Counts
Lowest Intensity in Range 221 - 350 nm	5000 Counts
Lowest Intensity in Range 351 - 500 nm	16150 Counts
Lowest Intensity in Range 351 - 500 nm	2000 Counts
Lowest Intensity in Range 501 - 950 nm	13102 Counts
Lowest Intensity in Range 501 - 950 nm	2000 Counts
Highest Intensity in Range 190 - 350 nm	81934 Counts
Highest Intensity in Range 190 - 350 nm	2000 Counts
Highest Intensity in Range 700 - 950 nm	62919 Counts
Highest Intensity in Range 700 - 950 nm	2000 Counts
Highest Intensity for D2 Alpha Line (600 - 700 nm)	157676 Counts
Highest Intensity for D2 Alpha Line	2000 Counts
Spectrum Integral	31863112
UV Integral (190 - 349 nm)	7384053

# UV Lamp Tests

Diode array and multiple wavelength

Intensity test



The profile of the intensity scan changes as a lamp ages

# UV Lamp Tests

Diode array and multiple wavelength

ASTM drift and noise

ASTM Drift and Noise Test 1260 HPLC » 1260 HPLC » G7115A:DEAC600377

General Limits Signals

**Test Name** ASTM Drift and Noise Test      **Description** The test performs ASTM drift and noise evaluation without reference.  
**Module** G7115A:DEAC600377 (1260 DAD WR)  
**Status** **Passed**  
**Start Time** 2/21/2019 4:56:26 PM  
**Stop Time** 2/21/2019 5:16:47 PM


Test Procedure

- ✓ 1. Check Prerequisites...
- ✓ 2. Remove Flow Cell.
- ✓ 3. Measure Noise...
- ✓ 4. Evaluate Data...

Result

Name	Value
Accumulated UV Lamp Burn Time	127.94 h
UV Lamp On-Time	1.60 h
Minimum Lamp On-Time	1 h
Accumulated Vis Lamp Burn Time	264.03 h
Vis Lamp On-Time	1.61 h
Minimum Lamp On-Time	1 h
Signal Drift value at 254 nm (UV)	-0.138 mAU/h
Maximum Allowed Drift	-1 ... 1 mAU/h
Signal Noise value at 254 nm (UV)	0.008 mAU
Maximum Allowed Noise	0 ... 0.02 mAU
Signal Drift value at 750 nm (Vis)	-0.415 mAU/h
Maximum Allowed Drift	-1 ... 1 mAU/h
Signal Noise value at 750 nm (Vis)	0.013 mAU
Maximum Allowed Noise	0 ... 0.02 mAU

ASTM Drift and Noise Test



Remove the flow cell from the lightpath of the detector.

OK Cancel

# UV Lamp Tests

Diode array and multiple wavelength


ASTM drift and noise

Run on a monthly basis, this test can help track the natural decline of the lamp and perhaps raise awareness of a dirty cell.

Name	Value	Description
Minimum Lamp On-Time	1 h	The minimum lamp on-time to perform a noise check.

Name	Lower limit	Upper limit	Description
Maximum Allowed Noise	0 mAU	0.02 mAU	The maximum allowed Signal noise in mAU.
Maximum Allowed Drift	-1 mAU/h	1 mAU/h	The maximum allowed Signal drift in mAU.





# UV Lamp Tests

Diode array and multiple wavelength

Cell test

Cell Test 1260 HPLC » 1260 HPLC » G7115A:DEAC600377

General Limits Signals

**Test Name** Cell Test **Description** The test compares the lamp intensity with and without the flow cell installed. The intensity ratio is an indicator of the amount of light absorbed by the flow cell.

**Module** G7115A:DEAC600377 (1260 DAD WR)

**Status** **Passed**

**Start Time** 2/21/2019 3:57:24 PM

**Stop Time** 2/21/2019 3:59:37 PM

Test Procedure

- ✓ 1. Check Prerequisites...
- ✓ 2. Remove Flow Cell.
- ✓ 3. Scan Intensity Spectrum...
- ✓ 4. Insert Flow Cell.
- ✓ 5. Scan Intensity Spectrum...
- ✓ 6. Evaluate Data...

Result


Name	Value
Accumulated UV Lamp Bum Time	126.95 h
UV Lamp On-Time	0.62 h
Accumulated Vis Lamp Bum Time	263.05 h
Vis Lamp On-Time	0.62 h
Intensity Integral without Flow Cell	32,088,720
Intensity Integral with Flow Cell	19,830,098
Intensity Ratio	0.62
Minimum Intensity Ratio	0.3

Diode array detectors with the fiber optic style flow cell require a Max Light test cell for this test (part number G4212-60011).

# UV Lamp Tests

Diode array and multiple wavelength

Cell test




Example of scans with and without the cell installed.

# How Do I Know if My UV Lamp is Good?

Diode array and  
multiple wavelength

Baseline inspection



Various factors contribute to the specific amplitude and pattern of baseline noise and drift, including the specific wavelength, mobile phase, room temperature, and data rate.

# How Do I Know if My UV Lamp Is Good?

Diode array and multiple wavelength

## Baseline inspection

### Long cycle wave

This is a rhythmic change in the baseline where the periodicity may be hours.

- Environmental influences

### Short cycle wave

This is a rhythmic change in the baseline where the periodicity may be seconds or minutes.

- Solvent mixing noise
- Mechanical issue in pump

If the cycle of the wave does not appear to be mixing noise, evaluate the health of the lamps through Lab Advisor intensity, noise and drift tests. Also, the cleanliness of the flow cell **can be evaluated** through the cell test.


## Excessive drift

In a UV baseline, light scattering shows up as drift. If the baseline is drifting more than expected, empty and rinse the solvent bottles, refilling with fresh solvent. Perform a cell test to check the cleanliness of the flow cell.

# Useful Parts

Parts that address potential issues and help to ease your daily tasks

Part Description	Information	Part number
InfinityLab Stay Safe caps	Prevents solvent evaporation; changes in mobile phase concentration	Various <a href="http://www.agilent.com/chem/staysafecaps">www.agilent.com/chem/staysafecaps</a>
InfinityLab Quick Connect and Quick Turn fittings	With <b>spring-loaded</b> function for optimized dead volume reduction	Various <a href="http://www.agilent.com/chem/infinitylabfittings">www.agilent.com/chem/infinitylabfittings</a>
Blank nut, long, 10-32	Blank nut, PEEK with steel core; for system diagnostic tests; finger tight up to 1300 bar, easy to use and gentle on receiving port	5043-0277
Agilent Captiva syringe filters	Solve issues like inlet clogging, increased backpressure, and retention time shift by filtering your samples	Various <a href="http://www.agilent.com/chem/filtration">www.agilent.com/chem/filtration</a>
InfinityLab Poroshell 120 columns	High efficiency and high resolution; available in 18 chemistries	Various <a href="http://www.agilent.com/chem/discoverporoshell">www.agilent.com/chem/discoverporoshell</a>



InfinityLab Poroshell 120 columns



InfinityLab Stay Safe cap on solvent bottle



InfinityLab Quick Connect fitting



InfinityLab Quick Turn fitting



Blank nut, 5043-0277

# LC Troubleshooting Poster Available

## LC Troubleshooting Guide

Your guide to solving common problems and staying productive

Agilent  
InfinityLab

### Places to Start

#### Solvents

- Use brown borosilicate bottles to avoid algae growth
- Prepare solvent volume to be used up within 1 to 2 days
- Use only HPLC-grade solvents filtered through 0.2 µm filters

#### Preparing and powering up the pump

- Inspect solvent bottles and inlet filters for damage or coloring
- Always use seal wash when installed and purge the pump
- Use the appropriate system conditioning method

#### Daily tasks

- Replace aqueous and organic mobile phases every second day
- Check seal wash solvent
- Flush the system with the composition of your application

#### Weekly tasks


- Change seal wash solvent and bottle and inspect solvent filters
- Check system backpressure and change filters if necessary

#### Pump shutdown

- Flush all channels to remove salt deposits and particulate matter
- Flush the system with appropriate storage solvent and power down the system

#### Handling of acetonitrile

- If possible, use 5 to 10% of water in your mobile phase
- Be sure to avoid ACN evaporation
- Don't leave ACN on the system for more than 2 to 3 days
- Perform a periodic warm water wash (60 to 70 °C) if you face problems



### Maintenance

Agilent Lab Advisor software helps you manage your Agilent LC instruments to achieve high-quality chromatographic results in the most efficient way by ensuring high instrument performance, productivity, and reliability. It is available free-of-charge.

- Diagnostic tests to evaluate performance
- Easier maintenance of all Agilent LC modules
- Comprehensive reports generated to ease communication with Agilent service

#### Retention Time Drift



Possible Cause	Solution
Inconsistent online mobile phase mixing	Ensure gradient system delivers constant composition; compare with manual preparation of mobile phase
Variation in column temperature	Thermostat or insulate column; ensure constant lab temperature
Insufficient equilibration time with gradient run or change in isocratic mobile phase	Make sure at least 10 column volumes pass through column after sample run
Selective evaporation of mobile phase component	Less vigorous helium sparging; keep solvent reservoirs covered; prepare fresh mobile phase
Contamination buildup	Occasionally flush column with strong solvent
Column overloaded with sample	Decrease injection volume or concentration

#### Pressure Fluctuation



Possible Cause	Solution
Leak in the system	Identify the channel and clean or replace check valve; replace pump seals
Buildup of particulates	Filter sample and mobile phase
Bubble in pump	Perform solvent degassing; sparge solvent with helium

#### Pressure Increase



Possible Cause	Solution
System blockage	Check flowpath (needle seat, capillaries, filter and fits)
Water/organic systems: buffer precipitation	Test buffer-organic mixtures to ensure compatibility
High Column Backpressure	<ul style="list-style-type: none"> <li>Column blockage: Buffer sample cleanup; use guard column</li> <li>Mobile phase viscosity too high: Use lower viscosity solvents or higher temperature</li> <li>Particle size too small: Use larger d<sub>p</sub> packing</li> <li>Plugged inlet fit: Replace column</li> </ul>

#### Drifting Baseline



Possible Cause	Solution
Positive/negative direction: contaminant buildup/dilution	Flush column; clean up sample; use pure solvents
Positive/negative: difference in refractive index of injection solvent	Use mobile phase for sample solvent
Temperature changes	Insulate and thermostat column and tubing

#### Noisy Baseline



Possible Cause	Solution
Contamination	Use degassed HPLC-grade solvents; flush system; clean up sample
Detector problems	Check number of hours of UV lamp; replace UV lamp or flow cell

#### Ghost Peaks



Possible Cause	Solution
Peaks from previous injection	Flush column to remove contaminants; check with blank injection
Contamination; unknown interferences in samples	Proper sample cleanup
Ion pair: disequilibrium	Prepare sample in actual mobile phase to minimize disturbance
Contaminated mobile phase	Check your mobile phase
Bubbles in solvent	Check and degas your solvents

#### Peak Tailing



Possible Cause	Solution
Unwashed dead volumes	Minimize number of connections; ensure injector seal is tight; ensure fittings are properly seated
Column performance	Change mobile phase; replace column
Silica-based: column degradation	Use specialty, polymeric, or sterically protected column
Silica-based: basic interactions with stationary phase	Use stronger mobile phase or add appropriate base (e.g., TEA)

#### Peak Broadening



Possible Cause	Solution
Injection volume too large	Decrease injection volume or solvent strength of injection solvent; use gradient methods
Low sampling rate of data system	Increase data rate
Detector cell volume too large	Use smallest possible cell volume
Injection volume too large	Decrease injection volume

#### Sensitivity Problems



Possible Cause	Solution
Peaks are outside of sensitivity range of detector	Dilute/concentrate sample to bring into linear region
Sample-related losses during preparation	Use internal standard during sample preparation; optimize sample preparation method

#### Leaks



Possible Cause	Solution
White powder at fittings/ loose fitting	Tighten fittings; replace capillaries
System leak	Identify location checking leak sensors/sensors; check flow cell

Discover more best practices for using an Agilent LC system:  
<https://www.agilent.com/chem/lc-best-practices>



Training courses are available at:  
<https://www.agilent.com/crosslab/university>



Get answers. Share insights. Join the Agilent Community at:  
<https://community.agilent.com>



For Lab Advisor software, please visit:  
<https://www.agilent.com/chem/lab-advisor>



Request yours today at  
[www.agilent.com/chem/troubleshootLC](https://www.agilent.com/chem/troubleshootLC)

# Resources for Support

- **New!** HPLC Advisor App: [HPLC Advisor app | Agilent](#)
- LC Troubleshooting poster: 5994-0709EN
- Resource page: <http://www.agilent.com/chem/agilentresources>
  - Quick reference guides
  - Catalogs, column user guides
  - Online selection tools, how-to videos
- InfinityLab Supplies catalog: [5991-8031EN](#)
- LC handbook: [5990-7595EN](#)
- YouTube – [Agilent channel](#) (maintenance videos)
- Agilent service contracts



AGILENT INFINITYLAB SUPPLIES FOR THE ANALYZER 1200 INFINITY HPLC WITH MULTISAMPLE  
Quick reference guide

Agilent supplies for Agilent instruments

Agilent Technologies is committed to ensuring your laboratory productivity, so we have produced this list of the most commonly ordered supplies and parts for the 1200 Infinity HPLC with Multisample.

This list of items is an abstract profile of LC materials, columns, and supplies designed to work together seamlessly to improve efficiency and performance.

Part Number	Description
0101-0001	Agilent Technologies
0101-0002	Agilent Technologies
0101-0003	Agilent Technologies
0101-0004	Agilent Technologies
0101-0005	Agilent Technologies
0101-0006	Agilent Technologies
0101-0007	Agilent Technologies
0101-0008	Agilent Technologies
0101-0009	Agilent Technologies
0101-0010	Agilent Technologies
0101-0011	Agilent Technologies
0101-0012	Agilent Technologies
0101-0013	Agilent Technologies
0101-0014	Agilent Technologies
0101-0015	Agilent Technologies
0101-0016	Agilent Technologies
0101-0017	Agilent Technologies
0101-0018	Agilent Technologies
0101-0019	Agilent Technologies
0101-0020	Agilent Technologies
0101-0021	Agilent Technologies
0101-0022	Agilent Technologies
0101-0023	Agilent Technologies
0101-0024	Agilent Technologies
0101-0025	Agilent Technologies
0101-0026	Agilent Technologies
0101-0027	Agilent Technologies
0101-0028	Agilent Technologies
0101-0029	Agilent Technologies
0101-0030	Agilent Technologies
0101-0031	Agilent Technologies
0101-0032	Agilent Technologies
0101-0033	Agilent Technologies
0101-0034	Agilent Technologies
0101-0035	Agilent Technologies
0101-0036	Agilent Technologies
0101-0037	Agilent Technologies
0101-0038	Agilent Technologies
0101-0039	Agilent Technologies
0101-0040	Agilent Technologies
0101-0041	Agilent Technologies
0101-0042	Agilent Technologies
0101-0043	Agilent Technologies
0101-0044	Agilent Technologies
0101-0045	Agilent Technologies
0101-0046	Agilent Technologies
0101-0047	Agilent Technologies
0101-0048	Agilent Technologies
0101-0049	Agilent Technologies
0101-0050	Agilent Technologies

For more information, visit: [www.agilent.com/chem/aten](http://www.agilent.com/chem/aten)

Agilent Technologies





# Contact Agilent Chemistries and Supplies Technical Support



1-800-227-9770 option 3, option 3:

Option 1 for GC and GC/MS columns and supplies

Option 2 for LC and LC/MS columns and supplies

Option 3 for sample preparation, filtration and QuEChERS

Option 4 for spectroscopy supplies

Option 5 for chemical standards

**Available in the U.S. and Canada 8-5 all time zones**

[gc-column-support@agilent.com](mailto:gc-column-support@agilent.com)

[lc-column-support@agilent.com](mailto:lc-column-support@agilent.com)

[spp-support@agilent.com](mailto:spp-support@agilent.com)

[spectro-supplies-support@agilent.com](mailto:spectro-supplies-support@agilent.com)

[chem-standards-support@agilent.com](mailto:chem-standards-support@agilent.com)