



### Concentrations of NDMA in malt and beer

Traditionally malt kilns (Figure 2) were direct fired i.e. the products of combustion from the fuel used to heat the air for drying passed directly through the grain bed. NDMA formation was therefore likely to have been a normal part of the process of malting, but its presence was probably not known about and so not a cause for concern. However, with increasing awareness of potentially harmful substances in the food chain and improvements in analytical equipment and techniques, NDMA in malt was identified as a problem. Although palliative treatments like reducing the pH at the surface of the grain during kilning by burning sulphur in the air stream offered a short term solution, indirect firing was seen as the ultimate answer. In indirect firing the hot air carrying the products of combustion from the heating fuel passes to exhaust through a heat exchanger which heats incoming, clean, ambient air to dry the malt. Major investment in indirect fired kilns (Figure 3) in malting companies around the world was set in train and now the ability to obtain NDMA levels in malt of around 1 µg/kg is commonplace; previously some levels had been recorded in the low hundreds. In the few remaining malting plants where direct fired kilns are still in operation the installation of low NO<sub>x</sub> burners and/or palliative sulphur burning in the early stages of kilning provides almost equivalent control of NDMA formation, thus enabling maltsters to continue to meet the very tight specifications set by brewers and distillers.

### Control of NDMA

Very low levels of NDMA in malt are now being consistently achieved and these levels are diluted a further 10-fold during the brewing of beer. Nevertheless, there remains a duty of care to all those who drink malt based beverages to ensure that the potentially harmful effects of NDMA are kept to an absolute minimum. To that end the monitoring of

NDMA levels in malt continues to be necessary.

### Analysis of NDMA at Crisp Malting Group

Previously a Thermal Energy Analyser was used to determine levels of NDMA in malt samples, but replacement with GCMS-QP2010 in early 2004 has proved very successful. Since that time hundreds of samples have been routinely analysed without problems and excellent correlation with laboratories using Thermal Energy Analysers has been achieved. Figure 4 shows the chromatogram of a malt sample with a concentration of 0.75 ppb (µg/kg). Details of the method are given below.

### NDMA in malt: method details

Malted barley is blended with water in a homogeniser/blender and subsequently filtered through

a Whatman No. 54 filter paper. NDPA internal standard and sodium chloride are added to the filtrate. After liquid/liquid extraction with DCM (dichloromethane) the DCM phase is dried with anhydrous sodium sulphate, filtered and evaporated in a water bath at 55 °C to about 1 mL. The cooled DCM phase is transferred to a GCMS vial and ready for injection.

### GC-MS conditions:

Column: Stabilwax, 30 metres, 0.25 mm ID, 0.5 µm film thickness  
Oven temperature: 50 °C for 1 min, then 20 °C/min to 175 °C holding at 175 °C for 1 min, then 40 °C/min to 250 °C holding at 250 °C for 10 min.

### Injection port temperature:

240 °C

### Interface temperature:

250 °C

### Ion source temperature:

200 °C

### References:

- [1] Nitrosamines in Malt and Beer; Wainwright, T., J. Inst. Brew., 1986, Vol 92, pp 73-80
- [2] The Chemistry of Nitrosamine Formation: Relevance to Malting and Brewing; Wainwright, T., J. Inst. Brew., 1986, Vol. 92, pp 49-64




Figure 4: Chromatogram of NDMA