

Fast screening of pesticides in
foods and agricultural products
with probe electrospray ionization
(PESI) tandem mass spectrometry

ASMS 2019 ThP-043

Zhenhe Chen, Satoshi Yamaki, Jing Dong, Yuki Hashi,
Naoki Hamada
Shimadzu China Mass Center,
Shimadzu China,
Beijing, 100020.

Introduction

Pesticides screening in food and agricultural products have been intensively emphasized due to the increasing food safety requirement. Recently, high resolution mass has been thought to a promising screening method, but it is very costly. Pesticides usually have broad polar range and both GCMS and LCMS methods should be used for their determination. Probe electro spray ionization (PESI) source, as an ambient ionization technique, which could be used for the determination of both polar and semi-polar compounds. Meanwhile, as an in situ ionization, PESI

source usually needs less sample pretreatments. In this study, 80 pesticides are selected as research targets, which are commonly monitored in vegetables every year in China. A PESI-MS/MS approach is developed aiming at both high throughput and quantification using MRM acquisition.




Figure 1 Sample pretreatment flow and schematic diagram of PESI

Material and Methods

Commercially available pesticides were used for this experiment. Standards of pesticides were diluted with water / isopropanol = 1/1 to an appropriate concentration and then analysed by DPiMS-8060. The water / isopropanol (1/1) solution consisted of 2 mM ammonium acetate buffer and 0.2% formic acid. The commercially available pesticides were used for the confirmation of MS spectrum and the optimization of MRM conditions.

Automated MRM parameter optimization with both APCI and ESI source were carried out by flow injection analysis

of authentic standards with a function of the LabSolutions LCMS control software. The DPiMS-8060 consisted of PESI source (Shimadzu) and a triple quadrupole mass spectrometer (LCMS-8060, Shimadzu).

Fast screening of pesticides in foods and agricultural products with probe electrospray ionization (PESI) tandem mass spectrometry

For sample pretreatment, 1.5 kg of fresh cucumber was pulverized with a fruit and vegetable processing machine (AUX-20B, Aux) for 3 min, and stored in 50 ml polyethylene centrifuge tube for use. 10 (± 0.01) g sample were put in a 50 ml centrifuge tube, 5ml ultrapure water and 10 ml of acetonitrile (5% acetic acid) were added and then vortexed for 3min. Afterwards, 4g Mg_2SO_4 and 1g NaCl were added and vortexed for another 3min. Finally, the sample was centrifuged at 4200 rpm for 5 min and supernatant was stored for use.

For sample analysis, 10 μ l of the supernatant was loaded in the sample plate and after evaporation at room temperature (ca.10mins), 10 μ l of the water / isopropanol was added in the same sample plate.

Table1 Instruments Parameters for DPiMS-8060

PESI	LCMS-8060
Ionization Outage time: 200 msec	Desolvent line: 200 °C
Take sample time: 50 msec	Heat block: 30 °C
Discharge voltage: 2.3 kV(+)/-3.0 kV(-)	Sample amount: 10 μ L
Cleaning time: 0.05 min(+)/0.05 min(-)	MRM method: 0.075 min/5 pesticides
Cycle time of probe: 2.78 Hz	Switch time: 1 msec
Ionization time: 250 msec	Residence time: 2 msec





Figure 2 Injection mode for DPiMS-8060

Results

Ionization confirmation and sensitivity check

Automated MRM parameter optimization with ESI or APCI source were carried out by flow injection analysis with an embedded function of Labsolutions. Afterwards, ionization confirmation of 80 pesticides was applied using DPiMS-8060. Fortunately, among the 80 pesticides, 74 compounds can be confirmed by DPiMS-8060 in

water/isopropanol (1/1) solution. The unconfirmed compounds include dicofol, fugiclor and hexacyclohexane (4 isomers), which are usually detected with GCMS.




Figure 3 Typical mass chromatograms of isocarbofos, carbofuran, and forchlorfenuron (1ppm for ionization confirmation and 1ppb for sensitivity check)

Single pesticide standard solution was prepared in 1, 10, and 100ppb for sensitivity check on DPiMS-8060. Among the 74 compounds, 27 compounds are traditionally detected with GCMS method and 47 compounds are with LCMS methods.

Generally, over 90% pesticides can be detected at lower than 10ppb, which demonstrated the high sensitivity of DPiMS-8060 (Table 2). 1ppb responses of typical pesticides are shown in Figure 3.

Fast screening of pesticides in foods and agricultural products with probe electrospray ionization (PESI) tandem mass spectrometry

Table 2 Linearity of typical pesticides

Traditional Method	Total	1 ppb	10 ppb	100 ppb	>100 ppb
GCMS	27	12	8	3	4
LCMS	47	45	2	0	0

Quantitative Analysis of real samples

Blank matrix was obtained the same way as normal samples. 1ppm stock solution was diluted to 5, 10, 20, 50, and 100ppb with water / isopropanol. The calibration curves were treated the same way as samples but standard solutions were added. Among the 74 compounds, the linear correlation coefficient of 36 compounds are over 0.98 (Table 3). The mass chromatogram of cucumber blank matrix with 100ppb standard are show in Figure 4. DPIMS-8060 can test 74 pesticides within 1.5 min, which is

suitable for fast screening of pesticides in vegetables. For most compounds traditionally detected with GCMS, they show poor or no linearity under PESI detection. This is probably due to the matrix component strongly affect the ionization of these compounds. Especially for pyrethroids, their precursors are usually ammonium adducts. Under current pretreatment procedure, none of them could be detected.

Table 3 Linearity of typical pesticides

No.	Name	Equations	R ²
1	Methamidophos	$y=10271x-3973$	0.997
2	Acephate	$y=242.1x+121.7$	0.988
3	Omethoate	$y=1324x-1757$	0.995
5	Sulfone aldicarb	$y=597.4x+1027.3$	0.997
7	Thiamethoxam	$y=436.4x+1740.4$	0.998
8	Imidacloprid	$y=347.2x-57.9$	0.998
9	Carbofuran-3-hydroxy	$y=1371x+1863$	0.998
10	Dimethoate	$y=503.7x+2466$	0.996
11	Acetamiprid	$y=335.8x+282.2$	0.998
13	Aldicarb	$y=271.5x+17.6$	0.992
14	Carbofuran	$y=545.7x+3104$	0.984
15	Carbaryl	$y=78.14x-327.39$	0.999
16	Phosemet	$y=63.11x+218.85$	0.995
17	Azoxystrobin	$y=858.1x+6042.2$	0.999
18	Malathion	$y=99.25x+10389.78$	0.99
19	Dimethomorph	$y=787.7x+5140.4$	0.992
20	Triadimefon	$y=58.36x+2903.33$	0.989
21	Triazophos	$y=794.9x-1005.5$	0.997
24	Diazinon	$y=12809x+3795$	0.999
26	Phoxim	$y=72.14x+3293.62$	0.992
29	Difenoconazole	$y=214.5x+3198.5$	0.998
30	Profenofox	$Y=146.2x+1198.6$	0.99

Fast screening of pesticides in foods and agricultural products with probe electrospray ionization (PESI) tandem mass spectrometry

31	Tridemorph	$y=911.2x+4730.8$	0.998
32	Chlorpyrifos	$y=47.29x+2052.55$	0.998
35	Pyridaben	$y=487.3x+14122.1$	0.982
37	Cyromazine	$y=830.1x+100260.9$	0.999
40	Pacloubtrazol	$y=188.1x+587.0$	0.997
43	Chlorantraniliprole	$y=94.29x+306.3$	0.997
45	Phorate	$y=43.58x+343.69$	0.993
46	Isofenphos-methyl	$y=164.7x+868.7$	0.992
60	Pyrimethanil	$y=1535x+7138$	0.999
64	Etofenprox	$y=133.8x+6467.5$	0.998
68	Fipronil	$y=708.3x+13232.6$	0.986
69	Fipronil Desulfanyl	$y=369.8x+21422.8$	0.987
70	Fipronil sulfide	$y=1083x+7228$	0.993
71	Fipronil Sulfone	$y=1562x+13836$	0.996




Figure 4 Typical mass chromatograms for 74 pesticides (100ppb)

10 µg/kg and 100 µg/kg spike recovery experiments were performed in cucumber. Table 4 shows the spike recovery rate for some pesticides, which demonstrate the DPiMS-8060 can be used for the fast

screening analysis. For the other pesticides, the recovery rate is not within the range of 60-120%. Another easy and simple extraction method should be developed for them.

Fast screening of pesticides in foods and agricultural products with probe electrospray ionization (PESI) tandem mass spectrometry

Table 4 Pesticides spike recovery in cucumber

No.	Name	10(µg/kg) spike	100(µg/kg) spike	MRL(µg/kg)*
1	Methamidophos	122%	66.6%	50
2	Acephate	85.4%	74.5%	300
3	Omethoate	99.5%	67.4%	20
5	Sulfone aldicarb	85.6%	85.0%	20
7	Thiamethoxam	85.0%	113%	500
8	Imidacloprid	70.6%	111%	200
9	Carbofuran-3-hydroxy	73.0%	93.0%	20
10	Dimethoate	67.1%	116%	200
14	Carbofuran	89.9%	124%	20
15	Carbaryl	93.7%	102%	1000
25	Emamectine	66.0%	86.8%	20
29	Difenoconazole	79.8%	130%	200
40	Paclobutrazol	63.9%	123%	50
43	Chlorantraniliprole	95.4%	147%	20
60	Pyrimethanil	69.1%	127%	2000
64	Etofenprox	124%	81.7%	500

Conclusion

- On DPiMS-8060, the analytical work flow is presented as high-throughput and robust measurement for pesticides in agriculture products.
- MRM methods of 74 pesticides have been confirmed on DPiMS-8060, which enables fast screening analysis in short time.
- Matrix calibration curves are necessary to obtain satisfactory recovery.

References

- 1) K. Hiraoka et al. Rapid Commun. Mass Spectrom. 2007, 21, 3139
- 2) K. Zaitzu et al. Anal. Chem. 2016, 88, 3556
- 3) K. Zaitzu et al. Anal. Chem. 2018, 90, 4695

First Edition: December, 2019



For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.