

Determination of 55 Residual Pesticides in Onion and Leek Samples Using GC-MS/MS

ASMS 2013 ThP-189

Feifei Tian, Furong Wang, Jiangtao Xing, Peng Gao,
Jun Fan, Shin-ichi Kawano
Shimadzu Global COE, Shimadzu (China) Co., Ltd.,
China

Determination of 55 Residual Pesticides in Onion and Leek Samples Using GC-MS/MS

1. Introduction

Determination of residual pesticides in vegetables is an important topic in food safety. Most commonly, pesticide residues in the matrix have been analyzed by GC/MS. But applying GC/MS method to onion and leek samples is often difficult, since many of the pesticide peaks are overlapped by matrix compounds. Triple quadrupole GC-MS/MS provides excellent sensitivity and selectivity in analyzing

compounds in complex matrices. In this report, a novel approach for determination of 55 residual pesticides in onion and leek samples using a GC-MS/MS in multiple reaction monitoring (MRM) mode has been developed. The established method was sensitive and reliable for analysis of the 55 pesticides in onion and leek samples.

2. Methods and Materials

Sample Preparation

The samples were extracted by acetonitrile and purified by Carb/PSA solid-phase extraction (SPE) cartridges. Internal standard substance (heptachlor-endo-epoxide, 0.028 mg/L) was used to overcome matrix effect and achieve

quantitative determination.

Organic onion and leek samples were purchased from local food market.

The samples were homogenized in a food processor.




Fig. 1 Sample pretreatment step

Determination of 55 Residual Pesticides in Onion and Leek Samples Using GC-MS/MS

Table 1 List of pesticides

NO	NAME	CAS	Quantitative ion-pair	Collision voltage(V)	Qualitative ion-pair	Collision voltage(V)	Ret. Time
1	Methamidophos	10265 - 92 - 6	141.0>95.0	8	141.0>126.0	4	5.913
2	Dichlorvos	62 - 73 - 7	185.0>93.0	14	185.0>109.0	14	6.033
3	Omethoate	1113 - 02 - 6	156.0>110.0	8	156.0>141.0	4	9.686
4	Sulfotep	3689 - 24 - 5	322.0>202.0	10	322.0>294.0	4	10.559
5	Phorate	298 - 02 - 2	260.0>75.0	8	260.0>231.0	4	10.799
6	alpha-HCH	319 - 84 - 6	218.9>182.9	8	218.9>144.9	20	10.916
7	Dimethoate	60 - 51 - 5	125.0>79.0	8	125.0>47.0	14	11.157
8	beta-HCH	319 - 85 - 7	218.9>182.9	8	218.9>144.9	20	11.447
9	Quintozene	82 - 68 - 8	294.8>236.8	16	294.8>264.8	12	11.502
10	gamma-HCH	58 - 89 - 9	218.9>182.9	8	218.9>144.9	20	11.598
11	Terbufos	13071 - 79 - 9	231.0>174.9	14	231.0>128.9	26	11.672
12	Fonofos	944 - 22 - 9	246.0>109.1	18	246.0>137.1	6	11.761
13	Diazinon	333 - 41 - 5	304.1>179.1	10	304.1>162.1	8	11.797
14	Pyrimethanil	53112 - 28 - 0	198.1>183.1	14	198.1>158.1	18	11.908
15	delta-HCH	319 - 86 - 8	218.9>182.9	10	218.9>144.9	20	12.171
16	Vinclozolin	50471 - 44 - 8	285.0>212.0	12	285.0>178.0	14	12.828
17	Parathion-methyl	298 - 00 - 0	263.0>109.0	14	263.0>136.0	8	12.869
18	Fenitrothion	122 - 14 - 5	277.0>260.0	6	277.0>109.1	14	13.365
19	Malathion	121 - 75 - 5	173.1>99.0	14	173.1>127.0	6	13.535
20	Chlorpyrifos	2921 - 88 - 2	313.9>257.9	14	313.9>285.9	8	13.683
21	Fenthion	55 - 38 - 9	278.0>109.0	20	278.0>125.0	20	13.76
22	Parathion	56 - 38 - 2	291.1>109.0	14	291.1>137.0	6	13.828
23	Triadimefon	43121 - 43 - 3	208.1>181.0	10	208.1>127.0	14	13.895
24	Isocarbophos	24353 - 61 - 5	289.1>136.0	14	289.1>113.0	6	13.909
25	Dicofol	115 - 32 - 2	250.0>139.0	14	250.0>215.0	8	13.988
26	Isufenphos-methyl	83733-82-8	199.0>121.0	14	241.1>121.1	22	14.217
27	Fipronil	120068 - 37 - 3	366.9>212.9	30	366.9>254.9	22	14.414
28	Phosfolan	947-02-4	255.0>227.0	6	255.0>140.0	22	14.524
IS	Heptachlor-epoxide	1024 - 57 - 3	352.8>262.9	14	352.8>281.9	12	14.536
29	Phenthoate	2597 - 03 - 7	273.9>125.0	20	273.9>246.0	6	14.62
30	Quinalphos	13593 - 03 - 8	157.1>129.0	14	157.1>93.0	10	14.643
31	Procymidone	32809 - 16 - 8	283.0>96.0	10	283.0>255.0	12	14.714
32	Methidathion	950 - 37 - 8	145.0>85.0	8	145.0>58.0	14	14.915
33	Profenofos	41198 - 08 - 7	336.9>266.9	14	336.9>308.9	6	15.549
34	p,p'-DDE	72 - 55 - 9	246.0>176.0	30	246.0>211.0	22	15.643
35	p,p'-DDD	72 - 54 - 8	235.0>165.0	24	235.0>199.0	14	16.436
36	o,p'-DDT	789 - 02 - 6	235.0>165.0	24	235.0>199.0	16	16.491
37	Triazophos	24017 - 47 - 8	257.0>162.0	8	257.0>134.0	22	16.685
38	p,p'-DDT	50 - 29 - 3	235.0>165.0	24	235.0>199.0	16	17.146
39	Iprodione	36734 - 19 - 7	314.0>245.0	12	314.0>56.0	22	17.853
40	Phosmet	732 - 11 - 6	160.0>133.0	14	160.0>77.0	24	17.989
41	Bifenthrin	82657 - 04 - 3	181.1>166.1	12	181.1>153.1	8	18.007
42	Fenprothrin	39515 - 41 - 8	265.1>210.1	12	265.1>172.1	14	18.192
43	Phosalone	2310 - 17 - 0	182.0>111.0	14	182.0>138.0	8	18.664
44	Cyhalothrin-1	68085 - 85 - 8	197.0>161.0	8	197.0>141.0	12	18.802
45	Cyhalothrin-2	68085 - 85 - 8	197.0>161.0	8	197.0>141.0	12	18.982
46	Coumaphos	56-72-4	362.0>109.0	16	362.0>226.0	14	19.929
47	Pyridaben	96489 - 71 - 3	147.1>117.1	22	147.1>132.1	14	19.969
48	Cyfluthrin-1	68359 - 37 - 5	226.1>206.1	14	226.1>199.1	6	20.331
49	Cyfluthrin-2	68359 - 37 - 5	226.1>206.1	14	226.1>199.1	6	20.428
50	Cyfluthrin-3,4	68359 - 37 - 5	226.1>206.1	14	226.1>199.1	6	20.488
51	Cypermethrin-1	52315 - 07 - 8	163.1>127.1	6	163.1>91.0	14	20.655
52	Cypermethrin-2	52315 - 07 - 8	163.1>127.1	6	163.1>91.0	14	20.756
53	Cypermethrin-3,4	52315 - 07 - 8	163.1>127.1	6	163.1>91.0	14	20.813
54	Flucythrinate-1	70124 - 77 - 5	199.1>157.1	10	199.1>107.1	22	20.815
55	Flucythrinate-2	70124 - 77 - 5	199.1>157.1	10	199.1>107.1	22	21.007
56	Fenvalerate-1	51630 - 58 - 1	419.1>225.1	6	419.1>167.1	12	21.546
57	Fluvalinate-1	69409 - 94 - 5	250.1>55.0	20	250.1>200.0	20	21.66
58	Fenvalerate-2	51630 - 58 - 1	419.1>225.1	6	419.1>167.1	12	21.724
59	Fluvalinate-2	69409 - 94 - 5	250.1>55.0	20	250.1>200.0	20	21.755
60	Difenoconazole-1	119446 - 68 - 3	323.0>265.0	14	323.0>202.0	28	22.02
61	Difenoconazole-2	119446 - 68 - 3	323.0>265.0	14	323.0>202.0	28	22.089
62	Deltamethrin-1	52918 - 63 - 5	252.9>93.0	20	252.9>171.9	8	22.126
63	Deltamethrin-2	52918 - 63 - 5	252.9>93.0	20	252.9>171.9	8	22.345

Determination of 55 Residual Pesticides in Onion and Leek Samples Using GC-MS/MS

GC/MS/MS Analysis

Treated samples were analyzed in MRM mode using a gas chromatograph coupled with a triple quadrupole mass spectrometer (GCMS-TQ8030, Shimadzu Corporation, Japan). The MRM transitions and collision energies for

every compound were acquired from the pesticide MRM database provided by Shimadzu Corporation which contains 440 pesticides.

Analytical Conditions

GC

Carrier gas	: He
Linear velocity	: 47.2 cm/sec
Injection mode	: splitless (1 min)
Injection port temperature	: 250°C
Column	: Rxi-5 Sil ms, 30 m × 0.25 mm, 0.25 μm
Temperature program	: 50°C (1 min)-25°C/min-125°C-10°C/min-300°C (15 min)

MS/MS

Ionization	: EI
Collision gas	: argon
Solvent cutting time	: 4 min
Ion source temperature	: 200°C
Interface temperature	: 250°C
Detector voltage	: tuning voltage + 0.6 kV
Monitoring mode	: MRM

3. Results and Discussion




Fig. 2 SIM, MRM chromatograms of phenthoate and o, p-DDT spiked in onion samples (0.005 mg/L)

Determination of 55 Residual Pesticides in Onion and Leek Samples Using GC-MS/MS

Pesticide analysis by GC/MS has mostly been performed in SIM (Selected Ion Monitoring) mode, but especially for complex matrices the analytical selectivity achieved was often not sufficient. Fig. 2 shows chromatograms phenthoate and o, p-DDT spiked in onion matrix (0.005 mg/L) in both SIM and MRM modes. The results showed a clear improvement in sensitivity and selectivity when using MRM.

In order to assess the method linearity, calibration curves were constructed for the 55 pesticides spiked in the sample

matrix, using heptachlor-endo-epoxide as the internal standard. The linear relation between peak area ratio and concentration ratio of target and internal standard substance were good from 0.001 mg/L to 0.1 mg/L ($r^2 > 0.99$). The limits of detection (LODs) of all the pesticides studied were 0.0001 mg/kg. The average recoveries were 70%-120% of target compounds and the relative standard deviations (%RSDs, n=6) were less than 11% in two spiked levels at 0.001 mg/kg and 0.005 mg/kg.




Fig. 3 Recovery of pesticides at 0.001mg/L spiked in onion sample




Fig. 4 Recovery of pesticides at 0.001mg/L spiked in leek sample

4. Conclusions

- The developed method is sensitive, rapid, precise, and widely linear, therefore it is suitable for the determination of pesticides.
- The recoveries were acceptable for multi-residue pesticide determination in onion and leek samples. GC/MS/MS technique (MRM) compared to SIM provides excellent sensitivity and selectivity for analysis of pesticides in complex matrices.