

How Do I Determine the Best Chromatography for my Sample

Jean Lane
Application Engineer
LC columns and Consumables Technical Support
January 20, 2021




Outline for Talk

Getting started

- Sample type and define the objective
- Particle options
 - Totally porous vs superficially porous
- Agilent column chemistries
 - Column chemistries C18 phases and more
- Alternate Selectivity
 - Effect on resolution
- Hilic, CS C18, and chiral chemistries
- Role of instruments
 - Instrument and column compatibility
 - Other instrument considerations




Pharma



BioPharma / Biotech



Fine Chemicals



Universities



Agricultural Biochemistry



Food & Health Supplements

Small molecule separations

Examples:

- Environmental pollutants
- Contaminants in food (pesticides)
- Food toxins
- Drugs of abuse

Large molecule separations

Examples:

- Biologics
 - Therapeutic proteins
 - mAbs (monoclonal antibodies)
- Synthetic polymers
 - Polystyrene
 - Polyethylene
- Natural polymers
 - Starch
 - Dextran

Types of Liquid Chromatography

Separation achieved using a physical or chemical characteristic of molecule in solution

LC Modes	Mechanism	Small Molecules	Large Molecules
Normal phase (NP)	Polar particle	✓	
Reversed phase (RP)	Nonpolar particle	✓✓✓✓	✓
Ion exchange (IEX)	Charged particles, +ve or -ve	✓	✓
Size exclusion gel permeation (SEC/GPC)	Porous particles with defined pore size	✓	✓
Affinity	Exploits biological interactions – antigen-antibody		✓
Chiral	Shape recognition	✓	
Hydrophilic interaction (HILIC)	Similar to normal phases but uses aqueous/organic mobile phases	✓	✓
Hydrophobic interaction (HIC)	Similar to RP but uses aqueous mobile phases		✓

Reversed phase is the primary technique for small molecule LC
Multiple techniques are required to characterize a biomolecule

Define the Objective


- What is the sample type?
- How complex is the sample?
- What are the goals of the separation?
- Is high efficiency important?
- Is speed important?
- Any instrument considerations or limitations?



Particle Options

Comparison of two porous particles

Novel superficially porous particle versus traditional totally porous particle



ZORBAX – Fully Porous Particles

Fully Porous Particle Column (ZORBAX)	Superficially Porous Particle Column (Poroshell)
Methods that require high sample loading	Faster, more efficient separations
Direct use in legacy fully porous particle methods	Newly developed or transferred methods
Direct scalability to prep-scale chemistries	Screening or method development
Some chemistries are unique to ZORBAX	High throughput methods

- Rapid trend to **superficially porous particles** due to their robustness and high efficiency
- More efficient Poroshell particles are typically the first choice for new methods
- In some cases, fully porous particles may still be considered

Why InfinityLab Poroshell?

Performance

Poroshell columns deliver improved performance compared to totally porous materials with the same particle size, or similar performance at significantly lower backpressures.

Compatibility

Poroshell 120 columns are ideal for the Agilent InfinityLab LC and LC/MSD instruments. Column ID tags enable convenient column tracking and reporting.

Scalability


Poroshell columns are available in 1.9, 2.7, and 4 μm pore size and simplify method transfer from HPLC to UHPLC to ultra-low dispersion UHPLC.

Robustness


Poroshell columns show excellent lifetime and are less demanding in sample preparation than fully-porous material with similar performance. This is due to bigger particle sizes and frits.

Selectivity


The stationary phase is the dominant factor for selectivity and, therefore, resolution. Poroshell offers a wide range of innovative chemistries to achieve selectivity differences.



Poroshell Technology – What Makes it Better?




Poroshell is made of a solid core with a porous outer layer



- Analytes travel through the particle more efficiently
- High efficiency allows you to use a larger SPP (2.7 μm) for nearly equivalent performance to a smaller sub-2 μm (STM) TPP column
- Using a larger particle allows for lower backpressure compared to similar, efficient totally porous STM columns and flexible use on HPLC or UHPLC systems

Which Particle Size is Best for Your Needs?


UHPLC
1.9 μm



Highest UHPLC performance

- Pressure rating: 1300 bar
- Typical pressure: similar to sub-2 μm totally porous
- Efficiency: ~120% of sub-2 μm totally porous
- Ideal for: Agilent 1290 Infinity II LC System

UHPLC/
HPLC




2.7 μm

UHPLC performance at lower pressures

- Pressure rating: 600 bar
- Typical pressure: 50% of sub-2 μm totally porous
- Efficiency: ~90% of sub-2 μm totally porous
- Ideal for: Agilent 1260 Infinity II LC System or Agilent 1260 Infinity II Prime LC System

HPLC



4.0 μm

Improved HPLC performance

- Pressure rating: 600 bar
- Typical pressure: often <200 bar
- Efficiency: ~200% of 5 μm totally porous
- Ideal for: Agilent 1220 Infinity II LC System

Which Particle Size is Best for Your Needs?

UHPLC




Highest UHPLC performance

- Pressure rating: 1300 bar
- Typical pressure: similar to sub-2 µm totally porous
- Efficiency: ~120% of sub-2 µm totally porous
- Ideal for: Agilent 1290 Infinity II LC System

Most flexible/common
particle size


UHPLC/
HPLC



UHPLC performance at lower pressures

- Pressure rating: 600 bar
- Typical pressure: 50% of sub-2 µm totally porous
- Efficiency: ~90% of sub-2 µm totally porous
- Ideal for: Agilent 1260 Infinity II LC System or Agilent 1260 Infinity II Prime LC System


HPLC



Improved HPLC performance

- Pressure rating: 600 bar
- Typical pressure: often <200 bar
- Efficiency: ~200% of 5 µm totally porous
- Ideal for: Agilent 1220 Infinity II LC System


Column Choices – Which Particle Type to Choose?



Totally porous particle

ZORBAX Eclipse Plus C18
4.6 x 250 mm, 5 μ m


Runtime: 35 minutes



Superficially porous particle

InfinityLab Poroshell 120 EC-C18
4.6 x 100 mm, 2.7 μ m

Runtime: 9 minutes



A: 0.1% formic acid in water, B: ACN
Gradient: 8-33% ACN in 30 or 8 min
1 or 2 mL/min, 25 $^{\circ}$ C, 254 nm
Agilent application note: 5990-5572EN

When to choose which product family

Special phases		
HPLC	UHPLC	LD-UHPLC
5 μm , 3 μm	---	---
Features	Phases	
High carbon load columns	Pursuit XRs, Pursuit XRs Ultra	
Analytical to prep	Pursuit, Polaris	
Alternative selectivity for polar and nonpolar	Polaris C18-Ether, C18 Amide, NH2	
<p>Unique chemistries that help to solve nonstandard applications from HPLC to prep</p>		

Agilent's Column Families for Small Molecule LC Columns

When to choose which product family

ZORBAX

HPLC

5 μ m, 3.5 μ m

UHPLC

1.8 μ m
(RRHT)

LD-UHPLC

1.8 μ m
(RRHD)

Features

Traditional, reliable columns that offer a vast amount of unique chemistries

Scalable phases that range from UHPLC to HPLC to Prep

Higher overall retention, accepts larger amounts of strong solvent during injection and often shows higher sample capacity.

Scalable from UHPLC to HPLC to Prep-LC with high sample capacity.

A proven and reliable portfolio of totally porous HPLC columns

The Agilent ZORBAX family offers all advantages of totally porous particle columns such as increased retention, loadability and resistance to sample solvents. Easily scale your methods all the way from UHPLC to preparative LC.



Agilent ZORBAX	Chemistry	Particle Size	Pore Size (Å)	Temperature Limit	pH Range	Endcapped	Carbon Load (%)	Surface Area	USP Designation	Benefits and Applications
EdgePlus C18		1.8, 3.5, 5	95	60 °C	2-9	Double	9	160 m ² /g	L1	General purpose Starting Point for LC method development
EdgePlus C8		1.8, 3.5, 5	95	60 °C	2-9	Double	7	160 m ² /g	L7	General purpose Lower retention of hydrophobic analytes vs. C18
EdgePlus PhenylHexyl		1.8, 3.5, 5	95	60 °C	2-8	Double	9	160 m ² /g	L11	Alternative selectivity for aromatic compounds Enhanced pi-pi interactions when using methanol
EdgePlus PMH	Polymethac C18	1.8, 3.5, 5	95	60 °C	2-9	Double	14	160 m ² /g	L1	Application-specific Designed for the separation of PAHs in LC
Edge XDB C18		1.8, 3.5, 5	80	60 °C	2-9	Double	10	180 m ² /g	L1	General purpose, higher carbon load Higher hydrophobicity with alternative selectivity for lipophilic analytes
Edge XDB C8		1.8 (RRHT) 3.5, 5, 7	80	60 °C	2-9	Double	7.6	180 m ² /g	L7	General purpose, higher carbon load Higher hydrophobicity with alternative selectivity for lipophilic analytes but reduced retention vs. XDB-C18
Edge XDB Phenyl		3.5, 5	80	60 °C	2-9	Double	7.2	180 m ² /g	L11	Alternative selectivity for aromatic compounds Enhanced pi-pi interactions when using methanol
Edge XDB CN		3.5, 5	80	60 °C	2-9	Double	4.2	180 m ² /g	L10	Polar analytes in RP low bleed Excellent peak shape of polar and mid-polar compounds
StableBond C18		1.8, 3.5, 5, 7	80	90 °C	0.8-8	No	10	180 m ² /g	L1	Low pH and high temperature Excellent stability and peak shape at highly acidic conditions
StableBond C8		1.8, 3.5, 5, 7	80	80 °C	1-8	No	5.5	180 m ² /g	L7	Low pH and high temperature Lower retention of hydrophobic analytes vs. C18
StableBond C2		1.8, 3.5, 5	80	80 °C	1-8	No	4	180 m ² /g	L56	Low pH and high temperature Reduced retention of hydrophobic analytes
StableBond AQ		1.8, 3.5, 5, 7	80	80 °C	1-8	No	Proprietary	180 m ² /g	L36	Polar analytes in RP Excellent peak shape and retention of polar compounds using reversed-phase LC, stable at 100% aqueous mobile phases
StableBond Phenyl		1.8, 3.5, 5, 7	80	80 °C	1-8	No	5.5	180 m ² /g	L11	Alternative selectivity for aromatic compounds Enhanced pi-pi interactions when using methanol
StableBond CN		1.8, 3.5, 5, 7	80	80 °C	1-8	No	4	180 m ² /g	L10	Polar molecules at low pH or high temperature, low bleed Excellent peak shape of polar and mid-polar compounds
Extend C18		1.8, 3.5, 5, 7	80	60 °C	2-11.5	Double	4	180 m ² /g	L1	High pH applications Robust performance and long lifetime under high pH
Bonus-8P		1.8, 3.5, 5, 7	80	60 °C	2-9	Triple	9.5	180 m ² /g	L60	Alternative Selectivity to C18 Improved peak shape for basic compounds, stable in 100% aqueous conditions
HILIC Plus		1.8, 3.5	95	Only mobile phase limits apply	1-8	No	0	180 m ² /g	L3	Polar analytes in HILIC mode Excellent retention of polar compounds by HILIC
Rx C18		3.5, 5, 7	80	60 °C	2-8	No	12	180 m ² /g	L1	General purpose High carbon load for increased retention
Rx C8		3.5, 5	80	80 °C	1-8	No	5.5	180 m ² /g	L7	General purpose
Rx C2		1.8 (RRHT) 5, 7	80	Only mobile phase limits apply	0.8-8	No	0	180 m ² /g	L3	Polar compounds in HILIC, NPLC and SFC mode Good starting point for method development

Which particle is best for my method?

	1.8 μ m ZORBAX (RRHT) highest UHPLC performance Maximum pressure: 1200 bar Ideal for: 1250 Infinity II LC or 1250 Infinity I Prime LC
	3.5 μ m ZORBAX (RRHT) ultra-fast chromatography at up to 600 bar Maximum pressure: 600 bar Ideal for: 1250 Infinity II LC
	5 μ m ZORBAX (RR) higher resolution of HPLC methods Maximum pressure: 400 bar Upgrade of traditional methods on general HPLC instruments
	8 μ m ZORBAX: Proven and reliable for HPLC methods Maximum pressure: 400 bar Used for traditional methods on general HPLC instruments and in preparative LC

What column ID and length should I choose?

Format	Comment
Column ID	4.6 mm for legacy methods 3.0 mm for lower solvent use than 4.6 mm 2.1 mm for lowest solvent use and MS applications
Column length	Shorter 50 to 100 mm for fastest separations Longer: 150 to 250 mm for increased resolution

Interested in modernizing your LC methods?

InfinityLab Poroshell chemistries are aligned with traditional ZORBAX chemistries—making it easy to transfer your methods from fully porous to super-fine porous particle columns.

ZORBAX Chemistry	InfinityLab Poroshell 120 Chemistry
ZORBAX EdgePlus C18	InfinityLab Poroshell 120 EC-C18
ZORBAX EdgePlus EC-C8	InfinityLab Poroshell 120 EC-C8
ZORBAX EdgePlus PhenylHexyl	InfinityLab Poroshell 120 PhenylHexyl
ZORBAX StableBond (RR)C18	InfinityLab Poroshell 120 (RR)C18
ZORBAX StableBond (RR)C8	InfinityLab Poroshell 120 (RR)C8
ZORBAX Bonus-8P	InfinityLab Poroshell 120 Bonus-8P
ZORBAX StableBond (RR)AQ	InfinityLab Poroshell 120 (RR)AQ
ZORBAX EdgePlus XDB-CN	InfinityLab Poroshell 120 EC-CN
ZORBAX HILIC Plus	InfinityLab Poroshell 120 HILIC

Agilent InfinityLab is an optimized portfolio of LC instruments, columns, and supplies that work together seamlessly for maximum efficiency and performance—regardless of application area. More information at www.agilent.com/chem/infinitylab

For more information about ZORBAX columns, go to www.agilent.com/chem/ZORBAX

Zorbax Selectivity Overview
Pub No. 5994-2212EN

Agilent's Column Families for Small Molecule LC Columns

When to choose which product family

InfinityLab Poroshell 120

HPLC

UHPLC

LD-UHPLC

4 μm

2.7 μm

1.9 μm

Features

Modern column technology that offers higher performance at similar backpressure

or comparable performance at reduced backpressure

Designed in with Agilent LC instruments and supplies

Universal column platform with offerings for all separation modes, i.e. RP, NP, HILIC, SFC as well as chiral separations

Modern, high-performance HPLC and UHPLC columns designed in for state-of-the-art instruments.

More Chemistries, More Choices For Solving Your Toughest Separation Challenges

The InfinityLab Poroshell 120 family has grown to include 3 particle sizes and 18 chemistries—including new phases for chiral and HILIC separations. So, you can efficiently separate the widest variety of compounds.



InfinityLab Poroshell 120	Chemistry	Particle Size	Flow Size	Temperature Limit	pH Range	Endcapped	Carbon Load	Surface Area	USP Designation	Benefits and Applications
EC-C18		1.9 μm, 2.7 μm, 4 μm	120 Å	60 °C	2.0-8.0	Yes	10%	130 m ² /g	L1	General purpose Excellent peak shape and efficiency for acids, bases, and neutrals
EC-C8		1.9 μm, 2.7 μm, 4 μm	120 Å	60 °C	2.0-8.0	Yes	9%	130 m ² /g	L7	General purpose Lower retention of hydrophobic analytes vs. C18
SB-C18		2.7 μm	120 Å	90 °C	1.0-8.0	No	9%	130 m ² /g	L1	Low pH Excellent stability and peak shape in highly acidic conditions
SB-C8		2.7 μm	120 Å	80 °C	1.0-8.0	No	5.5%	130 m ² /g	L7	Low pH Excellent stability at low pH Lower retention of hydrophobic analytes vs. C18
HPH-C18		1.9 μm, 2.7 μm, 4 μm	100 Å	60 °C	3.0-11.0	Yes	Proprietary	95 m ² /g	L1	High pH capable Robust performance and long lifetime Improved retention, resolution, and peak shape of basic compounds
HPH-C8		2.7 μm, 4 μm	100 Å	60 °C	3.0-11.0	Yes	Proprietary	95 m ² /g	L7	High pH capable Robust performance and long lifetime Lower retention of hydrophobic analytes vs. C18
Bonus-RP		2.7 μm	120 Å	60 °C	2.0-8.0	Yes	9.5%	130 m ² /g	L60	Alternate selectivity to C18 Improved peak shape for basic compounds, stable in 100% aqueous conditions
PPP		1.9 μm, 2.7 μm, 4 μm	120 Å	60 °C	2.0-8.0	Yes	5.1%	130 m ² /g	L43	Alternate selectivity Excellent peak shape for polar and nonpolar analytes Unique selectivity for aromatic and hydrogenated compounds
Phenyl-Hexyl		1.9 μm, 2.7 μm, 4 μm	120 Å	60 °C	2.0-8.0	Yes	9%	130 m ² /g	L11	Alternate selectivity with aromatic groups Highly nonpolar bonded phase takes advantage of π-π interactions
SB-Aq		2.7 μm	120 Å	80 °C	1.0-8.0	No	Proprietary	130 m ² /g	L96	Alternate selectivity Excellent peak shape and retention of polar compounds using reversed-phase LC Exceptional stability under high-pressure conditions, including 100% water
EC-CN		2.7 μm	120 Å	60 °C	2.0-8.0	Yes	3.5%	130 m ² /g	L10	Alternate selectivity Use in reversed-phase for alternate selectivity of polar and mid-polar compounds Use in normal-phase for excellent peak shape and retention of nonpolar analytes
HILIC-Z		2.7 μm	100 Å	80 °C	3.0-11.0	No	Proprietary	95 m ² /g	L114	Polar analysis Excellent retention of highly polar or charged compounds by HILIC Rugged performance at high pH or high temperature
HILIC		1.9 μm, 2.7 μm, 4 μm	120 Å	60 °C	0.0-8.0	No	NA	130 m ² /g	L3	Polar analysis Excellent retention of polar compounds by HILIC
HILIC-ONS		2.7 μm	120 Å	45 °C	1.0-7.0	Proprietary	Proprietary	130 m ² /g	L86	Polar analysis Fluorinated bonded phase offers alternate selectivity to other HILIC phases
Chiral-V		2.7 μm	120 Å	45 °C	2.5-7.0	Proprietary	Proprietary	130 m ² /g	L88	Chiral separations Amines, proteins, and complex basic and neutral compounds Reversed-phase, polar basic, normal phase, or polar organic modes
Chiral-T		2.7 μm	120 Å	45 °C	2.5-7.0	Proprietary	Proprietary	130 m ² /g	L63	Chiral separations Beta blockers, hydroxy acids, amino acids, proteins, benzodiazepines, and hydantins Reversed-phase, polar basic, normal phase, or polar organic modes
Chiral-CD		2.7 μm	120 Å	45 °C	3.0-7.0	Proprietary	Proprietary	130 m ² /g	L45	Chiral separations Sildenafil, fungicides, and protected amino acids Reversed phase or polar organic modes
Chiral-CF		2.7 μm	120 Å	45 °C	3.0-7.0	Proprietary	Proprietary	130 m ² /g	NA	Chiral separations Polar analysis Polar organic or normal phase modes

© 2018 Agilent Technologies, Inc. All rights reserved. See www.agilent.com for more information.

Note: HILIC-ONS and all four Chiral phases have a pressure limit of 400 bar.

Which particle is best for my method?

1.9 μm: Highest UHPLC performance
 - Maximum pressure: 1300 bar
 - Ideal for Agilent 1200 Infinity ILC

2.7 μm: UHPLC performance at lower pressures
 - Maximum pressure: 600 bar (unless otherwise noted)
 - Ideal for Agilent 1200 Infinity ILC or Agilent 1200 Infinity II Prime ILC

4 μm: Improved HPLC performance
 - Maximum pressure: 600 bar
 - Ideal for Agilent 1200 Infinity ILC

What column ID and length should I choose?

Format	Comment
Column ID	4.6 mm for legacy methods 3.0 mm for lower solvent use than 4.6 mm 2.1 mm for lowest solvent use and MS applications
Column length	Shorter 30 to 100 mm for fastest separations Longer 150 to 250 mm for increased resolution

Still using a legacy method?

InfinityLab Poroshell chemistries are aligned with traditional ZORBAX chemistries—making it easy to transfer your methods from fully porous to superficially porous particle columns.

InfinityLab Poroshell Chemistry	Aligned Chemistry
InfinityLab Poroshell 120 EC-C18	ZORBAX Eclipse Plus C18
InfinityLab Poroshell 120 EC-C8	ZORBAX Eclipse Plus EC-C8
InfinityLab Poroshell 120 Phenyl-Hexyl	ZORBAX Eclipse Plus Phenyl-Hexyl
InfinityLab Poroshell 120 SB-C18	ZORBAX StableBond SB-C18
InfinityLab Poroshell 120 SB-C8	ZORBAX StableBond SB-C8
InfinityLab Poroshell 120 Bonus-RP	ZORBAX Bonus-RP
InfinityLab Poroshell 120 SB-Aq	ZORBAX StableBond SB-Aq
InfinityLab Poroshell 120 EC-CN	ZORBAX Eclipse Plus EC-CN
InfinityLab Poroshell 120 HILIC	ZORBAX HILIC Plus

Agilent InfinityLab is an optimized portfolio of LC instruments, columns, and supplies that work together seamlessly for maximum efficiency and performance—regardless of application area. More information at: www.agilent.com/chem/infinitylab

For more information about InfinityLab Poroshell 120 Columns, go to www.agilent.com/chem/poroshell-120

Poroshell Selectivity Overview
Pub No. 5991-9013EN

Method transferability across product families

Traditional ZORBAX chemistries are aligned with InfinityLab Poroshell chemistries to offer simplified method transfer from fully porous particles to superficially porous particle columns.

InfinityLab Poroshell Chemistries

InfinityLab Poroshell 120 EC-C18

InfinityLab Poroshell 120 EC-C8

InfinityLab Poroshell 120 Phenyl-Hexyl

InfinityLab Poroshell 120 SB-C18

InfinityLab Poroshell 120 SB-C8

InfinityLab Poroshell 120 SB-Aq

InfinityLab Poroshell 120 Bonus-RP

InfinityLab Poroshell 120 EC-CN

InfinityLab Poroshell 120 HILIC



Aligned Chemistry

ZORBAX Eclipse Plus C18

ZORBAX Eclipse Plus C8

ZORBAX Eclipse Plus Phenyl-Hexyl

ZORBAX StableBond SB-C18

ZORBAX StableBond SB-C8

ZORBAX StableBond SB-Aq

ZORBAX Bonus-RP

ZORBAX Eclipse XDB-CN

ZORBAX HILIC-Plus

For more information on method transfer:
technical overview 5990-6588EN

The Poroshell 120 Family Offers Chemistries with Unique Selectivity

InfinityLab Poroshell 120 offers a broad portfolio to suit your needs

For Nonpolar Analytes					
Best All Around	Best for Low pH Mobile Phases	Best for High pH Mobile Phases	Best for Alternative Selectivity	Best for More Polar Analytes	Chiral
EC-C18 1.9 μm, 2.7 μm, 4 μm	SB-C18 1.9 μm, 2.7 μm, 4 μm	HPH-C18 1.9 μm, 2.7 μm, 4 μm	Bonus-RP 2.7 μm	SB-Aq 1.9 μm, 2.7 μm, 4 μm	Chiral-V 2.7 μm
EC-C8 1.9 μm, 2.7 μm, 4 μm	SB-C8 2.7 μm	HPH-C8 2.7 μm, 4 μm	PFP 1.9 μm, 2.7 μm, 4 μm	EC-CN 2.7 μm	Chiral-T 2.7 μm
Phenyl-Hexyl 1.9 μm, 2.7 μm, 4 μm		CS-C18 2.7 μm		HILIC 1.9 μm, 2.7 μm, 4 μm	Chiral-CD 2.7 μm
				HILIC-Z 1.9 μm, 2.7 μm, 4 μm	Chiral-CF 2.7 μm
				HILIC-OH5 2.7 μm	

The Poroshell 120 Family Offers Chemistries with Unique Selectivity

InfinityLab Poroshell 120 offers a broad portfolio to suit your needs

For Nonpolar Analytes					
Best All Around	Best for Low pH Mobile Phases	Best for High pH Mobile Phases	Best for Alternative Selectivity	Best for More Polar Analytes	Chiral
EC-C18 1.9 μm, 2.7 μm, 4 μm	SB-C18 1.9 μm, 2.7 μm, 4 μm	HPH-C18 1.9 μm, 2.7 μm, 4 μm	Bonus-RP 2.7 μm	SB-Aq 1.9 μm, 2.7 μm, 4 μm	Chiral-V 2.7 μm
EC-C8 1.9 μm, 2.7 μm, 4 μm	SB-C8 2.7 μm	HPH-C8 2.7 μm, 4 μm	PFP 1.9 μm, 2.7 μm, 4 μm	EC-CN 2.7 μm	Chiral-T 2.7 μm
Phenyl-Hexyl 1.9 μm, 2.7 μm, 4 μm		CS-C18 2.7 μm		HILIC 1.9 μm, 2.7 μm, 4 μm	Chiral-CD 2.7 μm
				HILIC-Z 1.9 μm, 2.7 μm, 4 μm	Chiral-CF 2.7 μm
				HILIC-OH5 2.7 μm	

The Poroshell 120 Family Offers Chemistries with Unique Selectivity


InfinityLab Poroshell 120 offers a broad portfolio to suit your needs

For Nonpolar Analytes					
Best All Around	Best for Low pH Mobile Phases	Best for High pH Mobile Phases	Best for Alternative Selectivity	Best for More Polar Analytes	Chiral
EC-C18 1.9 μm, 2.7 μm, 4 μm	SB-C18 1.9 μm, 2.7 μm, 4 μm	HPH-C18 1.9 μm, 2.7 μm, 4 μm	Bonus-RP 2.7 μm	SB-Aq 1.9 μm, 2.7 μm, 4 μm	Chiral-V 2.7 μm
EC-C8 1.9 μm, 2.7 μm, 4 μm	SB-C8 2.7 μm	HPH-C8 2.7 μm, 4 μm	PFP 1.9 μm, 2.7 μm, 4 μm	EC-CN 2.7 μm	Chiral-T 2.7 μm
Phenyl-Hexyl 1.9 μm, 2.7 μm, 4 μm	For RP start here	CS-C18 2.7 μm		HILIC 1.9 μm, 2.7 μm, 4 μm	Chiral-CD 2.7 μm
				HILIC-Z 1.9 μm, 2.7 μm, 4 μm	Chiral-CF 2.7 μm
				HILIC-OH5 2.7 μm	

Best All Round

Poroshell 120 EC C18 or EC C8


InfinityLab Poroshell 120 EC C18, 3.0 x 50 mm, 2.7 μ m, p/n 699975-302



Conditions

Mobile phase:	60% ACN, 40% water
Flow rate:	0.85 mL/min
Temperature:	26 °C
Detection:	254 nm
Sample:	2 μ L of RRLC checkout sample (p/n 5188-6529), alkylphenones

InfinityLab Poroshell 120 EC C8, 3.0 x 50 mm, 2.7 μ m, p/n 699975-306



Need For Use Over a Wide pH Range

Poroshell HPH C18 and HPH C8

Conditions

Flow rate: 2.0 mL/min

Temperature: 26 °C

UV detection: 254 nm




Gradient:

Time	%Buffer	%ACN
0	10	90
5	90	10
7	10	90


Sample:

1. Procainamide
2. Caffeine
3. Acetylsalicylic acid
4. Hexanophenone deg.
5. Dipyrimadole
6. Diltiazem
7. Diflunisal
8. Hexanophenone

InfinityLab Poroshell HPH-C18, 4.6 x 50 mm, 2.7 μm, p/n 699975-702




Polar Embedded Phase for Alternate Selectivity



Beta blockers

1. Atenolol
2. Pindolol
3. Naldolol
4. Metoprolol
5. Acebutolol
6. Propranolol
7. Alprenolol



**Polar embedded
group provides
unique selectivity**


0.35 mL/min, 10 mM pH 3.8 ammonium formate buffer and methanol, 10-70% methanol/12 min,
2.1 x 100 mm, 40 °C , DAD 260 nm

Selectivity impacts the resolution most

$$R_s = \frac{\sqrt{N}}{4} \left(\frac{\alpha - 1}{\alpha} \right) \cdot \frac{k'}{k' + 1}$$

Chemistries with Unique Selectivity


Selectivity impacts the resolution most



$$R_s = \frac{\sqrt{N}}{4} \left(\frac{\alpha - 1}{\alpha} \right) \cdot \frac{k'}{k' + 1}$$

Alpha	1.10	1.35	1.60	1.85	2.1
plates	5,000	10,000	15,000	20,000	25,000
k'	2.0	4.5	7.0	9.5	12.0

Selectivity impacts the resolution most



$$R_s = \frac{\sqrt{N}}{4} \left(\frac{\alpha - 1}{\alpha} \right) \cdot \frac{k'}{k' + 1}$$

Selectivity impacts resolution

- Stationary and mobile phase
- Temperature
- N is strongly influenced by alpha

Alpha	1.10	1.35	1.60	1.85	2.1
plates	5,000	10,000	15,000	20,000	25,000
k'	2.0	4.5	7.0	9.5	12.0

Why Not to Use C18 chemistry

Some reasons for trying an alternative

C18 columns are a very common 'go to' column for HPLC methods and may be a suitable column choice for simple methods. **But** the C18 chemistry may not be the optimal choice.

For analysis of polar analytes, it is typically suggested that these sample types can be better separated on chemistries that have a greater polarity than C18.

Reasons to **try** another chemistry:

- Too much retention or selectivity with C18 for desired analysis time
- Polar analytes are not well retained with low or no organic modifier
- Polar analytes not well resolved even if retained
- A C18 method already in use is not rugged enough (revalidate)

- Screening different column chemistries is commonly advised when sample mixtures are complex

Why is Changing the Bonded Phase Effective?

- Differences in interactions between polar and nonpolar compounds
- Other types of interactions with a bonded phase can be exploited (for example, pi-pi interactions)
- These all change with the bonded phase
- Changing the bonded phase can improve selectivity/resolution
- Reduce analysis time

Orthogonal: Orthogonality in chromatography refers to alternative selectivity between separations.

12 Reverse Phase LC Chemistries

InfinityLab Poroshell 120 offers a broad portfolio to suit your needs

Best All Around	Best for Low pH Mobile Phases	Best for High pH Mobile Phases	Best for Alternative Selectivity	Best for More Polar Analytes	Chiral
EC-C18 1.9 µm, 2.7 µm, 4 µm	SB-C18 1.9 µm, 2.7 µm, 4 µm	HPH-C18 1.9 µm, 2.7 µm, 4 µm	Bonus-RP 2.7 µm	SB-Aq 1.9 µm, 2.7 µm, 4 µm	Chiral-V 2.7 µm
EC-C8 1.9 µm, 2.7 µm, 4 µm	SB-C8 2.7 µm	HPH-C8 2.7 µm, 4 µm	PFP 1.9 µm, 2.7 µm, 4 µm	EC-CN 2.7 µm	Chiral-T 2.7 µm
Phenyl-Hexyl 1.9 µm, 2.7 µm, 4 µm		CS-C18 2.7 µm		HILIC 1.9µm, 2.7 µm, 4 µm	Chiral-CD 2.7 µm
				HILIC-Z 1.9 µm, 2.7 µm, 4 µm	Chiral-CF 2.7 µm
				HILIC-OH5 2.7 µm	

12 Reverse Phase LC Chemistries

InfinityLab Poroshell 120 offers a broad portfolio to suit your needs

Best All Around	Best for Low pH Mobile Phases	Best for High pH Mobile Phases	Best for Alternative Selectivity	Best for More Polar Analytes	Chiral
EC-C18 1.9 µm, 2.7 µm, 4 µm	SB-C18 1.9 µm, 2.7 µm, 4 µm	HPH-C18 1.9 µm, 2.7 µm, 4 µm	Bonus-RP 2.7 µm	SB-Aq 1.9 µm, 2.7 µm, 4 µm	Chiral-V 2.7 µm
EC-C8 1.9 µm, 2.7 µm, 4 µm	SB-C8 2.7 µm	HPH-C8 2.7 µm, 4 µm	PFP 1.9 µm, 2.7 µm, 4 µm	EC-CN 2.7 µm	Chiral-T 2.7 µm
Phenyl-Hexyl 1.9 µm, 2.7 µm, 4 µm			CS-C18 2.7 µm	Very orthogonal 1.9 µm, 2.7 µm, 4 µm	Chiral-CD 2.7 µm
				HILIC-Z 1.9 µm, 2.7 µm, 4 µm	Chiral-CF 2.7 µm
				HILIC-OH5 2.7 µm	

Control Selectivity with Choice of Phases

Poroshell 120 columns





InfinityLab Poroshell 120 columns, 2.1 x 100 mm, 2.7 μm

Conditions

Flow rate: 0.4 mL/min
Mobile phase: A: 0.1% formic acid in H₂O
B: MeOH + 0.1% formic acid
Temperature: 25 °C
Detection: 260 nm
Gradient: 40-80 and MeOH/14 min





Sample: 8 steroids

1. Hydrocortisone
2. B Estradiole
3. Androstadiene 3,17 dione
4. Testosterone
5. Etheestradione
6. Estrone



Agilent publication number: 5990-5951EN

The influence of stationary phase on selectivity and resolution is dominant




Time	% Organic
0	8
6	100
7	100
8	8

2mL/min 254 nm


Compounds

1. APAP
2. Phenacetin
3. Piroxicam
4. Tolmetin
5. Ketoprofen
6. Naproxen
7. Sulindac
8. Diclofenac
9. Diflunisal

Importance of Alternate Selectivity Chemistries



- Three compounds
 - Same molecular weight
 - Only differ by positional location of the functionality



InfinityLab Poroshell 120 columns, 4.6 x 50 mm, 2.7 μ m
70:30 – MeOH/H₂O, 1.5 mL/min, 40 °C, 254 nm

The Poroshell 120 Family of Chemistries

InfinityLab Poroshell 120 offers a broad portfolio to suit your needs

Best All Around	Best for Low pH Mobile Phases	Best for High pH Mobile Phases	Best for Alternative Selectivity	Best for More Polar Analytes	Chiral
EC-C18 1.9 µm, 2.7 µm, 4 µm	SB-C18 1.9 µm, 2.7 µm, 4 µm	HPH-C18 1.9 µm, 2.7 µm, 4 µm	Bonus-RP 2.7 µm	SB-Aq 1.9 µm, 2.7 µm, 4 µm	Chiral-V 2.7 µm
EC-C8 1.9 µm, 2.7 µm, 4 µm	SB-C8 2.7 µm	HPH-C8 2.7 µm, 4 µm	PFP 1.9 µm, 2.7 µm, 4 µm	EC-CN 2.7 µm	Chiral-T 2.7 µm
Phenyl-Hexyl 1.9 µm, 2.7 µm, 4 µm	Charged surface C18		CS-C18 2.7 µm	HILIC 1.9µm, 2.7 µm, 4 µm	Chiral-CD 2.7 µm
			HILIC-Z 1.9 µm, 2.7 µm, 4 µm	Chiral-CF 2.7 µm	
			HILIC-OH5 2.7 µm		
			HILIC chemistries		

What is HILIC and When Should it be Considered?

HILIC complements RPLC

Reversed-Phase LC	
	Nonpolar stationary phase (for example, C18)
	Polar mobile phase H ₂ O/CH ₃ OH, H ₂ O/CH ₃ CN
	Decrease retention by decreasing polarity of mobile phase CH ₃ CN ↑ = retention ↓
	Polar to nonpolar

Polarity

Mobile phase

Gradient

Elution order

Hydrophilic Interaction LC (HILIC)	
	Polar stationary phase (for example, silica)
	Polar mobile phase CH ₃ CN/H ₂ O
	Decrease retention by increasing polarity of mobile phase ddH ₂ O ↑ = retention ↓
	Nonpolar to polar

InfinityLab Poroshell 120 HILIC column options

Best for More Polar Analytes

HILIC

1.9 μm , 2.7 μm , 4 μm

HILIC-Z

1.9 μm , 2.7 μm , 4 μm

HILIC-OH5

2.7 μm

For more information on method development for HILIC,
see Agilent technical note 5991-9271EN


InfinityLab Poroshell 120 HILIC column options

Best for More Polar Analytes

HILIC
1.9 μm, 2.7 μm, 4 μm

HILIC-Z
1.9 μm, 2.7 μm, 4 μm

HILIC-OH5
2.7 μm



HILIC

- Bare silica chemistry
- For very simple mixtures

For more information on method development for HILIC,
see Agilent technical note 5991-9271EN


InfinityLab Poroshell 120 HILIC column options

Best for More Polar Analytes

HILIC
1.9 μm, 2.7 μm, 4 μm

HILIC-Z
1.9 μm, 2.7 μm, 4 μm

HILIC-OH5
2.7 μm



HILIC

- Bare silica chemistry
- For very simple mixtures

HILIC-Z

- Proprietary zwitterionic chemistry, high pH stable
- **The most modern and robust column – start method development here**
- PEEK-lined version available

For more information on method development for HILIC,
see Agilent technical note 5991-9271EN


InfinityLab Poroshell 120 HILIC column options

Best for More Polar Analytes

HILIC
1.9 μm, 2.7 μm, 4 μm


HILIC-Z
1.9 μm, 2.7 μm, 4 μm

HILIC-OH5
2.7 μm




HILIC

- Bare silica chemistry
- For very simple mixtures



HILIC-Z

- Proprietary zwitterionic chemistry, high pH stable
- **The most modern and robust column – start method development here**
- PEEK-lined version available




HILIC-OH5


- Brushed fructan chemistry
- Alternative selectivity

For more information on method development for HILIC, see Agilent technical note 5991-9271EN

Separation of 11 Sugars on Poroshell 120 HILIC-Z




1. Xylose
2. Arabinose
3. Fructose
4. Mannose
5. Glucose
6. Galactose
7. Sucrose
8. Maltose
9. Lactose
10. Maltotriose
11. Raffinose



Agilent publication
number: 5991-8984EN

Selectivity Options with Poroshell 120 HILIC Columns




Water soluble vitamins



Columns used: 2.1 x 100 mm, 2.7 μ m; A: 100 mM ammonium acetate + 0.5% acetic acid (pH ~4.6) in H₂O, B: CH₃CN, 0.5 mL/min, 87% B for 1 min, 87-50% B in 4 min, 3 min re-equilibration, 1 μ L injection of individual vitamin standards (0.1-0.4 mg/mL each), 40 $^{\circ}$ C, 260 nm, 80 Hz

InfinityLab Poroshell 120 CS-C18


Built using Agilent innovation



- High pH stable
- Alternate C18 selectivity


- Better peak shape for basic compounds
- Formic acid compatibility
- Reduced operating pressures
- Increased speed of analysis

- Column dimensions**
- 2.1, 3, 4.6 mm id x 50, 100, 150 mm length
 - PEEK-lined options
 - Corresponding Fast Guards

 PEEK-lined column options are rare in the reversed phase column market and help with challenging metal sensitive compounds.

Agilent InfinityLab Poroshell 120 CS-C18

Better peak shape than traditional C18 for analytes with FA mobile phase




A: 0.1% formic acid or 0.2% trifluoroacetic acid in water; B: acetonitrile; 0.4 mL/min; isocratic: %B varies; 2.1 x 100 mm columns, 1 µL injection, 30 °C, LC/MS: ESI+, dMRM; sample: 5 µg/mL of doxepin, desipramine, imipramine, nortriptyline, amitriptyline, trimipramine

Agilent publication number: 5994-2095EN

Agilent InfinityLab Poroshell 120 CS-C18

Better peak shape than traditional C18 for analytes with FA mobile phase





A: 0.1% formic acid or 0.2% trifluoroacetic acid in water; B: acetonitrile; 0.4 mL/min; isocratic: %B varies; 2.1 x 100 mm columns, 1 μ L injection, 30 $^{\circ}$ C, LC/MS: ESI+, dMRM; sample: 5 μ g/mL of doxepin, desipramine, imipramine, nortriptyline, amitriptyline, trimipramine

Agilent publication number: 5994-2095EN

Agilent InfinityLab Poroshell 120 CS-C18





Better peak shape than traditional C18 for analytes with FA mobile phase



A: 0.1% formic acid or 0.2% trifluoroacetic acid in water; B: acetonitrile; 0.4 mL/min; isocratic: %B varies; 2.1 x 100 mm columns, 1 µL injection, 30 °C, LC/MS: ESI+, dMRM; sample: 5 µg/mL of doxepin, desipramine, imipramine, nortriptyline, amitriptyline, trimipramine

Agilent publication number: 5994-2095EN

Choosing Between C18s

InfinityLab Poroshell 120	Chemistry	Pore Size	Endcapped	Carbon Load	Surface Area	Best For
EC-C18 1.9 μm , 2.7 μm , 4 μm		120 Å	Yes	10%	130 m ² /g	General purpose Excellent peak shape and efficiency for acids, bases, neutrals
SB-C18 1.9 μm , 2.7 μm , 4 μm		120 Å	No	9%	130 m ² /g	Low pH Excellent stability and peak shape in highly acidic conditions
HPH-C18 1.9 μm , 2.7 μm , 4 μm		100 Å	Yes	Proprietary	95 m ² /g	Capable of high pH Robust performance and long lifetimes
CS-C18 2.7 μm		100 Å	Yes	Proprietary	95 m ² /g	Alternate selectivity Improved peak shape and sample capacity for basic compounds with low ionic strength mobile phases Capable of high pH

The Poroshell 120 Family

19 chemistries with unique selectivity

InfinityLab Poroshell 120 offers a broad portfolio to suit your needs

Best All Around	Best for Low pH Mobile Phases	Best for High pH Mobile Phases	Best for Alternative Selectivity	Best for More Polar Analytes	Chiral
EC-C18 1.9 µm, 2.7 µm, 4 µm	SB-C18 1.9 µm, 2.7 µm, 4 µm	HPH-C18 1.9 µm, 2.7 µm, 4 µm	Bonus-RP 2.7 µm	SB-Aq 1.9 µm, 2.7 µm, 4 µm	Chiral-V 2.7 µm
EC-C8 1.9 µm, 2.7 µm, 4 µm	SB-C8 2.7 µm	HPH-C8 2.7 µm, 4 µm	PFP 1.9 µm, 2.7 µm, 4 µm	EC-CN 2.7 µm	Chiral-T 2.7 µm
Phenyl-Hexyl 1.9 µm, 2.7 µm, 4 µm		CS-C18 2.7 µm		HILIC 1.9µm, 2.7 µm, 4 µm	Chiral-CD 2.7 µm
				HILIC-Z 1.9 µm, 2.7 µm, 4 µm	Chiral-CF 2.7 µm
				HILIC-OH5 2.7 µm	

Why Do Chiral Separations

- Most small molecule drugs on the market today are either racemates or enantiomerically pure.
- Enantiomers: Same chemical and physical properties, but can have **very different behavioral properties.**
- It is important to characterize each enantiomer.

InfinityLab Poroshell 120 Chiral Chemistries


Column Chemistry	Chiral Selector (bonded chemistry)	Typical LC Mode	Typical Applications
InfinityLab Poroshell 120 Chiral-CF	Derivatized cyclofructan (CF6)	Polar Organic (PO)	Primary amines
		Normal Phase (NP)	Primary amines
InfinityLab Poroshell 120 Chiral-CD	Hydroxypropylated- β -cyclodextrin	Reversed Phase (RP)	Stimulants, fungicides, t-boc amino acids
		Polar Organic (PO)	Complex molecules
InfinityLab Poroshell 120 Chiral-V	Vancomycin (macrolide antibiotic)	Polar Ionic (PI)	Basic pharmaceuticals (various)
		Reversed Phase (RP)	Amines, profens
		Polar Organic (PO)	Complex neutral molecules
InfinityLab Poroshell 120 Chiral-T	Teicoplanin (macrolide antibiotic)	Polar Ionic (PI)	Beta blockers, hydroxyl acids
		Reversed Phase (RP)	Amino acids, hydroxyl acids, profens
		Polar Organic (PO)	Hydantoins, benzodiazepines

Application Examples

InfinityLab Poroshell chiral columns

1-benzyl-2,2-diphenylethylamine

Retention (min): 1.98/3.15 (Peak1/Peak 2)
Resolution: 9.47




Method Conditions

Column: InfinityLab Poroshell 120 Chiral-V (10 cm x 4.6 mm, 2.7 μ m)
Mobile phase: 100/0.1 wt %: Methanol/Ammonium Trifluoroacetate
Flow Rate: 1.0 mL/min
Temperature: Ambient (23 $^{\circ}$ C)
Injection Volume: 1.0 μ L
Detection: UV 220 nm

Benalaxyl – fungicide

Retention (min): 1.64/1.87 (Peak1/Peak 2)
Resolution: 3.13



Method Conditions



Column: InfinityLab Poroshell 120 Chiral-CD (10 cm x 4.6 mm, 2.7 μ m)
Mobile phase: 30/70: Acetonitrile/15 mM Ammonium Formate (pH 3.6)
Flow Rate: 1.0 mL/min
Temperature: Ambient (23 $^{\circ}$ C)
Injection Volume: 1.0 μ L
Detection: UV 220 nm

InfinityLab Chiral Applications Compendium

Publication number: 5991-8450EN

Fungicides

Haloxypop





Retention (min): 1.12/1.37 (Peak1/Peak 2)
Resolution: 4.48

Method Conditions

Column: InfinityLab Poroshell 120 Chiral-T (10 cm x 4.6 mm, 2.7 µm)
Mobile phase: 100/0.3 wt %: Methanol/Ammonium Trifluoroacetate
Flow Rate: 1.0 mL/min
Temperature: Ambient (23 °C)
Injection Volume: 1.0 µL
Detection: UV 220 nm

Mandipropamid





Retention (min): 7.65/8.22 (Peak1/Peak 2)
Resolution: 1.67

Method Conditions

Column: InfinityLab Poroshell 120 Chiral-V (10 cm x 4.6 mm, 2.7 µm)
Mobile phase: 30/70: Methanol/15 mM Ammonium Formate (pH 3.6)
Flow Rate: 0.5 mL/min
Temperature: 45 °C
Injection Volume: 1.0 µL
Detection: UV 230 nm

Mecoprop methyl ester





Retention (min): 12.73/13.39 (Peak1/Peak 2)
Resolution: 1.74

Method Conditions

Column: InfinityLab Poroshell 120 Chiral-T (15 cm x 4.6 mm, 2.7 µm)
Mobile phase: 30/70: Methanol/50 mM Ammonium Formate (pH 3.6)
Flow Rate: 0.5 mL/min
Temperature: 45 °C
Injection Volume: 1.0 µL
Detection: UV 230 nm

Mecoprop





Retention (min): 1.14/1.31 (Peak1/Peak 2)
Resolution: 3.59

Method Conditions

Column: InfinityLab Poroshell 120 Chiral-T (10 cm x 4.6 mm, 2.7 µm)
Mobile phase: 100/0.3 wt %: Methanol/Ammonium Trifluoroacetate
Flow Rate: 1.0 mL/min
Temperature: Ambient (23 °C)
Injection Volume: 1.0 µL
Detection: UV 220 nm

Amines

1-benzyl-2,2-diphenylethylamine





Retention (min): 2.61/2.89 (Peak1/Peak 2)
Resolution: 3.51

Method Conditions

Column: InfinityLab Poroshell 120 Chiral-CD (15 cm x 4.6 mm, 2.7 µm)
Mobile phase: 30/70: Acetonitrile/15 mM Ammonium Formate (pH 3.6)
Flow Rate: 1.0 mL/min
Temperature: Ambient (23 °C)
Injection Volume: 1.0 µL
Detection: UV 220 nm

1-benzyl-2,2-diphenylethylamine





Retention (min): 1.98/3.15 (Peak1/Peak 2)
Resolution: 9.47

Method Conditions

Column: InfinityLab Poroshell 120 Chiral-V (10 cm x 4.6 mm, 2.7 µm)
Mobile phase: 100/0.1 wt %: Methanol/Ammonium Trifluoroacetate
Flow Rate: 1.0 mL/min
Temperature: Ambient (23 °C)
Injection Volume: 1.0 µL
Detection: UV 220 nm

1,1'-binaphthyl-2,2'-diamine





Retention (min): 3.44/3.77 (Peak1/Peak 2)
Resolution: 2.11

Method Conditions

Column: InfinityLab Poroshell 120 Chiral-CD (10 cm x 4.6 mm, 2.7 µm)
Mobile phase: 90/10/0.3/0.2: Acetonitrile/Methanol/Trifluoroacetate/TEA
Flow Rate: 1.0 mL/min
Temperature: Ambient (23 °C)
Injection Volume: 1.0 µL
Detection: UV 254 nm

1,1'-binaphthyl-2,2'-diamine



Retention (min): 2.64/2.79 (Peak1/Peak 2)
Resolution: 1.31

Method Conditions

Column: InfinityLab Poroshell 120 Chiral-V (10 cm x 4.6 mm, 2.7 µm)
Mobile phase: 90/10: Methanol/15 mM Ammonium Formate (pH 3.6)
Flow Rate: 1.0 mL/min
Temperature: Ambient (23 °C)
Injection Volume: 1.0 µL
Detection: UV 220 nm



Role of the LC Instrument

Instrument considerations

Column compatibility with LC instruments

Extra column volume

Proper connections and fittings

Importance of cell volume


Data collection rate

Designed along with your LC instruments for highest performance


SPP Particle	For	Maximum Pressure	Typical Pressure	Efficiency	Target System
1.9 μm	Highest UHPLC performance	1300 bar	Similar to sub-2 μm totally porous	~120% of sub-2 μm totally porous	1290 Infinity II
2.7 μm	UHPLC performance at lower pressures	600 bar/ 1000 bar	50% of sub-2 μm totally porous	~90% of sub-2 μm totally porous	1290 Infinity II 1260 Infinity II
4 μm	Improved HPLC performance	600 bar	Typically <200 bar	~200% of 5 μm totally porous	1260 Infinity II VL 1220 Infinity II (VL)

Particle Size	ID	Optimum Flow	Column length	Recommended Use
1.9 μm	2.1 mm	0.4 – 0.5 mL/min	50	High speed
	3.0 mm	0.8 – 1 mL/min	100	High resolution
2.7 μm	2.1 mm	0.4 – 0.5 mL/min	>=150	Ultra-high resolution
	3.0 mm	0.8 – 1 mL/min		
	4.6 mm	1.5 – 2 mL/min		
4 μm	3.0 mm	0.5 – 0.75 mL/min		
	4.6 mm	1 – 1.25 mL/min		


Full Compatibility




Full Compatibility





Full Compatibility



Full Compatibility






Poor fitting connections



- Will broaden or split peaks, or cause tailing
- Will typically affect all peaks, but especially early eluting peaks
- Can cause carryover

Good


Properly fitted tubing, no dead volume




Importance of Correct Connections



Quick Turn



Quick Connect




Correct connection every time

Compatible to 1300 bar



Agilent technical note: Agilent InfinityLab UHPLC fittings

Publication number: 5991-5525EN

Optimizing Connecting Tubing Volume for UHPLC Columns



Length	10 mm	50 mm	100 mm	150 mm
Tubing id	Volume	Volume	Volume	Volume
0.17 mm (green)	0.227 μ L	1.1 μ L	2.27 μ L	3.3 μ L
0.12 mm (red)	0.113 μ L	0.55 μ L	1.13 μ L	1.65 μ L



Resolution and Importance of Flow Cell Volume

Differences in detector flow cell volume can affect N and R_s

Scenario: Agilent ZORBAX Rapid Resolution column: 75 mm, 3.5 μm ; flow rate: 1mL/min; $k = 3$


Flow Cell Volume	Band Broadening* (4.6 mm)	Band Broadening* (2.1 mm**)
1.7 μL	0.3%	6%
8 μL	6%	138%
14 μL	19%	423%

*Versus 8571 theoretical plates (HPLC Calculations Assistant, Version 2.1, Savant Audiovisuals)

**Flow rate, 0.2 mL/min

Maintaining Resolution at High Analysis Speed

Importance of data collection rate for narrow peaks



80Hz versus 10Hz (20Hz) Data Rate

- Peak Width: **- 55%** (- 30%)
- Resolution: **+ 90%** (+ 30%)
- Peak Capacity: **+ 120%** (+ 40%)
- App. Column Eff.: **+ 260%** (+ 70%)

Data Rate	Peak Width	Resolution	Peak Capacity
80 Hz	0.300	2.25	60
40 Hz	0.329	2.05	55
20 Hz	0.416	1.71	45
10 Hz	0.666	1.17	29
5 Hz	1.236	0.67	16

Sample: Phenones Test Mix
Column: Zorbax SB-C18, 4.6x30, 1.8um
Gradient: 50-100%ACN in 0.3min
Flow Rate: 5ml/min

Full Compatibility

Designed to seamlessly integrate into the InfinityLab family

Part of the
InfinityLab
family

Agilent **InfinityLab** products are designed to provide the best efficiency in your liquid chromatography workflow, regardless of application area. When relying on Agilent **InfinityLab instruments**, **LC columns** and **supplies** be assured that every part works together seamlessly.

InfinityLab LC columns



InfinityLab supplies and parts



InfinityLab LC series



InfinityLab StaySafe caps



InfinityLab accessories



Agilent Resources for Support

- Resource page <http://www.agilent.com/chem/agilentresources>
 - InfinityLab Poroshell 120 brochure – https://www.agilent.com/cs/library/brochures/5991-8750EN_InfinityLab_Poroshell120_brochure.pdf
 - Quick reference guides, product catalogs
 - Online selection tools, “How-to” videos
 - Column user guides - <https://www.agilent.com/en-us/support/liquid-chromatography/kb005965>
 - Tech support: <http://www.agilent.com/chem/techsupport>
- InfinityLab LC Supplies catalog ([5991-8031EN](#))
- Agilent University <http://www.agilent.com/crosslab/university>
- YouTube – [Agilent Channel](#)
- Your local product specialists
- Subscribe to Agilent Peak Tales podcasts at peaktales.libsyn.com



Contact Agilent Chemistries and Supplies Technical Support



Available in the USA and Canada 8-5 all time zones

1-800-227-9770 option 3, option 3:

Option 1 for GC and GC/MS columns and supplies

Option 2 for LC and LC/MS columns and supplies

Option 3 for sample preparation, filtration and QuEChERS

Option 4 for spectroscopy supplies

Option 5 for chemical standards



gc-column-support@agilent.com

lc-column-support@agilent.com

spp-support@agilent.com

spectro-supplies-support@agilent.com

chem-standards-support@agilent.com

Thank you for attending



Any questions?

Brochures



Poroshell 120
Portfolio brochure
5991-8750EN



Poroshell 120
ordering guide
5991-9123EN


More Chemistries, More Choices for Your Toughest Separation Challenges
Agilent InfinityLab Poroshell 120 columns



Efficiently separate the widest variety of columns
Based on superficially porous particle technology, Agilent InfinityLab Poroshell 120 columns feature a solid silica core and a porous outer layer. Compared to traditional fully porous particles of the same particle size, Poroshell 120 offers higher chromatographic efficiency and faster flow, high-resolution separations. The Poroshell 120 family includes 14 chemistries including new phases for chiral and HILIC separations. This ordering guide will help you choose the right particle size and bonded phase to achieve optimal separation for your samples.



Posters




Poroshell 120 selectivity overview
5991-9013EN



ZORBAX selectivity overview
5994-2212EN

Application Notes



Agilent Application Finder
>200 Poroshell 120 application notes


Selection Guide



Agilent InfinityLab Poroshell 120 LC Columns:
Novel Superficially Porous Columns


Comparison of Two Porous Particles

Novel Superficially Porous Particle vs. Traditional Totally Porous Particle



Poroshell 120
general selection
guide
5994-1982EN

HPLC Chemistry Selection: Poroshell 120




Legend

- Decision
- 1st choice
- Alternative

¹ C8 phases have nearly identical selectivity compared to C18. Suitable for faster analysis times when sufficient resolution and retention is achieved.

² Bare silica column. Shipped in RP solvents. For NP use, please properly flush out shipping solvent before.

Agilent InfinityLab Poroshell 120 LC columns
 General selection guide
 Publication number 5994-1982EN




Legend

- Decision
- 1st choice
- Alternative

¹ C8 phases have nearly identical selectivity compared to C18. Suitable for faster analysis times when sufficient resolution and retention is achieved.

Agilent InfinityLab Poroshell 120 LC columns
 General selection guide
 Publication number 5994-1982EN

Particle Size and Dimension: Poroshell 120



Particle Size	ID	Optimum Flow
1.9 μm	2.1 mm	0.4 – 0.5 mL/min
	3.0 mm	0.8 – 1 mL/min
	4.6 mm	1.5 – 2 mL/min
2.7 μm	2.1 mm	0.4 – 0.5 mL/min
	3.0 mm	0.8 – 1 mL/min
	4.6 mm	1.5 – 2 mL/min
4 μm	3.0 mm	0.5 – 0.75 mL/min
	4.6 mm	1 – 1.25 mL/min

Legend

- Decision** (Dark Blue Box)
- 1st choice** (Light Green Box)
- Alternative** (Blue Box)

- 1 ULD kit recommended (p/n 5067-5963)
- 2 Not possible with regulated gradient methods, not recommended for lab technicians that lack experience with UHPLC
- 3 Delay and dispersion volume. For example, 0.17 mm id tubing or bigger + 10 mm classic flow cell, valves, long tubing connections, old mixer design

Column Length	Recommended Use
50	High speed
100	High resolution
>=150	Ultrahigh resolution