

Simultaneous analysis for forensic drugs in human blood and urine using ultra-high speed LC-MS/MS

ASMS 2014 ThP-592

Toshikazu Minohata¹, Keiko Kudo², Kiyotaka Usui³,
Noriaki Shima⁴, Munehiro Katagi⁴, Hitoshi Tsuchihashi⁵,
Koichi Suzuki⁵, Noriaki Ikeda²

¹Shimadzu Corporation, Kyoto, Japan

²Kyushu University, Fukuoka, Japan

³Tohoku University Graduate School of Medicine, Sendai, Japan

⁴Osaka Prefectural Police, Osaka, Japan

⁵Osaka Medical Collage, Takatsuki, Japan

Simultaneous analysis for forensic drugs in human blood and urine using ultra-high speed LC-MS/MS

Introduction

In Forensic Toxicology, LC/MS/MS has become a preferred method for the routine quantitative and qualitative analysis of drugs of abuse. LC/MS/MS allows for the simultaneous analysis of multiple compounds in a single run, thus enabling a fast and high throughput analysis. In this study, we report a developed analytical system using ultra-high

speed triple quadrupole mass spectrometry with a new extraction method for pretreatment in forensic analysis. The system has a sample preparation utilizing modified QuEChERS extraction combined with a short chromatography column that results in a rapid run time making it suitable for routine use.

Methods and Materials

Sample Preparation

Whole blood sample preparation was carried out by the modified QuEChERS extraction method ⁽¹⁾ using Q-sep™ QuEChERS Sample Prep Packets purchased from RESTEK (Bellefonte, PA).

- 1) Add 0.5 mL of blood and 1 mL of distilled water into the 15 mL centrifugal tube and agitate the mixture using a vortex mixer.
- 2) Add two 4 mm stainless steel beads, 1.5 mL of acetonitrile and 100 µL of acetonitrile solution containing 1 ng/µL of Diazepam-d5. Then agitate using the vortex mixer.
- 3) Add 0.5 g of the filler of the Q-sep™ QuEChERS Extraction Salts Packet.

- 4) Vigorously shake the tube by hand several times, agitate well using the vortex mixer for approximately 20 seconds. Then centrifuge the tube for 10 minutes at 3000 rpm.
- 5) Move the supernatant to a different 15 mL centrifugal tube and add 100 µL of 0.1 % TFA acetonitrile solution. Then, dry using a nitrogen-gas-spray concentration and drying unit or a similar unit.
- 6) Reconstitute with 200 µL of methanol using the vortex mixer. Then move it to a microtube, and centrifuge for 5 minutes at 10,000 rpm.
- 7) Transfer 150 µL of the supernatant to a 1.5 mL vial for HPLC provided with a small-volume insert.

[ref.] (1) Usui K et al, Legal Medicine 14 (2012), 286-296




Figure 1 Scheme of the modified QuEChERS procedure

Simultaneous analysis for forensic drugs in human blood and urine using ultra-high speed LC-MS/MS

LC-MS/MS Analysis

Treated samples were analyzed using a Nexera UHPLC system coupled to a LCMS-8050 triple quadrupole mass spectrometer (Shimadzu Corporation, Japan) with LC/MS/MS Rapid Tox. Screening Database. The Database contains product ion scan spectra for 106 forensic and toxicology-related compounds of Abused drugs, Psychotropic drugs and Hypnotic drugs etc (Table 1) and

provides Synchronized Survey Scan[®] parameters (product ion spectral data acquisition parameters based on the MRM intensity as threshold) optimized for screening analysis.

Samples were separated on a YMC Triart C18 column. A flow rate of 0.3 mL/min was used together with a gradient elution.

Analytical Conditions

HPLC (Nexera UHPLC system)

Column	: YMC Triart C18 (100x2mm, 1.9µm)
Mobile Phase A	: 10 mM Ammonium formate - water
Mobile Phase B	: Methanol
Gradient Program	: 5%B (0 min) - 95%B (10 min - 13min) - 5%B (13.1 min - 20 min)
Flow Rate	: 0.3 mL / min
Column Temperature	: 40 °C
Injection Volume	: 5 µL

Mass (LCMS-8050 triple quadrupole mass spectrometry)

Ionization	: heated ESI
Polarity	: Positive & Negative
Probe Voltage	: +4.5 kV (ESI-Positive mode); -3.5 kV (ESI-Negative mode)
Nebulizing Gas Flow	: 3 L / min
Drying Gas Pressure	: 10 L / min
Heating gas flow	: 10 L / min
DL Temperature	: 250 °C
BH Temperature	: 400 °C
MRM parameter	:

Analytes	Ret. Time	Q1 m/z	Q3 m/z	Collision Energy
Diazepam-d5	9.338	290.15	154.05	-27
		290.15	198.20	-34
Alprazolam	8.646	309.10	281.10	-24
		309.10	205.10	-41
Atropine	5.378	290.15	124.15	-23
		290.15	93.20	-30
Estazolam	8.408	295.05	267.15	-24
		295.05	205.25	-37
Ethyl loflazepate	9.350	361.15	259.10	-30
		361.15	287.15	-19
Etizolam	8.786	343.05	314.10	-24
		343.05	138.15	-36
Haloperidol	8.253	376.15	165.15	-24
		376.15	123.10	-39

Analytes	Ret. Time	Q1 m/z	Q3 m/z	Collision Energy
Risperidone	7.993	411.20	191.05	-28
		411.20	69.05	-55
Triazolam	8.573	343.05	315.00	-27
		343.05	308.20	-25
Amobarbital (neg)	8.093	225.15	42.00	25
		225.15	182.00	14
Barbital (neg)	5.243	183.10	42.10	21
		183.10	140.10	15
Phenobarbital (neg)	6.762	231.10	42.20	19
		231.10	85.10	14
Thiamylal (neg)	8.883	253.00	58.10	23
		253.00	101.00	16


Simultaneous analysis for forensic drugs in human blood and urine using ultra-high speed LC-MS/MS

Type	Event#	+/-	Compound Name	m/z	Time (4.236 min - 10.415 min)
MRM	1	+	Atropine	290.15>124.15, 290.15>93.20	
MRM	2	+	Risperidone	411.20>191.05, 411.20>69.05	
MRM	3	+	Haloperidol	376.15>165.15, 376.15>123.10	
MRM	4	+	Estazolam	295.05>267.15, 295.05>205.25	
MRM	5	+	Triazolam	343.05>315.00, 343.05>308.20	
MRM	6	+	Alprazolam	309.10>281.10, 309.10>205.10	
MRM	7	+	Etizolam	343.05>314.10, 343.05>138.15	
MRM	8	+	0 Diazepam-d5	290.15>154.05, 290.15>198.20	
MRM	9	+	Ethyl loflazepate	361.15>259.10, 361.15>287.15	
MRM	10	-	Barbital (neg)	183.10>42.10, 183.10>140.10	
MRM	11	-	Phenobarbital (neg)	231.10>42.20, 231.10>85.10	
MRM	12	-	Amobarbital (neg)	225.15>42.00, 225.15>182.00	
MRM	13	-	Thiamylal (neg)	253.00>58.10, 253.00>101.00	



Figure 2 LCMS-8050 triple quadrupole mass spectrometer

Results and Discussion



Conc.	Area	Accuracy	%RSD	Conc.	Area	Accuracy	%RSD	Conc.	Area	Accuracy	%RSD	Conc.	Area	Accuracy	%RSD
0.01	9,004	112.1		0.01	4,865	114.4		0.01	29,832	108.4		0.01	3,047	107.0	
	8,288	105.1	6.57		5,109	119.9	8.71		32,436	116.7	5.14		3,064	109.2	5.63
	9,519	119.3			4,321	105.7			30,461	110.8			3,356	118.5	
0.1	75,236	89.6		0.1	48,038	84.0		0.1	335,202	91.3		0.1	27,991	94.8	
	75,983	89.6	6.04		49,152	85.1	1.82		309,273	83.7	4.74		25,542	85.7	7.83
	74,023	80.6			54,497	87.0			343,172	85.6			26,317	81.5	
1	829,519	99.9		1	604,640	103.7		1	3,826,373	102.8		1	288,776	99.0	
	831,098	99.6	2.53		581,207	99.2	2.22		3,718,854	99.4	1.66		297,332	101.5	1.96
	849,597	104.2			579,390	101.2			3,705,165	101.4			294,788	102.9	

Simultaneous analysis for forensic drugs in human blood and urine using ultra-high speed LC-MS/MS




Figure 3 Results of 8 drugs spiked in human whole blood using LCMS-8050

In this experiment, two different matrices consisting of human whole blood and urine were prepared and 18 drugs were spiked into extract solution. Calibration curves constructed in the range from 0.01 to 1 ng/mL for 12 drugs (Alprazolam, Aripiprazole, Atropine, Brotizolam, Estazolam, Ethyl loflazepate, Etizolam, Flunitrazepam,

Haloperidol, Nimetazepam, Risperidone and Triazolam) and from 1 to 100 ng/mL for 6 drugs (Bromovalerylurea, Amobarbital, Barbitol, Loxoprofen, Phenobarbital and Thiamylal). All calibration curves displayed linearity with an $R^2 > 0.997$ and excellent reproducibility was observed for all compounds ($CV < 12\%$) at low concentration level.

Simultaneous analysis for forensic drugs in human blood and urine using ultra-high speed LC-MS/MS




Figure 4 Results of 4 drugs spiked in human urine using LCMS-8050

Conclusions

- The validated sample preparation protocol can get adequate recoveries in quantitative works for all compounds ranging from acidic to basic.
- The combination of the modified QuEChERS extraction method and high-speed triple quadrupole LC/MS/MS with a simple quantitative method enable to acquire reliable data easily.