

# Application Data Sheet

## No.33

### GCMS

Gas Chromatograph Mass Spectrometer

## Analysis of Polycyclic Aromatic Hydrocarbons Using GC-MS

### Introduction

Polycyclic aromatic hydrocarbons (PAHs) are byproducts of burning fuels, such as mineral oils, and are known to be carcinogenic, mutagenic, and teratogenic. These compounds are regulated in many countries as harmful pollutants in environmental water and the atmosphere. In Japan, the "Measuring method manual for toxic air pollutants" and "Measuring method manual for polycyclic aromatic hydrocarbons (PAHs) in exhaust gas emissions" were revised in March 2011. In addition, environmental standards, including fine particle regulations (PM<sub>2.5</sub>) and fine particles, can include PAHs. Consequently, there is an increasing need to measure PAHs. This application datasheet introduces an example of measuring 36 PAHs.

### Analysis Conditions and Results

Table 1 shows the analysis conditions and Fig. 1 shows the total ion current chromatogram. Figure 2 (Page 2) shows enlarged chromatogram sections for PAHs with 4, 5, and 6 rings.

Table 1: Analysis Conditions

GC-MS	: GCMS-QP2010 Ultra	[MS]	
[GC]		Interface temperature	: 300°C
Column	: Rtx-35 (length 30 m, 0.32 mm I.D., df=0.25 μm)	Ion source temperature	: 230°C
Inlet mode	: Splitless	Solvent elution time	: 3.5min
Vaporizing chamber temp	: 300°C	Data sampling time	: 4.5 – 60min
Column oven temperature	: 90°C (2min) → (5°C/ min) → 320°C (12min)	Measurement mode	: Scan
Carrier gas	: Helium	Mass range	: m/z 45-450
Control mode	: Constant linear velocity (43.7cm/sec)	Event time	: 0.3sec
High pressure injection	: 150kPa (1.5min)		
Purge flow rate	: 3 mL/min		
Injection rate	: 1.0 μL		




Fig. 1: Total Ion Current Chromatogram of 36 PAHs

1.Naphthalene, 2.Biphenyl 3.Acenaphthylene 4.Acenaphthene 5.Flourene 6.Dibenzothiophene 7.Phenanthrene 8.Anthracene 9.4H-Cyclopenta[def]phenanthrene 10.Flouanthrene 11.Pyrene 12.Benzo[ghi]fluoranthene 13.Benzo[c]phenanthrene 14.Benzo[a]anthracene 15.Cyclopenta[cd]pyrene 16,17.Chrysene & Triphenylene 18.Benzo[b]fluoranthene 19.Benzo[k]fluoranthene 20.Benzo[j]fluoranthene 21.Benzo[a]fluoranthene 22.Benzo[e]pyrene 23.Benzo[a]pyrene 24.Perylene 25.Dibenz[a,j]anthracene 26.Dibenz[a,c]anthracene 27.Indeno[1,2,3-cd]pyrene 28.Dibenz[a,h]anthracene 29.Benzo[b]chrysene 30.Picene 31.Benzo[ghi]perylene 32.Anthanthrene 33.Dibenzo[b,k]fluoranthene 34.Dibenzo[a,h]pyrene 35.Coronene 36.Dibenzo[a,e]pyrene






Fig. 2: Enlarged Chromatograms  
Upper: 4-Ring Section (retention time 33.5 to 38.0 minutes);  
Middle: 5-Ring Section (retention time 40.5 to 44.5 minutes);  
Lower: 6-Ring Section (retention time 45.5 to 49.5 minutes)

## Summary

Because there are so many PAH isomers, separation by column is more important than separation by mass. Of the 36 polycyclic aromatic hydrocarbons, the Rtx-35 column can separate all except chrysene and triphenylene. Therefore, it enables separating PAHs with a wide range of boiling points.

