# **Packing Instructions for the Devices of the Thermo Scientific Dionex UltiMate 3000 Series**

For instructions in English, see page 1.

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To prevent the devices of the Thermo Scientific Dionex™ UltiMate™ 3000 series from being damaged during transit when they need to be returned to the factory, always use the original shipping container and observe the packing instructions in this document. We would also like to point out that inappropriate packaging automatically voids the product warranty. If the original packaging is no longer available, appropriate shipping containers and packing material can be ordered from the Thermo Fisher Scientific sales organization for Dionex HPLC products.

Follow the steps for the device that you want to ship:

|           |  |           |
|-----------|--|-----------|
| <b>1</b>  | <b>Autosamplers</b> .....  | <b>2</b>  |
| <b>2</b>  | <b>Flow Managers and Thermostatted Column Compartments</b> ..... | <b>3</b>  |
| <b>3</b>  | <b>Pumps</b> .....   | <b>4</b>  |
| <b>4</b>  | <b>Detectors</b> .....   | <b>4</b>  |
| <b>5</b>  | <b>Nano/Cap System</b> .....                                     | <b>7</b>  |
| <b>6</b>  | <b>Nano/Cap Pump</b> .....                                       | <b>9</b>  |
| <b>7</b>  | <b>Packing UltiMate 3000 Devices</b> .....                       | <b>10</b> |
|           | 7.1 Packaging with Foam Inserts .....                            | 11        |
|           | 7.2 Packaging with Foam Braces .....                             | 12        |
| <b>8</b>  | <b>Solvent Racks</b> .....                                       | <b>15</b> |
| <b>9</b>  | <b>Automated Fraction Collector</b> .....                        | <b>16</b> |
| <b>10</b> | <b>XRS Open Autosamplers</b> .....                               | <b>17</b> |

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## 1 Autosamplers WPS-3000 (all models) and ACC-3000

*For packing instructions for the OAS-3x00 Open Autosamplers refer to page 17.*

1. Shut down the autosampler as described in the Operating Instructions.

*This is a must! Verify that*


- ◆ All movable parts of the carousel (vial trays (racks), vials and well plates, and the 5-position vial holders for the 10-mL vials\* if applicable) and also the wash liquid reservoir if installed have been removed. If you have to return these parts, ship them in the accessories box. Make sure that no vials are present in the trays.  
\* For autosamplers with serial number 8013423 and higher, the 5-position holder may remain in the autosampler.

- ◆ The needle and needle arm are in the appropriate position for shipment (that is the needle arm is on the right and the needle is out of the needle seat and wash port). If they are not, move them into the correct position by either pressing the Standby button on the front of the autosampler or performing the Standby command from Chromeleon™.

2. Push the foam insert that secures the needle arm during transit into the sample compartment as far as it goes in.

*Note:*

Autosamplers with serial number 8010458 and higher are shipped with a black foam insert that secures also the carousel cover. Move the cover to the left before you push the foam insert into the sample compartment.

*Fig. 1: Foam insert securing the needle arm (left: WPS-3000PL; right: WPS-3000RS)*

3. Tilt the front panel upward and secure the panel as shown in the picture by using the two foam inserts.

Fig. 2: Foam insert for the front panel (here in a WPS-3000PL)

4. To proceed continue with the steps in Packing UltiMate 3000 Devices (→ page 10).

## 2 Flow Managers and Thermostatted Column Compartments FLM-3x00 and TCC-3x00

1. Shut down the device as described in the Operating Instructions.
2. *Only FLM-3x00 and TCC-3x00 (not TCC-3000SD and TCC-3000RS)*

Secure the front panel door. To do so, press the release button for the door and push the foam insert into the opening as far as it goes in.


Fig. 3: Securing the front panel door

3. *Only TCC-3000SD and TCC-3000RS*


*This is a must!* If the TCC is fitted with one or two column switching valves, you *have to* remove the entire valve unit (drive and drip tray) before you place the TCC in the shipping container. You may reinstall the long side panel that was installed when the TCC was shipped. If a valve unit needs to be returned, use the packaging in which the valve unit was shipped.

4. To proceed continue with the steps in Packing UltiMate 3000 Devices (→ page 10).

### 3 Pumps

#### All pumps of the UltiMate 3000 Pump Series

1. Shut down the pump as described in the Operating Instructions.
2. Tilt the front panel door upward and secure the door as shown in the picture by using the two foam inserts.

*Fig. 4: Securing the pump door (here insert installed on the left)*

3. To proceed continue with the steps in Packing UltiMate 3000 Devices (→ page 10).

### 4 Detectors

#### VWD-3x00, DAD-3000, MWD-3000, FLD-3x00, PDA-3000, UVD-3000, ECD-3000RS, Corona Veo (RS) (all models)

1. Shut down the detector as described in the Operating Instructions.
2. *Only VWD-3x00, DAD-3000, and MWD-3000 (all models)*
  - ◆ Verify that a flow cell (or dummy cell) is installed. This is to prevent that dust particles cause damage to the detector optics.
  - ◆ Secure the optics for shipment:
    - a) Push the cardboard with the information about the shipping locks under the orange shipping locks on the detector bottom as shown in the picture (→ page 5).  
(The cardboard was installed when the detector was shipped and is included in the shipping container for detectors, which can be ordered from the Thermo Fisher Scientific sales organization for Dionex HPLC products.)

- b) Tighten the shipping locks.
- c) Wrap the cardboard around the front panel door.
- d) Use tape to secure the cardboard to the detector top.


Fig. 5: Securing the detector optics (left: VWD-3x00; right: DAD/MWD-3x00)


Fig. 6: Cardboard with notice at detector front


### 3. Only FLD-3x00

Ship the FLD *without* a flow cell being installed. Instead, reinstall the flow cell cover that was installed when the detector was shipped.

Close the flow cell inlet and outlet with the plugs with which the flow cell was shipped. If a flow cell needs to be returned, use the packaging in which the flow cell was shipped.

4. *Only ECD-3000RS*

- ◆ Verify that no cell and potentiostat module is installed. This is to prevent that they damage the detector electronics during shipment.
- ◆ Tilt the front panel door of the detector upward and secure the door on the left and on the right as shown in the picture by using the two foam inserts.


*Fig. 7: Securing the detector door (here insert installed on the left)*

5. *Only Corona Veo (RS)*

- ◆ Make sure that you reinstall the protective tubing on the nebulizer that was installed when the nebulizer was shipped. Store and ship the nebulizer in its original foam packaging.
  - ◆ Reinstall the white cover on the right front side of the detector.
6. To proceed continue with the steps in Packing UltiMate 3000 Devices (→ page 10).


## 5 Nano/Cap System NCS-3500RS

1. Shut down the NCS-3500RS as described in the Operating Instructions.
2. Tilt the front panel door of the pump module upward and secure the door on the left and on the right as shown in the picture by using the two foam inserts.

*Fig. 8: Securing the pump door (here insert installed on the left)*

3. *This is a must!* If the NCS-3500RS is fitted with one or two column switching valves, you *have to* remove the valve before you place the NCS in the shipping container. To protect the valve drives, reinstall the foam inserts that were installed when the NCS was shipped. Push the side with the hole onto the drive; the plain side faces the column compartment door.

*Fig. 9: Protecting the valve drive*


If a valve needs to be returned, use the packaging in which the valve was shipped.

4. The NCS-3500RS packaging consists of a three-piece shipping container (with container bottom, spacer box, and container cover) and the actual packaging box.

*Always ship the NCS-3500RS on a pallet. This is a must!*

- a) First, place the container bottom on the pallet and place the spacer box in the bottom.

*Fig. 10: Shipping container with bottom and spacer box*

- b) Place one of the foam corners in each corner of the spacer box.

*Fig. 11: Spacer box with foam corners*

- c) Place the empty packaging box between the corners and continue with the steps in Packaging with Foam Inserts (→ page 11).
- d) Place the remaining four foam corners on the packaging box.

*Fig. 12: Packaging box with foam corners*

e) Close the shipping container with the container cover


*Fig. 13: Shipping container closed with cover*

f) Attach the shipping container on the pallet by using appropriate strapping bands.

## **6 Nano/Cap Pump NCP-3200RS**

1. Shut down the NCP-3200RS as described in the Operating Instructions.
2. Tilt the front panel door of the pump upward and secure the door on the left and on the right as shown in the picture by using the two foam inserts.

*Fig. 14: Securing the pump door (here insert installed on the left)*

3. To proceed continue with the steps in Packing UltiMate 3000 Devices (→ page 10).

## 7 Packing UltiMate 3000 Devices

*All UltiMate 3000 Devices except for Solvent Racks, Fraction Collectors, and Open Autosamplers. For information about how to pack these devices, see pages 15 and following pages.*

The shipping containers for some UltiMate 3000 devices have been changed. Follow the steps that are appropriate for your type of packaging.

- Newer shipping containers have two *foam inserts*, one of which covers the entire bottom of the container and the other one protects the device from above:

*Fig. 15: Foam insert, here covering the bottom of the shipping container*

Follow the steps in Packaging with Foam Inserts (→ page 11).

- Previous shipping containers have an outer box, *foam braces*, and an inner box.

*Fig. 16: Outer box with foam braces*

Follow the steps in Packaging with Foam Braces (→ page 12).

## 7.1 Packaging with Foam Inserts

*Note for Corona Veo detector packaging: The foam inserts of the packaging for the Corona Veo detector differ from the foam inserts shown in the pictures below. However, the packaging procedure for the detector corresponds to the procedure described in this section.*

1. Place the box on the floor, and then place one foam insert into the box.

*Fig. 17: Bottom foam insert*

2. After having shut down and prepared the device for shipment as described further up in this document, place the device on the foam insert and place the second foam insert on top of the device.

*Fig. 18: Device with foam insert on top*

3. *Only RS, SD, and BM pumps, NCP-3200RS, ECD-3000RS detectors and all VWD-3x00, DAD-3000, and MWD-3000 detectors*  
Place the additional spacer (box or foam insert) on the foam insert on top of the device.

4. Place the accessories box (with or without accessories, depending on whether you need to return them, too) on the upper foam insert and fill the free space by using the two cardboard stripes.

*Fig. 19: Accessories box*

5. Close and seal the box by using adhesive tape.

## 7.2 Packaging with Foam Braces

1. Place the larger box on the floor, place the foam braces along the left and right side in the box, and place the inner box between the braces.

*Fig. 20: Outer box with foam braces*

2. Place the foam insert on the bottom of the inner box.

*Fig. 21: Inner box with foam insert*

3. After having shut down and prepared the device for shipment as described further up in this document, place the device on the foam insert and surround the device with foam bars. The device must be protected by foam bars to the device height.

*Fig. 22: Placing the device in the box*

*Only autosamplers*

The autosampler is higher than other devices of the UltiMate 3000 series. Therefore, two foam inserts are required on each side.

*Fig. 23: Autosampler placed and protected in the box*

4. Place another foam insert on top of the device.

*Fig. 24: Top foam insert*

5. Close and seal the inner box by using adhesive tape. Place two foam braces along the edges of the inner box.

*Fig. 25: Outer box with completely packed inner box*

6. Place the accessories box (with or without accessories, depending on whether you need to return them, too) between the braces.

Accessories box

*Fig. 26: Outer box completely packed with inner and accessories boxes*

7. Close and seal the outer box by using adhesive tape.

## 8 Solvent Racks SR-3000 and SRD-3x00

1. Shut down the solvent rack as described in the Operating Instructions.
2. Push the foam braces onto the four corners of the solvent rack, and then place the rack in the inner box.

Fig. 27: Solvent rack with foam braces

3. Depending on whether the Solvent Rack is shipped with or without solvent bottles:

- ◆ Shipping the Solvent Rack *with* solvent bottles

Place the foam insert for the bottles (1) in the bottle tray of the Solvent Rack, place the bottles in the recesses provided in the insert, and then secure the bottle necks in the second foam insert (2).

Fig. 28: Securing the bottles (here for an SRD-3600)

- ◆ Shipping the Solvent Rack *without* solvent bottles

Place the foam insert for the bottles (1) in the bottle tray of the Solvent Rack, place the second foam insert (2) onto the first one.

4. Place the accessories box (with or without accessories, depending on whether you need to return them, too) on the foam insert. If necessary, fill the empty space between the Solvent Rack and the top of the box with packing material to prevent the Solvent Rack from moving inside the box during transit.

*Fig. 29: SRD packed with accessories box in the shipping container*

5. Close and seal the box by using adhesive tape.

## **9 Automated Fraction Collector AFC-3000**

1. Shut down the fraction collector as described in the Operating Instructions.
2. Before shipping the AFC, you have to remove the valve unit. If the valve unit needs to be returned, use the packaging in which the valve unit was shipped.
3. Surround the AFC with the protective wrap and place the device in the box.
4. Place the two foam spacers on the sides of the AFC as shown in the picture and fill the empty space between the fraction collection arm and the spacers with packing material to protect the arm during transit.

*Fig. 30: AFC placed in the shipping container*

5. Close and seal the box by using adhesive tape.


## 10 XRS Open Autosamplers OAS-3300TXRS and OAS-3600TXRS

1. Shut down the autosampler as described in the Operating Instructions.
2. Depends on whether you want to return the autosampler system or only a component or sub-assembly:
  - a) To return a component or sub-assembly, take the appropriate steps to remove the component or sub-assembly by following the related installation steps in the Operating Instructions in the reverse order.
  - b) To return the autosampler system, disassemble the autosampler by following the assembling steps in the Operating Instructions in the reverse order.
3. The pictures show how the main components should be arranged in the shipping boxes. Seal the boxes by using adhesive tape.

**Important:** If loose material (foam or similar items) is used to fill up a shipping box, it is mandatory to protect the components and sub-assemblies in a plastic bag. Removing the loose material and cleaning the components is very time consuming.

- ◆ For the basic autosampler system (that is, the x-, y-axes assembly, control terminal, injection unit (z-axis), basic autosampler power supply, safety guard, connecting cables, and miscellaneous accessory parts) - see further down.
- ◆ For the stack cooler - see page 19
- ◆ For the DLW option - see page 19
- ◆ For the table and other accessories - see page 20

*Shipping box — basic autosampler system (usually box 1)*

*Fig. 31: Packing the shipping box (basic autosampler) - 1*

| No. | Description                                     |
|-----|---|
| 1   | Foam inserts                                    |
| 2   | Box for power supply (basic autosampler system) |
| 3   | x-, y-axes assembly                             |


Fig. 32: Packing the shipping box (basic autosampler) - 2

| No. | Description  |
|-----|--|
| 1   | Injection unit (z-axis)  |
| 2   | Foam insert holding the injection unit   |
| 3   | Space for accessories box  |
| 3   | 3-valve housing with injection valve   |
| 4   | Space for miscellaneous connecting cables and additional autosampler leg (here hidden by the cables) |


Fig. 33: Packing the shipping box (basic autosampler) - 3

| No. | Description                           |
|-----|---------------------------------------|
| 1   | Control terminal                      |
| 2   | Mounting bracket for control terminal |
| 3   | Safety guard                          |
| 4   | Sample trays (VT54)                   |
| 5   | Reservoir tray                        |

*Shipping box — Stack cooler*


Fig. 34: Packing the shipping box (stack cooler)

| No. | Description                             |
|-----|---|
| 1   | Empty box (for packaging purposes only) |
| 2   | Foam inserts                            |
| 3   | Stack cooler                            |
| 4   | Box for power supply (stack cooler)     |
| 5   | Miscellaneous foam inserts              |

*Shipping box — DLW (Dynamic Load and Wash)*



Fig. 35: Packing the shipping box (DLW)

| No. | Description                   | No. | Description  |
|-----|-------------------------------|-----|--------------|
| 1   | Miscellaneous DLW accessories | 2   | Syringe      |
| 3   | Wash liquid reservoirs        | 4   | Wash station |
| 5   | Syringe holder assembly       | 6   | Pump module  |

*Shipping boxes — Table and Additional Accessories*

Use the shipping boxes in which the components were delivered. If loose material (foam or similar items) is used to fill up a shipping box, it is mandatory to protect the components and sub-assemblies in a plastic bag. Removing the loose material and cleaning the components is very time consuming.

*Fig. 36: Packing the shipping box (table)*

| No.   | Description   |
|-------|---|
| 1     | Shipping box (bottom part)  |
| 2+3   | Foam inserts to protect the table edges<br>The foam inserts are slightly different, depending on the table edge for which they are intended (front left and rear right edges (= nos. 2) and for front right and rear left edges (= nos. 3) of the assembled table, view from the front) |
| 4     | Tabletop (to be placed upside down onto the foam inserts)   |
| 5     | Foam insert to protect the table legs   |
| 6 + 7 | Table legs<br>(front left and rear right legs (= nos. 6) and front right and rear left legs (= nos. 7) of the assembled table, view from the front).  |
| 8     | Shipping box (top part)   |

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