

Method Development 101: From Beginner to Expert Part 1

From column selection to the first injection

Melissa Goodlad, Ph.D.
CSD Applications Engineer
February 20, 2024



Method Development 101

Agenda

Key Terminology

Method Goals

Column Selection

Mobile Phase Selection

Flow Rates & Injection Volumes

Scouting Gradient

Q&A





Method Development 101

Key Terminology



Key Terminology


Chromatographic process



Key Terminology


Chromatographic process

What makes one separation better than the other?



Good separation


How do we design methods that result in consistent, good separations?



Less good separation

Key Terminology


Peak anatomy



Key Terminology

Quantitative resolution

Resolution describes the ability of a column to separate the peaks of interest



$$R_s = \frac{t_{r2} - t_{r1}}{1/2 \cdot (W_{b2} + W_{b1})}$$

t_r Retention time compound


$W_{1/2}$ Peak width at half height

W_b Peak width at baseline

Key Terminology

Resolution

Resolution describes if baseline separation was achieved



- $R_s = 1$ separation of two peaks
- $R_s \geq 1.5$ Baseline separation and exact quantitation is possible
- $R_s \geq 1.7$ desirable for robust method

Key Terminology

Resolution: Influencing factors


$$R_s = \underbrace{\frac{1}{4} \sqrt{N}}_{\text{Efficiency}} \cdot \underbrace{\left(\frac{\alpha - 1}{\alpha} \right)}_{\text{Selectivity}} \cdot \underbrace{\left(\frac{k}{1 + k} \right)}_{\text{Retention}}$$

Improve resolution by improving any of these parameters:

- **Efficiency** describes the separation power of the column
- **Selectivity** has the highest influence on the resolution. Small changes in selectivity can lead to big changes in resolution
- **Retention** has only a significant influence at small k values

Key Terminology

Resolution: Influencing factors



Note:

Selectivity has highest impact on resolution

$$R_s = \underbrace{\frac{1}{4} \sqrt{N}}_{\text{Efficiency}} \cdot \underbrace{\left(\frac{\alpha - 1}{\alpha} \right)}_{\text{Selectivity}} \cdot \underbrace{\left(\frac{k}{1 + k} \right)}_{\text{Retention}}$$

[LC-Handbook-Complete-2.pdf \(agilent.com\)](https://www.agilent.com/LC-Handbook-Complete-2.pdf)

Key Terminology

Resolution: Retention factor

$$k = \frac{(t_R - t_0)}{t_0}$$


t_R = retention time for sample peak

t_0 = retention time for unretained peak

The **retention factor** measures the period of time that the sample component resides in the stationary phase relative to the time it resides in the mobile phase.

Main parameter affecting retention:
Mobile phase

Influence of Mobile Phase:



Key Terminology

Resolution: Selectivity

$$\alpha = \frac{k_2}{k_1}$$

α = selectivity

k_1 = retention factor of 1st peak


k_2 = retention factor of 2nd peak

Selectivity is a measure of the time or distance between the maxima of two peaks. If $\alpha = 1$, the two peaks have the same retention time and co-elute.

Main parameters influencing selectivity:

- Stationary phase
- Mobile phase

Influence of Stationary Phase:



Key Terminology

Resolution: Efficiency

$$N = 16(t_R/w_b)^2$$

N = efficiency

t_R = retention time

w_b = peak width at base

Columns with high plate numbers are more efficient. A column with a high N will have a narrower peak at a given retention than a column with a lower N number.

$$N \propto \frac{L}{d_p}$$

N = efficiency

L = column length

d_p = particle size


Parameters influencing column efficiency:

- Column length (increasing column length increases efficiency)
- Particle size (decreasing particle size increases efficiency)

Key Terminology

Efficiency: Van Deemter equation

- The Van Deemter curve describes the relationship between the mean linear velocity u of the mobile phase and the plate height H or plate number N .
- Minimum (H_{\min}) of the curve represents the optimal linear velocity u_{opt} .
- Van-Deemter curve is always analyte and method specific.




Lower h (reduced plate height) = higher efficiency $h = L/N$

Key Terminology

Efficiency: Van Deemter equation

- **A term:** eddy diffusion and flow distribution
 - particle size & packing quality important
 - narrow particle size distribution
- **B term:** longitudinal diffusion
 - Diffusion in the mobile phase
- **C term:** mass transfer
 - shorter diffusion paths
 - better with superficially porous particles
 - more effect on large molecules
- u : linear velocity
 - velocity of mobile phase through column
 - $u = L/t_0$ in cm/sec




| H | A | B | C | u |
|-----|-----|--------------------|-----|------------------|
| m | m | $m^2 \cdot s^{-1}$ | s | $m \cdot s^{-1}$ |

Key Terminology

Efficiency: Van Deemter

Eddy Diffusion



Note:

By flowing around the packaging vortices (eddy) are created.


→ Analyte molecules are carried through the column in different ways.

→ Molecules will reach detector at different time points!

Key Terminology

Efficiency: Van Deemter

Longitudinal Diffusion



Note:


- Can be observed in any pipeline with a laminar flowing fluid!
- Flow velocity is higher in the center than at the edges.
- Molecules migrating in the center are eluted faster than those migrating at the edge.
- Diffusion of the sample zone!

Key Terminology

Efficiency: Van Deemter

Mass Transfer

Diffusion path differences



Note:


- Equilibrium/mass transfer between mobile and stationary phase!
- Mass exchange takes time, so it is favored by low flow rates..
- Band broadening due to low flow through pores.

Key Terminology

Efficiency: Van Deemter

Mass Transfer

- Small particles lead to lower plate heights and therefore higher separation efficiency
- For smaller particles, the separation efficiency suffers less when increasing the flow



Key Terminology

Summary

| Condition | Retention (κ) | Selectivity (α) | Efficiency (N) |
|--|------------------------|--------------------------|----------------|
| % B | ● ● | ● | — |
| B-solvent (acetonitrile, methanol, etc.) | ● | ● ● | — |
| Temperature | ● | ● | ● |
| Column type (C18, phenol, etc.) | ● | ● ● | — |
| Mobile phase pH | ● ● | ● ● | ● |
| Buffer concentration | ● | ● | — |
| Ion-pair-reagent concentration | ● ● | ● ● | ● |
| Column length | NA | NA | ● ● |
| Particle size | NA | NA | ● ● |
| Flow rate | NA | NA | ● |

| Symbol | Meaning |
|--------|--|
| ● ● | Major effect |
| ● | Minor effect |
| — | Relatively small effect |
| blue | Conditions that are primarily used to control variable |

Method Development 101: Establishing Method Goals



Defining Method Objectives

- Is the primary goal detection, characterization or isolation of purified material?
- Is it necessary to resolve all sample components?
- What levels of accuracy and precision are required.
- How many samples will be analyzed at one time (Throughput)
- What HPLC equipment and operator skills are present in the laboratory that will use the final method?



Defining Method Objectives

Examples of common separation goals and method performance criteria

Good System Suitability Parameters

- Resolution: ≥ 2
- Peak shape: USP T_f close to 1 (< 2)
- Injection Repeatability: areas, T_f , etc. (RSD 0.1 - 0.25%)
- Absolute retention factors: $1 < k < 10$
- Relative Retention: α or k_2/k_1
- Signal-to-Noise Ratio: > 10


Avoid these for system suitability criteria:

- Column efficiency (theoretical plates)
- Absolute retention time

Method Performance Criteria

- Accuracy
- Precision
 - Ruggedness
 - Robustness
- Selectivity/Specificity
- Linearity
- Range
- Quantitation Limit (LOQ, 10x S/N)
- Detection Limit (LOD, 3x S/N)


Method Development 101: Column Selection



Column Selection


Selectivity

- Bonded phase affects selectivity (α)
- Different interactions for polar and non-polar compounds.
- Exploit other interactions with bonded phase (e.g., pi-pi)
- Changing the bonded phase can improve selectivity/resolution, reduce analysis time
- *Evaluating different bonded phase chemistries early can save time in optimization and generate a more robust method*




Column Selection

Overview



Column Selection

Base chemistry: silica





Column Selection

Column lifetime: low pH

Low pH methods

- Breaking of siloxane bond reduces column lifetime, especially with short chain alkyl ligands





Kirkland, J.J., J.L. Glajch, and R.D. Farlee, *Analytical Chemistry* (1989), 61, 2.

Column Selection

Column lifetime: high pH





High pH methods

- High pH methods can lead to the dissolution of the silica stationary phase
- Double end-capping and/or bidentate bonding improves peak shape of basic compounds and high pH stability



Column Selection


Stationary phase chemistry



| Chemistry | Alkyl | Alkyl polar | Phenyl, diphenyl, phenyl hexyl | Pentafluoro phenyl (PFP) | Cyano |
|---------------------------|---|---|---|---|---|
| Example | SB-C18, EC-C18 | Bonus RP, Aq-C18 | Phenyl, Phenyl hexyl | PFP | EC-CN |
| Hydrophobicity | ●●●● | ●●●● | ●●●● | ●●●● | ●●● |
| π - π interaction | — | — | ●●● (donor) | ●●● (acceptor) | ● |
| Dipole-dipole | — | ●● | ● | ●●●● | ●●● |
| Hydrogen bonding | ● | ●●●● | ●● | ●●● | ●● |
| Applicability | <ul style="list-style-type: none"> General purpose | <ul style="list-style-type: none"> Enhanced retention of polar analytes while also separating non-polar analytes | <ul style="list-style-type: none"> Alternate selectivity with aromatic and moderately polar groups | <ul style="list-style-type: none"> Alternate selectivity for halogenated, polar, and isomeric analytes Excellent peak shape for polar and non-polar compounds | <ul style="list-style-type: none"> Alternate selectivity for polar and mid-polar compounds |

Column Selection




Stationary phase chemistry



Mobile Phase 40 % ACN 60 % 25 mM Sodium Phosphate Buffer pH= 2.4 Flow Rate= 1.5 ml/min 4.6 x 50mm UV 210 nm

Column Selection

Stationary phase chemistry




1. Aniline
2. o-toluidine
3. Anisidine
4. Chloroaniline
5. Benzidine
6. Naphthylamine
7. o-tolidine
8. Dimethoxybenzidine
9. Dichlorobenzidine

Mobile phase A: 10 mM ammonium acetate, pH 4.7, Mobile phase B: methanol. 25% to 90% MPB in 9 minutes


Column Selection

Stationary phase chemistry

Moderate Sample Load: 30 ng amitriptyline



High Sample Load: 500 ng amitriptyline



Improved Sample Loading and Peak Shape with InfinityLab Poroshell 120 CS-C18 Columns (agilent.com)

Column Selection

Pore size


As a general rule, the pore size should be 3X the hydrodynamic radius of your analyte

Small molecules

- 80 – 120 Å
- Maximizes loading and retention




Peptides, proteins, other large biomolecules

- 120 Å (Peptides, small oligonucleotides)
- 300 Å to 450 Å (Proteins, mAb)
- 1000 Å (Larger proteins, larger oligonucleotides)
- 4000 Å (mRNA, pDNA, VLP)
- Maintain high efficiency



Column Selection


Particle size



Increasing resolution, pressure

Column Selection

Particle size



Column Selection


Particle Size

Original method

ZORBAX LC column Extend-C18

4.6 x 150 mm, 5 μ m

3 μ L inj




Updated method

ZORBAX RRHT column Extend C18

4.6 x 50 mm, 1.8 μ m

1 μ L inj.



Mobile phase: (70:30) MeOH: 50 mM pyrrolidine buffer Flow = 1.0 mL/min, Temp. : ambient

Column Selection

Particle Size



< 400 bar

InfinityLab Poroshell 120 4 μm

ZORBAX 3.5 & 5 μm

HPLC

400-800 bar

InfinityLab Poroshell 120 2.7 μm

ZORBAX RRHT 1.8 μm

UHPLC

800-1300 bar

Infinity Lab Poroshell 120 1.9 μm

ZORBAX RRHD 1.8 μm

Low dispersion UHPLC

Pressure

low

Moderate

High

Sample prep

low

SPP low / FPP moderate

High

Loadability

High


SPP medium / FPP high

SPP Low / FPP medium

Column Selection


Totally porous particles (TPP) vs. superficially porous particles (SPP)

- Analytes travel through the particle more efficiently
- High efficiency allows you to use a larger SPP (i.e., 2.7 μ m) for nearly equivalent performance to a smaller sub-2 μ m (STM) TPP column
- Using a larger particle allows for lower backpressure than comparably efficient totally porous STM columns and flexible use on HPLC or UHPLC systems



Column Selection

TPP vs. SPP



Column Selection


TPP vs. SPP

| SPP particle | For | Maximum pressure | Typical pressure | Efficiency | Target system |
|-------------------|--------------------------------------|--------------------|---|---|--|
| 1.9 μm | Highest UHPLC performance | 1300 bar | Similar to sub-2 μm totally porous | ~120% of sub-2 μm totally porous | 1290 Infinity II |
| 2.7 μm | UHPLC performance at lower pressures | 600 bar / 1000 bar | 50% of sub-2 μm totally porous | ~90% of sub-2 μm totally porous | 1290 Infinity II 1260 Infinity II |
| 4 μm | Improved HPLC performance | 600 bar | Typically < 200 bar | ~200% of 5 μm totally porous | 1260 Infinity II VL 1220 Infinity II (VL) |

Highend LC-Systems


1290 Infinity II LC

High-Speed Flexible



Core LC Modular

1260 Infinity II LC



Core LC Integr.

1220 Infinity II LC

Gradient Isocratic



Bio-inert LC

1260 Infinity II Bio-inert
1260 Infinity II Bio-SEC




Column Selection


TPP vs. SPP

The InfinityLab Poroshell column has 90 percent the efficiency and half the pressure of the ZORBAX column

InfinityLab Poroshell 120 EC-C18
2.7 μm , 3.0 x 100mm



ZORBAX Eclipse Plus C18
1.8 μm , 3.0 x 100mm



Column Selection

Column dimensions

Inner diameter (ID)

| ID | Optimum Flow Rate | Recommended Use |
|-----------|--------------------------|-------------------------------------|
| 4.6 mm | 1.00-1.25 mL/min | Legacy methods |
| 3.0 mm | 0.8-1.0 mL/min | Lower solvent use |
| 2.1 mm | 0.4-0.5 mL/min | MS applications, lowest solvent use |

Column length

| Column Length | Recommended Use |
|----------------------|------------------------|
| 50 mm | High throughput |
| 100 mm | High resolution |
| ≥150 mm | Ultra-high resolution |



Column Selection

Method develop kits

Column kits for method development

Simplify method development by offering one of **26** predefined method development kits that make method development easy!

- Particle Technology: Poroshell / Zorbax
- Different selectivities (classic RP / aqueous)
- pH ranges
- Column size options: 2.1 x 50 mm, 4.6 x 50 mm, and 4.6 x 100 mm with 1.8, 3.5 or 5 μm particle sizes

[Method development column kit Webpage](#)

[Method development column kit flyer](#)



Method Development Kits for HPLC | Agilent

Column Selection

Method develop kits













Poroshell Method Development Column Kits

| Product number | Description | Particle size |
|----------------|--|---------------|
| 5190-6155 | Poroshell 120 Selectivity Method Development Kit, includes Poroshell 120 EC-C18, Phenyl-Hexyl, Bonus RP columns , 2.1 x 50 mm | 2.7 µm |
| 5190-6156 | Poroshell 120 Selectivity Method Development Kit, includes Poroshell 120 EC-C18, Phenyl-Hexyl, Bonus RP columns , 4.6 x 50 mm | 2.7 µm |
| 5190-6157 | Poroshell 120 Aqueous Method Development Kit, includes Poroshell 120 SB-Aq, Phenyl-Hexyl, Bonus RP columns , 2.1 x 50 mm | 2.7 µm |
| 5190-6158 | Poroshell 120 Aqueous Method Development Kit, includes Poroshell 120 SB-Aq, Phenyl-Hexyl, and Bonus RP columns , 4.6 x 50 mm | 2.7 µm |
| 5190-6159 | Poroshell 120 L1, L7, and L10 USP Method Development Kit, includes Poroshell 120 EC-C18, EC-C8, EC-CN columns , 4.6 x 100 mm | 2.7 µm |
| 5190-6160 | Poroshell 120 L1, L7, and L10 USP Method Development Kit, includes Poroshell 120 EC-C18, EC-C8, EC-CN columns , 3.0 x 100 mm | 2.7 µm |

Method Development Kits for HPLC | Agilent














Column Selection

Superficially porous column (SPP) selection poster

| InfinityLab Poroshell 120 | Chemistry | Particle Sizes | Pore Size | Temperature Limit | pH Range | Endcapped | Carbon Load | Surface Area | USP Designation | Benefits and Applications |
|---------------------------|---|----------------------|-----------|-------------------|----------|-----------|-------------|-----------------------|-----------------|---|
| EC-C18 |  | 1.9 μm, 2.7 μm, 4 μm | 120 Å | 60 °C | 2.0–8.0 | Yes | 10% | 130 m ² /g | L1 | General purpose Excellent peak shape and efficiency for acids, bases, and neutrals |
| EC-C8 |  | 1.9 μm, 2.7 μm, 4 μm | 120 Å | 60 °C | 2.0–8.0 | Yes | 5% | 130 m ² /g | L7 | General purpose Lower retention of hydrophobic analytes vs. C18 |
| Aq-C18 |  | 2.7 μm | 120 Å | 90 °C | 1.0–8.0 | Yes | Proprietary | 130 m ² /g | L1 | Enhanced retention for challenging polar compounds while also separating non-polar analytes 100% aqueous mobile phase compatibility and low pH stability |
| SB-C18 |  | 1.9 μm, 2.7 μm, 4 μm | 120 Å | 90 °C | 1.0–8.0 | No | 9% | 130 m ² /g | L1 | Excellent stability at low pH Great peak shape in highly acidic conditions |
| SB-C8 |  | 2.7 μm | 120 Å | 80 °C | 1.0–8.0 | No | 5.5% | 130 m ² /g | L7 | Excellent stability at low pH Lower retention of hydrophobic analytes vs. C18 |
| HPH-C18 |  | 1.9 μm, 2.7 μm, 4 μm | 100 Å | 60 °C | 2.0–11.0 | Yes | Proprietary | 95 m ² /g | L1 | High pH capability designed for longest lifetime, especially under high pH conditions Robust performance and long lifetimes Similar selectivity compared to EC-C18 |
| HPH-C8 |  | 2.7 μm, 4 μm | 100 Å | 60 °C | 2.0–11.0 | Yes | Proprietary | 95 m ² /g | L7 | High pH capability Robust performance and long lifetimes Lower retention of hydrophobic analytes vs. C18 |
| CS-C18 |  | 2.7 μm | 100 Å | 90 °C | 1.0–11.0 | Yes | Proprietary | 95 m ² /g | L1 | High pH capability with alternate selectivity Improved peak shape and sample capacity for basic compounds with low ionic strength mobile phases |
| Bonus-RP |  | 2.7 μm | 120 Å | 60 °C | 2.0–8.0 | Yes | 9.5% | 130 m ² /g | L60 | Alternate selectivity to C18 Unique selectivity due to a polar embedded group, stable in 100% aqueous |
| PFP |  | 1.9 μm, 2.7 μm, 4 μm | 120 Å | 60 °C | 2.0–8.0 | Yes | 5.1% | 130 m ² /g | L43 | Alternate selectivity Excellent peak shape for polar and nonpolar analytes Unique selectivity for aromatic and halogenated compounds |
| Phenyl-Hexyl |  | 1.9 μm, 2.7 μm, 4 μm | 120 Å | 60 °C | 2.0–8.0 | Yes | 9% | 130 m ² /g | L11 | Alternate selectivity with aromatic groups Highly nonpolar bonded phase takes advantage of pi-pi interactions |
| SB-Aq |  | 1.9 μm, 2.7 μm, 4 μm | 120 Å | 80 °C | 1.0–8.0 | No | Proprietary | 130 m ² /g | L96 | Alternate selectivity Excellent peak shape and retention of polar compounds using reversed-phase LC Exceptional stability under high-aqueous conditions, including 100% water |

Column Selection

Totally porous column (TPP) selection poster

| Agilent ZORBAX | Chemistry | Particle Sizes | Pore Size (Å) | Temperature Limit | pH Range | Endcapped | Carbon Load (%) | Surface Area | USP Designation | Benefits and Applications | |
|---------------------------|---|---|----------------|-------------------|----------|-----------|-----------------|-----------------------|-----------------------|---|---|
| Eclipse Plus C18 |  | 1.8, 3.5, 5 | 95 | 60 °C | 2–9 | Double | 9 | 160 m ² /g | L1 | General purpose Starting Point for LC method development | |
| Eclipse Plus C8 |  | 1.8, 3.5, 5 | 95 | 60 °C | 2–9 | Double | 7 | 160 m ² /g | L7 | General purpose Lower retention of hydrophobic analytes vs. C18 | |
| Eclipse Plus Phenyl-Hexyl |  | 1.8, 3.5, 5 | 95 | 60 °C | 2–8 | Double | 9 | 160 m ² /g | L11 | Alternative selectivity for aromatic compounds Enhanced pi-pi interactions when using methanol | |
| Eclipse Plus PAH | Polymeric C18 | 1.8, 3.5, 5 | 95 | 60 °C | 2–9 | Double | 14 | 160 m ² /g | L1 | Application-specific Designed for the separation of PAHs in LC | |
| Eclipse XDB C18 |  | 1.8, 3.5, 5 | 80 | 60 °C | 2–9 | Double | 10 | 180 m ² /g | L1 | General purpose, higher carbon load Higher hydrophobicity with alternative selectivity for lipophilic analytes | |
| Eclipse XDB C8 |  | 1.8 (RRHT) 3.5, 5, 7 | 80 | 60 °C | 2–9 | Double | 7.6 | 180 m ² /g | L7 | General purpose, higher carbon load Higher hydrophobicity with alternative selectivity for lipophilic analytes but reduced retention vs. XDB-C18 | |
| Eclipse XDB Phenyl |  | 3.5, 5 | 80 | 60 °C | 2–9 | Double | 7.2 | 180 m ² /g | L11 | Alternative selectivity for aromatic compounds Enhanced pi-pi interactions when using methanol | |
| Eclipse XDB CN |  | 3.5, 5 | 80 | 60 °C | 2–9 | Double | 4.2 | 180 m ² /g | L10 | Polar analytes in RP, low bleed Excellent peak shape of polar and mid-polar compounds | |
| StableBond C18 |  | R ₁ =C18 | 1.8, 3.5, 5, 7 | 80 | 90 °C | 0.8–8 | No | 10 | 180 m ² /g | L1 | Low pH and high temperature Excellent stability and peak shape at highly acidic conditions |
| StableBond C8 |  | R ₁ =C8 | 1.8, 3.5, 5, 7 | 80 | 80 °C | 1–8 | No | 5.5 | 180 m ² /g | L7 | Low pH and high temperature Lower retention of hydrophobic analytes vs. C18 |
| StableBond C3 |  | R ₁ =C3 | 1.8, 3.5, 5 | 80 | 80 °C | 1–8 | No | 4 | 180 m ² /g | L56 | Low pH and high temperature Reduced retention of hydrophobic analytes |
| StableBond Aq |  | Proprietary | 1.8, 3.5, 5, 7 | 80 | 80 °C | 1–8 | No | Proprietary | 180 m ² /g | L96 | Polar analytes in RP Excellent peak shape and retention of polar compounds using reversed-phase LC, stable at 100% aqueous mobile phases |
| StableBond Phenyl |  | R ₁ =Phenylethyl | 1.8, 3.5, 5, 7 | 80 | 80 °C | 1–8 | No | 5.5 | 180 m ² /g | L11 | Alternative selectivity for aromatic compounds Enhanced pi-pi interactions when using methanol |
| StableBond CN |  | R ₁ =(CH ₂) ₆ -CN | 1.8, 3.5, 5, 7 | 80 | 80 °C | 1–8 | No | 4 | 180 m ² /g | L10 | Polar molecules at low pH or high temperature, low bleed Excellent peak shape of polar and mid-polar compounds |

A proven and reliable portfolio of totally porous HPLC columns (agilent.com)

ZORBAX Reversed-Phase Columns | Agilent

Method Development 101: Mobile Phase Selection



Mobile Phase Selection


Retention and selectivity

| Condition | Retention (κ) | Selectivity (α) | Efficiency (N) |
|--|------------------------|--------------------------|----------------|
| % B | ●● | ● | — |
| B-solvent (acetonitrile, methanol, etc.) | ● | ●● | — |
| Temperature | ● | ● | ● |
| Column type (C18, phenol, etc.) | ● | ●● | — |
| Mobile phase pH | ●● | ●● | ● |
| Buffer concentration | ● | ● | — |
| Ion-pair-reagent concentration | ●● | ●● | ● |
| Column length | NA | NA | ●● |
| Particle size | NA | NA | ●● |
| Flow rate | NA | NA | ● |

| Symbol | Meaning |
|--------|--|
| ●● | Major effect |
| ● | Minor effect |
| — | Relatively small effect |
| blue | Conditions that are primarily used to control variable |

Mobile Phase Selection

Chromatographic process




| Reversed-Phase LC | |
|-------------------|---|
| Polarity | Non-polar stationary phase (e.g., C18) |
| Mobile Phase | Polar mobile phase: H ₂ O/CH ₃ OH, H ₂ O/CH ₃ CN |
| Gradient | Decrease retention by decreasing polarity of mobile phase H ₂ O ↓ = retention ↑ CH ₃ CN ↑ = retention ↓ |
| Elution Order | polar to non-polar |

Mobile Phase Selection

Reversed-phase mobile phases

- HPLC or LCMS grade solvents
- UV transparency
- Low viscosity
- High boiling point
- Sample Solubility
- Low cost, low toxicity, non-corrosive



[Maximize Your Efficiency With Precision Solvents \(agilent.com\)](https://www.agilent.com)


Mobile Phase Selection

Reversed-phase mobile phases

HPLC or LCMS grade solvents

- Lowest impurity levels, reducing ghost peaks in gradient runs
- 0.2 μm pre-filtering safeguards system from contaminants and clogging
- Highest lot-to-lot reproducibility

Water/Methanol Gradient Overlay at 210 nm, 225 nm, and 254 nm



Gradient from 5-95% ACN. Detection wavelengths 210 nm (blue), 225 nm (turquoise), and 254 nm (purple); Range: 0-100 mAU

[Maximize Your Efficiency With Precision Solvents \(agilent.com\)](https://www.agilent.com)

Mobile Phase Selection


Reversed-phase mobile phases

UV transparency

- Mobile phase will have an absorbance of $A < 0.2$ AU at the wavelength used for detection of the sample



The Effect of Changing Wavelength on Sensitivity



The effect of changing wavelength on sensitivity of the analysis of beta blocker molecules. At lower wavelengths, the background (baseline) signal rises as the UV cutoff point of the solvent is approached

[Maximize Your Efficiency With Precision Solvents \(agilent.com\)](https://www.agilent.com)

Mobile Phase Selection

Reversed-phase mobile phases

Low viscosity

- Pressure is directly proportional to the viscosity of solvents



| Solvent | Boiling Point (°C) | Viscosity (cP) | UV cutoff (nm) |
|--------------------|--------------------|----------------|----------------|
| N-Hexane | 69 | 0.31 | 190 |
| Toluene | 78 | 0.59 | 285 |
| Methylene chloride | 40 | 0.44 | 233 |
| Tetrahydrofuran | 66 | 0.55 | 212 |
| Acetonitrile | 82 | 0.30 | 190 |
| 2-Propanol | 82 | 2.30 | 205 |
| Methanol | 65 | 0.54 | 205 |
| Water | 100 | 1.00 | <190 |

[Maximize Your Efficiency With Precision Solvents \(agilent.com\)](https://www.agilent.com)

Mobile Phase Selection

Mobile phase modifiers



Buffers

- To stabilize pHs of mobile phase (Phosphate, acetate, citrate)

Acidifiers

- To suppress ionization of acidic analytes (TFA, FA)

Ion-pairing reagents

- For separation of ionic compounds with reversed-phase methods (HFIP)

Amine modifiers

- To reduce tailing of basic analytes with reversed-phase methods (TEA)
-

Mobile Phase Selection

Buffer selection

Desired properties

- pK_a and buffer capacity
- Solubility
- UV absorbance (UV detection)
- Volatility (MS or ELSD)
- Ion-pairing properties
- Stability and compatibility with equipment



Mobile Phase Selection

Buffer selection

Desired properties

- pK_a and buffer capacity
 - Mobile phase pH should be ± 1.0 units from buffer pK_a
 - Concentration typically falls within 5 to 25 mM
 - pH of buffer should be at least one unit above or below the pK_a of the sample

| Name of Buffer | Range of pH | MS Compatible |
|------------------------------|-------------|---------------|
| Phosphate: pK_1 | 1.1-3.1 | No |
| Phosphate: pK_2 | 6.2-8.2 | No |
| Phosphate: pK_3 | 11.3-13.3 | No |
| Sodium acetate | 3.8-5.8 | No |
| Ammonium acetate (< 50 nM) | 3.8-5.8 | Yes |
| Trifluoro acetic acid (0.1%) | 2.0 | Yes |
| Phosphoric acid (0.1%) | 2.0 | No |
| Formic acid (0.1%) | 2.7 | Yes |
| Ammonium formate (< 50 nM) | 2.7-4.7 | Yes |
| Ammonium bicarbonate | 6.6-8.6 | Yes |
| TRIS | 7.3-9.3 | Yes |


Mobile Phase Selection

pH selectivity




Agilent InfinityLab Poroshell HPH-C18 4.6 x 50 mm, 2.7 µm


pH 3
10 mM HCO₂NH₄



pH 4.8
10 mM NH₄HCO₃



pH 10
10 mM NH₄HCO₃



1. Procainamide
2. Caffeine
3. Acetyl Salicylic Acid
4. Hexanophenone Deg.
5. Dipyrimadole
6. Diltiazem
7. Diflunisal
8. Hexanophenone

→ Acids
→ Bases

| Time | % Buffer | % MeCN |
|----------|----------|--------|
| 0 | 10 | 90 |
| 5 | 90 | 10 |
| 7 | 10 | 90 |
| 2 ml/min | | 254 nm |


Method Development 101: Flow Rates and Injection Volumes



Flow Rates and Injection Volumes

Flow rates


- The optimal flow rate depends on column diameter and particle size



Flow Rates and Injection Volumes

Injection volumes

- Injection volumes contribute to overall system volume
- Keep injection volumes to a minimum, while retaining solubility



Flowrates and Injection Volumes

Injection volumes

$$V_m = \pi \cdot r^2 \cdot L \cdot \sim 0.6$$

Column volume is calculated as the volume of a cylinder less the space occupied by the packing material. As an example, Agilent ZORBAX Eclipse Plus C18 packing material occupies 40% of the column, the remaining 60% of the cylinder would be considered as column volume.



| Column Dimensions (d _c x L, mm) | V _m (mL) |
|---|---------------------|
| 2.1 x 50 | 0.12 |
| 2.1 x 100 | 0.24 |
| 2.1 x 150 | 0.37 |
| 2.1 x 250 | 0.61 |
| 3.0 x 150 | 0.85 |
| 4.6 x 100 | 1.16 |
| 4.6 x 150 | 1.75 |
| 4.6 x 250 | 2.90 |


Method Development 101: Scouting gradient



Method Development 101:

Scouting gradient

- A good starting point when developing a method is a scouting gradient.
- Recommended starting conditions are 5–95% MP B with a low pH
- Gradient length is dependent on the column length




“Making to the most of a Gradient Scouting Run” LCGC North America Vol. 31, Number 1, 2013.

Method Development 101:

Scouting gradient

Table 1: Column volumes for various column configurations and recommended scouting gradient times

| L (mm) | d _c (mm) | V _m (mL) | F (mL/min) | t _G (min) |
|--------|---------------------|---------------------|------------|----------------------|
| 250 | 4.6 | 2.5 | 1 | 60 |
| 150 | 4.6 | 1.6 | 2 | 20 |
| 150 | 2.1 | 0.33 | 0.4 | 20 |
| 100 | 4.6 | 1.0 | 2 | 10 |
| 100 | 2.1 | 0.22 | 0.4 | 10 |
| 50 | 2.1 | 0.11 | 0.5 | 5 |




“Making to the most of a Gradient Scouting Run” LCGC North America Vol. 31, Number 1, 2013.

Method Development 101:


Scouting gradient

- Quick evaluation: how much of the gradient is occupied
 - $\frac{\Delta t_G}{t_G} \leq 25\%$ isocratic is recommended
 - $\frac{\Delta t_G}{t_G} \geq 40\%$ gradient is recommended

Ideal for Gradient Method



Ideal for Isocratic Method



Method Development 102



Jean Lane
Application Engineer
Agilent Technologies, Inc.

Title: HPLC Method Development: From Beginner to Expert, Part 2

Date: Thursday, March 28, 2024

Time: 11:00 AM Eastern Daylight Time

Duration: 1 hour

Method Development 102 will review and expand upon some of the 101 fundamentals as we cover advanced topics such as how to best transfer a method from one column dimension to another. We will explore reasons for why some methods are quite difficult to transfer to a different HPLC system. In addition, we will look at gradient method development and how to efficiently use a scouting gradient to quickly develop a good HPLC method.



[HPLC Method Development: From Beginner to Expert, Part 2 \(on24.com\)](https://on24.com)

Method Development 101: Conclusion



Conclusion


Tips for a robust method

- Always start method development with a new column
- Select columns with robust properties at pH of method
- Choose a quality column with long lifetimes
- Consider batch-to-batch reproducibility
- Consider scalability of particle sizes and chemistries for downstream method transfer
- Make sure mobile phase preparation is documented and transferrable



Agilent Resources for Support

- Resource page <http://www.agilent.com/chem/agilentresources>
 - Quick reference guides, product catalogs
 - Online selection tools, “How-to” videos
 - Column user guides - <https://www.agilent.com/en-us/support/liquid-chromatography/kb005965>
- Tech support: <http://www.agilent.com/chem/techsupport>
- InfinityLab LC Supplies catalog ([5991-8031EN](#))
- Agilent University <http://www.agilent.com/crosslab/university>
- YouTube – [Agilent Channel](#)
- Your local product specialists



Contact Agilent Chemistries and Supplies Technical Support



Available in the USA and Canada 8-5 all time zones

1-800-227-9770 option 3, option 3:

Option 1 for GC and GC/MS columns and supplies

Option 2 for LC and LC/MS columns and supplies

Option 3 for sample preparation, filtration and QuEChERS

Option 4 for spectroscopy supplies

Option 5 for chemical standards

gc-column-support@agilent.com

lc-column-support@agilent.com

spp-support@agilent.com

spectro-supplies-support@agilent.com

chem-standards-support@agilent.com

Questions?





Agilent

Trusted Answers