

Analysis of Fullerenes using LC-MS

Fullerenes (C_{60}) have been attracting attention ever since their structure including five-membered and six-membered rings was proposed (Fig.1). Many applications are expected for them as semiconductive or superconductive materials, or as materials for new plastics, due to their high electrical conductivity and a high degree of superconductivity when internally doped with metal.

It is known that the MALDI-TOFMS is effective in the analysis of fullerenes. The mass spectrum of a fullerene mixture obtained using the Shimadzu/Kratos AXIMA-CFR is shown in Fig.2. Positive ions created through laser ionization have been detected for C_{60} , C_{70} , C_{78} , and C_{84} .

In LCMS measurements, the column performs a vital role since components are individually separated. Fig.3 shows

a comparison of chromatograms obtained using the Develosil RPFULLERENE (made by Nomura Chemical), a column specially developed for fullerene, and a commercially available ODS column. The sample is a mixture of C_{60} and C_{70} (200ng each). A mixed solvent of A: 2-propanol and B: toluene/2-propanol = 90/10 was used for the mobile phase. (1) and (3) show a comparison with the same mobile phase. C_{30} (triaconthyl group) is chemically bonded in the RPFULLERENE, strengthening the retention of the components. (1) and (2) show a comparison at almost the same retention time. These comparisons prove that the use of RPFULLERENE improves the shape of the peaks (the symmetry coefficient of C_{60} is 1.16 with RPFULLERENE while it is 1.23 with ODS (2).




Fig. 1 Structure of C_{60}




Fig. 2 Mass spectrum of fullerenes (MALDI-TOFMS, positive)




Fig. 3 UV chromatograms of fullerenes (335nm)

(1) RPFULLERENE, A/B = 50/50 (2) ODS, A/B = 65/35

(3) ODS, A/B = 50/50

A : 2-propanol, B : toluene/2-propanol = 90/10

The APCI ionization method was employed in the LC-MS measurement, and the fullerene was detected as a negative ion. C₆₀ and C₇₀ were detected as ions of m/z 720 and 840

respectively. Fig.4 shows the mass chromatogram and Fig.5 shows the mass spectra.




Fig.4 Mass chromatograms of fullerenes




Fig.5 Mass spectra of C₆₀ and C₇₀

Table 1 Analytical conditions for LC-MS

Column	: Develosil RPFULLERENE (2.0 mmI.D. x 150 mmL, 5µm, Nomura Chemical)	Column temperature	: 40 degree C
Mobile phase	: A : 2-propanol, B : toluene/2-propanol = 90/10 A/B = 50/50	Probe temperature	: 400 degree C
Flow rate	: 0.2 mL/min		
Injection volume	: 2 µL		
Probe voltage	: -3.0 kV (APCI-Negative mode)		
Nebulizing gas flow	: 2.5 L/min (Air)		

SHIMADZU CORPORATION International Marketing Division

3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan
 Phone: 81 (3) 3219-5641 Fax: 81 (3) 3219-5710
 Cable Add. SHIMADZU TOKYO