

Application

Note
72

Blood Alcohols Analysis by Packed Column GC

Despite many changes in gas chromatography over the last several decades, analysis of blood alcohols by packed column GC, using a packing Supelco specifically tests for this purpose, continues to be a simple, reliable procedure. Supelco packing 1-1766 enables you to simultaneously monitor ethanol, methanol, and isopropanol, and their metabolites, acetaldehyde and acetone, in either glass or stainless steel columns.

Key Words:

- blood alcohols • ethanol • methanol

Packed column GC is a specific, rapid, precise method for identifying and quantifying blood alcohols and related compounds. Supelco™ packing 1-1766 (Carbopack™ B/5% Carbowax® 20M) enables you to simultaneously monitor ethanol, methanol, and isopropanol, and their metabolites, acetaldehyde and acetone. This packing can be used in either glass or stainless steel columns.

GC analyses of blood alcohols are routinely done by one of two methods: direct aqueous injection or headspace analysis. There are advantages to each method, but packing 1-1766 improves the separation obtained by either.


In *direct aqueous injection*, the whole blood or serum sample is diluted with an aqueous solution of the internal standard, and an aliquot is injected directly onto the column. Many analysts prefer this method because the instrumentation required is simpler than that for headspace analysis.

Direct sample injection onto Carbopack C/0.2% Carbowax 1500, a general purpose industrial packing for C1-C5 alcohols, has been used to separate blood alcohols (1). This packing, however, does not resolve acetaldehyde, a product of ethanol metabolism, from methanol. Packing 1-1766 provides better resolution. Figure A shows the isothermal separation of the alcohols, their metabolites, and n-propanol, an internal standard, obtained by direct sample injection.

If you use direct injection, nonvolatile sample components will be deposited at the column inlet. Peaks may begin to tail after the column has been in use for some time. You can eliminate tailing from a glass column by replacing the silanized glass wool plug and the first few inches of packing in the column inlet. This process can be repeated several times before the column must be completely repacked or replaced. If you prefer, you can inject the samples into a glass lined inlet or a PureCol™ column inlet liner (see our catalog), rather than into the column inlet. It is easier to replace these disposable products, when necessary, than to repack a column inlet.


Figure A. Isothermal Separation of Blood Alcohols

Packing: 60/80 Carbopack B/5% Carbowax 20M
Cat. No.: 11766 (15g)
Column: 6' x 2mm ID glass
Oven: 85°C
Carrier: helium, 20mL/min
Det.: FID
Inj.: 1µL, 0.05-0.10% each analyte



795-0296

Figure B. Pre-Column Installation



795-0297

Alternatively, if you wish to use direct injections, we recommend a 6' x 1/8" stainless steel column connected to a 6" x 1/4" OD x 2mm ID glass precolumn through a 1/4" x 1/8" Swagelok® reducing union (Figure B). To maintain column performance, simply repack the precolumn with packing 1-1766 or replace it as necessary.

In the second approach, *headspace analysis*, the internal standard is added to the sample and the mixture is held at a fixed temperature to equilibrate the volatiles between the liquid and the headspace. Only the vapor above the sample is injected onto the column, so contamination from nonvolatiles is avoided. The procedure is well suited to automated analysis with specially designed injection equipment. Manual injection of the vapor via a syringe can yield erratic results because the components can condense in the syringe.

When headspace analysis is used, it is especially important that the column be capable of separating acetaldehyde and methanol. The vapor concentration of acetaldehyde in the headspace can approach that of a clinically significant amount of methanol. If the two compounds are not adequately separated, acetaldehyde might be incorrectly identified as methanol. Acetaldehyde levels in the aqueous phase usually are too low to cause this problem when direct aqueous injection is used.

Always order the blood alcohols packing by catalog number (Cat. No. 1-1766). We test each lot of packing 1-1766 with the mixture shown in Figure A and approve *only* this packing for blood alcohol analyses. Other packings prepared from Carbowax B coated with Carbowax 20M are for different applications and are sold under different numbers.

Reference

1. Manno, B.R., and J.E. Manno, *J. Anal. Toxicol.*, **2**: 257-261 (1978).

Trademarks

Carbopack, PureCol, Supelco – Supelco, Inc.
 Carbowax – Union Carbide Corp.
 Swagelok – Crawford Fitting Co.

Ordering Information:

Description	Cat. No.
60/80 Carbopack B/5% Carbowax 20M Packing 15g	11766 [□]
Glass precolumn packed with Cat. No. 1-1766 6" x 1/4" OD x 2mm ID, pk. of 6 (for use with SS columns only)	12494
Swagelok bored-through reducing union 1/4" x 1/8" SS	22061

For packed columns, please refer to our catalog or call our Order Processing personnel.

[□] **NOTE:** Supelco cannot guarantee results obtained with similar packings, including custom-made packings of identical composition.

Contact our Technical Service Department
 (phone 800-359-3041 or 814-359-3041, FAX 814-359-5468)
 for expert answers to your questions.

Note 72

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA • Sigma-Aldrich de Argentina, S.A. • Buenos Aires 1119 AUSTRALIA • Sigma-Aldrich Pty. Ltd. • Castle Hill NSW 2154 AUSTRIA • Sigma-Aldrich Handels GmbH • A-1110 Wien
 BELGIUM • Sigma-Aldrich N.V./S.A. • B-2880 Bornem BRAZIL • Sigma-Aldrich Quimica Brasil Ltda. • 01239-010 São Paulo, SP CANADA • Sigma-Aldrich Canada, Ltd. • 2149 Winston Park Dr., Oakville, ON L6H 6J8
 CZECH REPUBLIC • Sigma-Aldrich s.r.o. • 186 00 Praha 8 DENMARK • Sigma-Aldrich Denmark A/S • DK-2665 Vallensbaek Strand FINLAND • Sigma-Aldrich Finland/YA-Kemia Oy • FIN-00700 Helsinki
 FRANCE • Sigma-Aldrich Chimie • 38297 Saint-Quentin-Fallavier Cedex GERMANY • Sigma-Aldrich Chemie GmbH • D-82041 Deisenhofen GREECE • Sigma-Aldrich (o.m.) Ltd. • Ilioupoli 16346, Athens
 HUNGARY • Sigma-Aldrich Kft. • H-1067 Budapest INDIA • Sigma-Aldrich Co. • Bangalore 560 048 IRELAND • Sigma-Aldrich Ireland Ltd. • Dublin 24 ISRAEL • Sigma Israel Chemicals Ltd. • Rehovot 76100
 ITALY • Sigma-Aldrich s.r.l. • 20151 Milano JAPAN • Sigma-Aldrich Japan K.K. • Chuo-ku, Tokyo 103 KOREA • Sigma-Aldrich Korea • Seoul MALAYSIA • Sigma-Aldrich (M) Sdn. Bhd. • Selangor
 MEXICO • Sigma-Aldrich Química S.A. de C.V. • 50200 Toluca NETHERLANDS • Sigma-Aldrich Chemie BV • 3330 AA Zwijndrecht NORWAY • Sigma-Aldrich Norway • Torshov • N-0401 Oslo
 POLAND • Sigma-Aldrich Sp. z o.o. • 61-663 Poznań PORTUGAL • Sigma-Aldrich Quimica, S.A. • Sintra 2710 RUSSIA • Sigma-Aldrich Russia • Moscow 103062 SINGAPORE • Sigma-Aldrich Pte. Ltd.
 SOUTH AFRICA • Sigma-Aldrich (pty) Ltd. • Jet Park 1459 SPAIN • Sigma-Aldrich Quimica, S.A. • 28100 Alcobendas, Madrid SWEDEN • Sigma-Aldrich Sweden AB • 135 70 Stockholm
 SWITZERLAND • Supelco • CH-9471 Buchs UNITED KINGDOM • Sigma-Aldrich Company Ltd. • Poole, Dorset BH12 4QH
 UNITED STATES • Supelco • Supelco Park • Bellefonte, PA 16823-0048 • Phone 800-247-6628 or 814-359-3441 • Fax 800-447-3044 or 814-359-3044 • email: supelco@sial.com

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

H