

Determination of 17 Organotin Compounds in Beverages Using Triple Quadrupole GC-MS/MS System

Application Note

Authors

Wenwen Wang, Zhe Cao
Agilent Technologies, Inc.
No.3, Wang Jing Bei Lu,
Chao Yang District,
Beijing, 100102,
China

Abstract

A gas chromatography-triple quadrupole mass spectrometer (GC-MS/MS) method has been developed for the determination of 17 organotin compounds in beverages. The sample was derivatized with NaBEt_4 , and used liquid-liquid extraction with hexane. Data were acquired in the MRM mode and an external standard was used for quantitative determination. Calibration curves were generated for all of the derivatized compounds with concentrations for most of them ranging from 0.001 mg/L to 0.200 mg/L. Samples were concentrated by a factor of 10 during preparation and so the working range for the method was effectively 0.0001 mg/L to 0.0200 mg/L. Correlation coefficients for all of the calibration curves were higher than 0.995. Recoveries were determined for three replicate extractions at two spike levels (0.001 mg/L and 0.005 mg/L). The average recovery for the vast majority of the organotins was between 70.0% and 120.0% with the relative standard deviation less than 10.0%. The method was found to be simple, fast and sensitive. It can be used for the trace analysis of multiple organotin compounds in beverages for both qualitative and quantitative analysis.



Agilent Technologies

Introduction

Organotin compounds (OTCs) have been widely used as polymer stabilizers, insecticides and pesticides. The highest amounts of OTCs are found in coating materials for boats and ships. They are used to prevent crustacean attachment to hulls. OTCs are identified as endocrine disruptive chemicals (EDCs) which may have a negative impact on human health due to exposure through the food chain.

OTCs have been found in ocean water, sediments, textiles and human urine. However, there are fewer studies on OTCs in food. These compounds can be analyzed using GC-PFPD, GC/MS (single quadrupole), ICPMS, or LC-ICPMS. In China, food is referenced against GB method 5009.215-2008 for organotin compounds. This method uses GC-PFPD, and can typically be calibrated down to between 0.5 and 1 µg/kg (depending on target compound). We needed a method that could be calibrated down to 0.1 µg/kg (in the sample) for various beverage samples, and would cover a wider range of compounds (as shown). We wanted to achieve this performance with a sample size no larger than 10 mL (commensurate with what was available) and with a preparation procedure that could be carried out conveniently with small-scale glassware and minimal use of reagents.

An Agilent 7000 Series GC-MS/MS was used in this research work to develop a precise, rapid, and reliable method to study 17 OTCs in beverages. Eleven beverage samples were tested to determine if they contained any of these organotins.

Experimental

Samples and reagents

- **Samples:** Beer, carbonated soft drinks, energy drinks, and drinks containing metabolic stimulants
- **Acetate buffer:** 82 g/L sodium acetate dissolved in water, solution adjusted to pH 4.5 with acetic acid
- **Derivatization reagent:** 2g NaBEt₄ dissolved in 10 mL of ethanol. This solution should be freshly prepared.
- Ethanol, hexane, methanol, chromatography grade purity

Standards and derivatization method


Standards

The standards were dissolved in methanol at 1,000 ppm (1 mg/mL), then further diluted to the required concentration.

Derivatization method

To 1 mL of standard solution, 1 mL of acetate buffer and 50 µL of derivatization reagents were added. The solution was shaken and allowed to react for 30 minutes. After the addition of 5 mL of water, the derivatized compounds were extracted in 1 mL of hexane. The mixture was vortexed for 10 seconds, and the two phases were allowed to separate. The clear upper layer (nonpolar hexane layer) was transferred to an autosampler vial for analysis.

Sample preparation



Chromatographic parameters

| | |
|----------------------------|---|
| GC system | Agilent 7890A GC |
| Column | Agilent HP-5 MS UI capillary column, 30 m × 0.25 mm, 0.25 µm (p/n 19091S-433UI) |
| Oven temperature program | 50 °C hold 1.5 minutes, at 10 °C/min to 300 °C, hold 1 minute |
| Carrier gas | Helium |
| Flow rate | 1.1 mL/min |
| Injection port temperature | 280 °C |
| Injection volume | 2 µL |
| Injection mode | Splitless, purge on after 1 minute |

Mass spectrum parameters

| | |
|------------------------|----------------------|
| Mass system | Agilent 7000B MS/MS |
| Ion source | EI |
| Ionization voltage | 70 eV |
| Ion source temperature | 230 °C |
| Interface temperature | 280 °C |
| Collision gas | Nitrogen 1.50 mL/min |
| Quenching gas | Helium 2.25 mL/min |
| Solvent delay | 2.0 minutes |

MRM parameters are shown in Table 1.

Table 1. Retention Time and MRM Parameters of 17 Derivatives of Organotin Compounds

| No. | Compound name | RT (t/min) | Precursor ion (m/z) | Product ion (m/z) | Collision energy (eV) |
|-----|------------------|------------|---------------------|-------------------|-----------------------|
| 1 | Trimethyltin* | 2.72 | 165 | 135 | 15 |
| | | | 163 | 133 | 15 |
| 2 | Dimethyltin* | 3.86 | 179 | 151 | 5 |
| | | | 151 | 135 | 10 |
| 3 | Monomethyltin* | 5.35 | 193 | 165 | 5 |
| | | | 165 | 137 | 5 |
| 4 | Monobutyltin* | 9.45 | 179 | 151 | 5 |
| | | | 179 | 123 | 10 |
| 5 | Tripropyltin* | 10.38 | 193 | 151 | 5 |
| | | | 193 | 123 | 10 |
| 6 | Tetrapropyltin* | 11.35 | 207** | 165 | 5 |
| | | | 207** | 123 | 10 |
| 7 | Dibutyltin* | 11.66 | 179 | 151 | 5 |
| | | | 263 | 207 | 5 |
| 8 | Monophenyltin* | 12.92 | 255 | 199 | 15 |
| | | | 255 | 277 | 5 |
| 9 | Monoheptyltin* | 13.29 | 179 | 151 | 5 |
| | | | 179 | 123 | 10 |
| 10 | Tributyltin* | 13.54 | 291 | 179 | 10 |
| | | | 207** | 123 | 15 |
| 11 | Monooctyltin* | 14.49 | 179 | 151 | 5 |
| | | | 179 | 123 | 10 |
| 12 | Tetrabutyltin* | 15.16 | 235 | 179 | 5 |
| | | | 291 | 179 | 10 |
| 13 | Diphenyltin* | 17.78 | 303 | 275 | 5 |
| | | | 303 | 197 | 15 |
| 14 | Diheptyltin* | 18.00 | 249 | 151 | 5 |
| | | | 249 | 123 | 15 |
| 15 | Diocetyl tin* | 19.83 | 263 | 151 | 5 |
| | | | 263 | 123 | 15 |
| 16 | Tricyclhexyltin* | 21.80 | 351 | 197 | 20 |
| | | | 349 | 195 | 20 |
| 17 | Triphenyltin* | 21.80 | 233 | 151 | 5 |
| | | | 233 | 123 | 15 |

* Organotin compound after derivatization with NaBEt₄

**While 207 m/z shares signal from column bleed for a 5 phase GC column, it produces high signal and low noise when run through the transitions shown.

Results and Discussion

Chromatographic separation results

GC analysis was finished in 22 minutes with baseline separation for 15 of 17 compounds. Although the last two OTCs could not be separated by this method, the MRM function of the 7000B GC/MS/MS allowed the two coeluting compounds to be separated based on the transition ions. The results are shown in Figure 1. The order and retention time is shown in Table 1.

Calibration curve, linear fit, and recovery results

Most of the OTCs calibration curves were prepared at 1.0–200.0 $\mu\text{g/L}$ (1.0, 5.0, 10.0, 20.0, 50.0, 100.0, 200.0 $\mu\text{g/L}$, seven points), and the results are shown in Table 2. The MRM chromatograms for the analysis of the 17 organotin standard mixture after derivatization with NaBEt_4 (5.0 $\mu\text{g/L}$) are shown in Figure 2.




Figure 1. GC-MS/MS chromatogram for the analysis of 17 organotin standard mixture after derivatization with NaBEt_4 .




Figure 2. MRM chromatogram for analysis of 17 organotin standard mixture after derivatization with NaBEt_4 (5.0 $\mu\text{g/L}$).

Precision and recovery studies were performed in triplicate on beverage samples that were spiked at concentrations of 0.01 mg/L and 0.05 mg/L in hexane, which correspond to 0.001 mg/L and 0.005 mg/L in sample. The majority of RSDs were within 10%, and recoveries within 70–120%. The results are shown in Table 2.

Analysis results of 11 types of beverage

Eleven beverages were tested with above established method for the 17 kinds of OTCs. Dimethyltin and monobutyltin were detected in one sample with concentrations of 0.27 µg/L and 0.39 µg/L respectively.

Conclusion

The Agilent 7000 Series GC-MS/MS system in the MRM mode has the advantage of eliminating most background interferences resulting in high selectivity. This experiment demonstrates this analytical system's effectiveness in analyzing 17 organotins a variety of beverages.

Table 2. Linear Range, R², Recovery and RSD of 17 Organotin Compounds

| No. | Compound name | Linear range (mg/L) | R ² | 0.001 mg/L (n = 3) | | 0.005 mg/L (n = 3) | |
|-----|-------------------------|---------------------|----------------|--------------------|-------|--------------------|-------|
| | | | | % Recovery | % RSD | % Recovery | % RSD |
| 1 | Trimethyltin (TMT)* | 0.001–0.200 | 0.999 | 87.3 | 6.7 | 119.5 | 10.6 |
| 2 | Dimethyltin (DMT)* | 0.001–0.100 | 0.997 | 136.3 | 5.1 | 136.8 | 2.4 |
| 3 | Monomethyltin (MMT)* | 0.001–0.100 | 0.999 | 113.0 | 7.0 | 116.8 | 5.5 |
| 4 | Monobutyltin (MBT)* | 0.001–0.200 | 0.999 | 133.6 | 6.5 | 135.7 | 2.8 |
| 5 | Tripropyltin (TPhT)* | 0.001–0.200 | 0.996 | 104.9 | 4.7 | 114.3 | 2.4 |
| 6 | Tetrapropyltin (TrPhT)* | 0.001–0.200 | 0.996 | 71.0 | 3.5 | 87.7 | 6.1 |
| 7 | Dibutyltin (DBT)* | 0.001–0.200 | 0.999 | 108.9 | 3.2 | 114.9 | 2.4 |
| 8 | Monophenyltin (MPhT)* | 0.001–0.200 | 0.998 | 115.4 | 4.6 | 118.0 | 1.7 |
| 9 | Monoheptyl tin (MHT)* | 0.001–0.200 | 0.999 | 116.2 | 9.4 | 114.4 | 9.0 |
| 10 | Tributyltin (TBT)* | 0.001–0.200 | 0.999 | 117.9 | 7.8 | 112.1 | 8.9 |
| 11 | Monooctyltin (MOcT)* | 0.001–0.200 | 0.997 | 129.4 | 3.8 | 112.8 | 8.1 |
| 12 | Tetrabutyltin (TeBT)* | 0.001–0.200 | 0.998 | 63.4 | 8.2 | 77.9 | 11.3 |
| 13 | Diphenyltin (DPhT)* | 0.001–0.200 | 0.998 | 76.0 | 6.5 | 83.7 | 1.8 |
| 14 | Diheptyl tin (DHT)* | 0.001–0.200 | 0.999 | 102.1 | 5.6 | 96.7 | 2.3 |
| 15 | Diocetyl tin (DOcT)* | 0.001–0.200 | 0.999 | 67.5 | 6.3 | 75.5 | 0.8 |
| 16 | Tricyclhexyltin (TCyT)* | 0.001–0.200 | 0.998 | 107.4 | 6.1 | 99.8 | 2.2 |
| 17 | Triphenyltin (TPhT)* | 0.001–0.200 | 0.999 | 82.2 | 7.5 | 82.4 | 3.4 |

References

1. Wang Xiuping, *et al.* "Development of an analytical method for organotin compounds in fortified flour samples using microwave assisted extraction and normal-phase HPLC with UV detection [J]". *Journal of Chromatography B*, 2006, **843**: 268-274.
2. H. Julien, *et al.* "Analytical advances in butyl-,phenyl- and octyltin speciation analysis in soil by GC-PFPD [J]" *Talanta*, 2008, **75**: 486-493.
3. Y. Mino, *et al.* "Determination of organotin in human breast milk by gas chromatography with flame photometric detection[J]" *Journal of Health*, 2008, **54**(2): 224-228.
4. F. David, P. Sandra, P.I. Wylie "Improving the Analysis of Organotin Compounds Using Retention Time Locked Methods and Retention Time Databases"
Agilent publication 5988-9256E.

For More Information

For more information on our products and services visit our Website at www.agilent.com/chem.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2014
Published in the USA
April 15, 2014
5991-4434EN



Agilent Technologies