



Solvents

Separation of chlorinated solvents

Application Note

Environmental

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography with an Agilent CP-Select 624 CB column separates 19 chlorinated solvents in 12 minutes.




Agilent Technologies

Conditions

Technique : GC-wide-bore
Column : Agilent CP-Select 624 CB fused silica WCOT
30 m x 0.53 mm, fused silica WCOT (df = 3.0 µm)
(Part no. CP7416)
Temperature : 50 °C → 200 °C. 10 °C/min
Carrier Gas : N₂, 10 mL/min
Injector : Direct,
T = 250 °C
Detector : FID
T = 250 °C
Sample Size : 0.02 µL
Solvent Sample : solvents mixture

Peak identification

- 1,1-dichloroethylene
- impurity
- 1,2-dichloroethylene
- trichloromethane (chloroform, CFC 20)
- tetrachloromethane
- 1,2-dichloroethane
- trichloroethylene
- epichlorohydrin
- tetrachloroethylene
- chlorobenzene
- 2-chlorotoluene
- 3-chlorotoluene
- 1,3-dichlorobenzene
- 1,4-dichlorobenzene
- benzyl chloride
- 1,2-dichlorobenzene
- impurity
- 1,2,4-trichlorobenzene
- α,α-dichlorotoluene (benzal chloride)



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01415



Agilent Technologies