


The Era of Mass Surveillance

Identify Suspects, Confirm Targets and
Explore Unknowns with GC/Q-TOF

Nathan Eno
GC/Q-TOF Product Manager

Kai Chen, Ph.D.
GC/MS Application Scientist



Panoramic Technology for Full Spectrum Surveillance





High Resolution and Mass Accuracy


Reproducible Spectral Performance

Simultaneous
High Resolution and
Wide Dynamic Range

Flexible Agilent
7890B GC

Sensitive Low Energy EI





Target Quantitation

Linearity and Detection Limit Performance for Targeted Assays



File Home View Method Tools Help

Edit New Open Append Method Report Save Validate Save As Exit




Method Setup Tasks Advanced Tasks Duplicate Compounds Manual Setup Tasks Update Swap Qualifier With Quantifier Outlier Setup Tasks Calibration Curve

Library Method Setup Reference Library Setup Reference Pattern Library


Batch Table

Sample: 10 ppb Avocado Sample Type: <All> Compound: Azoxystrobin

Sample				Azoxystrobin Results								Qualifier (388.0928) Results				Qualifier (345.1060) Results			
Name	Type	Level	Exp. Conc.	RT	Resp.	MI	Calc. Conc.	Final Conc.	Accuracy	Mass Accuracy	Ratio	Mass Accuracy	S/N	MI	Ratio	Mass Accuracy	S/N	MI	Ratio
2 ppb Avocado	Cal	2	2.0000	18.415	11667		2.0805	2.0805	104.0	-1.9689	29.6	-0.7354	7.38		22.4	-0.7069	3.61		15.8
5 ppb Avocado	Cal	3	5.0000	18.417	32806		4.3995	4.3995	88.0	-2.5398	29.7	-1.8039	16.99		19.5	-3.0019	4.66		18.3
10 ppb Avocado	Cal	4	10.0000	18.415	67221		8.1748	8.1748	81.7	-1.4376	28.1	-1.2841	26.68		18.5	-0.0052	6.28		18.1
20 ppb Avocado	Cal	5	20.0000	18.417	143907		16.5874	16.5874	82.9	-0.0917	27.0	0.3633	50.71		19.0	0.6008	21.55		15.4
50 ppb Avocado	Cal	6	50.0000	18.418	424513		47.3703	47.3703	94.7	-1.7079	24.8	-1.7057	77.07		18.4	-1.7582	39.17		15.0
100 ppb Avocado	Cal	7	100.0000	18.418	858975		95.0315	95.0315	95.0	-0.8646	26.6	-0.6051	455.66		19.3	-1.0061	48.39		16.1
200 ppb Avocado	Cal	8	200.0000	18.419	1845131		203.2142	203.2142	101.6	-1.5004	26.5	-0.7792	1020.33		20.5	-1.3373	104.66		16.2
500 ppb Avocado	Cal	9	500.0000	18.420	4638419		509.6419	509.6419	101.9	-0.0104	27.6	0.8329	1648.86		21.2	-0.2021	257.53		16.2




6-Hydroxychrysene - 9 Levels, 9 Levels Used, 10 Points, 10 Points Used, 0 QCs




Compound name	LOD (ng/mL) per matrix	
	Urine	Plasma
Naphthalene	0.52	0.38
Acenaphthylene	0.42	0.74
Acenaphthene	0.98	1.28
1-Hydroxynaphthalene	0.48	0.47
2-Hydroxynaphthalene	0.28	0.20
Fluorene	1.12	0.26
Phenanthrene	0.35	0.18
Anthracene	0.36	0.26
1,6-Dihydroxynaphthalene	0.34	0.16
2,7-Dihydroxynaphthalene	0.13	0.10
3-Hydroxyfluorene	0.22	0.36
2-hydroxyfluorene	0.25	0.21
4-Hydroxyphenanthrene	0.39	0.42
Fluoranthene	0.30	0.19
3-Hydroxyphenanthrene	0.60	0.27
1-Hydroxyphenanthrene	0.34	0.30
9-Hydroxyphenanthrene	0.50	0.33
Pyrene	0.42	0.10
Benzo(a)anthracene	0.16	0.23
Chrysene	0.26	0.09
6-Hydroxychrysene	0.14	0.16
Benzo(b)fluoranthene	0.38	0.26
Benzo(k)fluoranthene	0.38	0.27
Benz(a)pyrene	0.61	0.30
Indeno(1,2,3-cd)pyrene	0.65	0.35
Dibenz(a,h)anthracene	0.51	0.09
9-Hydroxybenzo(a)pyrene	0.47	0.34
Benzo(ghi)perylene	0.89	0.36
7,8-Hydroxybenzo(a)pyrene	6.65	3.63


6-Hydroxychrysene - 9 Levels, 9 Levels Used, 9 Points, 9 Points Used, 0 QCs



SureMass Signal Processing



Conventional Deconvolution
(Nominal Mass)



SureMass Signal Processing
(Profile Accurate Mass)

Suspect Screening

Analytical Performance Factors for Confident Compound Detection




Fingerprint image adapted from
The Photographer - Own work, CC BY-SA 3.0,
<https://commons.wikimedia.org/w/index.php?curid=8251516>


Suspect Surveillance – Now with even broader scope

GC/Q-TOF Pesticides & Environmental Pollutants PCDL – now with *1000+ compounds*:


- High Resolution Spectra
- Expert, human curation
- Better compound alignment with US EPA 8270 targets and Agilent GC/TQ MRM database




PAHs




Phthalates




Organophosphates



Nitroanilines



Amines



Chloronitrobenzenes

File Launch Edit View Find Identify Method Configuration Tools Help

Compound List: 1018 found, 32 shown, filtered on Flags (Tgt), Fb Conf.



Compound Identification		General			Target/Suspect Screening						
Name	CAS	Formula	m/z	Area	RT	RT Diff (Tgt)	Score (Frag Coelution)	Score (Frag Ratio)	Flags (Tgt)	FIs Conf.	FragMassDiff(ppm)
Diisobutylphthalate	84-69-5	C16 H22 O4	149.0233	93135	8.832	0.015	98.63	97.94	Qualified	5	0.68
Dimethomorph (E)	110488-70-5	C21 H22 Cl N O4	387.1241	21210	18.48	0.029	99.75	99.17	Qualified	6	2.42
Dimethomorph (Z)	113210-98-3	C21 H22 Cl N O4	387.1242	16614	18.79	0.09	99.8	99.73	Qualified	6	2.66
Fludioxonil	131341-86-1	C12 H6 F2 N2 O2	248.0394	27773	11.55	0.024	99.35	99.5	Qualified	6	0.6
Isophorone	78-59-1	C9 H14 O	138.104	20234	4.166	0.037	89.3	92	Qualified	5	1.03
Mefenoxam	70630-17-0	C15 H21 N O4	279.1414	1119	9.342	0.007	89.32	97.21	Qualified	5	0.11
Metalaxyl	57837-19-1	C15 H21 N O4	206.1176	12522	9.342	0.014	87.35	99.25	Qualified	6	0.11
Pentachlorobenzoni	20925-85-3	C7 Cl5 N	274.8437	23680	8.27	0.008	99.17	99.98	Qualified	6	0.3
Phenanthrene-D10	1517-22-2	C14 D10	188.1404	66032	8.396	0.089	98.69	99.81	Qualified	6	0.28
Phenol	108-95-2	C6 H6 O	94.0413	10093	3.506	0.097	97.14	98.8	Qualified	6	0.67
Thiamethoxam	153719-23-4	C8 H10 Cl N5 O3	212.0488	8652	10.29	0.023	87.04	87.71	Qualified	5	0.57
TPPA / Triphenyl ph	115-86-6	C18 H15 O4 P	326.07	14125	13.37	0.025	99.85	99.51	Qualified	6	0.81
Triethylphosphate	78-40-0	C6 H15 O4 P	155.0468	19149	4.099	-0.033	95.97	99.5	Qualified	6	0.36
Triisobutylphosphat	126-71-6	C12 H27 O4 P	98.9842	60751	6.088	-0.019	95.64	99.66	Qualified	5	0.65


Compound Identification Results: Cpd 322: Fludioxonil; C12 H6 F2 N2 O2; 11.551

ID Techniques Applied



FBF-FragConfirm

Best	Name	Formula	m/z	Mass	Mass (Tgt)	Diff (ppm)	Score
127.0413	127.0417		2.6	Qualified	51.9	11.6	561
153.0448	153.0447		0.4	Qualified	16.7	16.4	1667
154.0525	154.0525		0.3	Qualified	39.8	40.1	4212
182.0475	182.0475		0	Qualified	30.6	92.1	3952
248.0393	248.0392		0.6	Reference ion	100	363.9	14577
249.0424	249.043		2.4	Qualified	14.6	39.2	1596









Clostebol
Metabolite
(TMS-derivatized)



Sample Name	File Name	Compou	Hits
S1_Spike_1ul...	S1_Spike_1...	5047	69
S1_MRPL_1...	S1_MRPL_1...	5159	37
S1_Blank_1ul...	S1_Blank_1...	5323	77
S2_Blank_1ul...	S2_Blank_1...	5462	63

Component RT	Compound Name	Match Factor	Best Hit	Form
2.2143	Crotetamide	80.0	<input checked="" type="checkbox"/>	C12H
2.2429	Benzyl piperazine, N-trimethylsilyl-	89.1	<input checked="" type="checkbox"/>	C14H
2.3630	(+/-)-BDB, N-trimethylsilyl-	80.5	<input checked="" type="checkbox"/>	C14H
2.3737	Crotonamide, N-[1-(dimethylcarbamoyl)...	85.9	<input checked="" type="checkbox"/>	C13H
2.5367	Lidocaine, TMS derivative	80.5	<input checked="" type="checkbox"/>	C17H
2.5762	Benzoic acid, 4-[(2,4-dimethoxy-6-pent...	82.0	<input checked="" type="checkbox"/>	C28H
2.6210	Caffeine	88.4	<input checked="" type="checkbox"/>	C8H1
2.9101	Myristic acid, TMS derivative	85.8	<input checked="" type="checkbox"/>	C17H
2.9160	Methyl 2,6-dihydroxybenzoate, 2TMS ...	75.8	<input checked="" type="checkbox"/>	C14H
2.9903	Pipecolic acid, N-ethoxycarbonyl-, dod...	78.1	<input checked="" type="checkbox"/>	C21H
3.1435	Isoxanthopterin, tris(trimethylsilyl) deriva...	78.3	<input checked="" type="checkbox"/>	C15H
3.2312	Hexadecanoic acid, methyl ester	85.7	<input checked="" type="checkbox"/>	C17H
3.2996	Theobromine, TMS derivative	92.7	<input checked="" type="checkbox"/>	C10H
3.4140	Clobenzorex	87.6	<input checked="" type="checkbox"/>	C16H
3.5060	N,N-Diethyl-3-methoxy-4-(trimethylsilyl)...	87.8	<input checked="" type="checkbox"/>	C15H
3.6399	Succinic acid, 3-methylbut-2-en-1-yl di...	76.8	<input checked="" type="checkbox"/>	C22H
3.9212	L-Leucine, N-ethoxycarbonyl-N-methyl-...	85.8	<input checked="" type="checkbox"/>	C24H
4.0582	Silane, dimethyl(4-(2-phenylprop-2-yl)ph...	84.3	<input checked="" type="checkbox"/>	C20H
4.0813	Methyl propyl ether	76.3	<input checked="" type="checkbox"/>	C4H1
4.0908	2-Butanone, 3-methoxy-3-methyl-	93.1	<input checked="" type="checkbox"/>	C6H1
4.1340	Acetic acid, TBDMS derivative	77.1	<input checked="" type="checkbox"/>	C8H1




Unknowns Identification


Extensive Functionality and Tools for Identification and Elucidation




+EI Scan (rt: 12.515-12.525 min, 4 scans) Frag=70.0V U-D_BlankN_70eV_5Hz.D Subtract 70 eV




Candidate:
Escitalopram




+EI Scan (rt: 12.523-12.533 min, 4 scans) Frag=15.0V U-D_BlankN_15eV_5Hz.D Subtract 15 eV





+EI Scan (rt: 12.515-12.525 min, 4 scans) Frag=12.0V U-D_BlankN_12eV_5Hz.D Subtract 12 eV







Proprietary Low Energy EI Source












Clenbuterol (335.0690), MRPL = 0.2 ng/ml




4-OH-Stanozolol (560.3644), MRPL = 5 ng/ml


TIC










Compound ID = ???











Candidate #1



Candidate #2





Sort by **Score** Show structures for **All Formulas**

Sort

C15H24O; 48451598

Scores
MFG=100.0 MSC=94.8

ChemSpider: [484515](#)

More Info...

Fragment

C15H24O; 48451631

Scores
MFG=100.0 MSC=94.8

ChemSpider: [484516](#)

More Info...

Structure #10 -- elucidated: 76.7% Display Filters

	Mass	Intensity	Weight(%)
▶ 1	135.0803	53260.97	32.3
2	191.1429	31287.67	38.0
3	149.0960	27310.47	20.2
4	121.0645	6632.19	3.2
5	107.0491	3087.48	1.2
6	135.1201	865.06	0.5

Penalty=1.5 dM=1.0ppm F.D.S.=99.9 1 Of 1 Pen C9H

C9H12O-H Score=99.1

1.0 ppm

Penalty=1.5 dM=0.7ppm F.D.S.=99.9 1 Of 1

C13H20O-H Score=99.1

0.7 ppm

Penalty=1.5 dM=0.6ppm F.D.S.=100.0 1 Of 1

C10H14O-H Score=99.1

0.6 ppm

Penalty=2.0 dM=2.4ppm F.D.S.=99.7 1 Of 2

C8H10O-H Score=98.2

2.4 ppm

Penalty=2.0 dM=0.4ppm F.D.S.=100.0 1 Of 1


C7H8O-H Score=98.5

0.4 ppm

Candidate: Benzenemethanol derivative

Case Study: Food Safety

Full Spectrum Accurate Mass Data for Targets, Suspects and Unknowns



Safeguarding food supply is challenging



- 1000+ pesticides in use or remain in environment
- Other environmental pollutants are also of concern
 - Thousands of possible contaminants/food pairs
- High sensitivity and selectivity needed to meet “Default” MRLs in “dirty” extracts
 - Growing interests in broadest scope and even untargeted screening for risk assessment

Challenges for our customers



How much of each calibrated target compound is present?

Dozens of Compounds
e.g., 100~200 Targets

Are other target residues present?

100's of Compounds
e.g., 800 Targets

Is there anything else in my sample?

1000's of Compounds

Workflow Strategy

Acquire full-spectrum data



Targeted Method

Yes


Calibrate?

No

Targeted Quantitative Screening

Suspect Screening


Untargeted Screening





NIST

Deconvolution followed by matching public libraries

Configuration Optimized with Backflushing



- ✓ Reduced run times
- ✓ Longer column lifetime
- ✓ Less ion source contamination




7250 GC/Q-TOF Method Parameters

Agilent 7890 GC

Parameter	Value
Inert flow path configuration	Mid-column backflush
Columns	Agilent HP-5ms UI, 15 m, 0.25 mm id, 0.25 µm film (two each)
Inlet	MMI, 4 mm UI liner single taper w wool
Injection volume	1~2 µL
Injection mode	Cold Splitless
	60 °C for 0.2 minutes 600 °C/min to 300 °C, hold 330 °C, post run
Inlet flow (column 1)	1.0 mL/min (Chlorpyrifos-methyl locked at 9.143 min)
PUU flow (column 2)	column 1 flow + 0.2 mL/min
Oven temperature program	60 °C (hold 1 min) then 40 °C/min to 170 °C, then 10 °C/min to 310 °C (hold 3 min)
	Run time 20.75 min
Transfer line	280 °C
Midcolumn Backflush	
Timing	5 min duration during post-run
Oven temperature	310 °C
Aux EPC pressure	~50 psi
Inlet pressure	~2 psi

Agilent 7250 Q-TOF





Parameter	Value
Source temperature	280 °C (70 eV), 250 (°C)
Quad temperatures	150 °C
Collision cell gas flows	1 mL/min N ₂
	4 mL/min He
Electron energy	70 eV (Standard EI), 15 eV (Low energy EI)
Acquisition mass range	45-550 m/z
Spectral acquisition rate	5 spectra/sec



Method Evaluation

**Targeted Quantitative
Screening**

Sample Preparation

Matrix	Extraction	dSPE Cleanup
	QuEChERS (EN)	EMR-Lipid
	QuEChERS (EN)	Pigment matrix
	QuEChERS (EN)	Fruit/Vegetable
	QuEChERS (EN)	Fruit/Vegetable






- Blank food matrices spiked with pesticides for method development and calibration
- A mixture of 140+ pesticides represent OCs, OPs, carbamates, triazoles and pyrethroids, etc.
- 6 replicates



Replicate Precision


Targeted Quantitative Screening



Matrix Matched Calibration

Targeted Quantitative Screening


α -BHC - 8 Levels, 8 Levels Used, 8 Points, 8 Points Used, 0 QCs




2-500 ng/mL



Carbofuran - 8 Levels, 8 Levels Used, 8 Points, 8 Points Used, 0 QCs




Malathion - 8 Levels, 8 Levels Used, 8 Points, 8 Points Used, 0 QCs



Malathion 5 ng/mL



Deltamethrin 5 ng/mL

Deltamethrin (Decamethrin) - 8 Levels, 8 Levels Used, 8 Points, 8 Points Used, 0 QCs



Confidence Enhanced by Mass Accuracy

Targeted Quantitative Screening



Identify Suspect and Confirm Target

**Targeted Quantitative
Screening**

+

Suspect Screening

Guideline to Identify Compounds

EU SANTE/11813/2017 Guidelines

MS detector/Characteristics		Acquisition	Requirements for identification	
Resolution	Typical systems (examples)		minimum number of ions	other
Unit mass resolution	Single MS quadrupole, ion trap, TOF	full scan, limited m/z range, SIM	3 ions	S/N $\geq 3^{d)}$ Analyte peaks from both product ions in the extracted ion chromatograms must fully overlap.
	MS/MS triple quadrupole, ion trap, Q-trap, Q-TOF, Q-Orbitrap	selected or multiple reaction monitoring (SRM, MRM), mass resolution for precursor-ion isolation equal to or better than unit mass resolution	2 product ions	Ion ratio from sample extracts should be within $\pm 30\%$ (relative) of average of calibration standards from same sequence
Accurate mass measurement	High resolution MS: (Q-)TOF (Q-)Orbitrap FT-ICR-MS sector MS	full scan, limited m/z range, SIM, fragmentation with or without precursor-ion selection, or combinations thereof	2 ions with mass accuracy ≤ 5 ppm ^{a, b, c)}	S/N $\geq 3^{d)}$ Analyte peaks from precursor and/or product ion(s) in the extracted ion chromatograms must fully overlap. Ion ratio: see D12

^{a)} preferably including the molecular ion, (de)protonated molecule or adduct ion

^{b)} including at least one fragment ion

^{c)} < 1 mDa for $m/z < 200$

^{d)} in case noise is absent, a signal should be present in at least 5 subsequent scans

How to Confirm Identification

Suspect Screening

Compound List: 1018 found, 32 shown, filtered on Flags (Tgt), Fb Conf.

Compound Identification		General				Target/Suspect Screening					
Name	CAS	Formula	m/z	Area	RT	RT Diff (Tgt)	Score (Frag Coelution)	Score (Frag Ratio)	Flags (Tgt)	Fb Conf.	FragMassDiff(ppm)
Diisobutylphthalate	84-69-5	C16 H22 O4	149.0233	93135	8.832	0.015	98.63	97.94	Qualified	5	0.68
Dimethomorph (E)	110488-70-5	C21 H22 Cl N O4	387.1241	21210	18.48	0.029	99.75	99.17	Qualified	6	2.42
Dimethomorph (Z)	113210-98-3	C21 H22 Cl N O4	387.1242	16614	18.79	0.09	99.8	99.73	Qualified	6	2.66
Fludioxonil	131341-86-1	C12 H6 F2 N2 O2	248.0394	27773	11.55	0.024	99.35	99.5	Qualified	6	0.6
Isophorone	78-59-1	C9 H14 O	138.104	20234	4.166	0.037	89.3	92	Qualified	5	1.03
Mefenoxam	70630-17-0	C15 H21 N O4	279.1414	1119	9.342	0.007	89.32	97.21	Qualified	5	0.11
Metalaxyl	57837-19-1	C15 H21 N O4	206.1176	12522	9.342	0.014	87.35	99.25	Qualified	6	0.11
Pentachlorobenzoni	20925-85-3	C7 Cl5 N	274.8437	23680	8.27	0.008	99.17	99.98	Qualified	6	0.3
Phenanthrene-D10	1517-22-2	C14 D10	188.1404	66032	8.396	0.089	98.69	99.81	Qualified	6	0.28
Phenol	108-95-2	C6 H6 O	94.0413	10093	3.506	0.097	97.14	98.8	Qualified	6	0.67
Thiamethoxam	153719-23-4	C8 H10 Cl N5 O3	212.0488	8652	10.29	0.023	87.04	87.71	Qualified	5	0.57
TPPA / Triphenyl ph	115-86-6	C18 H15 O4 P	326.07	14125	13.37	0.025	99.85	99.51	Qualified	6	0.81
Triethylphosphate	78-40-0	C6 H15 O4 P	155.0488	19149	4.099	-0.033	95.97	99.5	Qualified	6	0.36
Triisobutylphosphat	126-71-6	C12 H27 O4 P	98.9842	60751	6.088	-0.019	95.64	99.66	Qualified	5	0.65

Compound Identification Results: Cpd 322: Fludioxonil; C12 H6 F2 N2 O2; 11.551

Best	Name	Formula	m/z	Mass	Mass (Tgt)	Diff (ppm)	Score
127.0413	127.041	C12 H6 F2 N2 O2	248.0394	248.0399	248.0397	-0.72	
153.0448	153.044						
154.0525	154.052						
182.0475	182.047						
248.0393	248.039						
249.0424	249.04						

Compound Chromatogram Results (zoomed)

Compound Fragment Spectrum Results (zoomed)



Targeted Screening


Real samples from the Supermarket

Targeted Quantitative Screening

+

Suspect Screening

Sample	Compound name	RT diff	Score		Mass diff	Amount (ng/mL)
		(min)	Coelution	Frag ratio	(ppm)	
	2,6-Dimethylphthalene	-0.006	97.87	88.53	0.21	ID only
	Azobenzene	0.051	91.18	67.01	0.54	ID only
	Benzylbenzoate	0.014	77.62	62.55	0.59	ID only
	Biphenyl	-0.002	97.05	86.43	1.34	ID only
	Boscalid	0.026	93.56	97.75	1.00	ID only
	DPA / Diphenylamine	0.003	98.54	99.8	0.51	ID only
	Eugenol	-0.003	93.08	86.83	1.02	ID only
	Fludioxonil	0.022	99.82	91.88	0.41	ID only
	Flutriafol	0.016	85.66	80.75	3.62	ID only
	Phenanthrene	0.005	86.57	95.45	0.29	ID only
	Pyraclostrobin	0.02	98.41	93.29	0.8	ID only
	Pyrimethanil	0.007	99.32	97.68	4.41	ID only
	Thiabendazole	0.012	99.67	99.37	0.22	ID only
		2,6-Dimethylphthalene	-0.004	99.44	63.82	0.17
Anthraquinone		0.009	94.04	87.87	0.19	ID only
Imazalil		0.028	88.91	84.45	1.05	ID only
Propiconazole		0.023	97.29	99.45	0.24	ID only
Pyrimethanil		0.013	99.12	97.68	4.69	ID only
Thiabendazole		0.017	88.02	99.29	0.18	ID only
Thiamethoxam		0.024	86.58	88.41	1.61	2.46
Tri(2-chloroethyl)phosphate		-0.008	82.55	75.39	2.85	ID only

Sample	Compound name	RT diff	Score		Mass diff	Amount (ng/mL)
		(min)	Coelution	Frag ratio	(ppm)	
	1,2,3,5-tetrachlorobenzene	0.008	88.8	85.1	0.48	ID only
	1-aminonaphthalene	0.021	96.6	72.3	0.96	ID only
	Anthraquinone	0.011	78.2	84.7	0.52	ID only
	Azobenzene	0.049	86.1	64.1	0.51	ID only
	Azoxystrobin	0.028	99.5	99.6	2	878
	Benzyl benzoate	0.009	74.3	55.8	0.07	ID only
	Boscalid	0.027	90.7	99.5	0.7	ID only
	Cyfluthrin	0.051	84.3	79.9	0.49	13.6
	Cyhalothrin	0.025	96.1	84.2	1.68	43
	Chlorthal-dimethyl	0.013	93.6	99.9	0.84	ID only
	Dimethomorph (E & Z)	0.029	99.8	99.2	2.42	972
	Fludioxonil	0.024	99.4	99.5	0.6	ID only
	Metalaxyl	0.014	87.4	99.3	0.11	ID only
	p,p'-DDE	0.015	90.8	87.2	1.7	4.65
	Pentachlorobenzonitrile	0.008	99.2	100	0.3	ID only
	Permethrin	0.022	99.5	98.9	1.8	61.3
	Phenathrene	0.003	91.8	97	1.08	ID only
	Pyraclostrobin	0.018	97.8	92.8	0.38	ID only
	Thiabendazole	0.023	88.7	86.7	0.35	ID only
	Thiamethoxam	0.023	87	87.7	0.57	2.46
	Triphenylmethane	0.016	84.4	85.7	0.29	ID only




Explore Unknown

Untargeted Screening

Untargeted Screening (NIST)

Untargeted Screening

+ TIC Scan Broccoli non-org 70eV8uA280C 02.D



Elucidate structures.....


low energy EI, MS/MS and accurate mass

Untargeted Screening


Step 1: Confirm M⁺

Step 2: Confirm fragment ions

Step 3: Structure elucidation on candidate



Low energy EI (15 eV)



Sort by # Reference: Show structures for All Formulas

Structure #1 -- elucidated: 68.4% ions, 97	Mass	Intensity	Weight(%)	No. of candid.	Best score
1	194.9504	2989.09	40.2	1	91.2
2	158.9741	2017.41	18.1	2	68.7
3	168.9477	1060.99	10.7	4	53.1

Structure #1 -- elucidated: 68.4% ions, 97

Penalty=13.0 dM=2.1ppm F.D.S.=99.6 1 Of 6 Penalty=17.0

C7H3Cl2N-2H Score=53.1 C7H7Cl2N-6H

C8-H-Cl3-N2 turns out to be 2,4,5-Trichloroisophthalonitrile. That is not NIST. It is in fact a degradation product of Chlorothalonil

7250 GC/Q-TOF

- Mass Accuracy
- Isotopic Fidelity
- MS/MS Performance
- Low Energy EI Capability
- MassHunter Software

Identify. Quantify. Simplify

Thank You!



- MP 238 Analysis of polycyclic aromatic hydrocarbons (PAH) and hydroxylated PAH metabolites in plasma and urine
- MP 284 Screening of contaminants in food and natural products
- TP 180 Low Energy EI and High Resolving Power Instrumentation for the Analysis of Arson Samples
- TP 213 Targeted and Untargeted Analysis of Pesticides and Other Contaminants in Fruits and Vegetables
- WP 280 Analysis of doping and forensic drugs in urine
- ThP 286 Analysis of halogenated polycyclic aromatic hydrocarbons in atmosphere
- ThP 295 Combination of chemical ionization (CI) and low electron energy ionization capabilities