

ASMS 2013

MP-685

Accurate Mass Retention
Time Locked Flavor
Database by GC/Q-TOF

¹Stephan Baumann, ¹Tim Conjelko,
¹Sofia Aronova, ²Sean LaFond, ³Frank
David, and ²Susan E. Ebeler

¹Agilent Technologies, Santa Clara,
CA;

²UC Davis Department of Viticulture
and Enology, Davis, CA;

³Research Institute for
Chromatography, Kortrijk, Belgium

Introduction

Companies in the flavor and fragrance industry often use proprietary GC/MS databases based on historical methods for research and quality control. Fortunately, many of their MS based libraries are built upon common non-polar columns and so their retention index information can be converted into an absolute retention time.

Although quadrupole based GC/MS is a powerful technique in the flavor and fragrance industry, it is not optimal for identifying new flavorants. The use of an accurate mass instrument with MS/MS capabilities allows for structure elucidation of novel flavorants in the part per million concentration range without the need for isolating compounds.

Using retention time locking, it is no longer necessary to calculate the retention index, but the absolute retention time can be used as an identification tool. Of course, retention times are still dependent on operating conditions, but small differences in carrier gas velocity and column length are compensated by re-locking the GC method by adjustment of the column head pressure. After re-locking, elution temperatures of the solutes are also constant. Retention time locking and retention time databases are excellent tools in essential oil and in flavor QA/QC analysis. Incorporating both accurate mass information and absolute retention times are useful in the identification of isomeric flavorants.

Methods

Sample Preparation

This database was developed on analytical standards run either in split mode or using Solid Phase Micro Extraction (SPME). The standards were originally made in either ethanol or hexanes at 1 mg/mL concentration dependent on solubility. The standards were diluted down to 10 ppm concentration in hexanes. Further dilutions were also made dependent on ionization efficiency of the specific flavorant.

Volatile samples were prepared in 20-mL headspace vials. HS-SPME was performed using a conditioned 50/30 µm DVB/Carboxen/PDMS StableFlex SPME fiber. Semi-volatile compounds were injected using a split/splitless inlet run in split mode. The Retention Time Locked (RTL) GC method was utilized for separation of flavor compounds to construct a targeted flavor database. The spectra were acquired in full spectrum acquisition mode.

Methods

Analytical Conditions

This study was performed using an Agilent 7890 GC coupled to an Agilent 7200 series Quadrupole-Time-of-Flight (Figure 1). GC and MS conditions are described in Table 1. A CTC-PAL was installed for the headspace SPME samples.

The SPME samples were equilibrated at 40 °C for 3 min with 500 rpm agitation. The fiber was then lowered into the headspace and the samples were extracted for 3 min before being desorbed in the inlet for 3 min. After 1 minute the purge flow flushed out the inlet. The fiber was left in the inlet for 3 additional minutes until the gas saver was turned on and the fiber was retracted.

The split injections were done using the same inlet and oven conditions to minimize the impact of the different injection techniques.



Figure 1. 7200 series GC/Q-TOF system.

| GC and MS Conditions: | |
|-----------------------------------|---|
| Column | HP-5 MS, 30 meter, 0.25 mm ID, 0.25 µm film |
| Injection volume | 1:10 Split or SPME |
| Purge to Split Vent | 60 mL/min at 1 min |
| Split/Splitless inlet temperature | 300 °C |
| Oven temperature program | 40 °C for 10 min 4 °C/min to 140 °C, 0 min hold 10 °C/min to 200 °C, 0 min hold 50 °C/min to 300 °C, 2 min hold 50 °C/min to 325 °C, 0.5 min hold |
| Carrier gas | Helium at 1.2 mL/min constant flow |
| Transfer line temperature | 300 °C |
| Ionization mode | EI |
| Source temperature | 230 °C |
| Quadrupole temperature | 150 °C |
| Scan range | 40 to 800 m/z |
| Spectral acquisition rate | 5 Hz, collecting both in centroid and profile modes |

Table 1. GC-MS conditions used in the study.

Results and Discussion

Chromatogram Deconvolution

The data was processed by chromatogram deconvolution using the Find by Chromatographic Deconvolution tool of MassHunter Qualitative software. Chromatogram deconvolution was able to extract clean spectra from background noise based on both retention time and peak shape. The compounds could then be searched against a unit mass database, generating a Match Factor score comparable to the score generated by NIST.




Figure 2. Find by Chromatographic Deconvolution tool was used to extract spectra.

Find by Formula and Fragment Formula Annotation

The Find by Formula and Fragment Formula Annotation tools were then used to automatically identify fragments of the compound identified by library search. In Find by Formula, a match score was generated based on abundance ratios, isotope spacing and isotope m/z match.




Figure 3. Find By Formula results for 2-heptanone. The match score is 99.59; it is based on abundance ratios, isotope spacing and isotope m/z match. These ions are compatible with the formula C₇H₁₄O. Moreover, the mass difference is only 0.73 ppm.

Library Editor

The Library Editor can read commercial EI spectral libraries. Currently the editor considers the standard *.L format databases as compressed and cannot write to that format. However it is possible to import custom library data in JCAMP formats. GC/MS ChemStation can create the JCAMP format files and LIB2NIST program can create the HPJ format. There is also a script to automatically import structures (.MOL files) from the ChemSpider web site.




Figure 4. A fragment is annotated and colored green when the formula of the fragment is a subset of the molecular formula identified by a library search. This is useful for rapid confirmation of the compound. Both spectra in this figure are the same, the one on the right is expanded in order to show low-abundance fragments that are not compatible with the C₇H₁₄O structure. These ions can be manually filtered out in the Library Editor.




Figure 5. Library Editor saves data in either binary or XML formats. It has a full featured search function and can include alternative names, boiling point, melting point, and user defined features. The spectra shown above has not been filtered and still contains matrix ions.

It is also possible to import spectra from the Find by Chromatographic Deconvolution and the Find by Integration algorithms. It would be useful to include spectra generated by these algorithms in the database if these techniques are frequently used for library searching.




Figure 6. Library curation with relative height and peak filters.

