

Addressing the World Shortage of Helium

Introducing the
Programmable Helium
Conservation Module

GPD Solutions
June 2013



Market Situation

Unreliable supply of helium worldwide and increasing prices have customers seeking alternative carrier gas solutions

ODD SPENDING
Helium Prices Hit the Roof, Leaving Balloon Sales Deflated
 By BRAD TUTTLE @bradtuttle
 citi Like 99
 Log In with Facebook

JULY 26, 2012
C&EN
 CHEMICAL & ENGINEERING NEWS

BUSINESS

INFLATION
 As price of helium rises
 By David Knowlton Wednesday,

HELIUM SUPPLIES ARE SCANT, AGAIN

Global **PRODUCTION SHORTFALL** means the element will be scarce into 2013
 MARC S. REISCH, C&EN NORTHEAST NEWS BUREAU



PHOTO: iStockphoto.com

FOR THE SECOND TIME in five years amid helium consumers are having trouble buying the inert and buoyant element. The situation is dire for some users, and helium suppliers say it is going to get worse before it gets better.

Dean Olson, nuclear magnetic resonance...

BUSINESS

WHY HELIUM IS IN SHORT SUPPLY
 Maintenance disruptions, new PROJECT DELAYS mean shortages for at least two more years

WHEN HELIUM shortages hit the market last year, it was the party helium balloons that first felt the pinch. Since then, suppliers of helium have grown more tight-fisted, pushing out on special contracts for balloons for research, medical, and industrial uses.

THE LATEST CRIMP in the helium supply comes from a plant in Wyoming that is supposed to be producing helium for the first time this year. When that happens over a two-year period this month, it will help ease the global helium supply.

HEALTH CARE is one of the most demanding users of helium. It is used in MRI machines to cool the superconducting magnets. A result of a large-scale and multi-faceted investigation, helium is the most abundant gas in the atmosphere, but it is not as easy to produce as other gases.



U.S. Senate Committee on Energy and Natural Resources
 "The Helium Stewardship Act of 2012"
 Washington, D.C.
 May 10, 2012
 Statement Submitted for the Record
 Airgas, Inc.

The U.S. is experiencing a significant helium shortage that is having an adverse effect on the nation's economy and ability to create jobs. A sufficient amount of helium is available in the U.S. to satisfy domestic demand, however, it is not fairly priced, nor is it reaching the U.S. businesses that need it. At issue is the U.S taxpayer-owned Federal Helium Reserve and the 1996 Helium Reserve Privatization Act.






RESOLVE
YOUR SEARCH
FOR VALUE

Agilent's COMPREHENSIVE Solutions

Carrier Gas Decision Tree



GC and GC/MS Migration to H₂ Carrier Gas

Helping Your Facility Address Issues Related to Helium


Summary of pre-recorded webinar: (September 11, 2012 C&E News)

- **System setup:**
 - H₂ safety, H₂ source, gas connection, system clean up
- **Method Migration:**
 - Method transfer SW, method migration consideration, revalidation
- **Analytical Performance Expectation:**
 - Sensitivity impact, MS spectrum impact, analytes compatibility
- **For More Details:**
 - Listen to recorded session at <http://www.agilent.com/chem/heliumupdate>

Only Agilent provides you with comprehensive solutions

Agilent's **COMPREHENSIVE** Solutions

Carrier Gas Decision Tree



GC Method Migration to H₂ and N₂

Helping Your Facility Address Issues Related to Helium

Summary of pre-recorded webinar: (Feb 12th, 2013 C&E News)


- **Many HPI GC Methods Suited to Nitrogen:**
 - No safety issue, chromatographic resolution suitable for N₂ carrier
- **Method Migration:**
 - Method transfer SW, method migration consideration, revalidation
- **Regulatory Compliance:**
 - ASTM method trends, Agilent's leadership position
- **For More Details:**
 - Listen to recorded session at <http://www.agilent.com/chem/heliumupdate>

Only Agilent provides you with comprehensive solutions



Agilent's **COMPREHENSIVE** Solutions


Carrier Gas Decision Tree



Reducing Helium Use With Conservation


New 7890B Helium Conservation Module

- Automatically switches carrier gas supply to N₂ Standby during idle time
- Integrates into the new 7890B Sleep and Wake function
- Combined with Helium Gas Saver to **GREATLY** reduce helium consumption
- Better alternative to just “shutting off the GC”
 - No system contamination with ambient air exposure
 - Faster re-start of heated zones



Carrier Gas Switch EPC Module


Ensures Business Continuity



- Built on 5th generation EPC
- Fully integrated and controlled by OpenLAB CDS
- Purge channel prevents cross contamination of gases
- Precise pressure control between tank and GC
- In AWAKE mode, switches between gases within 15-30 min for most detectors including MSD


How Does It Work ?

Helium Conservation Mode (Sleep Mode with Nitrogen Carrier)



How Does It Works ?

Normal Operation Mode (Wake Mode with Helium Carrier)



OpenLAB CDS: Configuring Sleep/Wake Operation

Simple, Straight Forward Setup

Agilent 7890B Configuration: Instrument 1

Connection Configuration **Resource Conservation**

Reduce gas and power consumption by setting gas saver and instrument schedule options

Instrument Schedule

Select a schedule that best matches how you use this instrument:

Custom

Schedule

Day	Set Wake Method	Wake Time	Set Sleep Method	Sleep Time
Sunday	<input type="checkbox"/>		<input type="checkbox"/>	
Monday	<input type="checkbox"/>		<input type="checkbox"/>	
Tuesday	<input type="checkbox"/>		<input type="checkbox"/>	
Wednesday	<input type="checkbox"/>		<input type="checkbox"/>	
Thursday	<input type="checkbox"/>		<input type="checkbox"/>	
Friday	<input type="checkbox"/>		<input type="checkbox"/>	
Saturday	<input type="checkbox"/>		<input type="checkbox"/>	

Wake Method: Sleep Method:


Wake to last active method before sleep

Perform a conditioning run on Wake



Performance: No Change in Chromatography After N₂ Carrier Sleep Method (GC/FID)

Day 1 - Original He carrier gas run



Day 2 – First He carrier gas run after overnight N₂ Sleep.M method




Day 3 – First He carrier gas run after overnight N₂ Sleep.M method




14 16 18 20 22 24 Min.

Performance: Pass MS Tune within 15 min after Switching from N₂ to He Carrier (GC/MSD)



Time (min)	Counts of Nitrogen Ion			
	5 mL/min He	Relative to Saturation	2 mL/min He	Relative to Saturation
3	1735168	20.69%	8388096	100.00%
4	1033280	12.32%	4959232	59.12%
5	590080	7.03%	1618944	19.30%
6	354112	4.22%	722944	8.62%
7	228480	2.72%	333696	3.98%
10	56984	0.68%	102576	1.22%
15	9052	0.11%	17080	0.20%

Performance: Switch between H₂ and N₂ Reduces Safety Concerns During Idle Time



Result: MSD system performance **restored quickly** (within 15-20 min) after switching back to hydrogen from overnight nitrogen carrier flow.

Helium Savings Calculator – Single GC Channel

Extend helium supply and lower cost using conservation techniques



Agilent Technologies

Method: ASTM D4815 - Ethanol in Gasoline

Column: PDMS 30m x 0.53mm x 2.65um

GC Flow Conditions

He Carrier Flow (mL/min):	8
He Split flow (mL/min):	70
Gas Saver Flow (mL/min):	20
Gas Saver On (min):	3
Run Time(min.):	20
Gas Volume in Cylinder (L):	8000
Runs per Day:	30
He Cylinder Cost (\$):	300
N2 Cylinder Cost (\$):	60

Parameter	No Conservation	Helium Conservation
Daily He Usage (L)	112	21
He Cylinder Life (days)	71	376
Daily N ₂ Usage (L)	0	24
N ₂ Cylinder Life (days)	0	340
Yearly He Cost (\$)	\$1,537	\$292
Yearly N ₂ Cost (\$)	\$0	\$64
Yearly Total Gas Cost (\$)	\$1,537	\$356

Example

- ASTM Method D4815
 - Widely used to measure ethanol in gasoline
 - Helium cylinder last 2 months under normal operation
- Helium Conservation
 - Helium cylinder life extended to 12 months
 - 4x yearly gas costs per year



Benefits of Helium Conservation

User Programmable for Automated Control

- **Seamless integration**

Fully integrated with 7890B and CDS (OpenLAB, Mustang, Mass Hunter)
Method includes carrier gas ID and set points for compliance and transfer
Easily implemented using new OpenLAB Sleep/Awake functions

- **Greater reliability**

Based on proven 5th generation AUX EPC
Warning from 7890 if set points are not reached
For hydrogen users, nitrogen substitution when GC idle

- **Greater performance**

Purge channel prevents cross contamination of gases
Acts as an intermediate pressure regulator from the tank to inlet EPC to ensure greater analytical precision



Reduce concerns about He availability!
No change to method eliminates need for revalidation!



Programmable Helium Conservation Module

Ensuring Business Continuity





*Thank
You!*

