

A Multiresidue Method for Pesticide Analysis Using an Orbitrap Tribrid Mass Spectrometer and Automatic Background Exclusion

Anastasia Kalli, Seema Sharma, and Dipankar Ghosh, Thermo Fisher Scientific, 355 River Oaks Pkwy, San Jose, CA 95134, USA

ABSTRACT

Purpose: Develop a multiresidue method for pesticides quantitation and improved untargeted screening on the Thermo Scientific™ Orbitrap ID-X™ Tribrid™ mass spectrometer

Methods: Chromatographic separation was performed on a Thermo Scientific™ Vanquish™ UHPLC system using a Thermo Scientific™ Accucore™ aQ C18 column

Mass spectrometric analysis was performed on an Orbitrap ID-X Tribrid mass spectrometer using AcquireX workflow, for automated generation of background exclusion list

Data processing was performed with the Thermo Scientific™ TraceFinder™ software utilizing the mzVault library for pesticide screening

Results: Excellent limits of quantitation, reproducibility, linearity and accuracies were obtained

Automated generation of background exclusion list, using the AcquireX workflow, resulted in a significant increase in the number of identified pesticides, through libraries matches, compared to DDA, at low ppb concentration levels

INTRODUCTION

Pesticides are routinely applied to crops for preventing, destroying or controlling pest activity. In order to protect the consumers and ensure they are not being exposed to pesticide levels harmful for their health, pesticides are regulated and several countries have established maximum residue levels (MRLs). Given the large number of pesticides used and the globalization of the food industry, multiresidue methods offer a great advantage allowing analysis of hundreds of pesticides in a single run. We have implemented a multiresidue method for the analysis of 250 pesticides on an Orbitrap ID-X Tribrid mass spectrometer utilizing AcquireX for automated generation of background exclusion list.

MATERIALS AND METHODS

Sample Preparation

Strawberry samples were obtained from a local retail store. Following homogenization, strawberry samples were extracted using a QuEChERS extraction kit. Briefly 10 g of sample was weighed into the QuEChERS extraction tube and 10ml of ACN was added. Samples were shaken, centrifuged and the supernatant was collected. Matrix-matched standards were prepared by spiking the pesticide standards (250 pesticides) into the extracted matrices at concentration levels ranging from 0.05 to 200 ppb.

Liquid Chromatography

Chromatographic separation was performed on a Vanquish UHPLC system using an Accucore aQ C18 column (100 × 2.1 mm, 2.6 μm). Mobile phase A consisted of 98% water, 2% methanol containing 5mM ammonium formate and 0.1% formic acid. Mobile phase B consisted of 98% methanol, 2% water containing 5mM ammonium formate and 0.1% formic acid. The flow rate was 300 μl/min and the column temperature was set at 25 °C. Analysis time was 15 min including 3 min equilibration time. Injection sample volume was 1 ul.

Mass Spectrometry


Mass spectrometric analysis was performed on an Orbitrap ID-X Tribrid mass spectrometer using AcquireX workflow, for automated generation of background exclusion list, or data dependent acquisition (DDA). The exclusion list was automatically generated from the matrix blank and automatically implemented in the acquisition method for all subsequent analyses to reduce background ion fragmentation. The AGC target value was set at 2e5 for the full MS and at 1e4 for the MS/MS scans. The resolution was set to 60,000 for the full MS scan and at 15,000 for the MS/MS scan events. HCD was performed with stepped collision energy (%) of 20, 40, and 70.

Data Processing

Data processing was performed with the TraceFinder software. Quantitation was performed on the precursor ion and mass tolerance was set to 5 ppm.

RESULTS

Figure 1. Representative calibration curves




Calibration curves for methoprotrolyne, mexacarbamate, thiabendazole and ametryn. Calibration ranges were from 0.1 to 200 ppb for methoprotrolyne and from 0.5 to 200 ppb for mexacarbamate, thiabendazole and ametryn.

Table 1. Obtained % Diff, % RSD and % at LOQ levels for representative pesticide residues

Pesticide Residue	LOQ (ppb)	% Diff Injection 1	% Diff Injection 2	% Diff Injection 3	% RSD	% CV
Ametryn	0.5	-0.07	-6.27	-8.01	4.38	4.48
Carbaryl	1.0	4.14	-6.79	3.43	6.10	5.10
Chloridazon	1.0	-0.13	-3.33	8.17	5.84	5.86
Ciometazone	0.5	0.54	4.46	-9.50	7.31	7.24
Cyanazine	0.5	-4.15	-1.42	-1.71	1.54	1.53
Cyazofamid	0.5	-4.48	-4.40	-7.92	2.13	2.01
Diclotofopos	0.5	3.66	2.87	9.13	3.24	3.62
Dinotefuran	1.0	2.14	8.90	1.67	3.88	4.88
Fensulfotioh	0.5	1.08	-4.20	5.61	4.87	4.74
Fuberidazole	0.5	-1.24	7.81	1.97	4.46	4.47
Hexazinone	0.5	-5.17	-0.35	-2.01	2.51	2.55
Heptonophos	1.0	-0.37	3.72	-0.37	2.34	2.40
Methabenzthiazuron	0.5	3.38	6.83	2.54	2.18	2.13
Metosulam	1.0	-1.21	-1.81	-2.98	0.92	0.72
Ofurace	0.5	6.41	-1.51	-0.11	4.16	4.14
Tebufepyrad	1.0	1.39	-5.94	6.06	6.02	5.60
Thiabendazole	0.5	1.49	-2.71	-5.99	3.84	3.81
Tricyclazole	0.5	-1.04	-4.31	-2.49	1.68	1.73


The % Diff represents the difference between the calculated amount and the expected amount, divided by the expected amount. % RSD is based on the calculated amount and % CV is based on the peak areas.

Figure 2. Obtained LOD (A) and LOQ (B) levels per number of analytes




LOD levels were at or below 1 ppb for 86% of the pesticide examined. LOQ levels were at or below 1 ppb for 82% of the pesticide residues tested.

Figure 3. Obtained % CV at LOQ levels




Obtained % CV values at LOQ levels were within 10% for 94% of all pesticides tested

Figure 4. Obtained % RSD at LOQ levels



Obtained % RSD values at LOQ levels were within 10% for 94% of all pesticides tested

Figure 5. Obtained % CV at all calibration levels for all pesticides residues tested



Obtained % CV were within 5% for the vast majority of the pesticide residues tested

Figure 6. AcquireX Workflow for Background Exclusion






Figure 7. Representative Examples of mzCloud Spectral Library Matches




Spectra library matches for trifloxystrobin at a spiked concentration of 0.5 ppb, coumaphos at spiked concentration of 0.5 ppb, imazali at spiked concentration of 1 ppb and buprofezin at spiked concentration of 1 ppb.

Figure 8. Number of Identified Pesticides with Background Exclusion using Acquire X and Data Dependent Acquisition (DDA)



Automated generation of background exclusion list with AcquireX workflow significantly increases the number of identified pesticides, compared to DDA, at low ppb concentration levels

Figure 9. Results Reproducibility with AcquireX workflow



Number of identified pesticides obtained with automated generation of background exclusion list with AcquireX workflow are highly reproducible. In the data presented above matrix-matched standards were prepared and analyzed on two different days

CONCLUSIONS

The Orbitrap ID-X Tribrid mass spectrometer delivers confident, accurate and reproducible quantitation

LOQ levels for 86% of the pesticide residues tested were at or below 1 ppb

AcquireX workflow enables the generation of background exclusion list in an automated fashion to increase productivity, efficiency and quality of results

Automated generation of background exclusion list with AcquireX workflow significantly increases the number of identified pesticide residues at low ppb levels compared to DDA

TRADEMARKS/LICENSING

© 2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

PO65482-EN0419S