

Poster Reprint

ASMS 2023
Poster number WP 231

An Innovative Approach to Q-TOF High Resolution Accurate Mass Analyte Screening Using an Improved Software Algorithm and Screener Tool

Peter JW Stone

Agilent Technologies, Santa Clara, CA, United States

Introduction

High resolution, accurate mass screening (HRAM) for high numbers of targets using a Q-TOF Mass Spectrometer can be difficult to reproducibly achieve in a single analysis due to the variability of sample matrices, the concentration ranges and the nature of the data mining algorithms.


This poster illustrates a proof of principle ability to detect a panel of over 200+ targets and suspects reproducibly in a single 8-minute analysis at high and low concentrations using a newly improved software mining algorithm SureMass.

The results presented herein this poster demonstrate near 100% identification of a comprehensive panel reproducibly to below 5ppm mass accuracy over a wide range of typical working concentrations.


SureMass Data Mining

What is SureMass?

- A low-level data processing algorithm designed to extract maximum information from accurate mass profile data.



Ridges and Features



Attributes of a Feature:

- Scan space (Collision Energy, Fragmentor Voltage)
- Assigned flight time and Retention Time
- Abundance profile and Flight time profile
- Assigned m/z

SureMass Data Mining (cont'd.)

Why SureMass?

Both Profile and Centroid data formats have individual Issues.

- Speed** - Profile data analysis is painfully slow
- Mass Accuracy** - Centroid data analysis is fast, but gives poor mass accuracy at high abundance levels
- Saturation** - Centroid can't handle Saturation (Mass shifts, Split peaks, False negatives)

SureMass handles saturation gracefully (No high-abundance negatives)

Experimental

Outline.

Repeat injections (N=10) containing 200+ positive polarity targets were injected onto an 8-minute reverse phase chromatography gradient for separation from isobars and interferences that may be present. Column phase was **Poroshell 120 EC C18 100 x 2.1mm**. (Aqueous MP= 0.01% FA, 5mM Amm F).

Pure solvent standard mixes at high and low concentrations relative to minimal and high working concentrations for the respective target classes included in this study were analyzed using an Agilent 6546 Q-TOF LC/MS operating in MS-only or All Ions positive polarity mode to attain the identification results.

Data mining software SureMass was utilized to identify and report low false negative and positive percentage rates via a novel screener software tool.

HPLC Gradient Conditions.

0.00 min 10% Methanol; 1.00 min 15% Methanol; 4.00 min 50% Methanol; 6.00 min 95% Methanol; 8.00 min 95% Methanol. 8 min total run time.

MS Conditions (AJS-ESI Positive Polarity).

Data collection Rate:

- 3x data points/s (MS-only)
- 9x data points/s (All Ions, CE = 0, 20, 30 eV)

Source Settings at HPLC Flow Rate of 0.5 mL/min:

- Nozzle Voltage: 500V
- Sheath Gas Temp/Flow: 360°C; 11L/min
- Dry Gas Temp/Flow: 300°C; 10 L/min
- Nebulizer Pressure: 55 psi
- Capillary Voltage: 3750 V

Screener Tool.

The Screener Tool Graphical User Interface is an integral element of MassHunter Quantitative Analysis Software. It displays each sample in the batch interactively and simplifies the analyst review process by displaying the targets and suspects using a color-coded 'traffic light' system of green, amber and red. Green=positive, amber=review required, Red= negative. Filters are provided to simplify the review process further.




Figure 1. Screener Tool Graphical User Interface.

Mass Match Score %.

The Mass Match Score in Table 1 is derived from the contributions of three scoring elements:

- Monoisotopic Mass Accuracy Score (50%)
- Isotope Spacing Score (25%)
- Isotope Abundance Score compared with Theoretical (25%)
 - Retention time Delta window also used for Identification.
 - Coelution Score available when using All Ions mode of acquisition.

Conclusions

- It has been shown that Comprehensive and large panels of Targets and Suspects can be accurately and reliably screened using SureMass Data Mining combined with the Agilent 6546 Q-TOF LC/MS.
- Mass match Scores are consistently in the high 90% ranges for high and low concentration replicate injections.
- Mass Accuracy measurements for very high and very low concentration ranges within target classes are consistently below 5 ppm.
- Fast data processing speed for SureMass converted data files of large target and Suspect lists.
- Easy and intuitive Screener Tool for fast data review.
- MS-only and All Ions data formats compatible and reliable.

References

1/ Karen E. Yannell,, and Manuel Gomez, Agilent Technologies, Inc., Pub 5994-1744EN

Results

Mass Accuracy Results for 200+ Target Mixes.

- Of the 200+ Targets and Suspects analyzed as part of this proof of principle study, the vast majority exhibited a Mass Match Score of High 90% values for Low and High Concentration Controls (Averaged over 10x different batches.)
- Mass Accuracy measurements were less than 5ppm for All Targets at low and high concentrations.
- MS-only data results are illustrated due to space constraints, however All Ions data yielded comparable results.

Compound	Transition	Low Controls N=10		High Controls N=10		Compound	Transition	Low Controls N=10		High Controls N=10	
		Mass Accuracy Average (ppm)	Mass Match Score Average	Mass Accuracy Average (ppm)	Mass Match Score Average			Mass Accuracy Average (ppm)	Mass Match Score Average	Mass Accuracy Average (ppm)	Mass Match Score Average
Name	Transition					Name	Transition				
25B-NBOMe	380.0856	0.7303	99.4506	0.5289	99.0208	Isobutyrylfentanyl/beta-Methylfentanyl	0.0865	99.1605	0.7147	99.4432	
25C-NBOMe	336.1361	0.9524	98.8240	1.2707	97.6275	Ketamine	1.3008	99.2295	1.3093	98.9875	
25H-NBOMe	302.1751	0.9877	99.5339	0.7183	99.4547	Lacosamide	1.4773	99.1081	1.4790	99.2036	
25I-NBF	416.0517	0.9026	99.5882	0.8628	99.5358	Lamotrigine	0.5244	98.1497	1.2585	98.1517	
25I-NBMD	442.0510	0.8085	99.5329	0.3036	99.6687	Lidocaine	0.9896	99.2873	0.8662	95.3565	
25I-NBOH	414.0561	0.9066	99.4942	0.5042	99.5230	Loperamide	1.4817	95.4532	0.6926	96.0532	
25I-NBOMe	428.0717	0.9317	99.5551	0.4652	99.4947	Loratadine	1.5419	96.2688	1.0261	97.7322	
2-MAPB	190.1226	1.0261	99.6256	1.5285	99.3207	Loxapine	1.1375	99.1810	0.7724	98.6391	
3-Chloroethcathinone	212.0836	1.9641	99.4559	1.1749	99.4396	LSA	-0.1899	91.0280	1.2810	83.8689	
3-Chloromethcathinone	198.0680	1.5290	99.5521	1.4755	98.7582	Maprotiline	1.4441	98.9089	0.8646	99.3495	
3-Fluorofentanyl	355.2180	0.4399	99.5225	0.4362	95.1465	MDA	180.1019	1.4183	99.5409	1.4508	
3-Fluorophenmetrazine	196.1132	1.6006	99.4387	1.0015	99.5006	MDEA	208.1332	1.5558	99.3208	1.5388	
4/5-APDB	178.1226	1.1267	99.6459	1.0252	97.9841	MDMA/Methedrone	194.1176	1.4041	98.8112	1.3169	
4-Chloro-Alpha-PVP	266.1306	1.4243	98.4546	0.8829	99.0547	MDPV	276.1594	1.1180	99.6175	0.8843	
4-Chloromethcathinone	198.0680	1.5290	99.5521	1.4755	98.7582	Meperidine	248.1645	1.4982	99.4144	1.0344	
4-Methoxybutyrylfentanyl	381.2537	1.0189	99.3459	0.9730	97.9662	Mephedrone/6-APDB	178.1226	1.1267	99.6459	1.0252	
5-DBFPV	274.1802	1.1097	99.4470	1.0646	99.1194	Mepivacaine	247.1805	1.4236	99.3452	0.9966	
5-Fluoro-PB-22	377.1660	1.3315	98.9360	1.1144	98.9819	Methadone	310.2165	1.1909	99.6076	0.0431	
6/5-MAPB	190.1226	1.0261	99.6256	1.5285	99.3207	Methamphetamine	150.1277	1.5506	99.6178	1.3213	
6-Acetylmorphine	328.1543	0.4403	98.8440	0.7928	96.4197	Methcathinone	164.1070	1.1943	99.6563	0.9424	
7/4/3-MAPB	190.1226	1.0261	99.6256	1.5285	99.3207	Methiopropamine	156.0841	1.0051	99.4840	0.9507	
7-Aminoclonazepam	286.0742	1.1568	99.3415	0.8375	98.9995	Methoxetamine	248.1645	1.6191	99.4144	1.1553	
7-Aminoflunitrazepam	284.1194	2.0733	79.7982	0.6751	99.0319	Methoxyacetylfentanyl	353.2224	0.0704	98.6097	0.7161	
7-APDB	178.1226	1.1267	99.6459	1.0252	97.9841	Methoxyphenidine	296.2009	1.1864	99.5226	1.0278	
AB-PINACA	331.2129	1.2961	98.7074	0.3699	98.9559	Methylone	208.0968	1.1695	99.2013	1.1699	
Acetaminophen	152.0706	1.3173	96.1147	-0.0152	95.2489	Methylphenidate	234.1489	0.9005	99.5964	0.3965	
Acetylcodeine	342.1700	1.2909	99.3449	0.7533	99.2765	Metoclopramide	300.1473	1.5702	99.0772	0.9129	
Acetylfentanyl/Benzylfentanyl	323.2118	-0.2557	86.9129	-0.3048	76.5476	Mexedrone	208.1332	1.5558	99.3208	1.5388	
Acrylfentanyl	325.2118	0.0969	99.3141	0.6229	95.0780	Midazolam	326.0855	1.0985	99.2322	0.6859	
AH-7921	339.1182	0.3177	99.4097	0.6871	97.5881	Mirtazapine	266.1652	0.8574	99.5545	0.8258	
alpha-Hydroxylprazolam	325.0851	0.3184	99.1008	0.2475	95.9373	Mitragynine	399.2278	1.3999	99.2904	0.8436	
alpha-Hydroxymidazolam	342.0804	0.9431	99.3600	0.6245	98.4171	Modafinil	296.0716	1.0768	97.5970	1.0625	
alpha-Hydroxytriazolam	359.0461	0.3155	99.5744	0.3510	98.6131	Morphine	286.1438	0.2784	99.5352	0.9745	
alpha-Methylacetylfentanyl	337.2274	1.0205	99.3151	0.7602	93.4747	Naloxone	328.1543	0.4403	98.8440	0.7928	
alpha-Methylbutyrylfentanyl	365.2587	0.4324	99.6015	1.0453	97.2639	N-Desalkylflurazepam	289.0538	1.1746	99.6046	1.2005	
alpha-Methylfentanyl	351.2431	0.3290	99.7223	0.6991	98.7657	Nefazodone	470.2317	1.1216	96.8872	0.7430	
alpha-Methylthiofentanyl	357.1995	0.4682	95.2031	0.8197	87.9848	N-Ethylpentylone/N,N-Dimethylpentylone	250.1437	2.2353	97.0007	1.0438	
alpha-PVP	232.1696	2.5272	98.4295	0.6026	98.7350	N-Methyl-Norcarfentanil	305.1859	0.3910	99.5332	0.4108	
Alprazolam	309.0902	0.8933	99.4558	0.4208	99.1752	Nordiazepam	271.0633	1.1876	99.3032	0.7893	
Amantadine	152.1434	0.9194	99.7589	0.4772	98.3491	Nordoxepin	266.1539	1.4257	98.9884	1.2922	
Amitriptyline	278.1903	0.9830	99.5747	0.8207	99.5786	Norflouxetine	296.1257	0.7822	99.5385	0.9228	
Amoxapine	314.1055	-0.2683	99.1872	-1.6177	96.2973	Norketamine	224.0837	1.3322	99.3716	1.2579	
Amphetamine	136.1121	0.6967	99.7644	1.2288	99.3957	Normeperidine	234.1489	0.9005	99.5964	0.3965	
Atomoxetine	256.1696	1.3800	99.0820	0.7558	99.2106	Norpropoxyphene	326.2115	1.0246	99.2507	0.6611	
Benocyclidine (BTCP)	300.1780	1.2402	96.8765	1.3403	96.4169	Nortriptyline	264.1747	1.0025	99.4706	0.9553	
Benzodioxolefentanyl	429.2173	-0.1015	99.8230	0.4187	98.2064	Ocfentanil	371.2129	0.4411	99.1990	0.6794	
Benzoylcegonine	290.1387	1.4494	99.0252	1.0271	96.4646	Ondansetron	294.1601	0.9560	99.5904	0.5222	
Benzotropine	308.2009	0.7651	99.5661	0.6244	99.5250	Orphenadrine	270.1852	1.4838	99.2224	1.1835	
beta-Hydroxyfentanyl	353.2223	0.2969	98.6097	0.9426	98.3199	ortho-Fluoroacrylfentanyl	353.2024	0.1671	99.7847	0.5610	
beta-Hydroxythiofentanyl	359.1788	0.3147	95.9565	0.5345	91.8893	ortho-Fluorofentanyl	355.2180	0.1662	98.5178	0.9839	
Brompheniramine	319.0804	0.5644	94.2589	0.6758	99.3564	Oxazepam	287.0582	-0.0699	97.9141	0.8636	
Bupivacaine	289.2274	1.1032	99.4875	0.7464	99.6823	Oxycodone	316.1543	0.8103	99.7328	1.2426	
Buprenorphine	468.3108	-0.1996	98.6039	0.4300	95.7811	Oxymorphone	302.1387	0.0979	85.4136	0.9646	
Bupropion	240.1150	0.9645	99.4560	0.5109	98.3737	para/meta/ortho-Fluorobutyrylfentanyl	369.2337	0.0421	98.5065	0.7074	
Buspirone	386.2551	1.3212	94.8839	0.7919	95.0845	para/meta-Fluorofentanyl	355.2180	0.1662	98.5178	0.9839	
Butyrylfentanyl	351.2431	-0.0559	99.1605	0.5723	99.4432	para/meta-Fluoroisobutyrylfentanyl	369.2337	0.0421	98.5065	0.7270	
Carbamazepine	237.1022	1.1831	99.4187	1.3417	99.3098	para-Chloroisobutyrylfentanyl	385.2041	0.1884	98.3865	0.4612	
Carfentanyl	395.2329	0.2101	99.2474	0.9409	92.0454	para-Fluoroacrylfentanyl	353.2024	0.1671	99.7847	0.5610	
Cetirizine	389.1626	1.5910	98.0111	0.7802	97.8271	Paroxetine	330.1500	1.1085	99.0307	0.8900	
Chlordiazepoxide	300.0898	1.3581	98.0470	1.2255	98.5859	PB-22	359.1754	1.9540	98.5461	1.8107	
Chlorpheniramine	275.1310	0.5823	99.7127	0.7041	99.1425	PCP	244.2060	1.2132	99.4043	1.0048	
Chlorpromazine	319.1030	1.6810	98.9511	0.9923	98.5398	Pentazocine	286.2165	1.3908	99.4933	1.1017	
cis-3-Methylbutyrylfentanyl	365.2587	0.4324	99.6015	1.0453	97.2639	Pentylone	236.1281	1.4895	99.3200	1.1191	
cis-3-Methylthiofentanyl	357.1995	0.4682	95.2031	0.8197	87.9848	Phenazepam	350.9716	0.5269	99.2149	1.0192	
Citalopram	325.1711	1.1514	98.9534	0.6463	99.5175	Pheniramine	241.1699	1.0795	99.5354	0.6692	
Clomipramine	315.1623	1.3570	98.4549	1.1720	98.2892	Phentermine	150.1277	1.5506	99.6178	1.3213	
Clozapine	327.1371	1.0386	98.9943	0.5687	96.7123	Phenylfentanyl	385.2274	-0.0273	93.8598	0.6260	
Cocaine	318.1700	1.2913	99.3682	0.6681	98.9057	Pregabalin	160.1332	0.5483	99.5620	1.5299	
Codeine	304.1543	1.6622	98.7518	1.0749	99.1163	Propoxyphene	340.2271	1.2738	99.4069	0.5193	
Crotonylfentanyl	300.1594	1.0750	99.5366	1.0536	99.0615	Propriptyline	264.1747	1.0025	99.4706	0.9553	
Cyclobenzaprine	349.2274	0.1884	98.6962	0.6017	97.0149	Pseudoephedrine/Ephedrine	166.1226	1.1219	99.7599	0.7868	
Cyclobenzaprine	276.1747	1.0516	99.6086	0.7257	99.4621	PVB	260.2009	1.2423	99.4275	1.1393	
Cyclopropylfentanyl	377.2587	0.5706	99.5663	0.6944	98.1511	Quetiapine	384.1740	1.0194	95.5122	1.0945	
Cyclopropylfentanyl	349.2274	0.1884	98.6962	0.6017	97.0149	Remifentanyl	377.2071	0.6801	99.6402	0.9869	
Desipramine	267.1856	1.0131	99.5078	0.8047	99.3615	Risperidone	411.2191	-0.2300	94.9033	0.2532	
Despropionylfentanyl (4-ANPP)	281.2012	0.7178	99.6724	0.7512	94.8923	SDB-005	359.1754	1.8704	98.5461	1.7272	
Dextromethorphan	272.2009	1.1780	99.2860	1.0161	99.0647	Sertraline	306.0811	1.0026	99.5621	1.1621	
Diazepam	285.0789	1.5970	99.1514	0.9648	99.0689	Strychnine	335.1754	1.4218	98.7952	1.0414	
Dibutylone	236.1281	1.4048	99.3200	1.0344	99.0859	Suvorexant	451.1644	1.3924	97.6641	1.0314	
Dicyclomine	310.2741	1.3248	99.2659	1.0877	91.3359	Tapentadol	222.1852	1.1988	96.6794	0.7673	
Diltiazem	415.1686	1.1079	94.4260	0.8306	94.5139	Tenocyclidine (TCP)	250.1624	0.9556	99.5439	1.2069	
Diphenhydramine	256.1696	1.3800	99.0820	0.7558	99.2106	Terfenadine	472.3210	1.2938	98.8366	1.1623	
Dothiepin	296.										