

# Flexibility of HRIM-MS Analysis for Targeted Lipid Profiling: Choose Your Own Adventure!



Rachel A. Harris<sup>1</sup>, Kelly L. Wormwood Moser<sup>1</sup>, Frances Carroll<sup>1</sup> | <sup>1</sup>MOBILion Systems, Chadds Ford, PA

## Abstract

- Lipids are challenging to analyze due to the presence of numerous isomeric species, many of which cannot be distinguished via traditional LC-MS approaches.
- The introduction of HRIM provides the analytical separation power to begin to unravel very structurally similar lipid species, including isomers.
- Flow injection and liquid chromatography are complementary front-end sample introduction methods for HRIM-MS that provide high throughput analysis and deeper sample characterization, respectively.

## Methods

- Data acquisition was carried out on a commercial equivalent high resolution ion mobility (HRIM) platform, MOBI<sup>TM</sup> (MOBILion Systems) coupled to a 6545XT QTOF (Agilent Technologies) with a 1290 Infinity II Autosampler (Agilent Technologies) for sample introduction.
- The sample analyzed in these experiments was a total ganglioside extract from Avanti Polar lipids diluted to approximately ~100 ug/mL
- Source settings and HRIM method parameters conserved between LC and FIA workflows





Figure 1. Optimized gradient used in LC-HRIM-MS workflow with an Ascentis Si column.

Figure 2. Experimental system configuration and set-up.

## FIA-HRIM-MS

- 2 min high-throughput analysis
- Quantitation of most abundant species in sample
- Isomers separated via high resolution ion mobility
- Ability to convert arrival time to collision cross section (CCS)
- Fast, reproducible sample profiling

## LC-HRIM-MS

- 30 min HILIC gradient separation
- Higher sensitivity quantitation of low abundance species
- Multidimensional isomer separation via LC and HRIM
- Ability to convert arrival time to collision cross section (CCS)
- Deep sample characterization

Table 1. Qualitative comparison of FIA-HRIM-MS and LC-HRIM-MS workflows that users should consider when deciding whether a given workflow meets their analytical needs.

## Results – HRIM resolves isomeric lipid species and methods are complementary between FIA and LC workflows




Figure 3. FIA-HRIM-MS and LC-HRIM-MS extracted ion mobiligrams of select gangliosides in porcine ganglioside extract. Overlaid extracted ion mobiligrams of six FIA-HRIM-MS replicate injections and six LC-HRIM-MS replicate injections for the gangliosides (A) GD1 d36:1, (B) GT1 d36:1, and (C) GM1 d36:1. Insets are magnifications of the six LC-HRIM-MS replicate injections alone.

## Results – LC-HRIM-MS reduces ion suppression for improved quantitation of low abundance species

Addition of LC to workflow enabled:

- The quantitation of nearly 6 times as many isomeric species as FIA alone
- Lower peak area %RSD per species, on average
- The detection of nearly 4.5 times as many high-quality features in an untargeted feature finding assessment




Figure 4. FIA-HRIM-MS and LC-HRIM-MS detection of the most abundant lipid species in porcine ganglioside extract. (A) In the LC dataset, the % abundance of the most abundant gangliosides, sulfatides, and glycerophospholipids in each class is plotted relative to the most abundant species, PS. The FIA dataset was handled similarly, plotted relative to a sulfatide. Only species with a relative standard deviation (RSD) < 20% in the LC-HRIM-MS datasets are included. The peak areas were summed for species that had multiple detected isomeric forms. \*Species in the FIA-HRIM-MS dataset with a %RSD>20% are noted. (B) The % abundance of the major isomers of GD1, GT1, GM1, GD3-Ac, GD1-Ac, GT3-Ac, GT1-Ac, GQ1-Ac, GQ1.

## Results – Ganglioside CCS in Conformational Space

CCS values for gangliosides were determined as follows:

- MOBILion TuneMix data was acquired at identical traveling wave settings as the lipid data.
- TuneMix data for non-surfing ions was used to generate a third-order polynomial function for CCS calculation.
- Ganglioside arrival times were extracted from the FIA-HRIM-MS dataset and converted to CCS using the polynomial function.




Figure 5. Ganglioside conformational space plot. (A) The CCS-calibrated FIA-HRIM-MS data for the d36:1 and d38:1 species from 12 ganglioside classes are plotted against m/z. (B) RSD values were calculated for the CCS values across three replicates.

## Conclusions

- FIA-HRIM-MS provides high-throughput ( $\leq 2$  min) analysis of the most abundant lipids in a ganglioside extract.
- Coupling HRIM-MS with LC enables deeper characterization of the sample by reducing ion suppression.
- Both workflows utilize HRIM for highly reproducible gas phase separation.
- Most exciting results:
  - CCS values of 40 gangliosides from 12 different classes were determined using FIA-HRIM-MS, all with %RSD < 0.2%
  - Adding LC to the workflow enabled the quantitation of an additional 24 low abundance isomers not detected via FIA