

# Improved Low Mass Transmission Efficiency in High Resolution Ion Mobility (HRIM) – Mass Spectrometry (MS) for Expanded Application Profiles



Joshua K. McBee<sup>1</sup>, Jacob W. McCabe<sup>1</sup>, Kelly L. Wormwood Moser<sup>1</sup>, Bud Buttrill<sup>2</sup>, Nathan Roehr<sup>1</sup>, Daniel DeBord<sup>1</sup> | <sup>1</sup>MOBILion Systems, Chadds Ford, PA; <sup>2</sup>S.E. Buttrill Jr. Scientific Consulting, Palo Alto, CA, USA.

## Introduction

- Current MOBILion MOBIETM hardware provides HRIM across a mass range of 300-3000 m/z.
- Improving Low Mass Transmission (LMT) would expand the breadth of applications MOBIETM can be utilized for.
- We investigated modifications of a MOBIETM platform to enhance LMT performance.
- Utility of improved LMT is demonstrated via Credentialed E. Coli metabolomics LC-HRIM-MS workflow.

## Hardware Modifications

- The SLIM RF frequency was increased from the default of 800 kHz to 1.07 MHz using a MIPS Ultra High Q RF Head (GAA Custom Electronics)
- Simulations demonstrated improved low mass ion confinement with higher RF frequencies.
- Nitrogen drift gas was replaced with 5.0 grade Helium at 2.5 Torr.
- MS data was generated on an Agilent 6545B QTOF MS coupled to an Acquiris SA220 14-bit ADC.

## Methodology




Figure 1. Diagram of hardware modifications to MOBIETM system.

Figure 2. Diagram of SLIM ion path.




Figure 3. Positive and negative spectra of 20% Agilent tune solution showing low mass ion transmission in HRIM mode. The 6545B QTOF was operated in 1700 m/z mode.

## Results




Figure 4. Separation of leucine and isoleucine at 132 m/z. 100pg each of leucine and isoleucine were infused at 5  $\mu$ l/min in 50:50 Methanol:Water.





Figure 5. HRIM achieves baseline separation of QC reference compounds. TOP: EIC of Waters LCMS QCRM (part num: 186006963-1) of all 9 compounds relative to HPLC Retention Time. BOTTOM: Mobiligram showing all 9 compounds (149 m/z – 608 m/z) sorted by mass in HRIM. Verapamil is commercially available as a racemic pair and are baseline separated possibly through a proton-induced chirality mechanism<sup>2</sup>.

References:

- Zheng et al. *Chem. Sci.* 8(11), 2017 pp7724-7736.
- Leritano C. et al. *Angew. Chem. Int. Ed. Engl.* 61(9), 2022.
- Chytil L. et al. *J. Chromatography B* 878(30), 2010.
- Mahieu NG et al. *Anal. Chem.* 86, 2014, pp 9583-9589.
- Di Poto C. et al. *J. Am. Soc. Mass Spectrom.* 32(8) 2021, pp 2072-2080.

## Results




Figure 6. Mobility Aligned Fragmentation (MAF) of Verapamil at 30 eV showing identical fragmentation patterns of both mobility peaks. The MS/MS spectra match well with previously published Verapamil spectra<sup>3</sup>.




Figure 7. Baseline separation of Syn- and Anti- Adenosine conformers. Identification of Adenosine from HILIC-Z chromatography of credentialed E. Coli kit (Cambridge Isotope Labs, part # MSK-CRED-KIT). Standard was mixed in a 1:1 C<sup>12</sup>/C<sup>13</sup> ratio and analyzed via HILIC-Z HPLC chromatography<sup>4</sup>. Mobility alignment of the two peaks confirms isotopologues with 10 carbons each. Identification of Adenosine was confirmed by MAF. The two peaks in the mobiligram represent the Syn- and Anti- conformers of Adenosine. These conformers have been previously observed but are baseline separated for the first time here<sup>5</sup>.

## Conclusions

- The use of Helium and higher frequency RF on the SLIM boards dramatically improves low mass range down to ~100 m/z.
- We demonstrate separation of isomers, conformers and racemic pairs in the mobility domain.
- This improved mass range extends the range of applications for the MOBIETM platform.