

Shimadzu's Recommendations for Dealing with Current Helium Gas Shortage

Shimadzu Corporation
GC/TA Business Unit, Analysis & Measurement Division
Sep 2022

Introduction

- [Introduction](#)
- [Reducing He Consumption](#)
- [Changing the Carrier Gas](#)
- [Summary](#)


SHIMADZU Global Analytical and Measuring Instruments

HOME / Products / Gas Chromatography / Countermeasures and Solutions for Helium Gas Supply Shortages

Countermeasures and Solutions for Helium Gas Supply Shortages

Due to helium gas supply shortages, soaring prices, and other factors, reducing the quantity of helium gas consumed and considering alternative carrier gases have become urgent issues. The following describes functionality for reducing helium gas consumption rates and indicates precautions for switching to a different carrier gas.

GC



Reducing Helium Consumption


- [Reducing Consumption During Analysis](#)
- [Reducing Consumption After Analysis](#)
- [Shutting off the Gas Supply After Analysis](#)

Changing the Carrier Gas

- [Switching to Nitrogen](#)
- [Switching to Hydrogen](#)

Application

GCMS



Reducing Helium Consumption

- [Reducing Consumption During Analysis](#)
- [Reducing Consumption During Standby or After Analysis](#)
- [Estimating Reductions by Ecology simulation](#)

Changing the Carrier Gas

- [Switching to Hydrogen](#)
- [Switching to Nitrogen](#)

Application




Because of the challenges seen by our customers, we have created a dedicated page for [Countermeasures and Solutions for Helium Gas Supply Shortages](#)

Please scan the above QR code to be taken to the webpage seen here.

- Introduction
- Reducing He Consumption
- Changing the Carrier Gas
- Summary

Introduction



Introduction

Alternative carrier gas compatibility table

GC Main Unit	Nitrogen	Hydrogen
Nexis GC-2030	✓	✓
GC-2014	✓	✓
GC-2025	✓	✓
GC -2010 Series	✓	✓
GC -8 Series	✓	-
GC -14 Series	✓	-
GC -17 Series	✓	-

GC Options	Nitrogen	Hydrogen
AOC-30/20	✓	✓
AOC-6000/5000	✓	✓
HS -20 Series	✓	✓
HS-10	✓	✓
TD-30/20	✓	-
Aqua PT6000/PT7000	✓	-
PY-3030D/PY-2020iD	✓	-

GC Detector	Nitrogen	Hydrogen
FID (packed)	✓	✓
FID (Capillary)	✓	✓
TCD (packed/capillary)	✓	✓
FPD (packed)	✓	-
FPD (Capillary)	✓	✓
ECD (packed)	✓	-
ECD (Capillary)	✓	✓
FTD (packed/capillary)	-	-
BID	-	-
SCD	✓	-

In order to use an alternative carrier, your “GC Main unit”, “Options”, and “Detector” must all be compatible. There are also other precautions and limitations. Please refer to the "Changing the Carrier Gas" session for details.

Additional notes:



- GC-8, 14, and 17 are not designed to be safe for use with hydrogen.
- Ar carriers have been used as substitutes for nitrogen in TCD detectors. For other detectors, substitution with nitrogen and hydrogen is recommended.

Reducing Helium Consumption

Reducing Consumption During Analysis



Carrier gas saving function

Shimadzu GC is equipped with a carrier gas saver mode that changes the split ratio after introducing a sample to the column to reduce gas consumption. **This mode is standard on the Nexis GC-2030, GC-2014, GC-2025 and GC-2010 series.**



Backflush function

Backflush function reduces carrier gas consumption by significantly shorten the analysis time. **This is an option for Nexis GC-2030 and GC-2010 series.**






Reducing Helium Consumption

Reducing Consumption After Analysis

Gas Selector

This is a dedicated option for Nexis GC-2030 and GCMS-NX series.

If the Nexis GC-2030 system is equipped with a gas selector, it can automatically switch to an alternative gas after analysis to save a significant amount of He consumption. (A manual gas switching valve is also available for conventional models.)



Reducing Helium Consumption


Shutting off the Gas Supply After Analysis

Automatic start/stop

GC system can be scheduled to automatically startup and shutdown. Prevents waste of helium when not analyzing sample.

This feature is included with LabSolutions GC and GCsolution.

This feature is standard on Nexis GC-2030, GC-2014, GC-2025, and GC-2010 series.




Remote start/stop

The LabSolutions Direct feature allows you to start and stop equipment remotely.

This feature is included with LabSolutions GC and GCsolution.

This feature is standard on Nexis GC-2030, GC-2014, GC-2025, and GC-2010 series.





Reducing Helium Consumption

Utilization of correct purity is key

Shimadzu Super Clean Gas Filter

The Shimadzu Super Clean gas filter adsorbs and removes impurities from the feed gas and purifies it to equivalent of 99.9999% purity. It is more effective to use a gas filter for both carrier gas and detector gas. (The recommended replacement cycle is one year.)

This is an option for Nexis GC-2030, GC-2014, GC-2025, and GC-2010 series.



Using the Shimadzu Super-Clean Gas Filter results in significantly lower detector noise.

Helium Purifier

This is effective for detectors that require high purity helium, such as BID detectors. For example, if there are 10 ppm impurities in the carrier gas (helium purity 99.999%), the helium purifier can be used to remove the impurity concentration to 10 ppb or less.

This is an option for Nexis GC-2030, GC-2014, GC-2025, and GC-2010 series. A power supply is required for use.



Changing the Carrier Gas

- **There are alternate options for helium as a carrier gas**
 - ✓ Hydrogen
 - ✓ Nitrogen
 - ✓ Argon (in some situations)
- **Based on selected, there may be additional considerations that need to be taken prior to use.**
- **Keep in mind, your laboratory needs and safety requirements when choosing an alternative gas!**

Hydrogen is a flammable gas and is not right for every lab. Please remember to take precautions when utilizing hydrogen gas with any instrument. Please refer to the Safety Instructions which is shipped with the Instrument.

Changing the Carrier Gas

Comparison table of nitrogen and hydrogen

In the case of a device configuration in which both nitrogen and hydrogen can be used as alternative carrier gases, consider the following advantages and disadvantages when choosing a carrier gas.

The two gases differ mainly in terms of "separation characteristics" and "safety".

	Nitrogen	Hydrogen
Advantage	<ul style="list-style-type: none"> ◆ No safety concerns compared to hydrogen. ◆ It's often cheaper than hydrogen. 	<ul style="list-style-type: none"> ◆ As the separation ability is superior to that of nitrogen, the analysis speed can be increased. ◆ Because the separation characteristics are close to those of helium, equivalent separation can be obtained without method review in many cases.
Disadvantage	<ul style="list-style-type: none"> ◆ When using carrier gas flow setting from original helium method, it is likely the separation will be bad, and the method needs to be reexamined. ◆ The analysis time is often extended compared to helium to improve the separation. 	<ul style="list-style-type: none"> ◆ Safety considerations and additional measures may be necessary, such as proper exhaust and hydrogen leakage countermeasures. ◆ In the case of a detector using hydrogen, a detector sensitivity variation due to a flow rate variation of the carrier gas (hydrogen) may occur. ◆ Because of its reducibility, it may react with the measured compounds and affect the results.

Changing the Carrier Gas

Alternative carrier gas compatibility table: GC Detectors

✓ in the table means that the carrier can be used for analysis and does not guarantee the same performance as helium.

Please be sure to perform a method validation prior to changing to an alternative carrier gas.

	Nitrogen	Hydrogen	Cautionary Notes
FID (packed)	✓	✓*1	*1: Adjustment of detector H2 flow rate is required due to sensitivity variation. If the column flow rate is high, sensitivity degradation is expected.
FID (Capillary)	✓	✓*2	*2: Adjustment of detector H2 flow rate is required due to sensitivity variation.
TCD (packed/capillary)	✓*3	✓	*3: When changing from He, a significant decrease in sensitivity is expected.
FPD (packed)	✓	-	
FPD (Capillary)	✓	✓*4	*4: Due to sensitivity fluctuation, adjustment of detector H2 flow rate is required.
ECD (packed)	✓	-	
ECD (Capillary)	✓*5	✓*5	*5: When using H2, sensitivity fluctuations are expected; N2 is recommended.
FTD (packed/capillary)	-*6	-	*6: Although it is possible to use only N2 as a carrier gas, a significant decrease in sensitivity is expected. In addition, helium is required for the detector to operate, so we have determined that this is not possible.
BID	-*7	-*7	*7 — Not available in principle
SCD	✓	-	

Please also refer to the annex at the end of this document.

Changing the Carrier Gas

Alternative carrier gas compatibility table: Sample Introduction Options

✓ in the table indicates that the alternative carrier gas can be used for analysis. However, it does not guarantee the same analytical results as the helium carrier gas.

Please be sure to perform a method validation prior to changing to an alternative carrier gas.

	Nitrogen	Hydrogen	Precautions and Auxiliary Gas
AOC-30/20	✓	✓	
AOC-6000/5000	✓	✓	Purge gas for HS/SPME: Only nitrogen can be used.
HS -20 Series	✓	✓	Vial Pressurization/Purge Gas (Trap only): Only Nitrogen can be used.
HS-10	✓	✓	Vial Pressurization: Only Nitrogen Can Be Used
TD-30/20	✓	-	
Aqua PT6000/PT7000	✓	-	
PY-3030D/PY-2020iD	✓	-	Furnace cooling gas: Nitrogen, compressed air can be used. Sample recovery gas: Only nitrogen can be used.

- Many GC option devcie require gas other than carrier gas.
 - In particular, hydrogen can be used as a carrier gas, but in many cases, it cannot be used as an auxiliary gas like as purge or etc.
- The available carrier gas varies also depending on the detector used, so please refer to the detector’s table in next page.

Please also refer to the annex at the end of this document.

Changing the Carrier Gas


Example method adjustment based on alternative carrier gas

Method development of alternative carriers is simplified by using **constant linear velocity control**.



Constant linear velocity control [40 cm/sec]

Constant linear velocity control [25 cm/sec]



- 1: Acetaldehyde
- 2: Acetone
- 3: Ethyl acetate
- 4: Methanol
- Concentration: 1% each
- Solvent: Ethanol

Comparison of separation degree by changing linear velocity

Changing the Carrier Gas

Method development of alternative carriers

EZGC Method Translator from Restek Corporation

Carrier Gas	Original	Translation
Carrier Gas	Helium	Nitrogen

Column	Original	Translation
Length	25.00	25.00 m
Inner Diameter	0.32	0.32 mm
Film Thickness	1.00	1.00 μ m
Phase Ratio	80	80

Control Parameters	Original	Translation
Outlet Flow	1.64	0.58 mL/min
Average Velocity	40.00	18.11 cm/sec
Holdup Time	1.04	2.30 min
Inlet Pressure (gauge)	95.43	38.65 kPa
Outlet Pressure (abs)	101.33	101.33 kPa

Oven Program	Ramp (°C/min)	Temp (°C)	Hold (min)	Ramp (°C/min)	Temp (°C)	Hold (min)
Number of Ramps	250	10	250	22.9	22.1	
1 (1-4)	50	300	5	22.9	300	10.9

Control Method	Original	Translation
Control Method	Constant Linear Velocity	

Results	Original	Translation
Run Time	16.00	35.18 min
Speed		0.45 x

Input the conditions when using He in the "Original" column.

Select Nitrogen as the alternate carrier gas in the "Translation" column.

The optimum parameters are calculated automatically.

Select the priority of the analysis using the "Result" item. The parameters will be optimized for that condition.

The EZGC Method Translator from Restek is useful translating original settings to alternative carrier gases. (You can use it for free.)

After you enter the current method and select an alternate carrier gas type and preferred control method, the conversion method appears based on the result you are trying to meet.



See the Restek Corporation website for more information.
<http://www.restek.com/ezgc-mtfc>

Changing the Carrier Gas

Gas Chromatograph Hydrogen Gas Safety


Hydrogen gas is classified as a dangerous substance and is extremely flammable. When using hydrogen gas, please understand the properties of hydrogen gas and handle it correctly. Please review the below website as well as SDS from supplier.

<https://www.shimadzu.com/an/service-support/technical-support/handling-precautions/gas-chromatography/index.html>




Excerpts of caution

The room in which the instrument is installed should be well ventilated. If a large amount of hydrogen gas leaks into a poorly ventilated room, there is a risk of explosion.



Always have a hydrogen leak check device available and periodically check the hydrogen gas line for leaks. We recommend the use of a hydrogen gas alarm in the room where the instrument is installed.



In addition, the gas chromatograph itself has safety measures for use with hydrogen.


Changing the Carrier Gas

Hydrogen Safety Measures in GC Main unit


SUS piping with high resistance to hydrogen embrittlement is used for the carrier gas flow path.

Copper piping in flowpath has safety risk when using hydrogen as carrier gas since hydrogen embrittlement. In Shimadzu GC, stainless steel material is used for entire flow path, a hydrogen carrier can be used without additional modify to the GC main unit. This is a standard feature of Nexis GC-2030, GC-2014, GC-2010 series, and GC-2025.

all metal inlet parts are SUS



all metal detector parts are SUS





Changing the Carrier Gas

Hydrogen Safety Measures in GC Main unit


Carrier gas leak check function

Gas leakage can occur from the sample vaporization port due to septum deterioration, untightened columns and other reasons. The presence of such gas leakage can be diagnosed with a single tap on the screen.
(This function is only available on Nexis GC-2030)



Pressure/Flow monitoring function with electronically controlled flow controller

The flow controller for carrier gas installed in the GC constantly monitors the pressure flow rate. If excessive gas flow is detected due to column breakage, etc., the flow controller closes the valve to shut off the hydrogen supply.
(This function is standard on the Nexis GC-2030, GC-2014, GC-2025 and GC-2010 series.)




Changing the Carrier Gas

Hydrogen Safety Measures – Optional GC Accessory

Hydrogen sensor

Shimadzu GC's newly designed hydrogen sensor monitors the hydrogen concentration in the GC oven and can detect potential leaks early. When leaks has been detected, it lowers the temperature and automatically switches to a safe standby mode. If the hydrogen concentration rises continuously, the main power is turned off to prevent accidents. **(This is an option for Nexis GC-2030)**



H₂ sensor functions

- Instrument forcibly power off
- Error message shown
- Oven flap opens
- Heaters turn off

Hydrogen sensor


The hydrogen sensor.
Monitoring inside the GC oven

Changing the Carrier Gas

Hydrogen Safety Measures - Optional GC Accessory

Hydrogen Flow limiter


The hydrogen flow limiter controls the excessive flow of hydrogen into the oven in case a column is forgotten to be connected or the flow controller occur uncontrol error (This is an option for Nexis GC-2030, GC-2014, GC-2025, and GC-2010 series.)



Changing the Carrier Gas

Hydrogen Safety Measures - Optional GC Accessory

When hydrogen gas is used as a carrier gas, the room should be well ventilated and the exhaust from the split vent of the AFC and the vents of the TCD and ECD should be exhausted to a ventilator where no retention occurs. We recommend installing a hydrogen gas exhaust option (see below, PN: S221-77174-41) or other equipment that can be connected to the vent area.



If the silicone rubber tubing is not long enough, the customer can prepare tubing with an inner diameter of 6 mm and an outer diameter of 8 mm for connection (up to 20 m).

Changing the Carrier Gas

- Introduction
- Reducing He Consumption
- Changing the Carrier Gas
- Summary

Hydrogen Safety Measures – Gas Generators

A hydrogen gas generator can be a good choice when installing a hydrogen cylinder is difficult. See the link below for more information.

<https://www.peakscientific.com/products/hydrogen/>



Features

- Suitable for flame gas and carrier gas at trace detection limits
- 99.99999%* Purity
- Internal leak detection with automatic shutdown features
- Proven PEM technology to generate hydrogen safely and reliably
- Regenerative PSA dryers to ensure highest level of purity
- Automatic loading pump as standard
- Maintenance limited to replacing de-ionizer cartridge
- Compact, space-saving modular design
- Creates hydrogen on demand, minimal storage of hydrogen in the system
- Combine multiple units for higher flow requirements
- GC in-oven hydrogen leak detector available as an optional extra
- Peak offers a 3 year cell warranty with this generator as standard.

*based on O2 content independently verified by National Physical Laboratory, UK

Changing the Carrier Gas

Applications using Alternate Carrier Gas

Applications for using alternative carrier gas are posted on the Shimadzu website.
<https://www.shimadzu.com/an/service-support/technical-support/technical-information/gas-chromatography/top/index.html>

The screenshot shows the Shimadzu website interface. The top navigation bar includes 'SHIMADZU', 'PRODUCTS', 'INDUSTRIES', 'RESOURCES', 'CASE STUDIES', 'NEWS/EVENTS', and 'SERVICE/SUPPORT'. Below the navigation, there is a search bar and user icons. The main content area is titled 'Application' and lists several articles:

- Example of Use of Alternate Carrier Gas in GC Analysis: Change from He Gas to N2 Gas
- Comparison of Separation Performance with Various Carrier Gases and Introduction of Gas Selector
- Analysis of Residual Solvents in Pharmaceuticals by Water-Soluble Samples Using N2 Carrier (JP18 Supplement II, USP 467)
- Analysis of Residual Solvents in Pharmaceuticals by Water-Insoluble Samples Using N2 Carrier (JP18, USP 467)
- Analysis of Residual Solvents in Pharmaceuticals by Water-Soluble Samples Using H2 Carrier (USP 467)
- Analysis of Residual Solvents in Pharmaceuticals by Water-Insoluble Samples Using H2 Carrier (USP 467)

Below the list, there is a 'Movies' section with a video player showing a thumbnail for 'Converting GC Methods from Helium to Nitrogen Using the EZGC Method Translator'.



The screenshot shows the Shimadzu application news page for 'No. G329'. The page title is 'Example of Use of Alternate Carrier Gas in GC Analysis: Change from He Gas to N₂ Gas'. The main text discusses the use of Helium (He) as a carrier gas and the challenges associated with its price and availability. It introduces the EZGC[®] Method Translator software, which is used to calculate the optimum linear velocity for different carrier gases. The page includes a graph showing the relationship between constant linear velocity and HETP for He and N₂. It also features a screenshot of the EZGC[®] Method Translator software interface, which allows users to input conditions and select the priority of the analysis. The page concludes with a section on the 'Use of Gas Selector' and a reference to Application News No. G328.

Application News
No. G329

Example of Use of Alternate Carrier Gas in GC Analysis: Change from He Gas to N₂ Gas

Helium (He) has mainly been used as the carrier gas for gas chromatography (GC), but limited availability and a sharp rise in the price of He have become issues. In recent years, hydrogen (H₂) and nitrogen (N₂) have been actively used as alternate carrier gases for He. As an example of an analysis in which the carrier gas was changed from He to N₂, this article introduces an example of an analysis by the vegetable oil (olive oil) immersion method for plastic materials using GC-FID. In this analysis, we used the overall migration test referred to the European Committee for Standardization, EN1186 "Materials and articles in contact with foodstuffs - Plastics (No 10/2011)".

Points to Note when Changing Carrier Gases

Fig. 1 shows an example of the characteristics of HETP (height equivalent of one theoretical plate), which shows the separation efficiency of a column when using He, N₂, or H₂ as the carrier gas.

Calculation of Optimum Linear Velocity by EZGC[®]

In this experiment, the optimum conditions when using different carrier gases were studied using the method conversion program EZGC[®] Method Translator, which is available from Restek Corporation via the internet. Fig. 3 shows the screen of EZGC[®] Method Translator. When the conditions for He are input in the "Original" column and the alternate carrier gas type is selected under "Translation" column, the optimum parameters for the alternate carrier gas are displayed. It is also possible to select the priority of the analysis by using the "Result" item. For example, select "Speed" when you wish to shorten the analysis time, and select "Efficiency" when separation efficiency is the priority. Because the focus in this experiment was separation of C18:0 and C18:1, the method was translated by selecting "Efficiency," which prioritizes separation efficiency, and the program calculated the optimum constant linear velocity of 18.11 cm/s.

Use of Gas Selector

Gas selector (P/N: S221-84916-41), which is a dedicated option of Nexis[®] GC-2030, not only changes the gas flow lines, but also makes it possible to switch between two types of carrier gas from the software. Use of the gas selector in combination with EZGC[®] enables easy calculation and control of N₂ for low-cost analyses which do not require high resolution, and He or H₂ when high resolution, high speed analysis is necessary. For details concerning the gas selector, please refer to Application News No. G328.

Summary

Recommendation: Reducing helium consumption

- Try carrier gas save mode.
- Use gas selector.
- Use the automatic start/stop function.
- Be mindful using low purity helium gas.

Recommendation: Changing the carrier gas

- Options for alternative carrier gases are based on not only the main unit but detectors and sample introduction accessories.
- Consider the advantages and disadvantages of nitrogen and hydrogen before selecting an alternative carrier gas.
- Implement necessary safety measures when using hydrogen.



Please be sure to read the "Safe Use of Hydrogen Gas" on our website before using hydrogen.

<https://www.shimadzu.com/an/service-support/technical-support/handling-precautions/gas-chromatography/index.html>

Summary

Compatibility table for Shimadzu GC model with recommended functions

	Helium gas reduction						Operation of alternative carriers						
	Carrier gas save function, backflush (P5)	Gas Selector (P6)	Manual changeover valve *Similar to gas selector (P6)	automatic shutdown automatic startup (P7)	remote shutdown remote start (P7)	Gas filter and He Purifier (P8)	EZGC Method Translator (P15)	Stainless steel piping for main body (P17)	Carrier gas leak check function (P18)	pressure flow monitoring function (P18)	Hydrogen sensor (P19)	H2 Safety options (P20,P21)	Hydrogen generator and hydrogen leak detection function *Similar to hydrogen sensor (P 22)
Nexis GC-2030	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
GC-2014	✓	-	✓	✓	✓	✓	✓	✓	-	✓	-	✓	✓
GC-2025	✓	-	✓	✓	✓	✓	✓	✓	-	✓	-	✓	✓
GC-2010 Series	✓	-	✓	✓	✓	✓	✓	✓	-	✓	-	✓	✓
GC-17 Series	✓	-	✓	✓	-	✓	✓	✓	-	✓	GC -8, GC -14 and GC -17 do not have additional safety features and are not recommended to be used with hydrogen carrier.		
GC-14 Series	-	-	✓	-	-	✓	✓	✓	-	-			
GC-8 Series	-	-	✓	-	-	✓	✓	✓	-	-			

Summary

Compatibility table for software programs with recommended functions

	Helium gas reduction related functions				Alternative Carrier Related Functions		
	Carrier gas save function (P5)	Gas Selector (P6)	Auto Start automatic shutdown (P7)	remote start remote shutdown (P7)	Carrier gas leak check function (P18)	pressure flow monitoring function (P18)	Hydrogen sensor (P19)
LabSolutions	✓	✓ (Nexis GC -2030 only)	✓	✓	✓ (Since this function is part of the GC main unit, it is available regardless of software.)		
GCsolution	✓	-	✓	-			
OpenLab	✓	✓ (Nexis GC -2030 only) GC Driver Version 2.20 or later)	-	-			
Empower	✓	✓ (Nexis GC -2030 only, GC Driver Ver 3.20 or later)	(Automatic shutdown is not supported, but similar functionality can be achieved by using in the shutdown method at the end of a sample batch/list.)	-			
Chromeleon	✓	✓ (Nexis GC -2030 only, GC Driver Ver 2.20 or later)		-			

LabSolutions support models:
 GCsolutions support models:
 OpenLab supported models:
 Empower compatible models:
 Chromeleon support models:

Nexis GC -2030, GC -2014, GC -2010 series, GC -2025
 GC -2014, GC -2010 series, GC -2025, GC -17 series
 Nexis GC -2030, GC -2014, GC -2010 series
 Nexis GC -2030, GC -2014, GC -2010 series
 Nexis GC -2030, GC -2014, GC -2010 series

(GC -2014 and GC -2010 series can be controlled using the Thermo Fisher Scientific driver included with Chromeleon TM 7.)
 (GC -2030 is the only GC that the Chromeleon GC driver developed by Shimadzu can control.)

Thank you very much.

The content of this video is based on the latest information available at the time of publication.

ClickTek and Nexis are trademarks of Shimadzu Corporation or its affiliates in Japan and other countries.

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, Whether or not they are used with trademark symbol "TM" or "®".