

# Automating component detection of small molecules in complex mixtures using HRAM Q-TOF data

**ASMS 2019**

Simon Ashton<sup>1</sup>; Kirsten Hobby<sup>1</sup>; Alan Barnes<sup>1</sup>; Neil Loftus<sup>1</sup>  
<sup>1</sup>Shimadzu Corporation, Manchester, United Kingdom

## Automating component detection of small molecules in complex mixtures using HRAM Q-TOF data

### Overview

- Applying high resolution accurate mass (HRAM) QTOF analysis for suspect screening of small molecules in complex matrices and automating component detection.
- The algorithm automatically locates components that behave as chromatographic features.

### Introduction

To accelerate small molecule component detection in complex mixtures a novel algorithm has been developed for high resolution accurate mass Q-TOF data sets. As with other techniques, the algorithm locates ions that behave as a recognized chromatographic feature (ion intensities rise and fall in abundance in a covariant manner) but differs in its ability to qualify spectral interpretation by identifying isotopes of various charge states in complex interlaced

mass spectra. It can also identify and correlate mass spectrometric signals from adducts, neutral losses and alternate charge states to give a simplified “single component” output when multiple ionized species are present. This algorithm was evaluated using a series of HRAM Q-TOF data sets including food safety and toxicology screening.

### Materials and Methods

The algorithm for component detection was applied to complex mixture analysis in food safety and drugs of abuse. Both matrices included a QuEChERS extract and the algorithm was applied to a generic HRAM LC-MS/MS method (Nexera LC and LCMS-9030 QTOF system, Shimadzu Corporation) as shown in figure 1. Components were typically separated using a Restek Raptor ARC18 or Restek Biphenyl column depending on

the test panel. The component detection stage requires minimal parameters to function effectively, and definition of typical chromatographic peak width is the only key parameter needed. Two parameters define acceptance criteria for the ion intensity in a mass spectrum and the sum of the ions identified in a component.

## Automating component detection of small molecules in complex mixtures using HRAM Q-TOF data




Figure 1. The component detection algorithm was applied to analysis of a complex QuEChERS extract in food safety and drugs of abuse using data acquired by the LCMS-9030 QTOF system. The LC-MS/MS method used in both test panels included full scan MS followed by sequential MS/MS scans. The algorithm was applied to the MS full scan data. (The mass range used for food safety was 140-900 Da and 100-500 Da for drugs of abuse screening).

## Results

### Food safety screening HRAM DIA-MS/MS data acquisition

Untargeted MS and DIA-MS/MS in Food Safety Screening  
TIC for a food commodity matrix extract spiked with a panel of pesticides  
Restek ARC column  
QC sample 0.1 mg/kg

HRAM QTOF mass chromatograms  
Food Safety screening  
Pesticide test panel  
Mass chromatograms for each target pesticide




Figure 2. The component detection algorithm was applied to analysis of a complex QuEChERS extract of an apple matrix spiked with 212 pesticides at a concentration of 0.1 mg/kg corresponding the QC concentration. Default parameters for threshold and a peak width (FWHM of 0.05 minute) were applied and the algorithm detected 3100 grouped components (including adducts, neutral losses and alternate charge states) in the extract.

## Automating component detection of small molecules in complex mixtures using HRAM Q-TOF data

Find - MatrixCalibnStd\_0-100mg-kg - 3848 components found - 3100 groups - (Screening Mode - Find\_TargetsList\_Pesticides\_Biphenyl\_10ppm.xlsx)

#	RT	m/z	Response	Hit #	Target Name	Target Formula	Target m/z	Target RT
740	3.865	250.01640	555845	1	Clothianidin	C6H8CIN5O2S	250.01600	3.872
747	3.869	263.01815	1777891	1	Demeton-S-methyl sulpho...	C6H15O5PS2	263.01713	3.874
786	3.956	165.10260	753261	1	Fenuron	C9H12N2O	165.10224	3.963
791	3.967	238.08471	1322191	1	Dicrotophos	C8H16NO5P	238.08389	3.969
814	4.007	255.13418	193787	1	Carbofuran-3-hydroxy	C12H15NO4	255.13393	4.016
816	4.008	292.02693	805955	1	Thiamethoxam	C8H10CIN5O3S	292.02656	4.015
853	4.080	218.10419	971152	1	Pymetrozine	C10H11N5O	218.10364	4.082
863	4.104	192.07744	1255210	1	Carbendazim	C9H9N3O2	192.07675	4.105
883	4.118	230.00743	1028657	1	Dimethoate	C5H12NO3PS2	230.00690	4.122
897	4.151	278.05728	443304	1	Sulfoxaflor	C10H10F3N3O5	278.05694	4.158
903	4.168	265.04293	1118282	1	Ethidimuron	C7H12N4O3S2	265.04236	4.174

### Data workflow

#### Acquire data using MS and DIA-MS/MS method (Figure 1)

MS 140-900 Da (40 msec; mass range m/z 140-900) followed by 38 DIA-MS/MS sequential mass scans (each DIA-MS/MS mass scan 20 msec; MS/MS mass range m/z 65-900; precursor ion isolation 20Da);

#### Open the component detection research application

#### Select default intensity threshold and chromatographic peak width

#### Run the component detection algorithm for the QC sample

Using a low response threshold a high number of components were detected and grouped together (this will detect all components including matrix derived ions from the QuEChERS extract).

#### Search for a target list of pesticides


#### Open an Excel spreadsheet with a target list of pesticides (the search criteria include name, formula, accurate mass and if available the retention time).

From the list of targets submitted in the suspect screening search, 207 components were matched by accurate mass (within 5 ppm; below m/z 200 the mass tolerance was set to <1 mDa) and retention time (within 0.25 minute). 5 compounds were not detected based upon default chromatographic peak width and signal to noise setting.

## Automating component detection of small molecules in complex mixtures using HRAM Q-TOF data

### Toxicological screening HRAM DIA-MS/MS acquisition

Untargeted MS and DIA-MS/MS in Drugs of Abuse (DoA) Screening  
TIC for a plasma sample spiked with a panel of DoA compounds  
Restek Biphenyl column  
Plasma concentration 50 ng/mL



HRAM QTOF mass chromatograms  
Drugs of Abuse (DoA) screening  
DoA Test panel  
Mass chromatograms for each DoA

Mass chromatograms rescaled for  
Morphine  
Benzoylcegonine  
Norcocaine  
(Factor of 0.1)




Figure 3. TIC and extracted mass chromatograms for a panel of 45 drugs of abuse targets spiked at 50 ng/mL into a human plasma extract. To find a list of DoA targets, the component detection algorithm was used to deconvolute the data file using default search parameters.

## Automating component detection of small molecules in complex mixtures using HRAM Q-TOF data

Find - Plasma\_50ng-mL\_DIA-mtd-50Hz-P\_003 - 7267 components found - 6223 groups - (Screening Mode - Target-List\_DoA.xlsx)

Sample: 3  
MS Event: 1  
Peak width (FWHM): 0.05 min  
Charge state: 1 - 3  
Thresholds  
Spectral intensity: 300  
Response: Very low  
User defined response: 1000

#	RT	m/z	Response	Hit #	Target Name	Target Formula	Target m/z	Target RT
3369	5.990	324.20690	4170403	1	LSD	C20H25N3O	324.20704	5.982
3410	6.091	284.11924	4103837	1	7-aminoflunitrazepam	C16H14FN3O	284.11937	6.083
3827	6.533	310.14113	1734742	1	Fluoxetine	C17H18F3NO	310.14133	6.534
3869	6.615	327.13690	2793447	1	Clozapine	C18H19ClN4	327.13710	6.602
3971	6.700	310.14168	3082	1	Fluoxetine	C17H18F3NO	310.14133	6.534
3975	6.706	337.22732	2911350	1	Fentanyl	C22H28N2O	337.22744	6.707
4220	6.916	244.20591	1794470	1	Phencyclidine	C17H25N	244.20598	6.919
4251	6.958	388.15855	2902008	1	Flurazepam	C21H23ClFN3O	388.15864	6.959
4309	7.034	340.22705	1228430	1	Dextropropoxyphene	C22H29NO2	340.22711	7.035
4433	7.180	321.01934	168657	1	Lorazepam	C15H10Cl2N2O2	321.01921	7.173
4597	7.338	287.05804	220873	1	Oxazepam	C15H11ClN2O2	287.05818	7.328
4703	7.421	316.04811	292206	1	Clonazepam	C15H10ClN3O3	316.04835	7.410

### Data workflow

Acquire data using a MS and DIA-MS/MS method (as in Figure 1; changes applied to the mass range)

MS 100-500 Da (20 msec; mass range m/z 100-500) followed by 19 DIA-MS/MS sequential mass scans (each DIA-MS/MS mass scan 20 msec; MS/MS mass range m/z 40-500; precursor ion isolation 20Da);

Open the component detection research application

Select default intensity threshold and chromatographic peak width

Run the component detection algorithm for the 50 ng/mL spiked plasma sample

The plasma sample was prepared as a protein crash (using acetonitrile) with QuEChERS clean-up. The component detection algorithm found 6223 grouped compounds.

Apply a target list of DoA; the forensic toxicology screening library of 1278 compounds was filtered to search for classes of drugs of abuse (the library was filtered to 153 targets)

Open an Excel spreadsheet with a target list of drugs of abuse (the search criteria include name, formula, accurate mass and if available the retention time).

Using a filtered library of drugs of abuse compounds all 45 targets were matched by accurate mass (within 5 ppm; below m/z 200 the mass tolerance was set to <1 mDa) and retention time (within 0.2 minute of the target list entries in the spreadsheet data base).

## Automating component detection of small molecules in complex mixtures using HRAM Q-TOF data

### Suspect screening in complex mixtures

Using a series of HRAM Q-TOF data sets for food safety and drugs of abuse screening, the component detection algorithm was used to locate covariant ions using a two-step process, 'search' followed by 'validation'.

In the first step, a set of processed mass spectral data is created by using a 'moving FWHM window' based on the user defined chromatographic FWHM value. This creates a set of combined scans over the whole chromatogram. Masses from each combined scan are then correlated with equivalent masses from adjacent combined mass spectra to link data across the chromatographic space. A peak identification method is applied to locate the apex of each proposed chromatographic component using the spatially linked masses. A background ion detection algorithm is then applied to the chromatographic data so that only chromatographically relevant components reported further. In parallel to locating covariant ions, a spectral

interpretation process is applied where isotopic clusters are found. Finally, the search process reports a list of proposed monoisotopic masses and approximate peak apex times.

In the second step, the algorithm further validates the results of the search process to accurately locate the chromatographic component in the raw data. Spectra are combined over the identified chromatographic peak, and then subject to spectral interpretation, to give the most accurate interpretation of isotopic distribution for each component. The validation process improves the accuracy of the component detection process and eliminates false positive reporting. Finally, specific species such as adducts, neutral losses and alternate charge states are identified and correlated so that the monoisotopic mass of the protonated species is reported as the component mass.




Figure 4. Mass spectrum for fentanyl processed by the component detection algorithm and the corresponding unprocessed mass spectrum. The component detection algorithm is designed to report a simplified "single component" output when multiple ionized species are present and the most accurate interpretation of isotopic distribution for each component.

## Automating component detection of small molecules in complex mixtures using HRAM Q-TOF data

### Conclusions

- A novel component detection algorithm has been applied to accelerate small molecule detection in complex mixtures for high resolution accurate mass Q-TOF data sets.
- Using existing HRAM MS and DIA-MS/MS methods, accurate masses of the ions detected and identified by the component detection algorithm were compared to the exact masses of compounds in a pesticide database and in a forensic toxicology database that was filtered for drugs of abuse targets.
- The algorithm can be applied to targeted and untargeted analysis of complex mixtures.

First Edition: July, 2019



Shimadzu Corporation  
[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.