

## Molar Mass Determination of Collagen Peptides

Establishing a cost-effective, easy-to-use QC method for the accurate, reliable, and precise determination of collagen peptide molar mass distributions



### Scope of work

Collagen peptides (also referred to as hydrolyzed collagen or collagen hydrolysate) are produced from natural proteins by hydrolysis. They are soluble in cold water, are highly digestible, and are used in dietary supplements and functional foods.

Collagen peptides are relatively small proteins, with typical molar masses of less than 10,000 g/mol. They comprise at least two (and at most, 100) amino acids. Common to nearly all macromolecules, collagen peptides do not exhibit only one specific molar mass, but a molar mass distribution (MMD). This MMD influences their properties.

GPC/SEC is an established technique for measuring MMDs of macromolecules, such as collagen peptides, but it suffers from two limitations:

- It is a relative technique that allows determination of close-to-true molar masses only if:
  - reference materials chemically matching the analyte structure are available
  - or if advanced detectors (such as online viscometers or online light scattering detectors) are applied.
- It is a liquid chromatographic technique that requires highly qualified and trained users, especially when advanced detectors, such as online viscometers or online light scattering detectors, are used.

The purpose of this study was to develop a simple, cost-effective, and robust high-precision GPC/SEC method for QC, enabling the identification of close-to-true molar masses for collagen peptides. The method should be applicable in any laboratory, allowing reliable interlaboratory comparison. Therefore, the use of molar-mass-sensitive detectors should be avoided.


The single steps in the study were to:

- Identify collagen peptides that are suitable as reference materials, and determine their characteristics
- Establish a suitable calibration procedure
- Perform a round-robin test of the final approach

## Introduction to GPC/SEC

GPC/SEC is a chromatographic technique in which a sample is separated according to the size of the molecules. The separation takes place in columns filled with porous particles. Molecules larger than the pores of the particles are excluded from the pores, and remain in the flowing eluent stream, eluting first from the column.

Molecules smaller than the pores diffuse in and out of the pores, resulting in later elution with decreasing size.<sup>1</sup>



## Calibration in GPC/SEC

To determine the required molar mass information by GPC/SEC, a calibration curve can be applied to relate the elution volume to the molar mass of the eluting macromolecules. The most common approach to constructing the calibration curve is injecting a series of chemically identical, narrowly distributed standards of known molar masses. The elution volumes at peak maxima are plotted against the molar masses of the standards, resulting in the desired calibration curve.

However, GPC/SEC separates by the size of the molecules in solution. Because at a given molar mass, the sizes of macromolecules of different chemical structures differ from each other, the molar mass determined by GPC/SEC depends on the chemical structure of the standards applied for calibration. In order to derive the true molar masses, the chemical structure of the analyte and the calibration standards must be identical. If materials with narrow molar mass distributions are not available, broadly distributed reference materials of known molar mass having the same structure as the analyte can also be applied. Alternatively, molar-mass-sensitive detectors (such as light scattering detectors or viscometers) can be used.<sup>1,2</sup>

## Broad standard calibration in GPC/SEC

Benoit, Rempp, and Grubisic<sup>3</sup> empirically showed that polymers eluting at the same GPC/SEC elution volume have identical sizes (hydrodynamic volumes).<sup>2</sup> The relation between chemically different macromolecules eluting at the same elution volume can be expressed as:

**Equation 1.**

$$M_2 = A \times M_1^B$$

A and B are parameters that depend on the Mark-Houwink parameters of the two polymers in the eluent. Usually, A and B, which relate the molar mass of the narrowly distributed standards ( $M_1$ ) to the molar mass of the analyte at the same elution volume ( $M_2$ ) are unknown, and need to be determined.

For this purpose, a calibration curve is first established, using standards that might differ in chemical structure from the analyte (base calibration). In the next step, one or more broadly distributed reference samples, for which one or more molar mass averages are known ( $M_n$  or  $M_w$ , target values) and which are of the same chemical structure as the material to be analyzed, are run on the column set to be calibrated. For each chromatogram of the reference samples, the molar mass averages are calculated using a test set of A and B to convert the previously established base calibration. The relations between the test parameters A and B, the chromatogram ( $S(V)$ ), and the calculated average molar masses are given by:

**Equation 2.**

$$M_{n,i}(A,B) = \frac{\sum S_i(V)}{\sum S_i(V)/M_2} = \frac{\sum S_i(V)}{\sum S_i(V)/(A \times M_1^B)}$$

**Equation 3.**

$$M_{w,i}(A,B) = \frac{\sum S_i(V) \times M_2}{\sum S_i(V)} = \frac{\sum S_i(V) \times A \times M_1^B}{\sum S_i(V)}$$

Here,  $M_{n,i}(A,B)$  and  $M_{w,i}(A,B)$  are the number and weight average molar masses for chromatogram i, when applying the parameters A and B.  $M_1(V)$  is the molar mass of the base calibration at elution volume V. The calculated molar mass averages depend on the selection of parameters A and B. Agreement between the calculated molar mass averages

and the target values will result only if A and B are chosen correctly. Applying suitable optimization algorithms allows varying A and B, until the best agreement between the calculated molar mass averages for the chromatograms of the broadly distributed reference samples with the target values is achieved. The identified optimal set of A and B can then be applied to convert the base calibration into a calibration curve representing the chemical structure of the analyte.

## Molar mass determination by GPC/SEC with advanced detection

The broad standard calibration approach, as described, requires a set of broadly distributed reference samples of known true molar masses. These molar masses can be determined using GPC/SEC with either online viscometry or online light scattering detection.

### GPC/SEC with viscosity detection

Benoit, Rempp, and Grubisic empirically showed that polymers eluting at the same GPC/SEC elution volume have identical products of intrinsic viscosity,  $[\eta]$ , and molar mass ( $[\eta] \times M$ ).<sup>3</sup>

An online viscometry detector attached to a GPC/SEC instrument provides, in combination with a concentration detector, information about the intrinsic viscosity at each point of the elution curve.

By running a series of narrowly distributed standards of known molar masses, and determining their intrinsic viscosities using the viscosity detector, a universal calibration curve can be established by plotting the product  $[\eta] \times M$  as a function of elution volume. The chemical structure of the standards applied for setting up the universal calibration curve can differ from the analyte's chemical structure.

For determining the analyte's molar mass, its intrinsic viscosity at each elution volume is derived using the viscosity detector. By taking the ratio of  $[\eta] \times M$  from the calibration curve and the experimentally derived intrinsic viscosity of the eluting sample at each elution volume, the true molar mass of the eluting macromolecule is calculated, allowing determination of the true molar mass averages of the sample.

### GPC/SEC with light scattering detection

Other than GPC/SEC with viscosity detection, static light scattering is one of the few absolute techniques to measure true molar masses.<sup>4</sup> Attached to a GPC/SEC instrument, an online light scattering detector provides, in combination with a concentration detector, information about the sample molar mass at each point of the elution curve. Calibration with reference materials is thus not required. The determination of system-related constants and parameters alone are sufficient.

However, light scattering is a technique generally more applicable to molecules of higher molar masses. If molecules of lower molar mass are to be analyzed, the concentration needs to be enhanced to achieve suitable signal intensity. The enhanced concentrations might result in column overloading, which in turn will reduce chromatographic resolution. The reduced resolution affects determination of the molar mass distribution, resulting in incorrect determination of molar mass averages.

In addition, molar mass determination by light scattering requires the accurate knowledge of the specific refractive index increment,  $dn/dc$ , which depends on the composition. The applied  $dn/dc$  value needs to be constant over the eluting peak, and accurately known. Otherwise, severe errors might result for the molar masses determined.<sup>5</sup>

### Suitable options for collagen peptides

For the characterization of future collagen peptides, reference materials for online viscometry were preferred over light scattering detection for two reasons:

- Collagen peptides have relatively low molar masses, and while both detectors are better with higher molar masses, viscometers often have a better signal quality at low molar masses, specifically when working at typical low GPC/SEC concentrations.
- Collagen peptides are not uniform with respect to amino acid composition; thus  $dn/dc$  might not only vary between different collagen peptide samples, but also across the peak of the sample. While the accuracy of GPC/SEC-light scattering results is highly influenced by the accuracy of the  $dn/dc$ , this is not the case for GPC/SEC viscometry results.

## Experimental

### Samples and mobile phase

**Dextran/pullulan:** Dextran and pullulan standards, with a molar mass at the peak maximum of  $M_p = 180$  g/mol up to  $M_p = 106,000$  g/mol were applied to establish the universal calibration curve.

**Collagen fragments:** Collagen fragments (CNBr peptides) for establishing the base calibration for the broad calibration approach were purchased from the Research Institute of Leather and Plastic Sheetting (Freiberg, Germany). The standard mixture contains molar masses between 3,500 and 38,000 g/mol, certified by the manufacturer. These materials exhibit chromophores for UV-Vis detection and are generally suitable for GPC/SEC calibration. However, the lack of knowledge on the exact amounts of each component prevents preparing accurately known concentrations (as is required for determining intrinsic viscosities to establish a universal calibration curve).

**Collagen peptide samples:** Five different samples, 13R078, 13R079, 833593, 891600 and 890435, were obtained from different manufacturers and used as reference materials after molar mass determination. The true molar masses of these samples were determined by GPC/SEC with online viscosity detection.

Six additional samples (sample A to F) were used in a round robin test to identify the suitability of the developed calibration approach.

**Mobile phase:** Phosphate buffer pH 5.30 with 0.2 M NaCl. (13.27 g  $KH_2PO_4$ , 0.44 g  $Na_2HPO_4$ , 11.69 g NaCl in 1.0 L de-ionized water) was applied as GPC/SEC eluent. All salts were obtained from Sigma Aldrich in p.A. purity. Water was de-ionized in-house using an ion exchanger (Behr, Germany).

### Instrumentation and software

The Agilent GPC/SEC system was equipped with an Agilent isocratic pump, a four-channel online degassing device, a column thermostat, and an autosampler. The detector system in this study was a combination of an Agilent variable wavelength UV detector ( $\lambda = 214$  nm), a viscometer detector, and an Agilent refractive index detector (RI). Data acquisition and analysis was performed using Agilent WinGPC.

Only the viscosity and the RI detector were applied to determine the properties and characteristics of the collagen peptide reference materials, as the pullulan and dextran samples that are required for universal calibration are not UV-Vis active.

For broad standard calibration with the collagen peptide reference materials and the collagen fragments, only the UV-Vis detector was used.

All measurements were carried out at a temperature of 40 °C and a flow rate of 0.5 mL/minute. A small amount of benzoic acid (internal flow marker) was added to the solvent to prepare the sample solution. At least two concentrations were prepared for each sample. After dissolving overnight, each sample solution was analyzed in duplicate without further pretreatment.

**GPC/SEC viscometry:** Typical sample concentrations were 3 to 7 g/L, and an injection volume of 100  $\mu$ L was applied.

**Broad standard calibration:** Typical sample concentration was 2 g/L, and an injection volume of 20  $\mu$ L was applied.

### Separation column

Agilent PROTEEMA columns (5  $\mu$ m particle size, 100 Å porosity, 8.0 mm id  $\times$  300 mm) and Agilent precolumn (8.0 mm id  $\times$  50 mm), Mainz, Germany.

## Results and discussion

### Characterization of collagen peptides using GPC/SEC-viscometry-RI


Monodisperse monosaccharides and narrowly distributed pullulans and dextrans with given  $M_p$  values were used to establish the universal calibration curve. The intrinsic viscosity of each standard was determined online. The universal calibration curve was obtained by plotting the product  $\log([\eta] \times M)$  against  $V$ , as shown in Figure 1. The resulting data points were fitted with a third order polynomial, resulting in residuals of less than <5%.

To obtain true molar masses, the broadly distributed collagen peptides were analyzed using the same conditions as those applied for establishing the universal calibration curve.

Table 1 lists the weight average molar masses ( $M_w$ ) and the average intrinsic viscosities (IV) of the collagen peptides, including their standard deviations.

**Table 1.** Weight average molar masses,  $M_w$ , and intrinsic viscosity, IV, of the collagen peptide reference samples (average values and standard deviations from duplicate injections of four individual sample preparations).

Sample	$M_w$ (g/mol)	IV (mL/g)
13R078	1,960 $\pm$ 130	4.57 $\pm$ 0.07
833593	2,820 $\pm$ 270	5.56 $\pm$ 0.27
13R079	3,920 $\pm$ 200	6.60 $\pm$ 0.20
891600	4,480 $\pm$ 380	7.49 $\pm$ 0.12
890435	9,020 $\pm$ 680	9.58 $\pm$ 0.80



**Figure 1.** Universal calibration curve on Agilent PROTEEMA columns, constructed using monosaccharides and narrowly distributed pullulan and dextran standards.




For the further study, the weight average molar masses ( $M_w$ ) were preferred over the number average molar masses ( $M_n$ ), which have a higher experimental error.<sup>6</sup> With the  $M_w$  provided in Table 1, the collagen peptides can be used as broadly distributed reference materials to convert a base calibration to a collagen calibration using the broad standard calibration approach.

### Application of collagen peptides for broad calibration using GPC/SEC-UV-Vis

While the pullulan and monosaccharide/dextran standards cover a suitable molar mass separation range and have the required purity to establish accurate concentrations (as required for the viscometry measurements and thus for constructing a universal calibration curve), they lack chromophores for UV detection. Because a UV-based method was requested, the dextrans and pullulans had to be replaced by UV-active calibrants suitable for establishing the base calibration.


A sample of collagen fragments (CNBr peptides) was identified as a suitable candidate for the base calibration when using GPC/SEC systems with only UV-Vis detection.

A typical chromatogram of the CNBr peptides with indication of the molar masses at the peak maximum (as provided by the supplier) is shown in Figure 2.



**Figure 2.** Chromatogram of the CNBr peptide sample on the Agilent PROTEEMA column. The peaks with given molar mass information and dark red arrows were used to construct the base calibration. Please note that the peak with 55 kDa molar mass is not as nicely separated on all columns as shown in this example. Since the peak elutes close to the exclusion limit of the columns recommended for this application, a loss in column resolution can result in only five peaks being visible.

The molar masses and elution volumes of the five peaks below 38 kDa were used to construct the base calibration. To convert the calibration curve constructed using the CNBr peptides into a calibration curve for collagen peptides, the five collagen peptide reference materials were run using the same experimental conditions as applied to run the CNBr peptide sample (Figure 3).




**Figure 3.** Overlay of the chromatograms of the collagen peptide reference materials on the Agilent PROTEEMA column.

The resulting chromatograms of the five collagen peptides (Figure 3) and their  $M_w$  (as shown in Table 1) were then applied to adjust the parameters of the broad standard approach, thereby establishing a calibration curve for collagen peptides.

Figure 4 shows a comparison of the base calibration corresponding to the CNBr peptides and the collagen peptide calibration curve. A significant shift between both calibration curves is observed, indicating the necessity to apply the broad calibration correction to the base calibration in order to achieve true molar masses. The application of the CNBr peptide calibration only would result in an overestimation of the molar masses for collagen fragments.

Once the calibration curve is established, it can be used for molar mass determination of collagen peptides of unknown molar mass.

Table 2 compares the target values as obtained from universal calibration with the weight average molar masses resulting after adjusting the parameters A and B of the broad calibration approach.



**Figure 4.** Overlay of calibration curves constructed using CNBr peptides (base calibration, red) and the collagen peptide calibration curve obtained by broad standard calibration (blue) on the Agilent PROTEEMA column.

**Table 2.** Comparison of weight average molar masses from universal calibration (target values) with weight average molar masses derived based on converted CNBr peptide calibration.

Sample	$M_w$ (g/mol)	
	Values (recalculated adjusted calibration curve)	Target Values (universal calibration)
13R078	1,929	1,960
833593	2,776	2,820
13R079	4,025	3,920
891600	4,596	4,480
890435	8,806	9,020


### Round robin test

To test the validity of the calibration approach, a round robin test was conducted involving four different laboratories. The task was determination of weight average molar masses for six collagen peptides, which were not applied during calibration.

The eluent, the concentration, and the injection volume were specified as described in the experimental section. A flow rate of  $0.5 \text{ mL/min}^1$  was requested. The laboratories were free to choose from one of two types of columns.

For calibration, the developed procedure with the CNBr peptides and the five broad reference materials was requested. The CNBr peptide sample was applied to establish the base calibration, which was converted into a calibration for collagen peptides using the five samples as shown above.

Afterwards, eight additional collagen peptide samples (A to F) were evaluated using the established calibration curve. The results are summarized in Table 3 and shown in Figure 5.



**Figure 5.** Comparison of weight average molar masses of collagen peptide samples derived in different laboratories using the suggested approach with broad standard calibration.

On average, the value derived for the molar mass of a sample at a laboratory differed by less than 4% from the mean of the sample.

In addition, the data were compared to true weight average molar masses, independently derived by GPC/SEC-viscometry (column UC). The average deviation between the molar mass derived at a laboratory and the true value was less than 9%, which confirms that the above-established calibration approach does not only allow for reliable interlaboratory comparison, but also results in close-to-true molar masses.

**Table 3.** Average molar masses of collagen peptide samples derived in different laboratories using the suggested approach with broad standard calibration.

Sample	M <sub>w</sub> (g/mol) UC	M <sub>w</sub> (g/mol) Lab 1	M <sub>w</sub> (g/mol) Lab 2	M <sub>w</sub> (g/mol) Lab 3	M <sub>w</sub> (g/mol) Lab 4
Sample A	2,950	2,720	2,730	2,650	2,550
Sample B	4,230	3,780	3,640	3,680	3,660
Sample C	2,120	1,950	2,140	1,950	1,870
Sample D	3,960	3,530	3,390	3,300	3,280
Sample E	4,500	4,640	4,450	4,495	4,450
Sample F	2,860	2,830	2,890	2,920	2,740


### Procedure for broad calibration with collagen peptide reference materials using Agilent WinGPC with guided broad calibration

The CNBr peptides and the five collagen peptide references are commercially available for the construction of GPC/SEC calibration curves for collagen peptides using only a UV-Vis detector. The details are described in the Collagen Peptides Monograph.<sup>7</sup>

The applied broad standard calibration approach is often implemented in commercial GPC/SEC software packages. A user-friendly software wizard guides the user through the task in WinGPC.


The procedure includes the following steps:

1. Run the CNBr peptide sample and construct a five-point base calibration based on the measured elution volume and the certified molar masses.
2. Run the broadly distributed collagen peptide reference materials, and add the resulting chromatograms to an overlay.
3. Start the Guided Broad Calibration wizard.
4. Use the overlay and assign the weight average molar masses (as provided in the reference materials) to each of the chromatograms (Figure 6).



**Figure 6.** Step 2 of the guided setup, processing the overlay of chromatograms of the five collagen peptide reference samples with assigned weight average molar masses.

5. Allow WinGPC to automatically adjust the base calibration curve such that best agreement between target values and calculated molar masses is achieved for the chromatograms (Figure 7).



**Figure 7.** Step 3 of the guided setup, showing the comparison of target values and derived values after adjusting the fit parameters.

6. Save the calibration curve for future analysis of unknown collagen peptides.



## References

1. Striegel, A. M. *et al.* Modern Size Exclusion Liquid Chromatography; 2nd ed.; John Wiley and Sons Inc. Hoboken, New Jersey, **2009**.
2. GPC/SEC Theory and Background. *Agilent Technologies eBook*, publication number 5994-5935EN, **2018**.
3. Grubišić, Ž.; Rempp, P.; Benoit, H. A Universal Calibration for Gel Permeation Chromatography. *Materials Science Journal of Polymer Science Part B: Polymer Letters* **1967**, B5, 753, 1067.
4. Kratochvil, P. Classical Light Scattering from Polymer Solutions. Elsevier, New York **1987**.
5. GPC/SEC Detection, section 1.6 How to Choose a Static Light-Scattering Technique. *Agilent Technologies eBook*, publication number 5994-5937EN, **2019**.
6. Held, D. GPC/SEC-Light Scattering: Garbage In, Garbage out. *The Column* **2016**. <https://www.chromatographyonline.com/view/tips-tricks-gpcsec-light-scattering-garbage-garbage-out>
7. Kilz, P.; Gores, F. Accuracy and Precision in GPC/SEC. *The Column* **2008**.
8. Data Integrity in GPC/SEC. *LC/GC Europe* **2020**.
9. GME Collagen Peptides Monograph, Final Version 1 – October **2020** [https://www.gelatine.org/fileadmin/user\\_upload/gme\\_content/GME\\_Statements/GME\\_CP\\_Monograph-version\\_1\\_-\\_public.pdf](https://www.gelatine.org/fileadmin/user_upload/gme_content/GME_Statements/GME_CP_Monograph-version_1_-_public.pdf)

## Contact information

Learn more:

[www.agilent.com](http://www.agilent.com)

Buy online:

[www.agilent.com/chem/store](http://www.agilent.com/chem/store)

Get answers to your technical questions and access resources in the Agilent Community:

[community.agilent.com](http://community.agilent.com)

U.S. and Canada

**1-800-227-9770**

[agilent\\_inquiries@agilent.com](mailto:agilent_inquiries@agilent.com)

Europe

[info\\_agilent@agilent.com](mailto:info_agilent@agilent.com)

Asia Pacific

[inquiry\\_lsca@agilent.com](mailto:inquiry_lsca@agilent.com)

[www.agilent.com](http://www.agilent.com)

DE07652947

This information is subject to change without notice.

© Agilent Technologies, Inc. 2024  
Printed in the USA, February 19, 2024  
5994-6747EN

