

CLP and EPA methods for Pesticides in Water using Agilent J&W DB-CLP1 / DB-CLP2 Columns

Pat Sasso, Doris Smith, and Ken Lynam Agilent Technologies Inc.

PittConn 2013
Poster # 1480-10P



Abstract

The Agilent J&W DB-CLP1 and DB-CLP2 column pair has the selectivity to separate many different halogenated compounds by GC/ECD. Chlorinated pesticides, herbicides, organohalides, and PCBs in drinking water are well resolved by these columns with multiple Contract Laboratory Program (CLP) and United States Environmental Protection Agency (US EPA) methods. Some manufacturers offer different column sets for the different methods. However, the Agilent CLP single column set avoids the downtime of switching between different columns and can analyze multiple methods. This productivity increase can lead to more billable instrument-hours. Compared to other vendors' columns, the J&W DB-CLP1/DB-CLP2 provides better combinations of faster analysis times and fewer coelutions, for analysis of a wider range of pesticides and pollutants covered by US EPA methods all in a single operation.

Introduction

GC column manufacturers design columns specifically for pesticide residue analysis, with differing capability and efficiency. This type of measurement is typically done on 2 columns for simultaneous primary and confirmation analysis, using a retention gap and a γ -splitter to connect the columns. In this poster, we analyzed pesticides and other potential water contaminants according to 9 United States Environmental Protection Agency (US EPA) and Contract Laboratory Program (CLP) methods. We used Agilent J&W DB-CLP1 and DB-CLP2 columns, a pair of arylene backbone phases that handle high oven temperatures for rapid analysis of CLP chlorinated pesticides in drinking and ground water. We demonstrated the columns' suitability for phenoxy acids, haloacetic acids, polychlorinated biphenyls, and other pollutants, according to EPA methods for the analysis of water samples.

The Nine CLP and EPA methods investigated were:

- CLP Pesticides – 19 targeted organohalides
- EPA 504 – 1 (EDB), (DBCP), and (123TCP)
- EPA 505 - organohalides
- EPA 508.1 – chlorinated series
- EPA 551 – solvents, (THMs), and (DBPs)
- EPA 552.3 – haloacetic acids and dalapon
- EPA 8081B – (extended analyte list)
- EPA 8082A – (PCBs) and arachlors
- EPA 8151A – chlorophenoxyacid herbicides



Results of Dual Column Confirmation

CLP Pesticides

DB-CLP1 30m 0.32mm ID 0.25 μ m (cat.# 123-8232)
DB-CLP2 30m 0.32mm ID 0.5 μ m (cat.# 123-8336)
5m x 0.32mm ID deact. guard column
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split

Instrument: Agilent 7890 GC with dual μ ECD
Sampler: Agilent 7693
Sample: 50 ng/mL CLP Pesticides
Inj. Vol.: 1 μ L splitless
Liner: Ultra Inert liner, single taper splitless (cat.#5190-2292)
Inj. Temp.: 250°C
Oven Temp.: 150°C (hold 0.2 min),
45°C/min to 250°C, 18°C/min to 300°C,
30°C/min to 330°C, hold 2.5 min
Carrier Gas: Helium, constant flow 3.5 mL/min
Detector: μ ECD @ 340 °C



1. Tetrachloro-m-xylene*	9. γ -Chlordane	16. Endosulfan II
2. α -BHC	10. α -Chlordane	17. 4,4'-DDT
3. γ -BHC	11. Endosulfan I	18. Endrin aldehyde
4. β -BHC	12. 4,4'-DDE	19. Endosulfan sulfate
5. Heptachlor	13. Dieldrin	20. Methoxychlor
6. δ -BHC	14. Endrin	21. Endrin ketone
7. Aldrin	15. 4,4'-DDD	22. Decachlorobiphenyl*
8. Heptachlor epoxide		*surrogate standard



EPA 508.1

DB-CLP1 30m 0.32mm ID 0.25 μ m (cat.# 123-8232)
DB-CLP2 30m 0.32mm ID 0.5 μ m (cat.# 123-8336)
5m x 0.32mm ID deact. guard column
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split

Instrument: Agilent 7890 GC with dual μ ECD
Sampler: Agilent 7673B
Sample: 100 ng/mL EPA 508.1 analytes
100 ng/mL Pesticide Surrogate Mix
Inj. Vol.: 2 μ L splitless
Liner: Ultra Inert liner, single taper splitless (cat.#5190-2292)
Inj. Temp.: 250°C
Oven Temp.: 80°C (hold 0.5 min) to 175°C at 26°C/min,
0.5°C/min to 235°C, 15°C/min to 300°C,
hold 5 min
Carrier Gas: Helium at 35 cm/sec (constant flow)
Detector: μ ECD @ 340 °C





Results of Dual Column Confirmation

EPA 8081B

DB-CLP1 30m 0.32mm ID 0.25 μ m (cat.# 123-8232)
DB-CLP2 30m 0.32mm ID 0.5 μ m (cat.# 123-8336)
5m x 0.32mm ID deact. guard column
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split

Instrument: Agilent 7890 GC with dual μ ECD
Sampler: Agilent 7693
Sample: 50 ng/mL EPA 8081B analytes
Inj. Vol.: 2 μ L splitless
Liner: Ultra Inert liner, single taper splitless (cat.#5190-2292)
Inj. Temp.: 250°C
Oven Temp.: 80°C (hold 0.5 min) to 150°C at 20°C/min,
5°C/min to 235°C, 15°C/min to 300°C,
hold 5 min
Carrier Gas: Helium at 43.5 cm/sec (constant flow)
Detector: μ ECD @ 325 °C

1. 1,2-Dibromo-3-chloropropane	11. α -BHC	21. DCPA	31. Chlorobenzilate (250ng/mL)	41. Captafol
2. Heptachlorocyclopentadiene	12. Pentachloronitrobenzene	22. Isodrin	32. Perthane (250ng/mL)	42. Methoxychlor
3. 1-Bromo-2-nitrobenzene	13. γ -BHC	23. Kelthane	33. Chloropropylate (250ng/mL)	43. Endrin ketone
4. Etradiazole	14. β -BHC	24. Heptachlor epoxide	34. Endrin	44. Mirex
5. Chlorobenz	15. Heptachlor	25. γ -Chlordane	35. Nitrofen	45. cis-Permethrin
6. Trifluralin	16. Dieldrin	26. trans-Nonachlor	36. 4,4'-DDD	46. trans-Permethrin
7. TCMX*	17. Alachlor	27. α -Chlordane	37. Endosulfan II	47. Decachlorobiphenyl*
8. Propachlor	18. δ -BHC	28. Endosulfan I	38. 4,4'-DDT	*surrogate standard
9. Diallate isomers (250ng/mL)	19. Chlorobalonil	29. 4,4'-DDE	39. Endrin aldehyde	
10. Hexachlorobenzene	20. Aldrin	30. Dieldrin	40. Endosulfan sulfate	





EPA 8151A

DB-CLP1 30m 0.32mm ID 0.25 μ m (cat.# 123-8232)
DB-CLP2 30m 0.32mm ID 0.5 μ m (cat.# 123-8336)
5m x 0.32mm ID deact. guard column
Inert Tee CFT device (cat.# G3184-60065) 1:1 Split

Instrument: Agilent 7890 GC with dual μ ECD
Sampler: Agilent 7693
Sample: 50 ng/mL Chlorophenoxyacid Herbicides (methyl esters)
50 ng/mL DCAA methyl ester
50 ng/mL 4,4'-Dibromooctafluorobiphenyl
Inj. Vol.: 2 μ L splitless
Liner: Ultra Inert liner, single taper splitless (cat.#5190-2292)
Inj. Temp.: 250°C
Oven Temp.: 85°C (hold 0.25 min) to 135°C at 25°C/min,
11°C/min to 200°C, 20°C/min to 300°C,
hold 1.5 min
Carrier Gas: Helium, constant flow 3 mL/min
Detector: μ ECD @ 325 °C

1. Dalapon methyl ester	11. Pentachloroanisole
2. 2,5-dichlorobenzoic acid methyl ester	12. 2,4,5-TP methyl ester
3. 2-Nitroanisole	13. 2,4,5-T methyl ester
4. DCAA methyl ester (SS)	14. Chloramben methyl ester
5. Dicamba methyl ester	15. Dinoseb methyl ester
6. MCPP methyl ester	16. 2,4-DB methyl ester
7. MCPA methyl ester	17. Bentazon methyl ester
8. 4,4'-DBOB (IS)	18. DCPA methyl ester
9. Dichloroprop methyl ester	19. Picloram methyl ester
10. 2,4-D methyl ester	20. Acifluorfen methyl ester



Materials and Methods

The same column pair was used throughout these analyses. Other conditions varied and are shown with their respective chromatograms. Traditionally, replacement of the retention gap or analytical column involved replacing a quartz γ -splitter and re-establishing all column connections. To avoid the downtime and extra maintenance associated with γ -splitters, we also used an Agilent Capillary Flow Technology (CFT) reusable splitter. The CFT splitter uses column connections that are individually connected to the splitter, allowing inlet and column maintenance independent of the other analytical column connection.

Column1: Agilent J&W DB-CLP1 30m x 0.32mm, 0.25 μ m (p/n 123-8232)
Column 2: Agilent J&W DB-CLP2 30m x 0.32mm, 0.25 μ m (p/n 123-8336)
Gap: 5m x 0.32 mm, deactivated fused silica tubing
CFT Device: Inert-tee (p/n G3184-60065), split 1:1
Inlet liner: Agilent UI, single taper splitless (p/n 5190-2292) (Method 505 liner was dual taper (p/n 5190-3983))
Instrument: Agilent 7890 GC w/ dual μ ECD
Auto-sampler: Agilent 7693 Automated Liquid Sampler (Method 505 and 508.1 : 7873B Liquid Sampler)
Magnifier: 20x Coddington magnifier (p/n 430-1020)

References

Please see the application notes listed below additional details at:
<http://www.chem.agilent.com/chem>
Agilent Application Notes: 5991-0207EN, 5991-0541EN, 5991-0615EN

Conclusions

Agilent J&W DB-CLP1 and DB-CLP2 columns analyze 9 EPA methods along with CLP pesticides, without changing column sets. This capability increases sample throughput with faster analysis times for greater lab productivity. In addition, this column pair delivers fewer coelutions, allowing more accurate determination. This column pair is the pair of choice for your dual ECD methods.