



# Polytron PT 1300 D





2.136.0100

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Scope of delivery 2.136.0100

Qt.	Order no.	Description
1 PCS	1.136.0100	<b>Polytron PT 1300 D</b>  Polytron PT 1300 D - Metrohm version  Homogenizer that can be controlled directly by <b>tiamo™</b> or Touch Control.
		
1 PCS	1.136.0110	<b>Power supply unit for Polytron PT 1300 D</b>
		

## Optional accessories

Order no.	Description	
6.02008.010	<b>Polytron holder</b>  For fastening the Polytron in place, suitable for 601405040.	 The image shows a teal-colored C-shaped holder with two black cylindrical caps. One cap is shown inserted into the holder, while the other is shown separately.
6.02008.030	<b>Holder for homogenizer for Eco Titrations instruments</b>  Holder for homogenizer for fastening the Polytron in place, suitable for 601405060 titration vessel lid	 The image shows a teal-colored holder with a central rectangular cutout and two black cylindrical caps. One cap is shown inserted into the holder, and another is shown separately.
6.1462.250	<b>Swing arm for Polytron and 786 Swing Head, left swinging</b>  Swing arm for Polytron and 786 Swing Head, left swinging, with one buret tip for the addition of solvent or titrant and three spray nozzles for cleaning. The swing arm allows the usage of the Polytron for sample preparation on multi row racks.	 The image shows a black metal swing arm assembly. It features a central vertical rod with a horizontal arm extending to the left. A buret tip is attached to the end of the horizontal arm. Three spray nozzles are also attached to the assembly.
6.2047.110	<b>Holder for Polytron</b>  Holder to fix the Polytron at the stand rod.	 The image shows a teal-colored holder with a central rectangular cutout and two black cylindrical caps. One cap is shown inserted into the holder, and another is shown separately.

6.9012.000

**Dispersing aggregate for Polytron, 125 mm**

Standard dispersing aggregate for Polytron 1300 D



---

6.9012.010

**Dispersing aggregate for Polytron, 157 mm**

Special dispersing aggregate (Metrohm version) with protruding knives

