




860 KF Thermoprep

2.860.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.860.0010

Qt.	Order no.	Description	
1 PCS	1.860.0010	860 KF Thermoprep The 860 KF Thermoprep is designed for thermal sample preparation in Karl Fischer titration. Numerous substances cannot be analysed by direct Karl Fischer titration as they are not soluble, react with the Karl Fischer reagent or release their water only very slowly or not until heated to high temperatures. Using tightly sealed sample vials, the sample is inserted in the oven. The samples can be subsequently analysed by any volumetric or coulometric KF titrator.	
1 PCS	6.1446.170	Heating tube stopper for coulometric KF vessels in connection with heatable tubing connection	
1 PCS	6.1448.057	Aluminum septum caps / 100 pieces Including Al sealing.	

2 PCS

6.1602.145

Drying flask cap

Screw cap with GL 45 thread for the 6.1608.050 drying flask.



2 PCS

6.1608.050

Drying bottle / 100 mL / GL 45



1 PCS

6.1805.010

FEP tubing / M6 / 13 cm

With light and kink protection.



2 PCS

6.1805.080 FEP tubing / M6 / 25 cm

With light and kink protection.



1 PCS

6.1805.100 FEP tubing / M6 / 40 cm

With light and kink protection



1 PCS

6.1808.040 Thread adapter M6 outer / M8 inner

Outer thread M6, inner thread M8.



1 PCS

6.1808.050 Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



2 PCS

6.1821.040 Filter tube

Filter tube for 6.1608.050 Drying bottle. For Rancimats and Karl Fischer ovens.



1 PCS

6.1830.030 Heating tubing

Heating tubing with M6 thread for 860 KF Thermoprep, 885 Compact Oven Sample Processor and 874 USB Oven Sample Processor.



1 PCS

6.2048.030 Guidance rod

Guidance rod for 860 KF Thermoprep.



1 PCS

6.2049.040

Needle holder with Luer lock

Needle holder (58 mm) with Luer lock for the Karl Fischer ovens 860, 874 and 885.



1 PCS

6.2060.010

Cable holder

Used with 860 KF Thermoprep.



1 PCS

6.2063.010

Sample holder for 6.2419.007 Sample vials (6 mL)



1 PCS

6.2419.007

Sample glass, 6 mL, 100 pieces

Sample glass, 6 mL, suitable for 6.1448.057 septum, 100 pieces



1 PCS

6.2621.110 Crimping tongs

Crimping tongs for sealing 6 mL sample vials (6.2419.000, 6.2419.007).



1 PCS

6.2621.120 Hexagon key 1.5 mm



1 PCS

6.2621.130 Hexagon key 2 mm

2 mm.



1 PCS

6.2627.000 Oven insert

Used with 860 KF Thermoprep.



1 PCS

6.2724.010 Dust filter

Dust filter for Rancimats and Karl Fischer ovens.



1 PCS

6.2739.000 Wrench

For tightening connectors



1 PCS

6.2811.000 Molecular sieve

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



1 PCS

6.2816.070 Septum piercing needle

Used with Karl Fischer ovens.



1 PCS

6.2816.080

Outlet needle

Used with KF Thermoprep and Oven Sample Processor.

