


TRACE 1300/1310: Deans Switch Microfluidics Kit Installation Note

The Deans Switch Microfluidics Kit **P/N 19005580** includes all the components required for installing the Deans Switch Microfluidics device on your TRACE 1310/1310 GC.

The DVD **P/N 31709595** provided with the kit contains:

- **Thermo Scientific Deans Switch Calculator** program P/N 27410335. It must be installed on a computer for the control of the device.



- **Deans Switch Microfluidics Installation Guide** P/N 31709740 in PDF format. It provides the instruction for installing the device, and for installing and using the Thermo Scientific Deans Switch Calculator program.

