

TraceFinder 4.1 System Requirements

These are the hardware and software configurations required for TraceFinder 4.1 operation.

System	Requirements
Hardware	<ul style="list-style-type: none">• 2.33 GHz dual-core processor enabled• 8 GB RAM with system managed memory enabled• DVD drive• Resolution display 1280 ×1024 (SXGA)• 5 GB available on drive C• NTFS format
Software	<ul style="list-style-type: none">• Microsoft™ Windows™ 7 SP1 Professional (64-bit) English only• Microsoft .NET Framework 4.5.1• Adobe™ Reader™ 10.1• Thermo Scientific software:<ul style="list-style-type: none">– Xcalibur™ 4.0– Foundation™ platform 3.1 SP3– mzVault™ 1.0 SP1– FreeStyle™ 1.1 SP1– LC Devices 3.1 or later– Exactive™ Series 2.7 or later– TSQ Endura™/TSQ Quantiva™ 2.0 or later– TSQ™ 2.6 or later– Fusion™ 1.2 SP1 or later– Aria™ MX 2.2 or later– GC Devices 3.2 or later– ISQ™ 3.2 SP1 or later– SII 1.2 for Xcalibur or later– TSQ 8000 3.2 or later– LTQ™ Series 2.7 SP4 or later