

Selecting between RoHS Cup LN and Eco-Cup LF

Products affected RoHS Cup LN series
 Eco-Cup LF series



< Overview >

The table below shows how to select a sample cup between RoHS Cup LN and Eco-Cup LF when Frontier Labs Auto-Shot Sampler and Pyrolyzer are used.

Sample cup	Analytical target		Description
	Phthalate	Polymer/additive	
RoHS Cup LN	Yes	No	<ul style="list-style-type: none"> No deactivation Disposable use Recommended for phthalates analysis
Eco-Cup LF	Yes	Yes	<ul style="list-style-type: none"> Deactivated Can be reused Disposable use recommended for analysis of brominated flame retardants

- Both made of stainless steel, height 8 mm / o.d. 4 mm / volume 80 μ L
- Volume discount available. Refer to the table below.

< Product number and quantity >

Sample cup	Product name	Product number	Quantity (ea.)	Front view of bottle
RoHS Cup LN series	RoHS Cup LN	PY1-EC80N	100 (100 x 1 bottle)	 <p>Orange label</p>
	RoHS Cup LN-V5	PY1-EC80N-V5	500 (100 x 5 bottles)	
	RoHS Cup LN-V10	PY1-EC80N-V10	1,000 (100 x 10 bottles)	
Eco-Cup LF series	Eco-Cup LF	PY1-EC80F	100 (100 x 1 bottle)	 <p>White label</p>
	Eco-Cup LF-V4	PY1-EC80F-V4	400 (100 x 4 bottles)	