


Powder Residue Found in a Plastic Bag Analyzed by GC-TOFMS for Illegal Drugs

Life Science and Chemical Analysis Centre
LECO Corporation, Saint Joseph, Michigan USA



Analysis of powder residue showing positive identification of 4 illegal drugs of abuse.

Column:
30 m x 0.25 mm x 0.25 μ m Rxi-5ms

TOFMS:
45 to 450 m/z @ 5 spectra/s

For questions on this analysis e-mail us at:
life_science@leco.com

For a complete listing of snapshots and application notes visit us on the web at www.leco.com

