


## Analysis of pesticides (USEPA 608)

**BPX5**

Column part number	054142		
Phase	BPX5	Detector	ECD, 320°C
Column	30 m x 0.22 mm x 0.25 µm	Injection mode	Split (300:1)
Initial temperature	160°C	Carrier gas	He, 25 cm/sec
Rate	7°C/min	Injection volume	0.5 µL
Final temperature	290°C, 5 min		


 For more information visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)