


Chlorophenoxyacid Herbicides
US EPA Method 615
Rtx®-5



1. dalapon
2. DCAA (S)
3. dicamba
4. MCPP
5. MCPA
6. DBOB (IS)
7. dichlorprop
8. 2,4-D
9. 2,4,5-TP
10. 2,4,5-T
11. dinoseb
12. 2,4-DB

30m, 0.53mm ID, 0.50 μ m Rtx®-5 (cat.# 10240)
0.5 μ L direct injection of chlorophenoxy herbicides,
on-column concentration 10–10,000 μ g/mL

Oven temp.: 60°C to 150°C @ 8°C/min. (hold 5 min.),
to 210°C @ 4°C/min.
Inj./det. temp.: 250°C / 275°C
Carrier gas: helium
Linear velocity: 35cm/sec. set @ 60°C
Det.: ECD w/anode purge