


# Aroclor® 1254 on Rtx®-PCB



GC\_EV00753

Column: Rtx®-PCB 30m, 0.25mm ID, 0.25 $\mu$ m (cat.# 13223)  
Sample: 200ng/mL Aroclor® 1254 (cat.# 32011)  
Inj.: 1.0 $\mu$ L splitless (hold 0.75 min.), 3.5mm ID single goose-neck inlet liner (cat.# 20962)  
Inj. temp.: 250°C  
Carrier gas: hydrogen, constant pressure  
Linear velocity: 71cm/sec. @ 110°C  
Oven temp.: 100°C (hold 1.0 min.) to 300°C @ 10°C/min. (hold 4 min.)  
Det.: ECD @ 310°C