


FAMES Saw Palmetto Stabilwax®

(split injection)



50m, 0.25mm ID, 0.25 μ m Stabilwax® (cat.# 10626-105).

1 μ L split injection.

Oven temp.: 110°C (hold 1 min.) to 240°C @ 8°C/min. (hold 25 min.)

Inj./FID temp.: 230°C/250°C

Carrier gas: hydrogen @ 2.5mL/min.

Split flow: 37.5mL/min.

Septum purge: 3mL/min.

Chromatogram provided by The Institute for Nutraceutical Advancement (INA)

Restek Corporation 110 Benner Circle Bellefonte, PA 16823
814-353-1300 • 800-356-1688 • Fax: 814-353-1309 • www.restekcorp.com