


Analysis of triglyceride distribution in cocoa butter

BPX5

Column part number	054133		
Phase	BPX5	Final temperature	360°C, 5 min
Column	12 m x 0.53 x 0.25 µm	Carrier gas	He, 6 psi
Initial temperature	100°C, 0.5 min	Injection mode	On-column
Rate 1	50°C/min to 280°C	Injection volume	1 µL
Rate 2	10°C/min to 360°C		


 For more information visit www.trajanscimed.com or contact techsupport@trajanscimed.com